Vorderseite

Suprasorb® G

Amorphes Gel

Gel amorpho

Amorphous Gel

Gel amorfó

Amorfí Gel

Amorfinen geeli

Amorfín gél

Amorf gél

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en | Instructions for Use

Product Description:

Suprasorb® G Amorphous Gel is a transparent hydrogel consisting of CMC polymers, propylene glycol and water. Suprasorb® G Amorphous Gel is sterile, individually packed, and is sterilisable. The gel amorpho Suprasorb® G is for single use only. The product is sterile unless the package is opened or damaged.

Indications:

Suprasorb® G Amorphous Gel hydrates the wound facilitating the debridement of necrotic tissue. Slough is absorbed by the gel and removed from the wound with the dressing change. Suprasorb® G Amorphous Gel is indicated for the desloughing of necrotic and fibrous tissue, particularly in deep wounds. Before using the product, the package is opened or damaged.

Application:

- Clean wound as usual
- Open the package where indicated
- Remove sterile syringe and pull protective cap off the tip
- Apply gel directly into the wound
- Use a suitable dressing to cover Suprasorb® G Amorphous Gel (depending on the condition of the wound, for example with Suprasorb® H Hydrocolloid Wound Dressing or Suprasorb® F Film Wound Dressing)

Dressing Change:

Suprasorb® G Amorphous Gel may remain in the wound for up to 3 days depending on the amount of exudate. The dressing should be changed before the wound becomes saturated. Before applying a new dressing, the wound should be cleaned, e.g. with saline or Ringer solution.

Precautions:

• In the case of infected wounds it is the responsibility of the physician in charge to decide whether the product should be used.

• The product may not be used in the presence of known allergies to any of its components.

• Before removing the cap, carefully press down the piston to avoid contact with the gel. Move the piston lightly to ease the subsequent controlled administration of the gel.

• The gel must be used as directed and should not be applied.

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