Auto Suture*
ENDO LUNG*, ENDO MINI-SHEARS* SHORT and ENDO SCIZ* SHORT
DISPOSABLE INSTRUMENTS

IMPORTANT!

This booklet is designed to assist in using the Auto Suture* ENDO LUNG*, ENDO MINI-SHEARS* Short and ENDO SCIZ* Short Instruments. It is not a reference to surgical endoscopic techniques. For information on endoscopic procedures, techniques, complications and hazards, please see the books: *Surgical Laparoscopy* (Zucker KA, ed. Quality Medical Publishing, Inc., St. Louis, MO 1991) and *Endoscopic Surgery* (White RA, Klein SR, Mosby-Year Book, Inc., St. Louis, MO, 1991).

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

CAUTION: Federal (USA) law restricts these devices to sale, distribution, and use by, or on the order of, a physician.

INDICATIONS & EFFECTS

The Auto Suture* ENDO MINI-SHEARS* Short and ENDO SCIZ* Short instruments are packaged with a cautery adapter. Remove the adapter from the instrument if not using cautery.The cautery adapter is compatible with the following manufacturers' male and female cables.

ELECTROCAUTERY COMPATIBILITIES

Aspen Beldon Birtcher Cabot Conmed Erbe Jarit Siemens Snowden Storz Stryker Valleylab Weck Wisap Wolf Zink

The Auto Suture* ENDO LUNG* instrument has application in a variety of general, thoracic, gynecologic, urologic and endoscopic procedures for temporary grasping/clamping of tissue and small tubular structures. The device has a 12 mm diameter insulated shaft, 19,3 cm in length, and is designed for introduction and use through all appropriately sized Auto Suture* trocar sleeves or larger sized trocar sleeves with the use of a convert-er. The scalloped dial located on the handle rotates the shaft 180°, in either direction. The jaws of the ENDO LUNG* instrument are closed by compression of the ring handles. To open the jaws of the ENDO LUNG*, pull back on the release trigger.

The Auto Suture* ENDO MINI-SHEARS* Short instrument with unipolar cautery has application in a variety of general, thoracic, gynecologic, urologic and endoscopic procedures for mobilization and transection of tissue. The device has a 5 mm diameter insulated shaft, 19,3 cm in length, and is designed for introduction and use through all appropriately sized Auto Suture* trocar sleeves or larger sized trocar sleeves with the use of a converter. The scalloped dial located on the handle rotates the shaft 360°. The scissor blades are closed and opened to transect the tissue by compression and release of the ring handles. The scissor blades have a curved cutting edge 10 mm in length, with maximum blade opening of 8 mm. The female cautery connection located on top of the handle may be utilized for unipolar cautery when attached to standard cautery cables and their generators (a male adapter is provided).

IMPORTANT: BIPOLAR CABLES HAVE A DIFFERENT TYPE OF PLUG THAT CANNOT BE CONNECTED TO A UNIPOLAR CONNECTION.

The Auto Suture* ENDO SCIZ* Short instrument with unipolar cautery has application in a variety of general, thoracic, gynecologic, urologic and endoscopic procedures for mobilizing and transecting tissue. The device has a 5 mm diameter insulated shaft, 19,3 cm in length, and is designed for introduction and use through all appropriately sized Auto Suture* trocar sleeves or larger sized trocar sleeves with the use of a converter. The scalloped dial located on the handle rotates the shaft 360°. The hooked scissor blades are closed and opened to transect tissue by compression and release of the ring handles. The hooked scissor blades have a cutting edge, 10 mm in length, with a maximum blade opening of 7,5 mm. The female cautery connection located on top of the handle may be utilized for unipolar cautery cables and their generators (a male adapter is provided).

IMPORTANT: BIPOLAR CABLES HAVE A DIFFERENT TYPE OF PLUG THAT CANNOT BE CONNECTED TO A UNIPOLAR CONNECTION.

CAUTIONS

- Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to complications, hazards and techniques, prior to performing endoscopic procedures.
- A thorough understanding of the principles and techniques involved in laser and electrosrurgical procedures is necessary to avoid shock and burn hazard to both patient and operator(s). Verify compatibility of instrumentation, and ensure that electrical isolation or grounding is not compromised.
- These devices are provided STERILE and are intended for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT RESTERILIZE.

CONTRAINDICATION

The ENDO MINI-SHEARS* and ENDO SCIZ* Short instruments with unipolar cautery are NOT intended for contraceptive coagulation of fallopian tubes, but may be used to achieve hemostasis following transection of the fallopian tube.

DO NOT EXPOSE TO TEMPERATURES ABOVE 130° F (54° C).

*Trademark of United States Surgical Corporation.
Manufactured by: United States Surgical Corporation, Norwalk, Connecticut 06856
© 1995 United States Surgical Corporation.
All Rights Reserved. Made in the U.S.A.

Auto Suture*
ENDO LUNG*, ENDO MINI-SHEARS* SHORT et ENDO SCIZ* SHORT
INSTRUMENTS

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Auto Suture European Services Center, S.A.,78990 Elancourt, France.
*Marchio di United States Surgical Corporation
Prodotto da: United States Surgical Corporation, Norwalk, Connecticut 06856 USA.
*Marca de United States Surgical Corporation
Fabricado por: United States Surgical Corporation, Norwalk, Connecticut 06856 USA.