
ETHICON ECHELON FLEX

Powered Vascular Stapler ECHELON Flex 35 with
Advanced Placement Tip, shaft length 320 mm

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Please read all information carefully.











Failure to properly follow the instructions may lead to serious surgical consequences, such as leakage or disruption.









Important: This package insert is designed to provide instructions for use of the ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip, shaft length 320 mm (further instrument). It is not a reference to surgical techniques.

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Instruction for use

Designation

The ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip, shaft length 320 mm and its reloads are intended for transection and resection of tissue and vasculature.

Indications

The Instrument and its reloads are intended for transection and resection of tissue and vasculature. The instruments have application in multiple open or minimally invasive general, gynecologic, urologic, thoracic, and pediatric surgical procedures.

Area of application

It is used only in medical and preventive health care institutions.

Users

The instrument is intended for use only by highly qualified physicians with sufficient knowledge in the field of surgical techniques.

Contraindications

Do not use the instrument on the aorta.

Do not use the instrument on ischemic or necrotic tissue.

Do not use any linear cutter on major vessels without making provision for proximal and distal control.

Tissue thickness should be carefully evaluated before firing any stapler. Refer to the Stapler Reload Product Code Chart below for the tissue compression requirements (closed staple height). If tissue cannot comfortably compress to the closed staple height, or easily compresses to less than the closed staple height, the tissue is contraindicated as it may be too thick or too thin for the staple identified in the Reload Product Code Chart.

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This instrument is not intended for use when surgical stapling is contraindicated.

Do not use the instrument where use of buttress material is indicated.

Do not use the instrument if upon closing, the tissue extends beyond the cut line marking, which represents the minimum cut distance, such as in creating an anastomosis.

Do not use the instrument on the parenchyma of the liver, pancreas, kidney, or spleen, or other tissues which cannot comfortably compress to the closed staple height.

Undesirable Side Effects/Residual Risks

Undesirable side effects and risks associated with powered surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

Device Description

The ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip and ENDOPATH ECHELON Vascular White Reload for Advanced Placement Tip (35mm, 4 Row) reloads are sterile, single patient use medical devices that simultaneously cut and staple tissue. There are four staggered rows of staples, two on either side of the cut line. The ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip and reloads have a staple line that is

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approximately 35 mm long and a cut line that is approximately 30 mm long. The shaft can rotate freely in both directions and an articulation mechanism enables the distal portion of the shaft to pivot to facilitate lateral access to the operative site.

The instrument is packaged with a primary lithium battery pack that must be installed prior to use. There are specific requirements for disposing of the battery pack. Refer to the Battery Pack Disposal section.

The instrument is packaged without a reload and must be loaded prior to use. A staple retaining cap on the reload protects the staple leg points during shipping, transportation, and installation of the reload. The instrument's lockout feature is designed to prevent a used or improperly installed reload from being refired or an instrument from being fired without a reload.

Caution: Do not load the instrument more than 12 times for a maximum of 12 firings per instrument.

Instrument Product Code Chart

Instrument Code	Description	Shaft Length
PVE35A	ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip	320 mm

Reload Product Code Chart

Reload	Open Staple Height	Closed Staple Height	Staple Line Length	Reload Color	Number of Staples	Staple Rows
VASECR 35	2,5 mm	1,0 mm	35 mm	White	36	4

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WARNING: The ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip may only be used with the ENDOPATH ECHELON Vascular White Reload for Advanced Placement Tip, (35mm, 4 Row).

Clinical Section

Examples of procedures in which linear cutters and staplers may be used include:

General	Gynecologic	Urologic	Thoracic
Appendectomy	Oophorectomy	Nephrectomy	Lung Resection
Splenectomy		Radical Nephrectomy	

Illustrations and Nomenclature (Illustration 1)

1. Closing Trigger
2. Red Firing Trigger Lock
3. Firing Trigger
4. Anvil Release Buttons (on either side)
5. Battery Pack
6. Battery Pack Release Tab
7. Manual Override Access Panel
8. Knife Reverse Switch
9. Rotating Knob
10. Articulation Lever
11. Reload
12. Reload Alignment Tab
13. Reload Alignment Slot
14. Staple Retaining Cap
15. Anvil Jaw
16. Channel Jaw
17. Staple Line Length

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Instructions for use

Prior to Using Instrument

- 1 Verify that the reload size matches the instrument size to be used (e.g. use an ENDOPATH ECHELON Vascular White Reload for Advanced Placement Tip, (35mm, 4 Row) with an ECHELON FLEX Powered Vascular Stapler instrument.
- 2 Verify compatibility of all instruments and accessories prior to using the instrument (refer to Warnings and Precautions).

Preparing the Instrument for Use

- 3 Using sterile technique, remove the instrument, battery pack, and reload from their respective packages. To avoid damage, do not flip the instrument, battery pack, or reload into the sterile field.
- 4 Install the battery pack. The battery pack must be installed prior to use. Insert the battery pack by lining up the release tabs on the battery pack with the slots on the rear of the instrument. It can be inserted in either orientation; there is no up or down (Illustration 3). Ensure battery pack is fully inserted into the device. An audible click will be heard when the battery pack is fully inserted.
Note: The instrument must be used within 12 hours of inserting the battery pack. Refer to the Battery Pack Disposal section for battery pack disposal instructions.
- 5 Prior to inserting the reload, ensure the instrument is in the open position (Illustration 1).
- 6 Examine the reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

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Caution: Tissue thickness should be carefully evaluated prior to using the instrument. Refer to the Reload Product Code Chart for proper reload selection.

WARNING: Loading the incorrect reload for the instrument size or model may result in transecting the tissue without stapling (e.g. loading a 35 mm reload into a 45 mm instrument).

7 Insert the reload by sliding it against the bottom of the channel jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. (Illustration 4) The instrument is now loaded and ready for use. Caution: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any blue-colored drivers are visible. (If blue-colored drivers are visible, the reload may not contain staples.)

Using the Instrument

8 Close the jaws of the instrument by squeezing the closing trigger until it locks in place (Illustration 5). An audible click indicates that the closing trigger and the jaws are locked. When the jaws of the instrument are closed, the red firing trigger lock and firing trigger are exposed.

Note: Applying pressure to the articulation lever during closure may cause the device to articulate.

Note: Do not pull the red firing trigger lock or firing trigger at this time. Doing so could cause the instrument to be partially or completely fired and the instrument will need to be reloaded before using on tissue.

9 Visually inspect the stapler to ensure proper reload seating. Introduce the instrument into the body cavity through a trocar of the appropriate size or through an incision (Illustration 6). When using a trocar, the instrument jaws must be visible past the trocar sleeve before opening the jaws.

Caution: For insertion and removal of articulating instruments, the jaws of the instrument must be straight, parallel to the shaft of the instrument. Failure to have the instrument jaws in the straight position will result in difficult insertion or withdrawal of the instrument and may result in damage to the instrument.

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Caution: When placing the instrument through the trocar or incision, avoid inadvertently pulling the red firing trigger lock and the firing trigger. Doing so could cause the instrument to be partially or completely fired and the instrument will need to be reloaded before using on tissue. If the instrument is partially fired, remove the instrument and slide the knife reverse switch forward to return the knife to home position (Illustration 7). To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side

of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger (Illustration 8b). Remove the reload and insert a new one (see Reloading the Instrument). If the instrument is fired completely, the knife will return to home position automatically. Remove the instrument and insert a new reload (see Reloading the Instrument).

Note: At any time, if the knife reverse switch does not return the knife to home position (Illustration 11b or 11c) and the jaws will not open:

a. First, ensure the battery pack is securely installed and the instrument has power (Illustration 9); then, try the knife reverse switch again (Illustration 7).

b. If the knife still does not return (Illustration 11b or 11c), use the manual override. After the manual override system is used, the instrument is disabled and cannot be used for any subsequent firings. To use the manual override, remove the access panel labeled “Manual Override” on the top of the instrument handle. The manual override lever will be exposed. Move the lever forward and backward until it can no longer be moved (Illustration 10). The knife will now be in the home position. This can be verified by viewing the position of the knife in the knife lockout window (Illustration 11a). Discard the instrument and dispose of the battery pack.

Once in the cavity, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument

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(Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger to reopen the jaws and return the

closing trigger to its original position (Illustration 8b).

11 If needed, rotate the jaws by pushing on the fins of the rotating knob. (Illustration 12). The shaft will rotate optimally when the jaws of the device are open.

12 To articulate the jaws inside the body cavity, turn the articulation lever in either direction (Illustration 2) until one of the three pre-determined angles is reached. The predetermined angles are equal increments up to 50° from the straight configuration. Articulating while closed does not provide optimal feedback and control, and may result in instrument damage.

WARNING: Do not attempt to articulate by pressing the front of the jaws against the grounding surfaces as device damage or tissue trauma may occur.

Caution: The instrument can only achieve a maximum articulation of angle of 50°. When the maximum angle is reached, the articulation lever force will increase, indicating the maximum angle has been reached.

Caution: Attempting to articulate beyond the stop in either direction, or while the instrument jaws are within the trocar, may result in instrument damage.

13 Position the instrument around the tissue to be stapled, ensuring that the tissue is proximal to the cut line markings, which represent the minimum cut distance (Illustration 13).

Caution: Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload, particularly in the crotch of the jaws, may result in an incomplete staple line.

WARNING: Do not use the instrument on the parenchyma of the liver, pancreas, kidney, or spleen, or other tissues which cannot comfortably compress to the closed staple height.

WARNING: Do not use the instrument if, upon closing, the tissue extends beyond the cut line marking, which represents the minimum cut distance, such as in creating an anastomosis

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
(Illustration 13). If this occurs, open and reposition the instrument. Firing with tissue in the tip of the jaws may result in incomplete cutting action, and/or improperly formed staples.

Caution: When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.

Caution: If using other devices to guide instrument placement on the tissue, such as catheters that are attached to the anvil or reload, those other devices should be removed prior to clamping and firing. Any device or tissue that is clamped within the tip of the instrument can cause improperly formed staples, or be stapled to the targeted tissue.

Caution: Crossing of staple lines may shorten the life of the instrument.

14 After positioning the instrument jaws, close the jaws by squeezing the closing trigger until it locks (Illustration 5). An audible click indicates that the closing trigger and the jaws are locked. When the jaws of the instrument are closed, the red firing trigger lock and firing trigger are exposed.


Caution: Ensure tissue has not extended (extruded) proximal to the proximal black line on the instrument (Illustration 1). Tissue forced into the instrument proximal to the black line may be transected without staples.

Caution: If the closing trigger is difficult to lock, reposition the instrument and take a smaller amount of tissue. Ensure that the proper reload and instrument selection has been made. (Refer to the Reload Product Code Chart.)

Caution: If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use instrument.

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Note: When firing across thick tissue, holding the jaws in place for 15 seconds after closing and prior to firing may result in better compression and staple formation.

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15 Pull back the red firing-trigger lock to enable the firing trigger to be pulled (Illustration 14).

16 Fire the instrument by pulling the firing trigger; the motor will activate audibly (Illustration 15). Continue to depress the trigger until the motor stops (audible feedback), in order to reach the end of the transection.

Caution: Since the motor may stop if it stalls, it is important to do a visual check to ensure that the knife is visible in the knife blade windows, and has reached the end of the transection.

Caution: Attempting to force the device to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in motor stall and the knife will stop. If this happens, release the firing trigger, slide the knife reverse switch forward (Illustration 7), and remove and reload the instrument. Then, reposition the instrument and take a smaller amount of tissue. Ensure that the proper reload and instrument selection has been made. (Refer to the Reload Product Code Chart).

Note: The motor will slow down when device encounters more force.

17 To complete the firing sequence, release the firing trigger to activate the motor and automatically return the knife to home position where the motor will stop. In this position, the instrument is locked out until the jaws are opened and re-closed.

Note: If it is necessary to interrupt the firing sequence or it is interrupted inadvertently by releasing the trigger during the firing sequence, pull the firing trigger again to continue. The status of the transection can be determined by observing the knife blade window on the underside of the channel jaw throughout the firing action (Illustration 1). This can be done as often as necessary until the knife reaches its end where releasing the trigger automatically returns the knife to home position.

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Caution: If the instrument locks out, the motor will stop. Release the firing trigger and slide the knife reverse switch forward to return the knife to the home position (Illustration 7). In this position, the instrument should be removed, opened, and reloaded in order to continue. To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger (Illustration 8b). Follow the instructions for Reloading the Instrument.

Note: At any time, if the knife reverse switch does not return the knife to home position (Illustration 11b or 11c), and the jaws will not open:

a. First, ensure the battery pack is securely installed and the instrument has power (Illustration 9). Then, try the knife reverse switch again (Illustration 7).

b. If the knife still does not return (Illustration 11b or 11c), use the manual override.

To use the manual override, remove the access panel labeled “Manual Override” on the top of the instrument handle. The manual override lever will be exposed. Move the lever forward and backward until it can no longer be moved (Illustration 10). The knife will now be in the home position. This can be verified by viewing the position of the knife in the knife lockout window (Illustration 11a). Discard the instrument and dispose of the battery pack.

Caution: After the manual override system is used, the instrument is disabled and cannot be used for any subsequent firings. Caution: Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficulty removing the device. Caution: If the firing mechanism becomes inoperative, do not continue to use the instrument.

18 To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument. While pressure is still on the anvil release button, slowly release the closing trigger (Illustrations 8a and 8b).

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Caution: If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the home position. The position of the knife can be determined by observing the knife in the knife lockout window. (Illustration 11a). If the knife is not in the home position or the position of the knife cannot be determined, slide the knife reverse switch to activate the motor and return the knife to home position (Illustration 7). Try opening the jaws again using the anvil release button on either side of the instrument.

Note: At any time, if the knife reverse switch does not return the knife to home position (Illustration 11b or 11c), and the jaws will not open:

a. First, ensure the battery pack is securely installed and the instrument has power (Illustration 9). Then, try the knife reverse switch again (Illustration 7).

b. If the knife still does not return (Illustration 11b or 11c), use the manual override.

Caution: If the jaws do not open at this point, forcing the closing trigger upward (away from the handle) until the closing trigger returns to its original position can enable the opening of the jaws. When an instrument is forced open, discontinue use as this may expose the knife, or damage the instrument, increasing the risk of injury to user or patient.

19 Gently pull the instrument away from the transected tissue and ensure it is released from the jaws.

Caution: Examine the staple lines for hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.

20 Before removing articulating instruments, move the jaws away from any obstruction inside the body cavity while keeping the jaws open and within the field of view, and turn the articulation lever to return to the straight position.

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21 To remove the instrument from the cavity, squeeze the closing trigger until it locks, closing the jaws (Illustration 5). Completely withdraw the instrument in the closed position.

Reloading the Instrument

22 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

23 Prior to reloading, ensure the instrument is in the open position (Illustration 1).

24 Push upward (toward the anvil) to unsnap the reload from the channel jaw. Discard the used reload (Illustration 16). **WARNING:** Prior to reloading the instrument, hold the instrument in a vertical position, with anvil and channel jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the anvil and channel jaw to clear any unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and channel jaw.

25 Examine the new reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload. Caution: Tissue thickness should be carefully evaluated

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prior to using the instrument. Refer to the Reload Product Code Chart. **WARNING:** Loading the incorrect reload for the instrument size or model may result in transecting the tissue without stapling (e.g. loading a 35 mm reload into a 45 mm instrument).

26 Insert the new reload by sliding it against the bottom of the channel jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use (Illustration 4).

Caution: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any blue-colored drivers are visible. (If blue-colored drivers are visible, the reload may not contain staples.)

Battery Pack Disposal

The battery pack contains a built-in battery drain and must be installed into the instrument to trigger the automatic battery drain feature. It does not need to remain in the instrument – it will continue to drain after it has been removed from the instrument. Once removed, depending on your local regulations, it can be put directly into the battery recycling bin or normal waste stream.

Note: Batteries used in this device contain 1,2-dimethoxyethane.

Prior to Installation into the Instrument

If battery pack needs to be disposed of prior to installation into instrument (e.g. product is beyond expiration date indicated on the package, or battery pack is dropped), first install the battery pack into the instrument, then remove, to trigger built-in battery drain.

After Instrument Use

The battery pack must be removed from the instrument prior to disposal.

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If battery pack requires decontamination prior to disposal, follow hospital protocol or the Battery Pack Cleaning and Disinfection instructions below.

To remove the battery pack, squeeze the release tabs and pull the battery pack straight back (Illustration 17). Note: It is not necessary to disassemble the battery pack.

Battery Pack Cleaning and Disinfection

WARNING: Do not use autoclaves, ethylene oxide, or radiation to sterilize or disinfect the battery pack.

Manual Cleaning

1 Remove the battery pack from the instrument before cleaning.

Note: Battery pack should not be submerged in water or cleaning solutions.

2 Clean the battery pack surfaces with a neutral pH detergent or neutral pH enzymatic detergent, prepared according to the manufacturer's instructions.

3 Use soft bristle brush to manually clean the battery pack with the cleaning solution.

4 Ensure areas containing crevices are scrubbed thoroughly.

5 Wipe off detergent thoroughly with lukewarm tap water.

6 Perform visual inspection to determine if debris is removed.

7 Repeat cleaning as necessary to obtain a visually clean battery pack.

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Disinfectants should be prepared and used according to the manufacturer's recommendations. It is recommended that the chemical disinfectant be wiped off with tap water.

Simple Green D Pro 3
10% Bleach (sodium hypochlorite solution)
70% Isopropyl alcohol

Standard Conventions Used

Information relative to the completion of a task in a safe and thorough manner will be supplied in the form of a WARNING, a Caution, or a Note statement. These statements are found throughout the documentation.

These statements should be read before continuing to the next step in a procedure.

WARNING: A WARNING statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in personal injury or loss of life.

Caution: A Caution statement indicates an operating or maintenance procedure, practice, or condition that, if not strictly observed, could result in damage to or destruction of the equipment.

Note: A Note statement indicates an operating or maintenance problem, practice, or condition that is necessary to accomplish a task efficiently.

Warnings and Precautions

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.

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- After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any blue-colored drivers are visible. (If blue-colored drivers are visible, the reload may not contain staples.)
- For insertion and removal of articulating instruments, the jaws of the instrument must be straight, parallel to the shaft of the instrument. Failure to have the instrument jaws in the straight position will result in difficult insertion or withdrawal of the instrument and may result in damage to the instrument.
- When placing the instrument through the trocar or incision, avoid inadvertently pulling the red firing trigger lock and the firing trigger. The instrument may be partially or completely fired and will need to be reloaded before using on tissue. If the instrument is partially fired, remove the instrument and slide the knife reverse switch forward to return the knife to home position (Illustration 7). To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger (Illustration 8b). Remove the reload and insert a new one (see Reloading the Instrument). If the instrument is fired completely, the knife will return to home position automatically. Remove the instrument and insert a new reload (see Reloading the Instrument).
- Do not attempt to articulate by pressing the front of the jaws against the grounding surfaces as device damage or tissue trauma may occur.
- The instrument can only achieve a maximum articulation of angle of 50°. When the maximum angle is reached, the articulation lever force will increase, indicating the maximum angle has been reached.
- Attempting to articulate beyond the stop in either direction, or while the instrument jaws are within the trocar, may result in instrument damage.
- Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the reload,

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particularly in the crotch of the jaws, may result in an incomplete staple line.

- Do not use the instrument on the parenchyma of the liver, pancreas, kidney, or spleen, or other tissues which cannot comfortably compress to the closed staple height.

- Do not use the instrument if upon closing, the tissue extends beyond the cut line marking, which represents the minimum cut distance, such as in creating an anastomosis (Illustration 13). If this occurs, open and reposition the instrument. Firing with tissue in the tip of the jaws may result in incomplete cutting action, and/or improperly formed staples.

- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the instrument jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the instrument jaws.

- If using other devices to guide instrument placement on the tissue, such as catheters that are attached to the anvil or reload channel, those other devices should be removed prior to clamping and firing. Any device or tissue that is clamped within the tip of the instrument can cause improperly formed staples, or be stapled to the targeted tissue.

- Crossing of staple lines may shorten the life of the instrument.

- Attempting to force the device to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in motor stall and the knife will stop. If this happens, release the firing trigger, slide the knife reverse switch forward (Illustration 7), and remove and reload the instrument. Then, reposition the instrument and take a smaller amount of tissue. Ensure that the proper reload and instrument selection has been made. (Refer to the Reload Product Code Chart)

- If the closing trigger is difficult to lock, reposition the instrument and take a smaller amount of tissue. Ensure that the proper reload and instrument selection has been made. (Refer to the Reload Product Code Chart.)

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- If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the instrument. Remove and do not continue to use instrument.
- Since the motor may stop if it stalls, it is important to do a visual check to ensure that the knife is visible in the knife blade windows, has reached the end of the transection.
- If the instrument locks out, the motor will stop. Release the firing trigger and slide the knife reverse switch forward to return the knife to the home position (Illustration 7). In this position, the instrument should be removed, opened, and reloaded in order to continue. To open the jaws, squeeze the closing trigger, then simultaneously press the anvil release button on either side of the instrument (Illustration 8a). While pressure is still on the anvil release button, slowly release the closing trigger (Illustration 8b). Follow the instructions for Reloading the Instrument.
- After the manual override system is used, the instrument is disabled and cannot be used for any subsequent firings.
- Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficulty removing the device.
- If the firing mechanism becomes inoperative, do not continue to use the instrument.
- If the jaws do not automatically open after the anvil release button is pressed, first ensure that the knife is in the home position. The position of the knife can be determined by observing the knife in the knife lockout window. (Illustration 11a). If the knife is not in the home position or the position of the knife cannot be determined, slide the knife reverse switch to activate the motor and return the knife to home position (Illustration 7). Try opening the jaws again using the anvil release button on either side of the instrument.
- If the jaws do not open at this point, pushing the closing trigger upward (away from the handle) until the closing trigger returns to its original position can enable the opening of the jaws. When an instrument is forced open, discontinue use as this

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may expose the knife, or damage the instrument, increasing the risk of injury to user or patient.

- Examine the staple lines for hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.
- Prior to reloading the instrument, hold the instrument in a vertical position, with anvil and channel jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the anvil and channel jaw to clear any unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and channel jaw.
- Do not use autoclaves, ethylene oxide, or radiation to sterilize or disinfect the battery pack.
- Use of any other type of battery other than the battery supplied with the device may result in increased EMISSIONS or decreased IMMUNITY of the ECHELON FLEX Stapler.
- Portable radio and mobile RF communications equipment can affect Medical Electrical Equipment. Follow the directions of the following Electromagnetic Emissions and Electromagnetic Immunity Guidance tables when using the ECHELON FLEX Stapler.
- Avoid use of the ECHELON FLEX Stapler adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the ECHELON FLEX Stapler and the other equipment to ensure normal operation.
- Risk of non-ionizing radiation. Interference may occur in the vicinity of this equipment.
- The battery may present a fire hazard if mistreated. Do not disassemble, heat above 100 °C, autoclave, crush, puncture, short external contacts or recharge.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only. Multiple patient use may compromise the device integrity or create a risk of contamination that, in turn, may result in patient injury or illness.

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Specifications

The ECHELON FLEX Stapler has a power rating of 25 W.

The ECHELON FLEX Stapler is resistant to water ingress and is classified per IEC 60601-1 as IPX0.

The ECHELON FLEX Stapler needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

WARNING: Use of any other type of battery other than the battery supplied with the device may result in increased EMISSIONS or decreased IMMUNITY of the ECHELON FLEX Powered Vascular Stapler.

WARNING: Portable radio and mobile RF communications equipment can affect Medical Electrical Equipment. Follow the directions of the following Electromagnetic Emissions and Electromagnetic Immunity Guidance tables when using the ECHELON FLEX Stapler.

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WARNING: Avoid use of the ECHELON FLEX Stapler adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the ECHELON FLEX Powered Vascular Stapler and the other equipment to ensure normal operation.

Storage Conditions

Temperature: -20 °C — 30 °C
Relative Humidity: 10% — 65%

Store away from direct sunlight and heat.

Transport Conditions

Temperature: -20 °C — 60 °C
Relative Humidity: 10% — 65%

Operating Environment Conditions

Temperature: 10 °C – 40 °C
Relative Humidity: 30% – 75%

Disposal

After expiration date of instrument and reloads (period of preservation of sterility) and/or in damaged packaging it must be disposed of as medical waste of class A according to SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" in accordance with the rules of the medical institution.

After use, the product is a potential biological hazard from contamination with blood and other biological materials. Recycling should also be in accordance with accepted medical practice standards, and relevant local, state, and Federal

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legislative acts and instructions, as the waste classes B according to SanPiN 2.1.7.2790-10.

Sterilization

Sterilization is made by gamma irradiation.
Sterility term – 3 years

Shelf life

3 years

Standard configuration

The instruments of cutters are supplied sterile for single patient use. Discard after use.

The contents of cardboard packaging (box):

- instrument 1 piece.

Instruction manual (1 piece) is supplied together with group package.

Note: for an explanation of the symbols on the package see the end of this manual.

ETHICON ENDO-SURGERY, LLC
475 Calle C, Guaynabo, PR 00969 USA

Manufacturer's warranty

Manufacturer is not liable for any damage or accidents resulting from improper use, human error or external emergencies. Regarding the product quality, please contact an authorized manufacturer's representative.

Authorized manufacturer's representative

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LLC Johnson & Johnson
Krylatskaya 17, Building 2, 121614 Moscow, Russia
tel. (495)5807777, fax (495)5807878.

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Interpretation of symbols

	See Instructions For Use
	Follow Instructions for Use (refer to blue symbol on outer packaging).
	Type CF Applied Part
	Caution: The battery may present a fire hazard if mistreated. Do not disassemble, heat above 100 °C, autoclave, crush, puncture, short external contacts or recharge.
	WARNING: Risk of non-ionizing radiation. Interference may occur in the vicinity of this equipment.
	Battery pack. Return waste to a collection system or treatment and recycling facilities. Applicable in the EU. Follow decontamination instructions before returning waste.
	Single Patient Use
	Sterilized by Irradiation. Sterility Guaranteed Unless Package Opened or Damaged. Do Not Resterilize.

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	Temperature limitation
	Humidity limitation
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Manufacturer/Date of Manufacture
	Catalogue Number
	Lot
	Authorized Representative in the USA
	Authorized Representative in the European Community
	Product is certified by a Nationally Recognized Testing Laboratory.
	Use Until Date
	UL Recognized Component
	Compliance with requirements of European Directive 93/42/EES

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