# Sonicision™

Curved Jaw Cordless Ultrasonic Dissector

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















## Sonicision™

**REF** SCDA13 **Curved Jaw Cordless Ultrasonic Dissector**  
5 mm - 13 cm

**REF** SCDA26 **Curved Jaw Cordless Ultrasonic Dissector**  
5 mm - 26 cm

**REF** SCDA39 **Curved Jaw Cordless Ultrasonic Dissector**  
5 mm - 39 cm

**REF** SCDA48 **Curved Jaw Cordless Ultrasonic Dissector**  
5 mm - 48 cm

Not made with natural rubber latex

### Indications for Use

The Sonicision Curved Jaw Cordless Ultrasonic Dissection Device is indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The device can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, plastic, pediatric, gynecologic, urologic, exposure to orthopedic structures (such as spine and joint space) and other open and endoscopic procedures. The Sonicision Curved Jaw Cordless Ultrasonic Dissection Device can be used to coagulate isolated vessels up to and including 5 mm in diameter.

The Sonicision 13 cm device is also indicated for use in otorhinolaryngologic (ENT) procedures.

### Contraindications

- The device is not indicated for incising bone.
- The device is not intended for contraceptive tubal occlusion.


### Warning

Safe use of the Sonicision Curved Jaw Cordless Ultrasonic Dissector requires the user to read, understand, and follow the accompanying instructions.

The dissector cannot be adequately cleaned or sterilized for safe reuse and is, therefore, intended for single use. Attempts to clean and sterilize the dissector for reuse may result in infection or product-failure risks to the patient and operator.

Visually inspect all system components for breaks, cracks, nicks, or other damage. Do not use damaged components. Use of damaged components may result in injury to the patient or operator.

Contact between the active blade and other metal objects (hemostats, clips, staples, retractors, etc.) may result in product damage, such as a broken blade. Pieces of a broken blade may fall into the surgical cavity causing unintended tissue damage.

- ① Clamping jaw – Holds tissue against active blade
- ② Active curved blade – Delivers energy effect
- ③ Shaft
- ④ Shaft-rotation wheel (black)– Rotates shaft 360°
- ⑤ Speaker – Emits system status tones
- ⑥ Dual-mode energy button (blue)- Selects energy activation in minimum or maximum mode
- ⑦ Jaw lever (black) – Closes and opens clamping jaw
- ⑧ Shaft measurement marks - measures up to 10 cm from the distal end (not available on all models)
- ⑨ Battery door (black) – Closes and seals battery compartment
- ⑩ Battery door release button (black) – Opens the battery door

Do not use if package is opened or damaged

Defibrillation-proof type BF applied part

## Warning

Use the device with caution when operating in close proximity to external or implanted electronic medical devices. Electromagnetic interference produced by the proximity of the ultrasonic generator can potentially cause a device to enter an unsafe mode. There is also a possibility that a nearby electronic device could alter the performance of the Sonicision device. Consult the manufacturer of the external or implanted electronic device or responsible hospital department for further information when use is planned in close proximity to such devices.

**Infection hazard** Non-sterile components may contaminate the sterile field and place the patient at risk for infection.

- Do not use the dissector if the sterile package is opened, damaged, or has exceeded the expiration date.
- Generators and battery insertion guides are provided non-sterile and must be cleaned and sterilized prior to first use and for each reuse. Refer to the cleaning and sterilization instructions in the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide*.
- Battery packs are non-sterile and cannot be sterilized. The battery pack must be aseptically inserted into the dissector according to instructions. See *Inserting the Battery Pack* on page 7.

**Fire and Burn Hazard** The clamping jaw, active blade, and distal 6 cm of the shaft are hot during, and for a period of time after activating the device. Avoid the potential for fire and unintended burns by always keeping the distal end of the device in sight and away from direct contact with adjacent tissue, textiles, and other flammable materials during and for a period of time after activation.

## Precaution

Ensure that the battery pack is sufficiently charged before use. Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for charging instructions.

Use of a battery pack that is not full charged can result in a procedural delay.

Use only Sonicision components listed in the compatibility table on the front cover. Components from other manufacturers are not compatible with this system and may cause injury to the patient and operator.

## Important

Instructions provided in this document are limited to the initial assembly, change of a battery pack during a procedure, and disassembly of the device. Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* provided with the Sonicision Reusable Generator for complete information and important warnings and cautions.

Medical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide*.

Contact Covidien to obtain additional or replacement copies of the instructions for use and *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* if needed.

## Notice

This device is intended to be used with a 5 mm trocar when used endoscopically. Verify proper trocar size and compatibility prior to using the device in a procedure.

Ensure the battery pack connector is dry. Attaching a battery pack or inserting into the battery charger when the connector is wet may damage the components.

Battery packs are shipped with a low state of charge and must be charged upon receipt.

To avoid electrostatic-damage to a battery pack, do not touch the battery pack connector.

## Description

The Sonicision Curved Jaw Cordless Ultrasonic Dissector is a sterile, single-use component for use with the Sonicision Reusable Generator and Reusable Battery Pack. The dissector provides control for device functions such as power levels, blade placement and position, grasping, coagulating, and dissecting tissue. It is designed to be inserted and extracted through a compatible 5 mm trocar, when used endoscopically. Refer to the figure on page 5 for dissector features.

## Assemble the Sonicision Curved Jaw Cordless Ultrasonic Dissection Device

Complete assembly of a working Sonicision Curved Jaw Cordless Ultrasonic Dissection Device requires the following components:

- Single-use dissector
- Reusable generator
- Reusable battery pack
- Reusable battery insertion guide

## Attaching the Generator

1. **Circulating Person:** Transfer the dissector, battery insertion guide(s), and generator aseptically from each sterile barrier into the sterile field.
2. **Scrubbed Person:** Slide the generator into the opening until it makes contact with the dissector (step 1.1 on page 2). Hold the dissector shaft-rotation wheel. Turn the torque knob clockwise until it clicks twice (step 1.2 on page 2).

## Inserting the Battery Pack

### Important

Attach the generator to the dissector prior to inserting the battery pack.

1. **Scrubbed Person:** Hold the handle with the open battery compartment up. Open the battery compartment door.
2. **Scrubbed Person:** Place the sterile battery insertion guide onto the open dissector handle and hold firmly (step 1.3 on page 2).
3. **Circulating Person:** Orient the battery pack as shown with the battery insertion guide. Insert the battery pack through the battery insertion guide into the dissector handle (step 1.4 on page 2).
4. **Circulating Person:** After the battery pack is seated in the dissector handle, remove the contaminated battery insertion guide from the sterile field (step 1.5 on page 2).

5. **Scrubbed Person:** Do not touch the battery pack. Close the battery compartment door from the outside until there is tactile feedback and an audible click (step 1.6 on page 2).

A series of audible tones and an illuminated green LED on the generator indicate correct assembly.

6. In order to ensure the device meets required essential performance, it should be tested after assembly and each battery pack exchange as follows:

With the clamping jaws open, test the minimum-power mode (energy button depressed to first position) and maximum-power mode (energy button fully depressed). If the device generates different pulsating tones for each power mode and a green LED is present, assembly is complete and the device is ready for use.

If inaudible or undifferentiated tones or a non-green LED are present, refer to the Device Status and Troubleshooting section of the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide*.

## Changing a Battery Pack During a Procedure

### Warning

Do not use a contaminated battery insertion guide from the previous assembly. A new sterile battery insertion guide must be used for each procedural battery pack exchange. Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for cleaning and sterilization instructions.

1. **Scrubbed Person:** Hold the device with the battery compartment up. Press the battery door release button to open the compartment (step 2.1 on page 3).
2. **Scrubbed Person:** Place the new sterile battery insertion guide over the exposed battery pack and seat it onto the open dissector handle and hold firmly (step 2.2 on page 3).
3. **Circulating Person:** Grasp the battery pack handle and pull the battery pack through the insertion guide out of the dissector and into a clean, non-sterile environment (step 2.3 on page 3).
4. **Circulating Person:** Orient the new battery pack as shown with the battery insertion guide. Insert the battery pack through the battery insertion guide into the dissector handle (step 2.4 on page 3).
5. **Circulating Person:** After the battery pack is seated in the dissector handle, remove the contaminated battery insertion guide from the sterile field (step 2.5 on page 3).
6. **Scrubbed Person:** Do not touch the battery pack. Close the battery compartment door from the outside until there is tactile feedback and an audible click (step 2.6 on page 3).
7. With the clamping jaws open, test the minimum-power mode (energy button depressed to first position) and maximum-power mode (energy button fully depressed) as described in *Assemble the Sonicision Curved Jaw Cordless Ultrasonic Dissection Device* on page 6.

## Disassembly

### Warning

Disassembly of the device with these instructions should only occur after the procedure is complete and sterile field is broken.

## Notice

Take particular care to keep the battery pack clean.

### Remove the Battery Pack

1. Hold the dissector handle with the battery compartment up. Press the battery door release button (step 3.1).
2. Use the battery pack handle to pull the battery pack out into a clean environment (step 3.2).

Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for cleaning, disinfecting, and charging instructions.

### Remove the Generator

1. To loosen the torque knob, hold the shaft-rotation wheel firmly with one hand while turning the torque knob counterclockwise with the other hand (step 3.3 on page 4).
2. Continue to turn the torque knob counterclockwise until the generator can be removed from the dissector (step 3.4 on page 4).

### After Use

Reprocess the generator and battery insertion guides as soon as possible following use. If reprocessing cannot be performed immediately, cover the components with a moist towel.

### More Information

Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for information regarding the following elements:

- Using the device
- Other component functions
- Disposal
- Cleaning and sterilizing battery insertion guides and generators, and cleaning and disinfecting battery packs
- Troubleshooting
- Technical specifications

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© 2016 Covidien.
Covidien Ireland Limited,