



JETSTREAM™ Atherectomy System

CALCIUM. PLAQUE. THROMBUS.  
**TREAT IT ALL**



# CALCIUM. PLAQUE. THROMBUS.

Jetstream rotational atherectomy is engineered to predictably treat real-world lesions. Patients with PAD display a wide range of lesion characteristics such as long, diffuse disease and CTOs — which often include mixed morphologies like calcium, plaque and thrombus. **Jetstream is the only atherectomy device designed to *treat it all*.**



## **ROTATIONAL BLADES – Create Concentric Lumens**

Rotational blades spin at ~70,000 RPMs to create concentric lumens, optimizing balloon-to-wall apposition for DCB or other adjunctive therapies.

## **FRONT-CUTTING BLADES – Immediately Engage Lesions**

Five front-cutting blades immediately engage lesions and help enable the treatment of tight or occluded vessels.

## **EXPANDABLE BLADES – Provide Sizing Flexibility**

“Blades Down/Blades Up” technology enables maximum luminal gain while providing the flexibility to treat multiple vessel diameters with the same catheter.

## **ACTIVE ASPIRATION – Helps Reduce Embolization Risk**

Dynamic and continuous aspiration mechanically removes debris, helping to minimize the risk of distal embolization, and debulk the lesion.

## **DIFFERENTIAL CUTTING – Deflects Away from Healthy Tissue**

The mechanism of action allows the blades to cut the diseased, inelastic, tissue while deflecting away from the healthy, elastic, tissue.


# TREAT IT ALL




# JETSTREAM CASE EXAMPLES

## CASE 1:

**Chronic Total Occlusion of the Superficial Femoral Artery**



Hydrophilic 0.035" wire and support catheter used to cross SFA CTO




**Stand-alone Jetstream result**  
2 passes blades down,  
1 pass blades up with  
2.4/3.4 mm XC catheter



**Post DCB**  
Two 6.0 x 100 mm drug-coated balloons

## CASE 2:

**Adductor Canal Disease**




**Adductor Canal disease**

**Stand-alone Jetstream result**  
(IVUS image above)  
2.1/3.0 mm XC catheter

**Final angiogram**  
following  
5 x 80 balloon

## CASE 3:

**Common Femoral Artery Disease**



**Left Common Femoral disease**

**IVUS baseline of 4.4 mm<sup>2</sup> pre-treatment**  
(270 degree arc of calcium)

**Stand-alone Jetstream result,**  
revealed by angiogram  
(prior to PTA)

**IVUS imaging revealed lumen area of 12.5 mm<sup>2</sup> post-treatment**

# REAL-WORLD CLINICAL DATA

## JET REGISTRY<sup>1</sup> — Treatment effects of Jetstream Atherectomy System

The JET Registry demonstrated a high freedom from TLR rate at 12-months and low distal embolization rate in patients with long (16.4 cm), real-world lesions.

### Patient and Lesion Characteristics:

- 241 patients with 258 lesions
- 41% diabetic
- 16.4 cm lesion length
- 36.1% occluded
- 90% of lesions had visible calcium
- 47.7% Grade 3 and 4 calcium present

### Procedure Details:

- 22.4% of cases used embolic protection
- 4.7 minutes average Jetstream Runtime
- 1.4% distal embolization rate

### Key Clinical Results:

**77.2%**  
Patency\*

**81.7%**  
Freedom from TLR

**73.4%**  
Rutherford Category Improvement

**58.2%**  
ABI Improvement

Drug Coated Balloons were not used in this study

70% of patients had no or minimal symptoms (Rutherford Category 0-1)

**Post-Procedure:** 98.3% of patients had  $\leq 30\%$  residual diameter stenosis

\*Patency based on a DUS PSVR  $\leq 2.5$ ; Binary Restenosis was reported as 22.8%. The JET Registry had limited DUS follow-up at 12 months (57,241 patients)

## JET-SCE2 — Jetstream + DCB

In the JET-SCE, the TLR rate was significantly reduced with atherectomy and adjunctive DCB compared to atherectomy with adjunctive PTA at 18-months.

### Patient and Lesion Characteristics:

- 81 patients
- 53.1% Diabetic
- 25.9% CTOs
- 14.9 cm average lesion length in PTA cohort
- 12.0 cm average lesion length in DCB cohort

### Key Clinical Results:

At 18-months results demonstrated...

**91%**  
**JETSTREAM + DCB**  
Freedom from TLR\*

\*18-month Jetstream + PTA = 63.7% Freedom from TLR


## JETSTREAM CALCIUM STUDY<sup>3</sup>

The Jetstream Calcium study demonstrated Jetstream's ability to create statistically significant luminal gain in severe and moderate calcium as measured by IVUS.

### Patient and Lesion Characteristics:

- 55 patients treated with Jetstream
- 56% Diabetic
- 63.6% Severe Calcium
- $>90^\circ$  superficial calcium,  $>5$  mm in length

### Key Clinical Results:



# SPECIFICATIONS

Catheter Length	Min. Introducer Size	Max. Guidewire Diameter	Tip Diameter	Target Therapy Speed	GTIN	UPN/Order Code	Catalog Number	Unit	Qty
<b>Jetstream™ 2.4/3.4 mm XC Atherectomy Catheter</b>									
120 cm	7 F	0,014"	2,4 mm 3,4 mm	70K rpm	08714729889922	112266-001	PV41340	Each	1
<b>Jetstream™ 2.1/3.0 mm XC Atherectomy Catheter</b>									
135 cm	7 F	0,014"	2,1 mm 3,0 mm	70K rpm	08714729889892	112264-001	PV31300	Each	1
<b>Jetstream™ 1.85 mm SC Atherectomy Catheter</b>									
145 cm	7 F	0,014"	1,85 mm	73K rpm	08714729889861	112262-001	PV3118F	Each	1
<b>Jetstream™ 1.6 mm SC Atherectomy Catheter</b>									
145 cm	7 F	0,014"	1,6 mm	73K rpm	08714789889830	112260-001	PV3116F	Each	1

Consoles	UPN/Order Code
<b>Jetstream™ Console EU</b>	<b>050599-010</b>
<b>Jetstream™ Console UK</b>	<b>050500-020</b>

The C-Code used for the Jetstream Atherectomy System is C1724. C-Codes are used for hospital outpatient device reporting for Medicare and some private payers. Note: Boston Scientific Corporation is not responsible for correct use of codes on submitted claims; this information does not constitute reimbursement or legal advice.

- Garcia, L. (2017). Jetstream atherectomy in treating de novo or non-stent restenotic femoropopliteal disease: One-year results from the JET registry. Registry results presented at the Leipzig Interventional Course (LINC), Leipzig, Germany
- Shammas, N (2017). Long Term Outcomes with Jetstream Atherectomy System with or without Drug Coated Balloons in Treating Femoropopliteal Arteries: A Single Center Experience (JET-SCE). JET-SCE results presented as a poster at Society for Cardiovascular Angiography and Interventions (SCAI) Scientific Sessions, New Orleans, LA.
- Maehara A, Mintz G, Shimshak T, Ricotta J, Ramaiah V, Foster M, Davis T, Gray W. Intravascular ultrasound evaluation of JETSTREAM atherectomy removal of superficial calcium in peripheral arteries. EuroIntervention 2015;11:96-103



**Jetstream System Components**

## Jetstream™ Catheters combined with Console

CAUTION: All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.

PI-852208-AA Printed in Germany by medicalvision.

**Boston  
Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.