
Valleylab™

Reusable Pencils

REF E2100**REF** E2100E**en** Instructions for Use**fr** Mode d'emploi**pt** Instruções de uso**da** Brugervejledning**de** Gebrauchsanleitung**es** Instrucciones de uso**fi** Käyttöohjeet**it** Istruzioni per l'uso**nl** Gebruiksaanwijzing**no** Bruksanvisning**sv** Bruksanvisning**ru** Инструкция по применению**zh** 使用说明**Compatible with:**

| | |
|--------|--|
| E2100 | FORCETRIAD SW v3.2 or higher VLFT10GEN SW v1.0 or higher Force FX Series |
| E2100E | ERBE™* ICC Series ERBE™* ACC Series |

① Electroswitching Pencil

② Generator

REF E2100 **Rocker Switch Pencil, Reusable**
15' (4.6 m) Cord

REF E2100E **Rocker Switch Pencil, ERBE™*, Reusable**
15' (4.6 m) Cord

E2100E not for sale in U.S.

For use with a maximum peak voltage of 5600V.

Not made with natural rubber latex

Handswitching pencil

Do not use if package is opened or damaged

Type CF applied part

Indications for Use

This device is intended to be used in general surgical procedures where electroswitching cutting and/or coagulation is required.

Product Description

The E2100 is a reusable electroswitching instrument intended for use in general surgical procedures when monopolar electroswitching cutting and/or coagulation is required. It is for

use with electrodes with a 0.150 inch diameter molded insulator.

The E2100 series handswitching instruments are designed with a rocker style switch capable of CUT or COAG mode. The instruments are supplied non-sterile with a reusable electrode with locking hex.

Warning

Do not activate an E2100 Valleylab™ Rocker Switch Pencil while installing or removing the electrode. The user may also turn the generator off, place the generator in lowest power settings, or unplug the pencil from the generator during manipulation of the electrode.

Prior to use, ensure the pencil has been cleaned and sterilized in accordance with the following processing instructions.

Fire Hazard Always place the active electrode in a clean, dry, insulated safety holster when not in use.

- Electroswitching accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.
- Electroswitching accessories may cause fire or burns if placed close to or in contact with flammable materials, such as gauze or surgical drapes. Place longer electrodes, such as extended electrodes, away from the patient and drapes.

Warning

Fire/Explosion Hazard Sparking and heating associated with electrosurgery can provide an ignition source.

- Keep gauze and sponges wet.
- Keep electrosurgical electrodes away from flammable materials and environments enriched with oxygen (O₂).
- Use of electrosurgery in oxygen-enriched environments increases the risk of fire. Therefore, take measures to reduce the oxygen concentration at the surgical site. If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.
- Avoid enriched oxygen (O₂) and nitrous oxide (N₂O) atmospheres near the surgical site, especially during head and neck surgery. Both oxygen and nitrous oxide support combustion and may result in fires and burns to patients or surgical personnel.
- Prevent pooling of flammable fluids and the accumulation of flammable or oxidizing gases or vapors under surgical drapes or near the surgical site.
- Do not activate the generator until flammable vapors from skin prep solutions and tinctures have dissipated.
- Avoid the accumulation of naturally occurring flammable gases in body cavities such as the bowel.

The instrument is intended for use only with the equipment listed on the cover of this document and with ERBE™ ICC series and ACC series. Use of this instrument with other generators may not result in the desired tissue effect, may result in injury to the patient or surgical team, may cause damage to the instrument, or may result in electromagnetic interference.

Tissue buildup (eschar) on the tip of an active electrode may create embers that pose a fire hazard, especially in oxygen-enriched environments. Keep the electrode clean and free of all debris.

Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.

Warning**Fire Hazard During Oropharyngeal Surgery**

Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.

If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube, and be sure to keep sponges wet throughout the procedure.

Question the need for 100% oxygen during oropharyngeal or head and neck surgery.

If necessary, scavenge excess oxygen with separate suction.

Confirm proper power settings prior to and during a procedure. Use the lowest power settings that achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power-setting adjustments.

Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may carry electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids:

- Always keep the external surface of the active electrode away from adjacent tissue while activating the generator.
- Clear conductive fluid from the electrode before activating the electrosurgical device.

The electrode must fit completely and securely into the device. An incorrectly seated electrode may result in burns to the patient or surgical personnel.

Fire Hazard Do not modify or add to the insulation of active electrodes. If insulated electrodes are required, use an appropriate Covidien insulated electrode.

Warning

Some surgeons may elect to “buzz the hemostat” during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon's hands may result. To minimize the risk, take these precautions:

- Do not “buzz the hemostat” with a needle electrode.
- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate CUT rather than COAG. CUT has a lower voltage than COAG.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.
- “Buzz the hemostat” below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon's hands.
- When using a coated or nonstick blade electrode, place the edge of the electrode against the hemostat or other metal instrument.

Fire/Explosion Hazard Position surgical electrode cables to avoid contact with the patient or other leads.

Confirm proper electro-surgical generator setting before proceeding with surgery. Use the lowest power setting to achieve the desired effect.

Explosion Hazard Do not use electro-surgery in the presence of flammable anesthetics.

Electric Shock Hazard Do not connect wet accessories to the generator.

Do not use on patients who have electronic implants such as cardiac pacemakers without first consulting a qualified professional. A possible hazard exists because interference with the action of the electronic implant may occur, or the implant may be damaged.

Warning

Due to concerns about the carcinogenic and infectious potential of electro-surgical by products (such as tissue smoke plume and aerosols), protective eye wear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.

The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.

Precaution

Inspect accessories and cords for breaks, cracks, nicks, or other damage before every use. If damaged, do not use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team.

Visual inspection alone may not be sufficient to ensure that the insulation is intact.

Do not exceed maximum power limits as stated in these or other instructions for use. Exceeding recommended power settings may result in patient injury or product damage.

To prevent unintended injury to patient or operating room staff personnel, the conductive metal post of the electrode must be properly and completely installed into the nose of the pencil. Any exposed metal can result in arcing at the electrode and pencil connection.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Do not activate electrodes while in contact with or near other instruments, including cannulas. Localized burns to the patient or physician may result.

Do not activate the generator in an open circuit condition. To reduce the chances of unintended burns, activate the generator only when the active electrode is near or touching the target tissue.

Before use, examine the generator and accessories for defects. Do not use cables or accessories with damaged insulation or connectors.

Activate the generator only when ready to deliver electro-surgical current and when the active tip is in view.

Deactivate the generator before the tip leaves the surgical site.

Notice

When using an electrode, inspect the electrode frequently to verify that the coating is not damaged. If the coating becomes damaged, discard the electrode. Some of the factors that can increase the likelihood of coating damage are higher power settings, longer activation times, and more arcing.

Do not clean the electrode with a scratch pad or other abrasive object. Scraping with a sharp object or bending may damage the electrode. If damaged, discard the electrode.

Wipe the electrode often with moist gauze or other material.

Processing Instructions

This procedure has been developed in accordance with the professional recommendations of the Association for the Advancement of Medical Instrumentation and the Association of Operating Room Nurses.

The product life, cleaning efficacy, and sterility have been validated for the following cleaning and processing instructions when using the recommended detergents according to the manufacturer's instructions.

- Manual cleaning with Hexanios™ or Prolystica™.
- Manual cleaning followed by automatic cleaning in a Steris model Reliance™ 444 single chamber washer with Prolystica™, H-Klenz II™ (11 pH), citric acid, and Hinge-Free™ at 180°F (82°C).
- 50-use life validated with vacuum assisted steam up to 280°F (138°C) 20 minutes duration

Cleaning

Warning

Automated washing without manual cleaning (steps 1 - 4) is ineffective as a stand-alone method.


Notice

Do not use terpene or petroleum-based cleaners. Product damage will result.


Precaution

The electrode must be removed from the pencil for effective cleaning.


1. Remove the electrode from the pencil. Discard the tip protector if present.


- Discard single-use electrodes. Covidien single-use electrodes are manufactured with blue insulation. Covidien reusable electrodes are manufactured with black or gray insulation.
- Clean reusable electrodes with the E2401 Valleylab™ Tip Cleaner.

2. Clean the pencil with a soft brush in a detergent cleaning solution. Rinse the pencil thoroughly for five minutes.


3. Clean the electrode channel with a small brush and a detergent cleaning solution. Rinse the pencil thoroughly for five minutes.


Optional Step

- Place the electrode and the handswitching pencil with the cord loosely coiled in a basket. Process electrode, pencil, and cord in a mechanical washer.


- Inspect the pencil and cable for:
 - Cracks in the rocker cover

- Exposed wires or cracks in the cord

- Bent electrode contact

Discard if the pencil or cord is damaged.

Sterilization Recommendations

Precaution

The electrode must be removed from the pencil for effective sterilization.

Please contact the sterilizer manufacturer for instructions and recommendations specific to the sterilizer system in use.

Steam

Wrapped (Central Supply Room Wrap)

Vacuum Assisted, 270°-280°F (132°-138°C) for 4 minutes, dry for 20 minutes

Gravity Displacement, 270°-280°F (132°-138°C) for 10 minutes, dry for 30 minutes

Gravity Displacement, 250°-268°F (121°-131°C) for 20 minutes, dry for 30 minutes


Unwrapped ("flash" sterilization)

Vacuum Assisted, 270°-280°F (132°-138°C) for 4 minutes, dry for 1 minute



Gravity Displacement, 270°-280°F (132°-138°C) for 10 minutes, dry for 1 minute

Before Surgery


- Follow the instructions for applying the return electrode to the patient and connecting the return electrode to the generator.



- Remove the tip protector from the electrode. Connect the pencil cord to the generator.

- Attach the holster to the surgical drape by pulling the drape material through the slots (E2100 Accessory Holster sold separately).

Notice

Do not attach the grey pencil cord to the holster using the holster cord lock as this may result in damage to the cord insulation.

4. Test the pencil before surgery:

- Turn on the generator.
- Set the generator to 1 watt for both cut and coagulate modes.
- Press the yellow cut button on the pencil. Verify that the yellow cut indicator on the generator illuminates.
- Press the blue coag button on the pencil. Verify that the blue coag indicator on the generator illuminates.

Warning

Hold the pencil safely away from patients, personnel, and surgical drapes during testing.

Precaution

Always confirm correct power settings if the electrode in this device is replaced or exchanged.

The available high-power settings may cause patient injury or product damage. Verify that generator power settings are appropriate for the procedure, device, and return electrode being used.

During Surgery

Place the pencil in an insulated holster or away from the patient and drapes when not in use.

Cleaning the Instrument During Use

Wipe the electrode with a wet gauze pad as needed.


Notice

Do not activate the instrument while cleaning the electrode. Injury to operating room personnel may result.

After Surgery

Cover the instrument with a damp towel to prevent drying. Clean as soon as possible. Do not exceed twelve hours between surgery and cleaning.

Follow the processing instructions described above for proper cleaning and sterilization procedures.

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© 2011, 2021 Covidien.
COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, TM* brands are trademarks of their respective owner.
Medtronic B.V.,