

Y-Tubing 152cm with Single Check Valve

REF: 801106

Part Number: 223374













For use with MRI Injectors such as **OptiStar**® **LE / OptiStar**® **Elite**

Description: 152 cm Y-tubing with single check valve - 200 psi

Components / Materials

Luer Type

1 male luer + cap 2 female luer + 2 caps

Materials

Clear Polyvinyl Chloride (PVC), Polycarbonate (PC), Polypropylene (PP), High-Density Polyethylene (HDPE)

Packaging / Sterilization

Quantity per case - 100

Lid - Tyvek® Pouch

Sterilant - EtO

Shipper dimensions - 67.6 x 29.2 x 25.4cm

OptiStar® is a trademark of Guerbet or its affiliates. DuPont™ and Tyvek® are trademarks or registered trademarks of DuPont or its affiliates.

Edition date: May 2020 - P20000160

Exclusively distributed by:

Guerbet/Liebel-Flarsheim Company LLC
Made in Costa Rica.

Liebel-Flarsheim Company LLC.,

ECREP Guerbet, France

Regulation may differ from country to country. For more specific legal information please refer to your local Guerbet representative

Non contractual information.

Presentations may differ from country to country. Some features might be paid option.

For complete information please refer to the local Summary of Product Characteristics and instructions for use. Please refer to legal mentions available on our website: consumables.guerbet.com



