

Liebel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

P/N: 801106
Lengths: 150 cm
Injectors: Liebel-Flarsheim Optistar MR



Contents: Y-connector with
150 cm patient line and single check valve

Package: 100 units/box

Pressure Limit Syringe : 200 PSI / 14 bar

Materials: Polycarbonate
Polyvinyl Chloride
Silicone
Latex Free

Sterilisation: Ethylene Oxide (EO)

Shelf Life: 3 years from Date of Sterilisation

Use: Single Use – do not reuse

CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237

Country of Origin: USA



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

