

# Liebel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

---

**P/N:** FLEXxx  
**Lengths:** 50 cm / 75 cm / 120 cm  
**Description:** Flexible High Pressure Line  
for Angiography



---

Connectors: LL f/m  
Package: 10 units/box  
Pressure Limit Syringe : 1200 PSI / 84 bar

Materials: Polyurethane  
ABS  
**Latex Free**

Sterilisation: Ethylene Oxide (EO)  
Shelf Life: 5 years from Date of Sterilisation  
Use: Single Use – do not reuse  
CE: 0123

Manufacturer : Medtron AG,  
Hauptstraße 255, 66128 Saarbrücken  
Country of Origin: Germany

---



Tyco Healthcare Europe  
Liebel-Flarsheim Products  
Josef-Dietzgen-Str. 1-3  
53761 Hennef  
Germany

