

Masimo Product Catalog

2019 – Q3



This catalog contains Masimo proprietary material and is provided for information purposes only. No part of this catalog may be reproduced or retransmitted without Masimo's prior written consent.

The purpose of the product catalog is to provide a general listing of Masimo products, and it is not intended to represent products that are offered for sale. Please see local Masimo representatives for the commercially available products in specific regions/countries with the appropriate regulatory approvals/clearances.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Regulatory Notice

The following products, features and/or functions are not available for sale in the United States:

ATA (Adaptive Threshold Alarm); Eve™ (Newborn CCHD Screening Application); Halo Index; iSpO₂ Rx; ORi (Oxygen Reserved Index); Pronto with Next Generation SpHb; RD SedLine Pediatric Sensor; rainbow Super DCI-mini; RD rainbow Lite SET-1 sensors; and continuous RRp (Respiration Rate from Pleth).

The following products are not available for sale in Canada:

Doctella, EMS Cables; iSpO₂; ORi (Oxygen Reserved Index); Pronto with Next Generation SpHb; Rad-67; Rad-97 with NIBP; rainbow Super DCI-mini; RD rainbow Lite SET-1 sensors; and TIR-1.

Masimo Corporation
52 Discovery
Irvine, CA 92618, USA
Tel.: 949-297-7000
Fax.: 949-297-7001
www.masimo.com




Authorized representative in the European Community:



MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany

Patents: www.masimo.com/patents.htm

Masimo,  3D Alarm, Blue, Cleanshield, DCI, Doctella, EMMA, E1, FastSat, Iris, iSpO₂, Kite, LNCS, LNOP, LNOPv, M-LNCS, MOC-9, MX-1, NomoLine, O3, Phasein, Pronto, Pronto-7, Pulse CO-Oximeter, PVi, Rad-5, Rad-5v, Rad-57, Rad-67, Rad-8, Radical, Radical-7, Radius-7, rainbow, rainbow Acoustic Monitoring, rainbow 4D, RAM, RDS, RRa, RRp, Root, SatShare, SedLine, SET, SpCO, SpHb, SpMet, TFA-1, TF-I, and uSpO₂ are registered trademarks of Masimo Corporation.

Eve, MightySat, Ori, Oxygen Reserve Index, RAS-45, Rad-97, rainbow SET, RD SedLine, RD SET, RD rainbow SET, Sketch, Trace, and UniView are trademarks of Masimo Corporation.

Apple, iPad, iPhone, iPod, iPod touch, and Lightning are trademarks of Apple Inc.

"Made for iPod," "Made for iPhone," and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, or iPad may affect wireless performance.

The use of the trademark PATIENT SAFETYNET is under license from University HealthSystem Consortium.

All other trademarks and registered trademarks are property of their respective owners.

The devices provided in the catalog may not have been licensed in accordance with applicable Federal Laws.

© 2019 by Masimo Corporation.

CONFIDENTIAL AND PROPRIETARY – INTERNAL USE ONLY

Contents

- Monitoring Devices.....5
 - Masimo rainbow® Pulse CO-Oximeter® Devices.....5
 - SedLine® Brain Function Monitoring.....9
 - O3® Regional Oximeter Monitoring.....11
 - NomoLine® Capnography13
 - Noninvasive Blood Pressure (NIBP)22
 - Temperature28
 - Masimo SET® Pulse Oximeter Devices.....29
 - Masimo Patient SafetyNet™ System.....32
 - Replica™38
 - Doctella™40
 - Kite® and UniView™ Supplemental Display40
 - Philips® Compatible Devices42
 - GE® Healthcare Compatible Devices44
 - MightySat™ Rx.....45
 - iSpO2® Rx.....47
- Parameter Options and Upgrade Tools51
 - Masimo rainbow SET™ Parameter Options.....51
- Sensors and Cables.....55
 - RD rainbow SET™ Adhesive Sensors (for rainbow SET Devices with SpHb®, SpO2, and SpMet®)55
 - RD rainbow Lite SET™ Adhesive Sensors (for rainbow SET Devices with SpO2, and ORi).....56
 - rainbow® Adhesive Sensors R1 series (for rainbow SET Devices with SpHb, SpO2, and SpMet)57
 - rainbow® Adhesive Sensors (for rainbow SET Devices with SpO2, SpCO®, and SpMet)58
 - rainbow® Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO2, SpMet; SpHb required)59
 - rainbow® Reusable Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet)61
 - rainbow® Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet)62
 - rainbow® Patient Cables (for rainbow SET Devices)63
 - rainbow Acoustic Monitoring® Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRa®).....71
 - RD SET™ Disposable Sensors74
 - RD SET™ Cabled Sensors.....75
 - RD SET Specialty Sensors78
 - RD SET Reusable Sensors79
 - RD SET Patient Cables81
 - RD SET Adapter Cables85
 - LNCS® Adhesive Sensors86
 - LNCS Single Patient Use Sensors with Replaceable Tapes.....89
 - LNCS Specialty Sensors91
 - LNCS Sterile Sensors92
 - LNCS Reusable Sensors94
 - LNCS Red Patient Cables (for use with Masimo rainbow SET Devices)95
 - LNC Patient Cables (for use with Masimo SET® Devices)97
 - LNCS Adapter Cables.....102

M-LNC Patient Cables..... 103

M-LNCS Adapter Cables..... 107

Direct Connect Sensors (for use with Masimo SET® Devices) 109

Red Direct Connect Sensors (for Masimo rainbow SET Devices)..... 109

Mountable Patient Cable Extensions and Adapter..... 111

Accessories 113

 Training Kits 113

 Demo Tools 115

 Software..... 116

 Bedside Device Accessories 116

 Handheld Device Accessories 121

 Sensor Accessories..... 124




 SatShare® Interface Cables 128





 SatShare Extension Cable..... 131




Index 133




Monitoring Devices

Masimo rainbow® Pulse CO-Oximeter® Devices






PN	Description	UM	Product Image
9500	<p>Radical-7® Handheld</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (Pi) parameters as well as one RD rainbow SET, Red LNCS® or rainbow Patient Cable. Can be configured with optional rainbow parameters and features in cleared countries only. Check with customer service on availability:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb®) • Pleth Variability Index (PVi®) • Oxygen Reserve Index™ (ORi™) • Carboxyhemoglobin (SpCO®) • Methemoglobin (SpMet®) • Oxygen Content (SpOC™) • Acoustic Respiration (RRa®) • Respiration Rate from Pleth (RRp™)* • Adaptive Threshold Alarms (ATA) • Eve™ (Newborn CCHD Screening Application)* <p>This unit ships standard with a 802.11 wireless radio. *Available in countries where cleared.</p>	EA	
9515	<p>Root®</p> <p>For use with Radical-7, Radius-7® and MOC-9™ modules.</p> <ul style="list-style-type: none"> • Screen Size: 10.1 in (25.65 cm) • Resolution: 1280 x 800 Pixels • Color: 24-bit RGB • Display Type: Backlit Active Matrix TFT LCD • Touchscreen: Multi-Touch P-Cap <p>Provides Ethernet connection, nurse call connector, USB 2.0 Ports, MOC-9 ports, and Iris® Connectivity Ports. Comes with choice of power cords.</p>	EA	
9695	<p>Root</p> <p>With Noninvasive Blood Pressure and Temperature</p> <p>For use with Radical-7, Radius-7 and MOC-9 modules.</p> <ul style="list-style-type: none"> • Screen Size: 10.1 in (25.65 cm) • Resolution: 1280 x 800 Pixels • Color: 24-bit RGB • Display Type: Backlit Active Matrix TFT LCD • Touchscreen: Multi-Touch P-Cap • Integrated Blood Pressure Module • Integrated Temperature Module <p>Provides Ethernet connection, nurse call connector, USB 2.0 Ports, MOC-9 ports, and Iris Connectivity Ports. Comes with choice of power cords.</p>	EA	

PN	Description	UM	Product Image
9021	<p>Radical® Docking Station, RDS®-1</p> <p>Comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare cable connector.</p>	EA	
9023	<p>Radical Docking Station, RDS-3</p> <p>limited communication capabilities.</p> <p>Comes with choice of power cord and serial output connector and Analog output / nurse call connector.</p>	EA	
9738	<p>Kit, Rad-97™ Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi) and one M20 patient cable. Can be configured with optional rainbow parameters and features in cleared countries only. Check with customer service on availability:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Oxygen Reserve Index (ORi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) • Acoustic Respiration (RRa) • Respiration Rate from Pleth (RRp) • Eve™ (Newborn CCHD Screening Application) 	EA	
9797	<p>Kit, Rad-97 Pulse CO-Oximeter with Noninvasive Blood Pressure</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi), noninvasive blood pressure (NIBP) and one M20 patient cable. Can be configured with optional rainbow parameters and features in cleared countries only. Check with customer service on availability:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Oxygen Reserve Index (ORi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) • Acoustic Respiration (RRa) • Respiration Rate from Pleth (RRp) <p><i>Cuff must be ordered separately.</i></p>	EA	

PN	Description	UM	Product Image
9799	<p>Kit, Rad-97™ Pulse CO-Oximeter with NomoLine® Capnography</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi), integrated NomoLine Capnography module and one M20 patient cable. Can be configured with optional rainbow parameters and features in cleared countries only. Check with customer service on availability:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Oxygen Reserve Index (ORi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) • Acoustic Respiration (RRa) • Respiration Rate from Pleth (RRp) <p><i>NomoLine Sampling Lines must be ordered separately.</i></p>	EA	
9664	<p>Kit, Radius-7</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi) parameters, and Acoustic Respiration (RRa). Can be configured with optional rainbow parameters and features in cleared countries only. Check with customer service on availability:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) <p>Bluetooth radio allows for transfer of parameter data to the Root patient monitoring and connectivity platform.</p>	EA	
9794	<p>Kit, Rad-67™ Spot Check Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi), Total Hemoglobin (SpHb) and one rainbow M20 patient cable. Can be configured with optional rainbow parameters in cleared countries only. Check with customer service on availability:</p> <ul style="list-style-type: none"> • Pleth Variability Index (PVi)* • Carboxyhemoglobin (SpCO)* • Methemoglobin (SpMet)* • Oxygen Content (SpOC)* <p>rainbow Sensor must be ordered separately.</p> <p>Next Generation SpHb technology is only available with the DCI®-mini sensor.</p> <p>Shown with optional rainbow Super DCI-mini sensor.</p> <p><i>*Available in countries where cleared.</i></p>	EA	




PN	Description	UM	Product Image
9167	<p>Pronto® Spot Check Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi) and SpHb parameters.</p> <p>rainbow Sensor not included.</p> <p>rainbow Sensor must be ordered separately.</p>	EA	
9757	<p>Pronto Spot Check Pulse CO-Oximeter, with Next Generation SpHb</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi) and SpHb parameters.</p> <p>rainbow Sensor not included.</p> <p>rainbow Sensor must be ordered separately.</p> <p>Next Generation SpHb technology is only available with the DCI-mini sensor.</p> <p><i>Available in countries where cleared.</i></p>	EA	
9216-U 9216	<p>Rad-57® Configurable Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (Pi) parameters.</p> <p>The Rad-57 can ship with or be field upgraded with:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) <p>When shipped with rainbow parameters, rainbow Sensor must be purchased separately.</p> <p>rainbow Sensors not included.</p> <p>9216-U is for SpO₂ only Devices.</p> <p>Refer to Handheld Device accessories for available colors.</p> <p>Shown with rainbow Adult Reusable Sensor and in optional protective boot.</p>	EA	


SedLine® Brain Function Monitoring

PN	Description	UM	Product Image
9513	<p>SedLine Kit For use with Root and RD SedLine Sensors. Includes: 1 SedLine MOC-9 Module 1 RD SedLine Patient Cable</p>	EA	
9829	<p>SedLine Module only Includes: 1 SedLine MOC-9 Module</p>	EA	
4248	<p>RD SedLine™ Sensor For use with SedLine MOC-9 Module and RD SedLine Patient Cable. 25/ box, non-sterile Not made with natural rubber latex.</p>	Box	
4329	<p>RD SedLine™ Pediatric Sensor For use with SedLine MOC-9 and RD SedLine Patient Cable 25/box, non-sterile Not made with natural rubber latex Available in countries where cleared.</p>	Box	
4298	<p>RD SedLine Patient Cable For use with RD SedLine Sensor and SedLine MOC-9 Module.</p>	EA	

PN	Description	UM	Product Image
4316	<p>RD SedLine Training Kit For use with Root and SedLine MOC-9 Module</p> <p>5 RD SedLine Sensors 1 RD SedLine Patient Cable 1 DFU</p>	EA	
4549	<p>RD SedLine Pediatric Training Kit For use with Root and SedLine MOC-9 Module</p> <p>5 RD SedLine Pediatric Sensors 1 RD SedLine Patient Cable 1 DFU Available in countries where cleared.</p>	EA	
3637	<p>Replacement SedLine Patient Cable For use with SedLine MOC-9 and SedLine Sensors.</p>	EA	
3443	<p>Masimo Cable Tether For use with Masimo sensors and patient cables 20/ pack, non-sterile Not made with natural rubber latex</p>	Pack	



O3® Regional Oximeter Monitoring

PN	Description	UM	Product Image
4273	<p>O3 Adult Training Kit 6 O3 Adult Sensors 1 DFU</p>	EA	
4315	<p>O3 Pediatric Training Kit 6 O₃ Pediatric Sensors 1 DFU</p>	EA	
4616	<p>O3 Infant/Neonatal Training Kit 6 O₃ Infant/Neonatal Sensors 1 DFU</p>	EA	
9637	<p>O3 Regional Oximeter Module For use with O3 sensors.</p>	EA	





PN	Description	UM	Product Image
3756	<p>O3 Adult Sensor Single Patient Use 20/Box Weight ≥ 40 kg</p>	Box	
4235	<p>O3 Pediatric Sensor Single Patient Use 20/Box Weight: CE Weight < 40 kg FDA Weight ≥ 5 kg and < 40 kg</p>	Box	
4384	<p>O3 Infant and Neonatal Sensor Single Patient Use 20/Box Weight: FDA Weight < 10 kg Available in countries where cleared</p>	Box	
4206	<p>O3 Regional Oximeter Module Demo tools (2 /Packet) For use with an O3 Regional Oximeter Module (not included).</p>	EA	

NomoLine® Capnography





ISA™ CO2





PN	Description	UM	Product Image
9921	<p>NomoLine® ISA CO₂ For use with Root and NomoLine Sampling Lines. Includes:</p> <ul style="list-style-type: none"> • 1 NomoLine ISA CO₂ MOC-9 Module • 1 ISA CO₂ Bracket 	EA	
9685	<p>ISA OR+ Kit For use with Root and NomoLine Sampling Lines. Includes:</p> <ul style="list-style-type: none"> • 1 ISA OR+ MOC-9 Module • 1 ISA OR+ Bracket 	EA	





EMMA™





PN	Description	UM	Product Image
3639	<p>EMMA (mmHg) Capnograph Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 EMMA Capnograph mmHg • 1 EMMA Pouch • 1 EMMA Lanyard • 2 AAA Batteries 	EA	
3678	<p>EMMA (kPa) Capnograph Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 EMMA Capnograph kPa • 1 EMMA Pouch • 1 EMMA Lanyard • 2 AAA Batteries 	EA	
17448	<p>EMMA Adult / Pediatric Airway Adapter</p> <p>For use with EMMA Capnometer and EMMA Capnograph. Adult / Pediatric airway adapter 25 / Box, Non-sterile</p>	Box	
17449	<p>EMMA Infant Airway Adapter</p> <p>For use with EMMA Capnometer and EMMA Capnograph. Infant airway adapter 10 / Box, Non-sterile</p>	Box	





NomoLine Sampling Lines





PN	Description	UM	Product Image
3814	<p>NomoLine LH Adult / Pediatric Airway Adapter Set, 2m Adult/Pediatric Airway Adapter Set For short-term applications in low-humidity conditions For intubated patients 25 / Box</p>	Box	
3815	<p>NomoLine LH Adult / Pediatric Airway Adapter Set, 3m Adult/Pediatric Airway Adapter Set For short-term applications in low-humidity conditions For intubated patients 25 / Box</p>	Box	
3816	<p>NomoLine LH Infant Airway Adapter Set, 2m Infant Airway Adapter Set For short-term applications in low-humidity conditions For intubated patients 25 / Box</p>	Box	
4369	<p>NomoLine LH Infant / Neonatal Airway Adapter Set, 2m Infant/Neonatal Airway Adapter Set For short-term applications in low-humidity conditions For intubated patients 25 / Box</p>	Box	





PN	Description	UM	Product Image
3827	<p>NomoLine HH Adult / Pediatric Airway Adapter Set, 2m Adult/Pediatric Airway Adapter Set For long-term applications in high-humidity conditions For intubated patients 25 / Box</p>	Box	
3828	<p>NomoLine HH Adult / Pediatric Airway Adapter Set, 3m Adult/Pediatric Airway Adapter Set For long-term applications in high-humidity conditions For intubated patients 25 / Box</p>	Box	
3829	<p>NomoLine HH Infant Airway Adapter Set, 2m Infant Airway Adapter Set For long-term applications in high-humidity conditions For intubated patients 25 / Box</p>	Box	
4367	<p>NomoLine HH Infant / Neonatal Airway Adapter Set, 2m Infant/Neonatal Airway Adapter Set For long-term applications in high-humidity conditions For intubated patients 25 / Box</p>	Box	

PN	Description	UM	Product Image
3817	<p>NomoLine LH Adult Nasal CO₂ Cannula Adult Nasal Cannula For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3818	<p>NomoLine LH Pediatric Nasal CO₂ Cannula Pediatric Nasal Cannula For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3819	<p>NomoLine LH Infant Nasal CO₂ Cannula Infant Nasal Cannula For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3822	<p>NomoLine LH Adult Nasal / Oral CO₂ Cannula Adult Nasal/Oral Cannula For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	

PN	Description	UM	Product Image
3823	<p>NomoLine LH Pediatric Nasal / Oral CO₂ Cannula Pediatric Nasal/Oral Cannula For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3826	<p>NomoLine LH Adult Single Nasal Prong CO₂ Cannula Adult Single Nasal Prong For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3830	<p>NomoLine HH Adult Nasal CO₂ Cannula Adult Nasal Cannula For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3831	<p>NomoLine HH Pediatric Nasal CO₂ Cannula Pediatric Nasal Cannula For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	




PN	Description	UM	Product Image
3832	<p>NomoLine HH Infant Nasal CO₂ Cannula Infant Nasal Cannula For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3835	<p>NomoLine HH Adult Nasal / Oral CO₂ Cannula Adult Nasal/Oral Cannula For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3836	<p>NomoLine HH Pediatric Nasal / Oral CO₂ Cannula Pediatric Nasal/Oral Cannula For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3839	<p>NomoLine HH Adult Single Nasal Prong CO₂ Cannula Adult Single Nasal Prong For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	





PN	Description	UM	Product Image
3820	<p>NomoLine LH Adult Nasal CO₂ Cannula with O₂ Adult Nasal Cannula with O₂ Delivery For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3821	<p>NomoLine LH Pediatric Nasal CO₂ Cannula with O₂ Pediatric Nasal Cannula with O₂ Delivery For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3824	<p>NomoLine LH Adult Nasal/Oral CO₂ Cannula with O₂ Adult Nasal/Oral Cannula with O₂ Delivery For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3825	<p>NomoLine LH Pediatric Nasal/Oral CO₂ Cannula with O₂ Pediatric Nasal/Oral Cannula with O₂ Delivery For short-term applications in low-humidity conditions For non-intubated patients 25 / Box</p>	Box	

PN	Description	UM	Product Image
3833	<p>NomoLine HH Adult Nasal CO₂ Cannula with O₂ Adult Nasal Cannula with O₂ Delivery For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3834	<p>NomoLine HH Pediatric Nasal CO₂ Cannula with O₂ Pediatric Nasal Cannula with O₂ Delivery For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3837	<p>NomoLine HH Adult Nasal/Oral CO₂ Cannula with O₂ Adult Nasal/Oral Cannula with O₂ Delivery For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	
3838	<p>NomoLine HH Pediatric Nasal/Oral CO₂ Cannula with O₂ Pediatric Nasal/Oral Cannula with O₂ Delivery For long-term applications in high-humidity conditions For non-intubated patients 25 / Box</p>	Box	

Noninvasive Blood Pressure (NIBP)





Disposable Blood Pressure Cuffs



PN	Description	UM	Product Image
25242	<p>Disposable Neonate Cuff Size 1 Single Patient Use BP cuff with range of 3-6 cm Male Bayonet connection, for use with BP Hose with Female Bayonet connection 20 / box</p>	Box	 <p>A white disposable neonate blood pressure cuff with a range of 3-6 cm. It features a male bayonet connection on the left side and a female bayonet connection on the right side. The cuff is connected to a clear plastic tube. The cuff has a white label with orange and black text and graphics.</p>
25243	<p>Disposable Neonate Cuff Size 2 Single Patient Use BP cuff with range of 4-8 cm Male Bayonet connection, for use with BP Hose with Female Bayonet connection 20 / box</p>	Box	 <p>A white disposable neonate blood pressure cuff with a range of 4-8 cm. It features a male bayonet connection on the left side and a female bayonet connection on the right side. The cuff is connected to a clear plastic tube. The cuff has a white label with blue and black text and graphics.</p>
25244	<p>Disposable Neonate Cuff Size 3 Single Patient Use BP cuff with range of 6-11 cm Male Bayonet connection, for use with BP Hose with Female Bayonet connection 20 / box</p>	Box	 <p>A white disposable neonate blood pressure cuff with a range of 6-11 cm. It features a male bayonet connection on the left side and a female bayonet connection on the right side. The cuff is connected to a clear plastic tube. The cuff has a white label with green and black text and graphics.</p>

PN	Description	UM	Product Image
25245	<p>Disposable Neonate Cuff Size 4 Single Patient Use BP cuff with range of 7-13 cm Male Bayonet connection, for use with BP Hose with Female Bayonet connection 20 / box</p>	Box	
25246	<p>Disposable Neonate Cuff Size 5 Single Patient Use BP cuff with range of 8-15 cm Male Bayonet connection, for use with BP Hose with Female Bayonet connection 20 / box</p>	Box	
25296	<p>Disposable Infant Cuff Single Patient Use BP cuff with range of 8-13 cm Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 20 / box</p>	Box	
25297	<p>Disposable Child Cuff Single Patient Use BP cuff with range of 12-19 cm Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 20 / box</p>	Box	




PN	Description	UM	Product Image
25298	<p>Disposable Small Adult Cuff Single Patient Use BP cuff with range of 17-25 cm Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 20 / box</p>	Box	 <p>The image shows a white rectangular package for a SunTech Disposable Small Adult Cuff. The package features the SunTech logo at the top, a blue horizontal band with the text 'Small Adult' and '17-25 cm', and a graphic of three stylized human figures. A white plastic cuff with a female quick-connect port is attached to the bottom of the package.</p>
25299	<p>Disposable Adult Cuff Single Patient Use BP cuff with range of 23-33 cm Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 20 / box</p>	Box	 <p>The image shows a white rectangular package for a SunTech Disposable Adult Cuff. The package features the SunTech logo at the top, a blue horizontal band with the text 'Adult' and '23-33 cm', and a graphic of three stylized human figures. A white plastic cuff with a female quick-connect port is attached to the bottom of the package.</p>
25300	<p>Disposable Large Adult Cuff Single Patient Use BP cuff with range of 31-40 cm Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 20 / box</p>	Box	 <p>The image shows a white rectangular package for a SunTech Disposable Large Adult Cuff. The package features the SunTech logo at the top, a red horizontal band with the text 'Large Adult' and '31-40 cm', and a graphic of three stylized human figures. A white plastic cuff with a female quick-connect port is attached to the bottom of the package.</p>
25301	<p>Disposable Thigh Cuff Single Patient Use BP cuff with range of 38-50 cm Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 20 / box</p>	Box	 <p>The image shows a white rectangular package for a SunTech Disposable Thigh Cuff. The package features the SunTech logo at the top, a black horizontal band with the text 'Thigh' and '38-50 cm', and a graphic of three stylized human figures. A white plastic cuff with a female quick-connect port is attached to the bottom of the package.</p>

Reusable Blood Pressure Cuffs




PN	Description	UM	Product Image
25303	<p>Reusable Infant Cuff Reusable BP cuff with range of 8-13 cm. Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 5 / box</p>	Box	 An orange reusable blood pressure cuff with a black hose and a white connector. The cuff is shown in its packaging.
25260	<p>Reusable Child Cuff, Masimo Reusable BP cuff with range of 12-19 cm. Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 5 / box</p>	Box	 A green reusable blood pressure cuff with a black hose and a white connector. The cuff is shown in its packaging with the SunTech logo and 'One-Piece Cuff' text.
25261	<p>Reusable Small Adult Cuff, Masimo Reusable BP cuff with range of 17-25 cm. Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 5 / box</p>	Box	 A blue reusable blood pressure cuff with a black hose and a white connector. The cuff is shown in its packaging with the SunTech logo and 'Small Adult' text.
25262	<p>Reusable Adult Cuff, Masimo Reusable BP cuff with range of 23-33 cm. Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 5 / box</p>	Box	 A dark grey reusable blood pressure cuff with a black hose and a white connector. The cuff is shown in its packaging with the SunTech logo and 'Adult' text.


PN	Description	UM	Product Image
25263	<p>Reusable Large Adult Cuff, Masimo Reusable BP cuff with range of 31-40 cm. Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 5 / box</p>	Box	
25304	<p>Reusable Thigh Cuff Reusable BP cuff with range of 38-50 cm. Female Quick Connect connection, for use with BP Hose with Male Quick Connect connection 5 / box</p>	Box	

Blood Pressure Hoses and Adapters



PN	Description	UM	Product Image
25302	<p>NIBP Patient Hose, 3 Meters, Female Bayonet Metal Female Bayonet connection For use with BP cuffs with Male Bayonet connection.</p>	EA	
25247	<p>NIBP Patient Hose, 3 Meters, Male Quick Connect Plastic Male Quick Connect connection For use with BP cuffs with Female Quick Connect connection.</p>	EA	
26426	<p>NIBP Patient Hose Adapter, 10 cm, with Male Bayonet and Female Slip-Luer Male Slip-Luer BP cuffs to Female Bayonet NIBP Patient Hose</p>	EA	




Temperature




PN	Description	UM	Product Image
4301	<p>TIR-1 Thermometer Packaging includes:</p> <ul style="list-style-type: none"> • TIR-1 Unit • 2 AA Batteries • TIR-1 Operator's Manual 	EA	
25251	<p>Oral Temperature Probe and Well Kit, Welch Allyn Includes:</p> <ul style="list-style-type: none"> • Oral/axillary temperature probe • Oral/axillary temperature probe well 	EA	
25252	<p>Rectal Temperature Probe and Well Kit, Welch Allyn Includes:</p> <ul style="list-style-type: none"> • Rectal temperature probe • Rectal temperature probe well 	EA	

PN	Description	UM	Product Image
26289	<p>Temperature Probe Covers, Welch Allyn 7500 / case Case contains 30 dispensers, each with 10 packs of 25 probe covers</p>	EA	







Masimo SET® Pulse Oximeter Devices







PN	Description	UM	Product Image
9190	<p>Rad-8®, Horizontal Comes with either an RD SET or LNCS Patient Cable and a power cord.</p>	EA	
9191	<p>Rad-8, Vertical Comes with either an RD SET or LNCS Patient Cable and a power cord.</p>	EA	



PN	Description	UM	Product Image
9196	<p>Rad-5® Pulse Oximeter</p> <p>Comes with a Red DCI®-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of operator's manual.</p> <p>Shown in optional handheld boot and with DCSC Sensor.</p>	EA	
9198	<p>Rad-5 Pulse Oximeter with Sample Pack</p> <p>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of operator's manual.</p> <p>Shown in optional handheld boot and with DCSC Sensor.</p>	EA	
9214	<p>Rad-5 Pulse Oximeter (International Only)</p> <p>Comes with a choice of sensor, a choice of available colors of Handheld Boots and a choice of operator's manual.</p> <p>Shown in optional handheld boot and with DCSC Sensor.</p>	EA	






PN	Description	UM	Product Image
9197	<p>Rad-5v® Pulse Oximeter</p> <p>Comes with a choice of DCI sensor and choice of patient cable, a choice of available colors of Handheld Boots and a choice of operator's manual.</p> <p>Shown in optional handheld boot and with DCSC Sensor.</p>	EA	
9199	<p>Rad-5v Pulse Oximeter with Sample Pack</p> <p>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of operator's manual.</p> <p>Shown in optional handheld boot and with DCSC Sensor.</p>	EA	
9215	<p>Rad-5v Pulse Oximeter (International Only)</p> <p>Comes with a choice of sensor, a choice of available colors of Handheld Boots and a choice of operator's manual.</p> <p>Shown in optional handheld boot and with DCSC Sensor.</p>	EA	

Masimo Patient SafetyNet™ System









PN	Description	UM	Product Image
9095	Patient SafetyNet System	System	
25405	Patient SafetyNet View Station with Software For use with Patient SafetyNet	EA	
25134	Kit Patient SafetyNet Appliance 1UM For use with Patient SafetyNet	EA	
25151	Kit Patient SafetyNet Appliance 2UM For use with Patient SafetyNet	EA	
25133	Kit Patient SafetyNet Appliance 1UM and View For use with Patient SafetyNet	EA	
25150	Kit Patient SafetyNet Appliance 2UM and View For use with Patient SafetyNet	EA	

PN	Description	UM	Product Image
25353	Halo™ Appliance For use with Patient SafetyNet <i>Available in countries where cleared.</i>	EA	
9712	Iris Gateway™ Appliance	EA	
9920	AI Aql Mini-PC for Iris Gateway	EA	
9937	Iris Device Management System	EA	
2445	Patient SafetyNet Bed License For use with Patient SafetyNet	EA	No Image
2125	Halo Bed License For use with Patient SafetyNet <i>Available in countries where cleared.</i>	EA	No Image
2447	Notification License For use with Patient SafetyNet <i>Note: For non-proprietary notification interface</i>	EA	No Image
2146	Patient SafetyNet Compatible Radical Docking Station, RDS-1 For use with Patient SafetyNet Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	
2147	Patient SafetyNet Compatible Radical Docking Station, RDS-3 For use with Patient SafetyNet Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet® Compatible Radical Docking Stations.</i>	EA	

PN	Description	UM	Product Image
3963	HL7, ADT, Interface For use with Patient SafetyNet	EA	No Image
3964	HL7, Results Out, Interface For use with Patient SafetyNet	EA	No Image
15007	Colubris Workgroup Bridge For use with Patient SafetyNet One required per pulse oximeter; alternative for 3Com Workgroup Bridge	EA	
15275	NET232-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter For use with Patient SafetyNet	EA	
15274	NET232-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter For use with Patient SafetyNet	EA	No image
16190	NET232-POE-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter For use with Patient SafetyNet	EA	No image
16191	NET232-POE-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter For use with Patient SafetyNet	EA	No image
15719	NET232-U-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter For use with Patient SafetyNet	EA	No image
15720	NET232-U-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter For use with Patient SafetyNet	EA	No image
23416	Kit HP ProCurve Serial (Male-Male) For use with Patient SafetyNet	EA	No image
23415	Kit HP ProCurve Serial (Male-Female)	EA	No image




PN	Description	UM	Product Image
17100	<p>Cables to GO Patch Cable 3 ft. Blue For use with Patient SafetyNet</p>	EA	
15768	<p>Cables to GO Patch Cable 7ft. Blue For use with Patient SafetyNet</p>	EA	
15769	<p>Cables to GO Patch Cable 14 ft. Blue For use with Patient SafetyNet</p>	EA	
15770	<p>Cables to GO Patch Cable 25 ft. Blue For use with Patient SafetyNet</p>	EA	
17644	<p>HP V1910-24G-POE Port Switch For use with Patient SafetyNet</p>	EA	
14530	<p>Null Modem Cable, DB-9 Male to DB-9 Male, 1" Length For use with Patient SafetyNet</p>	EA	No Image

PN	Description	UM	Product Image
14857	XRackPro2 Rackmount 4U Cabinet, Noise Reduction Enclosure For use with Patient SafetyNet	EA	
14424	Paging Transmitter UHF For use with Patient SafetyNet Frequency range (450 to 470 MHz). One required per floor in hospital.	EA	
16365	Paging Transmitter VHF For use with Patient SafetyNet Frequency range (135 to 175 MHz). One required per floor in hospital.	EA	
16262	Pager 440- 444 MHz, Synthesized, POCSAG For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 440 to 444 MHz range.	EA	
16263	Pager 444 - 448 MHz, Synthesized, POCSAG For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 444 to 448 MHz range.	EA	
16158	Pager 448 - 452 MHz, Synthesized, POCSAG For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 448 to 452 MHz range.	EA	
16264	Pager 452 - 456 MHz, Synthesized, POCSAG For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 452 to 456 MHz range.	EA	

PN	Description	UM	Product Image
14253	<p>Pager 456 - 460 MHz, Synthesized, POCSAG</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 456 to 460 MHz range.</p>	EA	
14254	<p>Pager 460 - 464 MHz, Synthesized, POCSAG</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 460 to 464 MHz range.</p>	EA	
14255	<p>Pager 464 - 468 MHz, Synthesized, POCSAG</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 464 to 468 MHz range.</p>	EA	
16527	<p>Pager 135 – 142.999 MHz, Synthesized, POCSAG VHF Band A</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 135 to 142.999 MHz range.</p>	EA	
16528	<p>Pager 143 – 150.999 MHz, Synthesized, POCSAG VHF Band B</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 143 to 150.999 MHz range.</p>	EA	
16366	<p>Pager 151 – 158.999 MHz, Synthesized, POCSAG VHF Band C</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 151 to 158.999 MHz range.</p>	EA	
16529	<p>Pager 159 – 166.999 MHz, Synthesized, POCSAG VHF Band D</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 159 to 166.999 MHz range.</p>	EA	
16530	<p>Pager 167 – 174.999 MHz, Synthesized, POCSAG VHF Band E</p> <p>For use with Patient SafetyNet. Requires Paging Transmitter. Pager tunable in the 167 to 174.999 MHz range</p>	EA	


PN	Description	UM	Product Image
2456	Patient SafetyNet Installation Services For use with Patient SafetyNet.	EA	No image
2457	Masimo Wireless Site Survey For use with Patient SafetyNet.	EA	No image

Replica™


PN	Description	UM	Product Image
9863	Replica System	1	
26525	Replica application installed on Zebra TC51-HC phone, US	1	
26526	Replica application installed on Zebra TC51-HC phone, OUS	1	
4506	Replica application to be installed on customer's existing phone	1	No Image
19989	Zebra TC51-HC Power Supply Attaches to Zebra Charge Cable (19941) For OUS installations, requires Universal Power Supply (19755)	1	No Image
19990	Zebra TC51-HC Power Cord Attaches to Zebra Multi-Slot ShareCradles For OUS installations, requires Universal Power Supply (19755)	1	No Image






PN	Description	UM	Product Image
19941	Zebra TC-51 Charge Cable Charges one TC-51 phone Requires PN 19989 Power Supply	1	
19942	Zebra TC-51 HC Replacement Battery	1	
19943	Zebra Multi-Slot ShareCradle Requires Zebra TC51-HC Power Cord Charges 4 Zebra TC51-HC phones and 4 replacement batteries	1	
100188	Zebra Carrying Clip	1	
19755	Universal Power Supply	1	


Doctella™

PN	Description	UM	Product Image
4593	<p>Doctella, Data Sharing Subscription License, Consumer Devices Phone not included</p> <p>Doctella requires the Masimo Implementation Team to visit on-site for configuration of CarePrograms and the provider dashboard.</p>	EA	



Kite® and UniView™ Supplemental Display





PN	Description	UM	Product Image
9766	<p>Kite® Mini-PC Includes Masimo Sketch USB Installer (PN 9899)</p>	EA	
9856	<p>UniView™ Mini-PC</p>	EA	

PN	Description	UM	Product Image
19446	<p>32" Samsung TV (Optional) TV to be connected with the Kite or UniView Mini-PC to display</p>	EA	
19715	<p>Medical Grade Keyboard (Optional) Keyboard to be connected to Kite or UniView Mini-PC to navigate software</p>	EA	
17100	<p>Cables to GO Patch Cable 3 ft. Blue (Optional) For connecting the Kite or UniView Mini-PC or Masimo device to the network</p>	EA	
15768	<p>Cables to GO Patch Cable 7 ft. Blue (Optional) For connecting the Kite or UniView Mini-PC or Masimo device to the network</p>	EA	
15769	<p>Cables to GO Patch Cable 14 ft. Blue (Optional) For connecting the Kite or UniView Mini-PC or Masimo device to the network</p>	EA	

PN	Description	UM	Product Image
15770	<p>Cables to GO Patch Cable 25 ft. Blue (Optional) For connecting the Kite or UniView Mini-PC or Masimo device to the network</p>	EA	



Philips® Compatible Devices

PN	Description	UM	Product Image
3924	<p>Masimo rainbow SET IntelliVue® Module For use with IntelliVue Monitor. <i>Patient Cable must be ordered separately.</i> <i>Available in countries where cleared.</i></p>	EA	
3915	<p>Masimo SET® IntelliVue® Module For use with IntelliVue Monitor. <i>Patient Cable must be ordered separately.</i></p>	EA	




PN	Description	UM	Product Image
4319	<p>uSpO₂ Masimo SET Oximetry Cable, ODU Connector uSpO₂ Patient Cable with 5-pin ODU connector, 6 FT. For use with:</p> <ul style="list-style-type: none"> • Philips Respironics System One • Philips Respironics Alice 6 Lab Diagnostics System • Philips Respironics Devices Using the Broadband Oximetry Module <p>Compatible with Masimo LNCS oximetry sensors</p>	EA	 <p>A white 6-foot patient cable with a 5-pin ODU connector at one end and a Masimo SET connector at the other. An inset image shows a close-up of the ODU connector.</p>
4321	<p>uSpO₂ Masimo SET Oximetry Cable, NightOne Connector uSpO₂ Patient Cable with NightOne connector, 1 FT. For use with Philips Respironics Alice NightOne Device Compatible with Masimo LNCS oximetry sensors</p>	EA	 <p>A white 1-foot patient cable with a NightOne connector at one end and a Masimo SET connector at the other. An inset image shows a close-up of the NightOne connector.</p>
4323	<p>uSpO₂ Masimo SET Oximetry Cable, DreamStation uSpO₂ Patient Cable with DreamStation connector, 6 FT. For use with Philips DreamStation Devices Compatible with Masimo LNCS oximetry sensors</p>	EA	 <p>A white 6-foot patient cable with a DreamStation connector at one end and a Masimo SET connector at the other. An inset image shows a close-up of the DreamStation connector.</p>
18710	<p>USB-to-RS232 Converter Cable For connecting Root with IntelliBridge® (EC5 portion) and VueLink® modules.</p>	EA	 <p>A white USB-to-RS232 converter cable with a blue USB-A connector at one end and a white RS232 connector at the other.</p>




Philips®, IntelliVue® and VueLink® are registered trademarks of Philips Corporation.

GE® Healthcare Compatible Devices

PN	Description	UM	Product Image
9477	<p>uSpO₂ Pulse Oximeter Cable Kit</p> <p>For use with GE ApexPro FH Telemetry Transmitters.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNCS® Adtx-3 • 1 LNCS Pdtx-3 • 1 LNCS Inf-3 • 1 LNCS Neo-3 sensor <p><i>For US Markets Only</i></p> <p>May only be sold to End-Customers with pre-existing installed ApexPro Telemetry Systems. For more information, see GE Upgrade Profiles under OEM on RadNet.</p>	EA	
9476	<p>uSpO₂ Pulse Oximeter Cable Kit</p> <p>For use with GE ApexPro, ApexPro CH, and CARESCAPE T4 & T14 Telemetry Transmitters.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNCS Adtx-3 • 1 LNCS Pdtx-3 • 1 LNCS Inf-3 • 1 LNCS Neo-3 sensor <p><i>Available in countries where cleared.</i></p> <p>May only be sold to End-Customers with pre-existing installed ApexPro Telemetry Systems. For more information, see GE Upgrade Profiles under OEM on RadNet®.</p>	EA	



MightySat™ Rx



PN	Description	UM	Product Image
9707	<p>MightySat Rx Fingertip Pulse Oximeter (CE Markets Only)</p> <p>The MightySat Rx Fingertip Pulse Oximeter is intended as a noninvasive device that measures and displays arterial oxygen saturation (SpO₂), Pulse Rate (PR), and Perfusion Index (Pi).</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 MightySat Rx • 1 Carrying Case • 1 Lanyard • 2 AAA Batteries <p><i>Available in countries where cleared.</i></p>	EA	
9807	<p>MightySat Rx Fingertip Pulse Oximeter w/Bluetooth (CE Markets Only)</p> <p>The MightySat Rx Fingertip Pulse Oximeter is intended as a noninvasive device that measures and displays arterial oxygen saturation (SpO₂), Pulse Rate (PR), and Perfusion Index (Pi). Optional Bluetooth LE wireless technology for the wireless transfer of patient data to smart devices.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 MightySat Rx • 1 Carrying Case • 1 Lanyard • 2 AAA Batteries <p><i>See Masimo.co.uk for a list of compatible devices.</i></p> <p><i>Available in countries where cleared.</i></p>	EA	
9907	<p>MightySat Rx Fingertip Pulse Oximeter w/Bluetooth, PVi, RRp(CE Markets Only)</p> <p>The MightySat Rx Fingertip Pulse Oximeter is intended as a noninvasive device that measures and displays arterial oxygen saturation (SpO₂), Pulse Rate (PR), Perfusion Index (Pi), optional Pleth Variability Index (PVi), and Respiration Rate from the Pleth (RRp). Optional Bluetooth LE wireless technology for the wireless transfer of patient data to smart devices.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 MightySat Rx • 1 Carrying Case • 1 Lanyard • 2 AAA Batteries <p><i>See Masimo.co.uk for a list of compatible devices.</i></p> <p><i>Available in countries where cleared.</i></p>	EA	

PN	Description	UM	Product Image
9709	<p>MightySat Rx Fingertip Pulse Oximeter (US Markets Only)</p> <p>The MightySat Rx Fingertip Pulse Oximeter is intended as a noninvasive device that measures and displays arterial oxygen saturation (SpO₂), Pulse Rate (PR), and Perfusion Index (Pi).</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 MightySat Rx • 1 Carrying Case • 1 Lanyard • 2 AAA Batteries 	EA	
9809	<p>MightySat Rx Fingertip Pulse Oximeter w/Bluetooth (US Markets Only)</p> <p>The MightySat Rx Fingertip Pulse Oximeter is intended as a noninvasive device that measures and displays arterial oxygen saturation (SpO₂), Pulse Rate (PR), and Perfusion Index (Pi). Optional Bluetooth LE wireless technology for the wireless transfer of patient data to smart devices.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 MightySat Rx • 1 Carrying Case • 1 Lanyard • 2 AAA Batteries <p><i>See Masimo.com for a list of compatible devices.</i></p>	EA	
9909	<p>MightySat Rx Fingertip Pulse Oximeter w/Bluetooth, PVi, RRp (US Markets Only)</p> <p>The MightySat Rx Fingertip Pulse Oximeter is intended as a noninvasive device that measures and displays arterial oxygen saturation (SpO₂), Pulse Rate (PR), Perfusion Index (Pi), and optional Pleth Variability Index (PVi), and Respiration Rate from the Pleth (RRp). Optional Bluetooth LE wireless technology for the wireless transfer of patient data to smart devices.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 MightySat Rx • 1 Carrying Case • 1 Lanyard • 2 AAA Batteries <p><i>See Masimo.com for a list of compatible devices.</i></p>	EA	

iSpO₂® Rx




PN	Description	UM	Product Image
3588	<p>iSpO₂ Rx 30-Pin Series, Adult See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included.</i></p>	EA	
3590	<p>iSpO₂ Rx 30-Pin Series, Pediatric See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included.</i></p>	EA	
3721	<p>iSpO₂ Rx 30-Pin Series, M-LNCS Available for adults and pediatric patients only. See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</i></p>	EA	






PN	Description	UM	Product Image
9575	<p>iSpO₂ Rx Lightning Series, Adult See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included.</i></p>	EA	
9576	<p>iSpO₂ Rx Lightning Series, Pediatric See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included.</i></p>	EA	
9569	<p>iSpO₂ Rx Lightning Series, M-LNCS Available for adult and pediatric patients only. See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</i></p>	EA	
9578	<p>iSpO₂ Rx Micro USB Series, Adult See www.masimo.com for compatibility information. <i>Available in countries where cleared. Phone not included.</i></p>	EA	


PN	Description	UM	Product Image
9579	<p>iSpO₂ Rx Micro USB Series, Pediatric See www.masimo.com for compatibility information. Available in countries where cleared. Phone not included.</p>	EA	
9577	<p>iSpO₂ Rx Micro USB Series, M-LNCS Available for adult and pediatric patients only. See www.masimo.com for compatibility information. Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</p>	EA	

Parameter Options and Upgrade Tools

Masimo rainbow SET™ Parameter Options

PN	Description	UM	Product Image
9081	<p>Field Upgrader</p> <p>Used to load optional rainbow® parameters and features onto compatible monitors; available in countries where cleared.</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Oxygen Reserve Index™ (ORi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) • 3D Desat Index Alarm™ • 3D Pi Delta Alarm <p>Each optional parameter and feature sold separately. Maximum quantity of 25 of any individual parameter or feature per Field Upgrader.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 Field Upgrader • 1 to 5 operator's manuals • Parameter upgrade indicator sticker(s) • Postage-paid envelope to return Field Upgrader 	EA	
9792	<p>M20 Field Upgrade Tool</p> <p>Used to load optional rainbow parameters and features onto compatible monitors; available in countries where cleared.</p> <p>Available upgrade options include:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Pleth Variability Index (PVi) • Oxygen Reserve Index (ORi) • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Oxygen Content (SpOC) • Acoustic Respiration (RRa) • Respiration Rate from Pleth (RRp) <p><i>Each optional parameter and feature sold separately.</i></p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 M20 Field Upgrade Tool • 1 to 5 Operator's Manuals • Parameter upgrade indicator sticker(s) • Postage-paid envelope to return Field Upgrade Tool 	EA	
3688	<p>Root (Non-EWS), Software Upgrade Tool, US</p> <p>Includes:</p> <ul style="list-style-type: none"> • Upgrade tool • Upgrade tool directions for use • Postage-paid envelope to return upgrade tool 	EA	





PN	Description	UM	Product Image
3689	<p>Root (Non-EWS), Software Upgrade Tool, OUS</p> <p>Includes: Upgrade tool Upgrade tool directions for use Postage-paid envelope to return upgrade tool <i>Available in countries where cleared.</i></p>	EA	
4382	<p>Root with EWS, Software Upgrade Tool, US</p> <p>Includes: Upgrade tool Upgrade tool directions for use Postage-paid envelope to return upgrade tool</p>	EA	
4383	<p>Root with EWS, Software Upgrade Tool, OUS</p> <p>Includes: Upgrade tool Upgrade tool directions for use Postage-paid envelope to return upgrade tool <i>Available in countries where cleared.</i></p>	EA	
4397	<p>Rad-97, Software Upgrade Tool, US</p> <p>Includes: Upgrade tool Upgrade tool directions for use Postage-paid envelope to return upgrade tool</p>	EA	
4398	<p>Rad-97, Software Upgrade Tool, OUS</p> <p>Includes: Upgrade tool Upgrade tool directions for use Postage-paid envelope to return upgrade tool <i>Available in countries where cleared.</i></p>	EA	
3525	<p>Radical-7, Software Upgrade Tool, US</p> <p>Includes: Upgrade Tool Upgrade Tool Directions For Use Postage-paid envelope to return upgrade tool</p>	EA	No Image
3526	<p>Radical-7, Software Upgrade Tool, OUS</p> <p>Include: Upgrade Too Upgrade Tool Directions For Use Postage-paid envelope to return upgrade tool <i>Available in countries where cleared</i></p>	EA	No Image
3645	<p>Eve™ (Newborn CCHD Screening Mode)</p> <p>The Eve Newborn Screening Application on the Radical-7 Pulse CO-Oximeter automates the screening steps with animated instructions, including sensor application, measurement selection, and screening result determination.</p> <p>Includes: Upgrade card Upgrade card directions for use Radical-7 Operator's Manual supplement, CCHD Postage-paid envelope to return upgrade card <i>Available in countries where cleared.</i></p>	EA	No image

PN	Description	UM	Product Image
9848	<p>Eve™ (Newborn CCHD Screening Mode)</p> <p>The Eve Newborn Screening Application on the Rad-97 Pulse CO-Oximeter automates the screening steps with animated instructions, including sensor application, measurement selection, and screening result determination.</p> <p>Includes:</p> <ul style="list-style-type: none"> Upgrade tool Upgrade tool directions for use CCHD addendum to Rad-97 Operator's Manual Postage-paid envelope to return upgrade tool <p><i>Available in countries where cleared.</i></p>	EA	
3710	<p>Service Tool, MX-5, US</p> <p>Upgrades up to 100 instruments.</p>	EA	No image
4185	<p>Service Tool, MX-5, OUS</p> <p>Upgrades up to 100 instruments.</p> <p><i>Available in countries where cleared.</i></p>	EA	No image
4552	<p>Service Tool, MX-5, M20, US</p> <p>Upgrades up to 100 instruments.</p>	EA	No image
4553	<p>Service Tool, MX-5, M20, OUS</p> <p>Upgrades up to 100 instruments.</p> <p><i>Available in countries where cleared.</i></p>	EA	No image
2341	<p>Service Tool, MX-1®</p> <p>Upgrades up to 100 instruments.</p>	EA	No image
2897	<p>Service Tool, MX-3</p> <p>Upgrades up to 100 instruments.</p>	EA	No image
3351	<p>Service Tool, MX-3 with ATA</p> <p>Upgrades up to 100 instruments.</p> <p><i>Available in countries where cleared.</i></p>	EA	No image
The following parameters are applicable to both PN 9081 and PN 9792:			
2630	<p>Total Hemoglobin (SpHb) parameter</p> <p>Can be ordered on a Radical-7, Rad-97, Rad-57 or on an upgrade tool to upgrade a compatible monitor in the field</p>	EA	No image
2313	<p>Pleth Variability Index (PVi)</p> <p>Can be ordered on a Radical-7, Rad-97, Rad-67*, Rad-57 or on an upgrade tool to upgrade a compatible monitor in the field.</p> <p><i>* Available in countries where cleared</i></p>	EA	No image
3950	<p>Oxygen Reserve Index (ORi)</p> <p>Can be ordered on a Radical-7 or Rad-97 or on an upgrade tool to upgrade a compatible monitor in the field.</p> <p><i>Available in countries where cleared.</i></p>	EA	No image





PN	Description	UM	Product Image
2296	<p>Carboxyhemoglobin (SpCO) parameter</p> <p>Can be ordered on a Radical-7, Rad-97, Rad-67*, Rad-57 or on an upgrade tool to upgrade a compatible monitor in the field.</p> <p><i>* Available in countries where cleared</i></p>	EA	No image
2297	<p>Methemoglobin (SpMet) parameter</p> <p>Can be ordered on a Radical-7, Rad-97, Rad-67*, Rad-57 or on an upgrade tool to upgrade a compatible monitor in the field.</p> <p><i>* Available in countries where cleared</i></p>	EA	No image
2940	<p>Acoustic Respiration (RRa)</p> <p>Can be ordered on a Rad-97 or on an upgrade tool to upgrade a compatible monitor in the field.</p>	EA	No image
3518	<p>Respiration Rate from Pleth (RRp)</p> <p>Can be ordered on a Radical-7, Rad-97 or on an upgrade tool to upgrade a compatible monitor in the field.</p>	EA	No image
2636	<p>Pi Delta 3D Alarm®, Fixed Term License</p> <p>Can be ordered on a Radical-7, Rad-97 or on a Field Upgrader (PN 9081) to upgrade a compatible monitor in the field.</p>	EA	No image
2637	<p>Desat Index 3D Alarm, Fixed Term License</p> <p>Can be ordered on a Radical-7, Rad-97 or on a Field Upgrader (PN 9081) to upgrade a compatible monitor in the field.</p>	EA	No image

Sensors and Cables

RD rainbow SET™ Adhesive Sensors (for rainbow SET Devices with SpHb®, SpO₂, and SpMet®)

PN	Description	UM	Product Image
4026	<p>RD rainbow SET-2 Adt Adult Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i></p>	Box	
4027	<p>RD rainbow SET-2 Pdt Pediatric Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight 10-50 kg Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i></p>	Box	
4028	<p>RD rainbow SET-2 Inf Infant Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight 3-30 kg Use replacement tape 4439 Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i></p>	Box	
4029	<p>RD rainbow SET-2 Neo Adult / Neonatal Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight <3 kg*, > 30 kg Use replacement tape 4440 Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i> *SpO₂ and SpMet only</p>	Box	





RD rainbow Lite SET™ Adhesive Sensors (for rainbow SET Devices with SpO₂, and ORi)

PN	Description	UM	Product Image
4042	<p>RD rainbow Lite SET-1 Adt Adult Adhesive Sensors (SpO₂, ORi) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i> <i>Available in countries where cleared.</i></p>	Box	
4043	<p>RD rainbow Lite SET-1 Pdt Pediatric Adhesive Sensors (SpO₂, ORi) Single patient use 10 / box, non-sterile Weight 10-50 kg Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i> <i>Available in countries where cleared.</i></p>	Box	
4044	<p>RD rainbow Lite SET-1 Inf Infant Adhesive Sensors (SpO₂, ORi) Single patient use 10 / box, non-sterile Weight 3-30 kg Use replacement tape 4046 Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i> <i>Available in countries where cleared.</i></p>	Box	
4045	<p>RD rainbow Lite SET-1 Neo Adult / Neonatal Adhesive Sensors (SpO₂, ORi) Single patient use 10 / box, non-sterile Weight <3 kg*, > 30 kg Use replacement tape 4047 Not made with natural rubber latex. <i>Requires RD rainbow SET patient cable.</i> <i>Available in countries where cleared.</i> *SpO₂ only</p>	Box	






rainbow® Adhesive Sensors R1 series (for rainbow SET Devices with SpHb, SpO₂, and SpMet)



PN	Description	UM	Product Image
3792	<p>rainbow R1 25 Adult Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. Requires rainbow RC patient cable</p>	Box	
3793	<p>rainbow R1 20 Pediatric / Slender Digit Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. Requires rainbow RC patient cable</p>	Box	
2415	<p>rainbow R1 20L Infant Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement tape 2624. Requires rainbow RC patient cable</p>	Box	
2414	<p>rainbow R1 25L Adult / Neonatal Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight <3 kg* or >30 kg Not made with natural rubber latex. Use replacement tape 2623. Requires rainbow RC patient cable *SpO₂ And SpMet only</p>	Box	


rainbow® Adhesive Sensors (for rainbow SET Devices with SpO₂, SpCO®, and SpMet)

PN	Description	UM	Product Image
2221	<p>rainbow R25 Adult Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i></p>	Box	
2222	<p>rainbow R20 Pediatric Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i></p>	Box	
2220	<p>rainbow R20-L Infant Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement tape 2624. <i>Requires rainbow RC patient cable</i></p>	Box	
2219	<p>rainbow R25-L Adult / Neonatal Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight <3 kg* or >30 kg Not made with natural rubber latex. Use replacement tape 2623. <i>Requires rainbow RC patient cable</i> <i>*SpO₂ And SpMet only</i></p>	Box	

rainbow® Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO₂, SpMet; SpHb required)

PN	Description	UM	Product Image
3418	<p>rainbow DCI SC 200 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3419	<p>rainbow DCI-P SC 200 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3420	<p>rainbow DCI SC 400 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3421	<p>rainbow DCI-P SC 400 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3647	<p>rainbow DCI SC 1000 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	



PN	Description	UM	Product Image
3648	<p>rainbow DCI-P SC 1000 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3797	<p>rainbow DCI-mini SC-200 Adult/Pediatric/Infant reusable sensor (SpHb, SpO₂ and SpMet*) Weight ≥ 3 kg Not made with natural rubber latex 200 Spot Checks <i>Requires rainbow RC series patient cable</i> *Available in countries where cleared.</p>	EA	
3798	<p>rainbow DCI-mini SC-400 Adult/Pediatric/Infant reusable sensor (SpHb, SpO₂ and SpMet*) Weight ≥ 3 kg Not made with natural rubber latex 400 Spot Checks <i>Requires rainbow RC series patient cable</i> *Available in countries where cleared.</p>	EA	
4309	<p>rainbow Super DCI-mini SC-200 Adult/Pediatric/Infant reusable sensor (SpO₂, SpHb, SpCO, and SpMet) Weight ≥ 3 kg Not made with natural rubber latex 200 Spot Checks <i>Requires rainbow RC series patient cable</i> Available in countries where cleared.</p>	EA	
4310	<p>rainbow Super DCI-mini SC-400 Adult/Pediatric/Infant reusable sensor (SpO₂, SpHb, SpCO, and SpMet) Weight ≥ 3 kg Not made with natural rubber latex 400 Spot Checks <i>Requires rainbow RC series patient cable</i> Available in countries where cleared.</p>	EA	

PN	Description	UM	Product Image
4311	<p>rainbow Super DCI-mini SC-1000 Adult/Pediatric/Infant reusable sensor (SpO₂, SpHb, SpCO, and SpMet) Weight ≥ 3 kg Not made with natural rubber latex 1000 Spot Checks <i>Requires rainbow RC series patient cable</i> Available in countries where cleared.</p>	EA	






rainbow® Reusable Sensors (for rainbow SET Devices with SpCO, SpO₂, and SpMet)






PN	Description	UM	Product Image
2696	<p>rainbow DCI Adult Reusable Sensor, 3 ft. (SpCO, SpMet, SpO₂) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i></p>	EA	
2697	<p>rainbow DCI-P Pediatric Reusable Sensor, 3 ft. (SpCO, SpMet, SpO₂) 1 / box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i></p>	EA	
3657	<p>rainbow DCI-6 Adult Reusable Sensor, 6ft (SpCO, SpMet, SpO₂) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i></p>	EA	
3658	<p>rainbow DCI-P-6 Pediatric Reusable Sensor, 6ft (SpCO, SpMet, SpO₂) 1 box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex <i>Requires rainbow RC Cable</i></p>	EA	





rainbow® Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO₂, and SpMet)





PN	Description	UM	Product Image
2201	<p>rainbow DCI-dc3 Adult Reusable Direct Connect Sensor, 3 ft. (SpO₂, SpCO, SpMet) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2069	<p>rainbow DCI-P-dc3 Pediatric Reusable Direct Connect Sensor, 3 ft. (SpO₂, SpCO, SpMet) 1 / box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex.</p>	EA	





rainbow® Patient Cables (for rainbow SET Devices)





PN	Description	UM	Product Image
4071	<p>RD rainbow SET MD20-1.5 RD rainbow 20-pin Patient Cable, 1.5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors Not for use with Rad-8 or Rad-5/5v</p>	EA	
4072	<p>RD rainbow SET MD20-05 RD rainbow 20-pin Patient Cable, 5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors Not for use with Rad-8 or Rad-5/5v</p>	EA	
4109	<p>RD rainbow SET MD20-08 RD rainbow 20-pin Patient Cable, 8 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors Not for use with Rad-8 or Rad-5/5v</p>	EA	
4073	<p>RD rainbow SET MD20-12 RD rainbow 20-pin Patient Cable, 12 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors Not for use with Rad-8 or Rad-5/5v</p>	EA	
4255	<p>RD rainbow SET M20-1.5 RD rainbow M20 Patient Cable, 1.5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	





PN	Description	UM	Product Image
4256	<p>RD rainbow SET M20-05 RD rainbow M20 Patient Cable, 5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4257	<p>RD rainbow SET M20-12 RD rainbow M20 Patient Cable, 12 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4236	<p>rainbow M20-01 rainbow M20 Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p>	EA	
4237	<p>rainbow M20-04 rainbow M20 Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p>	EA	
4238	<p>rainbow M20-12 rainbow M20 Patient Cable, 12 ft For use with rainbow and M-LNCS sensors</p>	EA	





PN	Description	UM	Product Image
2405	<p>rainbow RC-1 rainbow 20-pin Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p>	EA	
2406	<p>rainbow RC-4 rainbow 20-pin Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p>	EA	
2404	<p>rainbow RC-12 rainbow 20-pin Patient Cable, 12 ft For use with rainbow and M-LNCS sensors</p>	EA	
3659	<p>rainbow RC-EXT-8 rainbow extension cable, 8 ft For use with rainbow and M-LNCS sensors</p>	EA	

PN	Description	UM	Product Image
4074	<p>RD rainbow SET R25-1.5 RD rainbow 25-pin Patient Cable, 1.5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4075	<p>RD rainbow SET RA25-1.5 RD rainbow 25-pin, right angle Patient Cable, 1.5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4076	<p>RD rainbow SET R25-05 RD rainbow 25-pin Patient Cable, 5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4077	<p>RD rainbow SET RA25-05 RD rainbow 25-pin, right angle Patient Cable, 5 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	





PN	Description	UM	Product Image
4106	<p>RD rainbow SET R25-08 RD rainbow 25-pin Patient Cable, 8 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4107	<p>RD rainbow SET RA25-08 RD rainbow 25-pin, right angle Patient Cable, 8 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4078	<p>RD rainbow SET R25-12 RD rainbow 25-pin Patient Cable, 12 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	
4079	<p>RD rainbow SET RA25-12 RD rainbow 25-pin, right angle Patient Cable, 12 ft For use with RD rainbow SET, RD rainbow Lite SET, or RD SET series sensors</p>	EA	





PN	Description	UM	Product Image
3933	<p>rainbow RC25-1 rainbow 25-pin, Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p>	EA	
3934	<p>rainbow RC25-1RA rainbow 25-pin, right angle Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p>	EA	
3628	<p>rainbow RC25-4 rainbow 25-pin, Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p>	EA	
3375	<p>rainbow RC25-4RA rainbow 25-pin, Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p>	EA	






PN	Description	UM	Product Image
3342	<p>rainbow RC25-12 rainbow 25-pin, Patient Cable, 12 ft For use with rainbow and M-LNCS sensors</p>	EA	
3343	<p>rainbow RC25-12RA rainbow 25-pin, right angle Patient Cable, 12 ft For use with rainbow and M-LNCS sensors</p>	EA	
4480	<p>RC-1 EMS, rainbow Patient Cable rainbow 20-pin Patient Cable, 1 ft For use with rainbow and M-LNCS sensors Note: EMS Cables have larger cable diameter compared to standard patient cables including Kevlar materials for enhanced strength</p>	EA	
4481	<p>RC-4 EMS, rainbow Patient Cable rainbow 20-pin Patient Cable, 4 ft For use with rainbow and M-LNCS sensors Note: EMS Cables have larger cable diameter compared to standard patient cables including Kevlar materials for enhanced strength</p>	EA	

PN	Description	UM	Product Image
4482	<p>RC25-1 EMS, rainbow Patient Cable rainbow 25-pin Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p> <p>Note: EMS Cables have larger cable diameter compared to standard patient cables including Kevlar materials for enhanced strength</p>	EA	
4483	<p>RC25-4 EMS, rainbow Patient Cable rainbow 25-pin Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p> <p>Note: EMS Cables have larger cable diameter compared to standard patient cables including Kevlar materials for enhanced strength</p>	EA	
4484	<p>RC25-1RA EMS, rainbow Patient Cable rainbow 25-pin, right angle Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p> <p>Note: EMS Cables have larger cable diameter compared to standard patient cables including Kevlar materials for enhanced strength</p>	EA	
4485	<p>RC25-4RA EMS, rainbow Patient Cable rainbow 25-pin, right angle Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p> <p>Note: EMS Cables have larger cable diameter compared to standard patient cables including Kevlar materials for enhanced strength</p>	EA	





rainbow Acoustic Monitoring® Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRa®)




PN	Description	UM	Product Image
4171	<p>RAS-45, Acoustic Respiration Sensor, 3 ft. Acoustic Respiration Sensor, 3 ft. Adult / Pediatric Adhesive Sensors (RRa) Single use 10/box, non-sterile Weight > 10 kg Not made with natural rubber latex</p>	Box	
4425	<p>RAS-45 Inf/Neo, Acoustic Respiration Sensor, 3 ft. Acoustic Respiration Sensor, 3 ft. Infant / Neonatal Adhesive Sensors (RRa) Single use 10/box, non-sterile Weight < 10 kg Not made with natural rubber latex For US markets Only</p>	Box	
3475	<p>RAS-125c, Acoustic Respiration Cloth Sensor Adult / Pediatric Adhesive Sensors (RRa) Single patient use 10 / box, non-sterile Weight > 10kg Not made with natural rubber latex.</p>	Box	
3483	<p>RAS-125c, Short Term Monitoring Acoustic Respiration Cloth Sensor Adult / Pediatric Adhesive Sensors (RRa) Single patient use 10 / box, non-sterile Weight > 10kg Not made with natural rubber latex.</p>	Box	

PN	Description	UM	Product Image
4279	<p>RD rainbow SET RAM Dual Cable MD20-12</p> <p>RD & RAS rainbow 20 pin dual patient cable, 12 ft</p> <p>For use with RD SpO2 series or rainbow® series sensors and RAS series sensors</p>	EA	
4512	<p>RD rainbow SET RAM Dual Cable M20-12</p> <p>RD & RAS rainbow M20 dual patient cable, 12 ft</p> <p>For use with RD SpO2 series or rainbow® series sensors and RAS series sensors</p>	EA	
4278	<p>RD rainbow SET RAM Dual Cable R25-12</p> <p>RD & RAS rainbow 25 pin dual patient cable, 12 ft</p> <p>For use with RD SpO2 series or rainbow® series sensors and RAS series sensors</p>	EA	
3661	<p>RAM Dual Cable LNC-10</p> <p>LNCS & RAS 20 pin dual patient cable, 10 ft</p> <p>For use with LNCS SpO2 series sensors and RAS series sensors</p>	EA	


PN	Description	UM	Product Image
4099	<p>RAM Dual Cable LNC M20-10 LNCS & RAS M20 dual patient cable, 10 ft For use with LNCS SpO2 series sensors and RAS series sensors</p>	EA	
3906	<p>RAM Dual Cable, Red 25 LNC-10 LNCS & RAS 25 pin dual patient cable, 10 ft For use with LNCS SpO2 series sensors and RAS series sensors</p>	EA	
3660	<p>RAM Dual Cable RC-10 M-LNCS & RAS rainbow 20 pin dual patient cable, 10 ft For use with M-LNCS SpO2 series or rainbow® series sensors and RAS series sensors</p>	EA	
4100	<p>RAM Dual Cable M20-10 M-LNCS & RAS rainbow M20 dual patient cable, 10 ft For use with M-LNCS SpO2 series or rainbow® series sensors and RAS series sensors</p>	EA	
3503	<p>RAM Dual Cable RC25-10 M-LNCS & RAS rainbow 25 pin dual patient cable, 10 ft For use with M-LNCS SpO2 series or rainbow® series sensors and RAS series sensors</p>	EA	






RD SET™ Single Patient Use Sensors






PN	Description	UM	Product Image
4000	<p>RD SET Adt Adult SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Available in countries where cleared</i></p>	Box	
4001	<p>RD SET Pdt Pediatric SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex <i>Available in countries where cleared</i></p>	Box	
4002	<p>RD SET Inf Infant SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight 3 kg – 20 kg Not made with natural rubber latex Use Replacement tape 4046 <i>Available in countries where cleared</i></p>	Box	
4003	<p>RD SET Neo Adult / Neonatal SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex Use Replacement tape 4047 <i>Available in countries where cleared</i></p>	Box	


PN	Description	UM	Product Image
4004	<p>RD SET NeoPt Neonatal SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. Use Replacement tape 4048 <i>Available in countries where cleared</i></p>	Box	
4005	<p>RD SET NeoPt-500 Neonatal SpO₂ Non-adhesive Sensor Single-patient-use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. Use Replacement tape 4094 <i>Available in countries where cleared</i></p>	Box	
4014	<p>RD SET Blue® Neonatal / Infant / Pediatric SpO₂ Adhesive Sensor Single patient use 10 / box, Non-Sterile Weight 2.5 – 30kg Not made with natural rubber latex Use replacement tape 4046</p>	Box	

RD SET™ Cabled Sensors

PN	Description	UM	Product Image
4470	<p>RD SET Adt CS-2 Adult SpO₂ Adhesive Sensors, 2 ft. Single-patient-use 20/Box, non-sterile Weight > 30 kg Not made with natural rubber latex</p>	Box	



PN	Description	UM	Product Image
4471	<p>RD SET Pdt CS-2 Pediatric SpO₂ Adhesive Sensors, 2 ft. Single-patient-use 20/Box, non-sterile Weight 10-50 kg Not made with natural rubber latex</p>	Box	
4472	<p>RD SET Inf CS-2 Infant SpO₂ Adhesive Sensors, 2 ft. Single-patient-use 20/Box, non-sterile Weight 3-20 kg Not made with natural rubber latex Use replacement tape 4046</p>	Box	
4473	<p>RD SET Neo CS-2 Adult / Neonatal SpO₂ Adhesive Sensors, 2 ft. Single-patient-use 20/Box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex Use replacement tape 4047</p>	Box	
4474	<p>RD SET NeoPt CS-2 Neonatal SpO₂ Adhesive Sensors, 2 ft. Single-patient-use 20/Box, non-sterile Weight < 1 kg Not made with natural rubber latex Use replacement tape 4048</p>	Box	
4475	<p>RD SET Adt CS-3 Adult SpO₂ Adhesive Sensors, 3 ft. Single-patient-use 20/Box, non-sterile Weight > 30 kg Not made with natural rubber latex</p>	Box	


PN	Description	UM	Product Image
4476	<p>RD SET Pdt CS-3 Pediatric SpO₂ Adhesive Sensors, 3 ft. Single-patient-use 20/Box, non-sterile Weight 10-50 kg Not made with natural rubber latex</p>	Box	
4477	<p>RD SET Neo CS-3 Adult / Neonatal SpO₂ Adhesive Sensors, 3 ft. Single-patient-use 20/Box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex Use replacement tape 4047</p>	Box	
4478	<p>RD SET Neo CS-3 Adult / Neonatal SpO₂ Adhesive Sensors, 3 ft. Single-patient-use 20/Box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex Use replacement tape 4047</p>	Box	
4479	<p>RD SET NeoPt CS-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single-patient-use 20/Box, non-sterile Weight < 1 kg Not made with natural rubber latex Use replacement tape 4048</p>	Box	
4015	<p>RD SET E1® Adult SpO₂ Disposable Ear Sensor, 3 ft. Single-patient-use 10 / box, non-sterile, non-adhesive Weight > 30 kg Not made with natural rubber latex. Requires Masimo SET® MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology</p>	Box	

PN	Description	UM	Product Image
4016	<p>RD SET TFA-1® Adult / Pediatric SpO₂ disposable tranreflectance forehead sensor, 3 ft. Single-patient-use 10 / box, non-sterile Weight > 10 kg Not made with natural rubber latex Requires Masimo SET® MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology</p>	Box	



RD SET Specialty Sensors




RD SET Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of averaging, sensitivity, or FastSat® setting.


PN	Description	UM	Product Image
4011	<p>RD SET Trauma Adult SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex Use Replacement tape 4048</p>	Box	
4012	<p>RD SET Newborn Infant/Pediatric Infant/Pediatric SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight 3 – 30 kg Not made with natural rubber latex Use Replacement tape 4094</p>	Box	

PN	Description	UM	Product Image
4013	<p>RD SET Newborn Neonatal Neonatal SpO₂ Adhesive Sensor Single-patient-use 20 / box, non-sterile Weight < 3 kg Not made with natural rubber latex Use Replacement tape 4048</p>	Box	



RD SET Reusable Sensors





PN	Description	UM	Product Image
4050	<p>RD SET DCI Adult Reusable Sensor, 3 ft 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex. <i>Available in countries where cleared.</i></p>	EA	
4051	<p>RD SET DCI-P Pediatric/Slender Digit Reusable Sensor, 3 ft 1 / box, non-sterile Weight 10 - 50 kg Not made with natural rubber latex. <i>Available in countries where cleared.</i></p>	EA	


PN	Description	UM	Product Image
4052	<p>RD SET DBI Adult Soft Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex</p>	EA	
4053	<p>RD SET TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex. <i>Available in countries where cleared.</i></p>	EA	
4054	<p>RD SET YI Multi-site Reusable Sensor, 3 FT 1/box, non-sterile Multiple foam and adhesive wraps Weight >1 kg Not made with natural rubber latex</p>	EA	





PN	Description	UM	Product Image
4055	<p>RD SET TF-I® Adult Transreflectance Forehead Sensor, 3 FT 1/box, non-sterile Weight >30 kg Not made with natural rubber latex</p>	EA	

RD SET Patient Cables





PN	Description	UM	Product Image
4102	<p>RD SET MD20 - 1.5 RD SET 20-pin SpO₂ Patient Cable</p>	EA	
4103	<p>RD SET MD20 - 05 RD SET 20-pin SpO₂ Patient Cable, 5 ft.</p>	EA	



PN	Description	UM	Product Image
4104	<p>RD SET MD20 - 12 RD SET 20-pin SpO₂ Patient Cable, 12 ft.</p>	EA	 <p>A 12-foot patient cable with a red connector at one end and a white connector at the other.</p>
4080	<p>RD SET MD14 - 05 RD SET Patient Cable, 5 ft For use with RD SET series sensors</p>	EA	 <p>A 5-foot patient cable with white connectors at both ends.</p>
4108	<p>RD SET MD14 - 08 RD SET Patient Cable, 8 ft. For use with RD SET series sensors</p>	EA	 <p>An 8-foot patient cable with white connectors at both ends.</p>
4081	<p>RD SET MD14 - 12 RD SET Patient Cable, 12 ft For use with RD SET series sensors</p>	EA	 <p>A 12-foot patient cable with white connectors at both ends.</p>

PN	Description	UM	Product Image
4082	<p>RD SET MP - 05 RD SET Dual Key patient cable, 5 ft RD SET to Intellivue Masimo SET® or Philips FAST SpO₂ For use with RD SET series sensors</p>	EA	
4083	<p>RD SET MP - 12 RD SET Dual Key patient cable, 12 ft RD SET to Intellivue Masimo SET or Philips FAST SpO₂ For use with RD SET series sensors</p>	EA	
4084	<p>RD SET GE - 05 RD SET GE Patient Cable, 5 ft For use with RD SET series sensors</p>	EA	
4085	<p>RD SET GE - 12 RD SET GE Patient Cable, 12 ft For use with RD SET series sensors</p>	EA	


PN	Description	UM	Product Image
4087	<p>RD SET SL - 12 RD SET SpaceLabs Patient Cable, 12 ft For use with RD SET series sensors</p>	EA	 A coiled white patient cable with a standard connector on one end and a SpaceLabs connector on the other.
4191	<p>RD SET NK - 05 RD SET Nihon Kohden Patient Cable, 5 ft For use with RD SET series sensors</p>	EA	 A coiled white patient cable with a standard connector on one end and a Nihon Kohden connector on the other.
4192	<p>RD SET NK - 12 RD SET Nihon Kohden Patient Cable, 12 ft For use with RD SET series sensors</p>	EA	 A coiled white patient cable with a standard connector on one end and a Nihon Kohden connector on the other.
4095	<p>RD SET EXT - 4 RD SET series extension cable, 4 ft</p>	EA	 A coiled white extension cable with a standard connector on one end and a sensor connector on the other.


RD SET Adapter Cables



PN	Description	UM	Product Image
4088	<p>RD to PC Adapter Cable RD SET Series sensors to LNOP PC patient cable, 1.5 ft</p>	EA	
4089	<p>RD to LNC Adapter Cable RD SET Series sensors to uSpO₂ cable, 3 ft</p>	EA	
4105	<p>RD to LNC Adapter Cable RD SET Series sensors to LNC, uSpO₂ cable, or LNC-MAC patient cable, 1.5 ft</p>	EA	
4090	<p>RD to M-LNC Adapter Cable RD SET Series sensors to M-LNC patient cable, 1.5 ft</p>	EA	


PN	Description	UM	Product Image
4091	<p>LNOP to RD Adapter Cable LNOP Series to RD patient cable, 1.5 ft.</p>	EA	
4092	<p>LNCS to RD Adapter Cable LNCS Series to RD patient cable, 1 ft.</p>	EA	

LNCS® Single Patient Use Sensors




PN	Description	UM	Product Image
1859	<p>LNCS Adtx Adult SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.</p>	Box	





PN	Description	UM	Product Image
2317	<p>LNCS Adtx-3 Adult SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.</p>	Box	
1860	<p>LNCS Pdtx Pediatric SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
2318	<p>LNCS Pdtx-3 Pediatric SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
1861	<p>LNCS Inf-L Infant SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex.</p>	Box	

PN	Description	UM	Product Image
1862	<p>LNCS Neo-L Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex.</p>	Box	
1901	<p>LNCS NeoPt-L Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex.</p>	Box	
2918	<p>LNCS E1® Adult SpO₂ disposable ear sensor, 3 ft. Single patient use 10 / box, non-sterile, non-adhesive Weight > 30 kg Not made with natural rubber latex. Requires Masimo SET® MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology</p>	Box	

PN	Description	UM	Product Image
3858	<p>LNCS TFA-1™ Adult / Pediatric SpO₂ disposable tranreflectance forehead sensor, 3ft Single patient use 10 / box non-sterile Weight > 10 kg Not made with natural rubber latex.</p>	Box	




LNCS Single Patient Use Sensors with Replaceable Tapes

PN	Description	UM	Product Image
2328	<p>LNCS Inf Infant SpO₂ Adhesive Sensors, 18 in. Single Patient Use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	
2319	<p>LNCS Inf-3 Infant SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	
2329	<p>LNCS Neo Neonatal SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	

PN	Description	UM	Product Image
2320	<p>LNCS Neo-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	
2330	<p>LNCS NeoPt Neonatal SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2321	<p>LNCS NeoPt-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2331	<p>LNCS NeoPt-500 Neonatal SpO₂ non-adhesive sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	




LNCS Specialty Sensors




LNCS Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.

PN	Description	UM	Product Image
2411	<p>LNCS Trauma Adult SpO₂ Adhesive Sensors Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2412	<p>LNCS Newborn Neonatal Neonatal SpO₂ Adhesive Sensors Single patient use 20 / box, sterilizable Weight < 3 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2413	<p>LNCS Newborn Infant / Pediatrics Infant / Pediatric SpO₂ Adhesive Sensor Single patient use 20 / box, sterilizable Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	



LNCS Sterile Sensors



Note: The sensors are not intended to be used for evaluations or carried as sales samples.

PN	Description	UM	Product Image
2601	<p>LNCS Amtx-3 Adult SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight > 30 kg Not made with natural rubber latex.</p>	Box	
2602	<p>LNCS Pmtx-3 Pediatric SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
2603	<p>LNCS Inf-3 Infant SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	



PN	Description	UM	Product Image
2604	<p>LNCS Neo-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	
2605	<p>LNCS NeoPt-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2606	<p>LNCS NeoPt-500 Neonatal SpO₂ non-adhesive sensor, 18 in. Single patient use 20 / box, Sterile Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	





LNCS Reusable Sensors



PN	Description	UM	Product Image
1863	<p>LNCS DCI Adult Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1864	<p>LNCS DCI-P Pediatric Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	EA	
2653	<p>LNCS DBI Adult Soft Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1895	<p>LNCS TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	

PN	Description	UM	Product Image
1896	<p>LNCS TF-I® Adult Transflectance Forehead, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2258	<p>LNCS YI Multisite Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.</p>	EA	



LNCS Red Patient Cables (for use with Masimo rainbow SET Devices)





PN	Description	UM	Product Image
4251	<p>Red LNC M20-01 LNCS M20 SpO₂ Patient Cable, 1 ft.</p>	EA	
4252	<p>Red LNC M20-04 LNCS M20 SpO₂ Patient Cable, 4 ft.</p>	EA	





PN	Description	UM	Product Image
4253	<p>Red LNC M20-10 LNCS M20 SpO₂ Patient Cable, 10 ft.</p>	EA	
4254	<p>Red LNC M20-14 LNCS 20-pin SpO₂ Patient Cable, 14 ft.</p>	EA	
2365	<p>Red LNC-01 LNCS 20-pin SpO₂ Spot Check Patient Cable.</p>	EA	
2055	<p>Red LNC-04 LNCS 20-pin SpO₂ Patient Cable, 4 ft.</p>	EA	





PN	Description	UM	Product Image
2056	Red LNC-10 LNCS 20-pin SpO ₂ Patient Cable, 10 ft.	EA	
2057	Red LNC-14 LNCS 20-pin SpO ₂ Patient Cable, 14 ft.	EA	





LNC Patient Cables (for use with Masimo SET® Devices)

PN	Description	UM	Product Image
2364	LNC-1 LNCS Series Patient Cable.	EA	
2017	LNC-4 LNCS Series Patient Cable, 4 ft.	EA	




PN	Description	UM	Product Image
1814	<p>LNC-10 LNCS Series Patient Cable, 10 ft.</p>	EA	
3935	<p>Red 25 LNC-1 LNCS Series 25-pin SpO₂ Patient Cable, 1 ft.</p>	EA	
3630	<p>Red 25 LNC-4 LNCS Series 25-pin SpO₂ Patient Cable, 4 ft</p>	EA	
3936	<p>Red 25 LNC-1RA LNCS Series 25-pin right angle SpO₂ Patient Cable, 1 ft</p>	EA	




PN	Description	UM	Product Image
3376	<p>Red 25 LNC-4RA Series LNCS Series 25-pin right angle SpO₂ Patient Cable, 4 ft</p>	EA	
3345	<p>Red 25 LNC-10 Series LNCS Series 25-pin SpO₂ Patient Cable, 10 ft</p>	EA	
2013	<p>LNC-14 LNCS Series Patient Cable, 14 ft.</p>	EA	
2021	<p>LNC-4-Ext LNC Extension Cable, 4 ft.</p>	EA	

PN	Description	UM	Product Image
2281	<p>LNC MP10 LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ LNCS Dual Key Patient Cable, 10 ft.</p>	EA	
3755	<p>LNC MP4 LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ LNCS Dual Key Patient Cable, 4 ft</p>	EA	
2430	<p>LNC-4-GE LNCS to GE Patient Cable, 4 ft.</p>	EA	
2016	<p>LNC-10-GE LNCS to GE Patient Cable, 10 ft.</p>	EA	


PN	Description	UM	Product Image
2432	<p>LNC-SL-10 LNCS to Spacelabs Patient Cable, 10 ft.</p>	EA	
2285	<p>LNC-04-DB9 LNCS to DB9 Connector Patient Cable, 4 ft.</p>	EA	
2286	<p>LNC-10-DB9 LNCS to DB9 Connector Patient Cable, 10 ft.</p>	EA	
3984	<p>LNC-10-NK LNCS to Nihon Kohden Patient Cable, 10 ft.</p>	EA	



LNCS Adapter Cables


PN	Description	UM	Product Image
1816	<p>LNCS to PC Adapter Cable LNCS Series to LNOP PC Adapter Cable, 18 in.</p>	EA	
2263	<p>LNC CMS LNCS Series to Philips 12-pin Instrument Adapter Cable, 3 m.</p>	EA	
2264	<p>LNC MAC-GE LNCS Series to GE® Adapter Cable, 3 m.</p>	EA	
2266	<p>LNC MAC-SL LNCS Series to Spacelabs Adapter Cable, 3 m.</p>	EA	




PN	Description	UM	Product Image
2336	LNC MAC-SL2 LNCS Series to Spacelabs Adapter Cable, 3 m.	EA	 A 3-meter cable with a blue connector on one end and a white connector on the other.
2268	LNC MAC-180 LNCS Series to N-180 Adapter Cable, 10 ft.	EA	 A 10-foot cable with a white connector on one end and a larger white connector on the other.
2269	LNC MAC-395 LNCS Series to N-395 Adapter Cable, 10 ft.	EA	 A 10-foot cable with a white connector on one end and a larger white connector on the other.

M-LNC Patient Cables




PN	Description	UM	Product Image
2523	M-LNC-1 M-LNCS Series Patient Cable, 1 ft.	EA	 A 1-foot patient cable with a white connector on one end and a red connector on the other.





PN	Description	UM	Product Image
2524	<p>M-LNC-4 M-LNCS Series Patient Cable, 4 ft.</p>	EA	
2525	<p>M-LNC-10 M-LNCS Series Patient Cable, 10 ft.</p>	EA	
2526	<p>M-LNC-14 M-LNCS Series Patient Cable, 14 ft.</p>	EA	
2831	<p>M-LNC-MP10 M-LNCS to IntelliVue Masimo SET® or IntelliVue Philips FAST-SpO₂ M-LNC Dual Key Patient Cable, 10 ft.</p>	EA	


PN	Description	UM	Product Image
3754	<p>M-LNC-MP4 M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ M-LNC Dual Key Patient Cable, 4 ft.</p>	EA	 A 4-foot patient cable with a blue connector on one end and a white connector with a red key on the other.
2655	<p>M-LNC-4-GE M-LNCS to GE Patient Cable, 4 ft.</p>	EA	 A 4-foot patient cable with a blue connector on one end and a white connector with a red key on the other.
2533	<p>M-LNC-10-GE M-LNCS to GE Patient Cable, 10 ft.</p>	EA	 A 10-foot patient cable with a blue connector on one end and a white connector with a red key on the other.
2714	<p>M-LNC-SL-10 M-LNCS to Spacelabs Patient Cable, 10 ft.</p>	EA	 A 10-foot patient cable with a blue connector on one end and a white connector with a red key on the other.

PN	Description	UM	Product Image
2534	<p>M-LNC-04-DB9 M-LNCS Series to DB9 Connector Patient Cable, 4 ft.</p>	EA	
2535	<p>M-LNC-10-DB9 M-LNCS Series to DB9 Connector Patient Cable, 10 ft.</p>	EA	
3985	<p>M-LNC-10-NK M-LNCS to Nihon Kohden Patient Cable, 10 ft.</p>	EA	


M-LNCS Adapter Cables

PN	Description	UM	Product Image
2542	<p>M-LNC MAC-SL2 M-LNC Series to Spacelabs Adapter Cable, 3m.</p>	EA	
2541	<p>M-LNC MAC-SL M-LNC Series to Spacelabs Adapter Cable, 3 m.</p>	EA	
2537	<p>M-LNC MAC-180 M-LNCS Series to N-180 Adapter Cable, 10 ft.</p>	EA	


PN	Description	UM	Product Image
2538	<p>M-LNC MAC-395 M-LNCS Series to N-395 Adapter Cable, 10 ft.</p>	EA	
2540	<p>M-LNC MAC-GE M-LNCS Series to GE Adapter Cable, 3 m.</p>	EA	
2539	<p>M-LNC MAC CMS M-LNCS Series to Phillips 12-pin Instruments Adapter Cable, 3 m.</p>	EA	
2654	<p>M-LNCS to LNC Adapter Cable M-LNCS Series to LNC Patient Cable, 1.5 ft.</p>	EA	




PN	Description	UM	Product Image
2750	<p>M-LNCS to LNC-MAC Adapter Cable M-LNCS Series to LNC-MAC Adapter Cable, 1.5 ft.</p>	EA	

Direct Connect Sensors (for use with Masimo SET® Devices)




PN	Description	UM	Product Image
1396	<p>DCSC Adult Reusable Direct Connect SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	

Red Direct Connect Sensors (for Masimo rainbow SET Devices)

PN	Description	UM	Product Image
2053	<p>Red DCI-dc3 Adult Reusable Direct Connect SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30kg Not made with natural rubber latex.</p>	EA	





PN	Description	UM	Product Image
2256	<p>Red DCIP-dc3 Pediatric Reusable Direct Connect SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight 10kg – 50kg Not made with natural rubber latex</p>	EA	
2643	<p>Red DBI-dc3 Adult Soft Reusable Direct Connect Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2644	<p>Red DBI-dc8 Adult Soft Reusable Direct Connect Sensor, 8 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	




Mountable Patient Cable Extensions and Adapter

PN	Description	UM	Product Image
2478	<p>14-Pin to 20-Pin Adapter Cable Connect PXC2 to Masimo SET® 20-pin devices, SpO₂ only.</p>	EA	
2371	<p>PXC2-16 Mountable Extension Cable, 16 ft. For use with Masimo RD SET MD14, LNC, M-LNC Patient Cables and DCSC, DBI-dc3 Direct Connect sensors or Masimo 14-pin to 20-pin Adapter Cables.</p>	EA	
2373	<p>PXC2-30 Mountable Extension Cable, 30 ft. For use with Masimo RD SET MD14, LNC, M-LNC Patient Cables and DCSC, DBI-dc3 Direct Connect sensors or Masimo 14-pin to 20-pin Adapter Cables.</p>	EA	




Accessories


Training Kits

PN	Description	UM	Product Image
4286	<p>RD rainbow SET-2 Training Kit (SpHb, SpMet, SpO₂)</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 RD rainbow SET-2 Adt • 1 RD rainbow SET-2 Pdt • 1 RD rainbow SET-2 Inf • 1 RD rainbow SET-2 Neo • Sensor application and troubleshooting card • Directions for use 	Box	
4303	<p>RD rainbow Lite SET-1 Training Kit (SpO₂ and ORi)</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 RD rainbow Lite SET-1 Adt • 1 RD rainbow Lite SET-1 Pdt • 1 RD rainbow Lite SET-1 Inf • 1 RD rainbow Lite SET-1 Neo • Sensor application and troubleshooting card • Directions for use 	Box	
4201	<p>RD SET Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 RD SET Adt • 1 RD SET Pdt • 1 RD SET Inf • 1 RD SET Neo • 1 RD SET NeoPt • 1 neonatal foot with application card • 1 RD SET sensor application card • Directions for use 	EA	
4617	<p>RD SET CS Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 RD SET Adt CS-2 • 1 RD SET Neo CS-2 • 1 RD SET Adt CS-3 • 1 RD SET Neo CS-3 • Directions for Use 	EA	



PN	Description	UM	Product Image
2345	<p>LNCS® Sensor Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNCS Amtx • 1 LNCS Pmtx • 1 LNCS Inf • 1 LNCS Neo • 1 LNCS NeoPt • 1 neonatal foot with application card • Directions for use 	EA	
4287	<p>RAS Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 2 RAS-125c • 2 RAS-45 • 1 Sensor Application and Troubleshooting card • Directions for use 	EA	
4533	<p>RAS-45 Inf/Neo Sensor Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 4 RAS-45 Inf/Neo Sensors • Directions for use <p>For US markets Only</p>	EA	No image
4381	<p>NomoLine Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 2 - 3827 NomoLine HH Adult/Pediatric Airway Adapter Set • 2 - 3837 NomoLine HH Adult Nasal/Oral CO2 Cannula with O2 • 2 - 3820 NomoLine LH Adult Nasal CO2 Cannula with O2 • 2 - 3814 NomoLine LH Adult/Pediatric Airway Adapter Set • Directions for Use 	EA	

Demo Tools


PN	Description	UM	Product Image
2389	<p>Mini-Demo Connector - 30 Min Used to demo product (best used to display trending information). Includes RAM waveform</p>	EA	
4186	<p>Mini-Demo Connector ORi - 30 Min Used to demo ORi (best used to display trending information). Not available for the US. For approved markets only.</p>	EA	
4266	<p>M20 Demo Connector w/ORi - 30 Min Used to demo ORi (best used to display trending information) Not available for the US. For approved markets only.</p>	EA	
4265	<p>M20 Demo Connector - 30 Min Used to demo product (best used to display trending information)</p>	EA	
4352	<p>25-Pin rainbow Demo Tool - 30 Min Used to demo product (best used to display trending information).</p>	EA	







PN	Description	UM	Product Image
4353	<p>25-Pin rainbow Demo Tool - 2 Min Used to demo product (best used for short demonstration).</p>	EA	






Software







PN	Description	UM	Product Image
9899	<p>Masimo Sketch™ Configuration Software Sketch is a network configuration software for Masimo Devices. Compatible with:</p> <ul style="list-style-type: none"> Radical-7 (with or without RDS) Root and Root NIBPT (with Radical-7) Rad-97 product family 	EA	
9853	<p>Masimo Trace™ Data Analysis and Reporting Tool (USB) Trace is a data download and reporting tool for Masimo devices. Compatible with:</p> <ul style="list-style-type: none"> Radical-7 (with or without RDS) Root and Root NIBPT (with Radical-7, Radius-7 and connected MOC-9 modules) Rad-97 product family <p>Licensing:</p> <ul style="list-style-type: none"> Requires license activation License duration is 1 year 	EA	
4330	<p>Masimo Trace Renewal License</p> <ul style="list-style-type: none"> Renewal license duration is 1 year 	EA	No image




Bedside Device Accessories

PN	Description	UM	Product Image
1315	<p>Radical Series Replacement Battery For use with Radical and Radical-7 Blue Screen Handheld.</p>	EA	

PN	Description	UM	Product Image
2386	<p>Radical-7 Color Screen Replacement Battery For use with Radical-7 Color Screen only.</p>	EA	
23794	<p>Radical-7 Touchscreen Replacement Battery, Rechargeable Lithium Ion</p>	EA	
1395	<p>Radical Series Handheld Lock To secure Radical and Radical-7 Blue Screen Handheld in Radical Docking Stations.</p>	EA	
30534	<p>Radical Series Handheld Lock Key To unlock Radical and Radical-7 Blue Screen Handhelds from Docking Stations.</p>	EA	
1584	<p>Radical Series Power Cord Lock To secure a power cord to Radical Docking Stations 5 / pack</p>	EA	
2228	<p>Rad-8 Power Cord Lock</p>	EA	






PN	Description	UM	Product Image
4249	<p>RD SET Tester Tester for devices containing Masimo SET or Masimo rainbow SET Connects to RD SET or RD rainbow SET patient cables Testers for SpO₂ and PR only</p>	EA	
4275	<p>Masimo rainbow M20 Tester Tester for devices containing Masimo SET or Masimo rainbow SET, M20 instrument connector Tests SpO₂, PR, SpCO, SpMet and SpHb</p>	EA	
3776	<p>Masimo SET® M15 Tester Tester for devices containing Masimo SET and connects to M-LNC or rainbow® RC patient cables Testers SpO₂ and PR only</p>	EA	
2367	<p>9-Pin Masimo SET Tester Tester for devices containing Masimo SET and connects to LNC patient cables Tests SpO₂ and PR only.</p>	EA	
1795	<p>Masimo SET LNOP Tester For Oximeter & Patient Cable SpO₂ and PR Testing. Tester for devices containing Masimo SET and connects to LNOP PC patient cables. Tests SpO₂ and PR only.</p>	EA	

PN	Description	UM	Product Image
2368	<p>Masimo rainbow SET Tester</p> <p>Tester for devices containing Masimo SET® or Masimo rainbow SET. Tests SpHb, SpCO, SpMet, PVi, Pi, SpO₂ and BPM.</p>	EA	
3502	<p>Masimo rainbow SET Tester</p> <p>Tester for devices containing Masimo SET or Masimo rainbow SET, 25-pin instrument connector. Tests SpO₂, PR, SpCO, SpMet and SpHb.</p>	EA	
1593	<p>14-Pin Masimo SET Tester</p> <p>For Oximeter SpO₂ and PR Testing</p> <p>Tester for devices containing Masimo SET, 14-pin instrument connector. Tests SpO₂ and PR only.</p>	EA	
2350	<p>Radical-7 Color Handheld Lock</p> <p>Locks Radical-7 Color Handheld to the Docking Station. Used with PN 2351.</p>	EA	
2351	<p>Radical-7 Color Handheld Lock Key</p> <p>Used for installation of Radical-7 Color Handheld Lock and unlocking handheld from docking station. Used with PN 2350</p>	EA	
38879	<p>Arm Band, Radius-7, Small, White</p> <p>20 / box</p>	Box	

PN	Description	UM	Product Image
38880	Arm Band, Radius-7, Medium, White 20 / box	Box	
38881	Arm Band, Radius-7, Large, White 20 / box	Box	
19118	Barcode Reader, Black	EA	
19356	Locking Power Cord, US, 10A, 125V, RoHS	EA	
18313	USB-Null Modem Cable For connecting Root and Rad-97 to Windows-based PC for downloading trend data with MICT.	EA	
17417	USB-to-Serial Adapter	EA	









Handheld Device Accessories



PN	Description	UM	Product Image
2208	<p>Water Resistant Protective Carrying Case for Battery Operated Portable Handheld, Red</p> <p>This can be used with Rad-57 and Rad-5/5v Includes rubber label parameter patches</p>	EA	
10368	<p>Radical-7 Ambulation Pouch</p> <p>Pack of 10</p>	EA	No image
1909	<p>Data Transfer Download Cable</p> <p>Used for connecting Rad-5 (with 14-pin connector) to Personal Computer to enable download of trending data.</p>	EA	
2063	<p>Data Transfer Download Cable, Red – 20-Pin</p> <p>Data Transfer Download Cable, 20-pin, for Rad-5, Rad-57 and Pronto Trend Data Download.</p>	EA	
4361	<p>Data Transfer Download Cable, Rad-67</p> <p>Connects to Rad-67 power charging connector (Masimo Proprietary Port) to transfer data (download only) to computer storage via USB connector, and to update software to the device via USB connector.</p>	EA	








PN	Description	UM	Product Image
1908	<p>TrendCom, Trend Download Software</p> <p>Used with appropriate Data Transfer Download cable to initiate download and de-compression of Rad-5 and Rad-57 trending data to a Personal Computer. Used with serial cable to initiate download and decompression of Radical-7 and Rad-8 trending data.</p>	EA	
1842	<p>Handheld Boot, Gray</p> <p>For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds</p>	EA	
1980	<p>Handheld Boot, Yellow</p> <p>For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds</p>	EA	
1981	<p>Handheld Boot, Red</p> <p>For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds</p>	EA	
1982	<p>Handheld Boot, Orange</p> <p>For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds</p>	EA	





PN	Description	UM	Product Image
2097	<p>Handheld Boot, Royal Blue For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds</p>	EA	
2098	<p>Handheld Boot, Light Blue For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds.</p>	EA	
2099	<p>Handheld Boot, Pink For use with Rad-5, Rad-5v, Pronto and all Rad-57 Handhelds</p>	EA	
2357	<p>Light Shield For use in strobing environment and high sunlight intensity areas. 5 / Pack</p>	EA	
2815	<p>Masimo Ambient Shield Accessory For use with rainbow® Sensors. 10 / pack</p>	EA	

Sensor Accessories

PN	Description	UM	Product Image
4439	RD rainbow SET Inf Replacement Tapes For use with RD rainbow SET Inf sensors 50 / box, non-sterile Not made with natural rubber latex	EA	
4440	RD rainbow SET Neo Replacement Tapes For use with RD rainbow SET Neo sensors 50 / box, non-sterile Not made with natural rubber latex	EA	
2624	Replacement Tapes rainbow® R1 20L and R20-L For use with rainbow R1 20L and R20-L 50 / box, non-sterile Not made with natural rubber latex	Box	
2623	Replacement Tapes for rainbow R1 25L and R25-L For use with rainbow R1 25L and R25-L 50 / box, non-sterile Not made with natural rubber latex	Box	
4046	RD SET Inf Replacement Tapes For use with RD SET Inf/Inf CS-2/Inf CS-3, RD SET Blue, and RD rainbow Lite SET Inf 102 / box, non-sterile Not made with natural rubber latex	Box	
4047	RD SET Neo Replacement Tapes For use with RD SET Neo/Neo CS-2/Neo CS-3, and RD rainbow Lite SET Neo 102 / box, non-sterile Not made with natural rubber latex	Box	
4048	RD SET SofTouch™ Replacement Wraps For use with RD SET NeoPt/NeoPt CS-2/NeoPt CS-3, RD SET Trauma, RD SET Newborn Neonatal 10 / pack, non-sterile Not made with natural rubber latex	Pack	
4094	RD SET NeoPt-500 SofTouch Replacement Wraps For use with RD SET NeoPt-500 and Newborn Infant/Pediatric 10 / pack, non-sterile Not made with natural rubber latex	Pack	






PN	Description	UM	Product Image
2307	<p>Replacement Tapes for Inf Series Sensors For use with LNCS Inf/Inf-3 102 / box, non-sterile Not made with natural rubber latex</p>	Box	
2308	<p>Replacement Tapes for Neo Series Sensors For use with LNCS Neo/Neo-3 102 / box, non-sterile Not made with natural rubber latex</p>	Box	
2309	<p>Replacement SofTouch Wraps for NeoPt, NeoPt-3, Trauma and Newborn Neonatal Sensors For use with LNCS NeoPt/NeoPt-3, Trauma and Newborn Neonatal 10 / pack, non-sterile Not made with natural rubber latex</p>	Pack	
2322	<p>Replacement SofTouch Wraps for NeoPt-500 and Newborn Infant / Pediatric Sensors For use with LNCS NeoPt-500 and LNCS Newborn Infant/Pediatric 10 / pack, non-sterile Not made with natural rubber latex</p>	Pack	
1596	<p>CleanShield® Multisite Wrap For use with YI Sensors Adult / Pediatric / Neonatal 100 / box, non-sterile Not made with natural rubber latex</p>	Box	
1597	<p>Standard Multisite Wrap For use with YI Sensors Adult / Pediatric / Neonatal 100 / box, non-sterile Not made with natural rubber latex</p>	Box	
1598	<p>Standard Petite Wrap For use with YI Sensors Adult slender digit / Pediatric 100 / box, non-sterile Not made with natural rubber latex</p>	Box	







PN	Description	UM	Product Image
4317	<p>Adhesive Tabs</p> <p>For use with Disposable Adhesive Sensors Adult/ Pediatric/ Infant/ Neonatal/ Neo-Pt 120 / pack, Non Sterile</p>	Pack	
1600	<p>Accessory Adhesive Squares</p> <p>For use with YI Sensors 144 / pack, non-sterile</p>	Pack	
1903	<p>Transflectance Sensor Adhesive Pad</p> <p>For use with TF-I Sensors 20 / pack, non-sterile</p>	Pack	
16179	<p>Tele Plast</p> <p>Philips SpO₂ Connector Plug Compatible with Philips MMS module</p>	EA	
1602	<p>Foam Wrap</p> <p>For use with YI Sensors 12 / pack, non-sterile Not made with natural rubber latex</p>	Pack	
1601	<p>Accessory Clothing Clips</p> <p>For securing reusable sensor cables to patient clothing 5 / pack, non-sterile</p>	Pack	
3444	<p>Masimo Soft Wrap Armband</p> <p>Used to secure the sensor cable along the patients arm. 20 / pack, non-sterile Not made with natural rubber latex</p>	Pack	







PN	Description	UM	Product Image
3929	<p>Masimo Forehead Sensor Headbands For use with forehead sensors 20 / box, non-sterile Not made with natural rubber latex</p>	Box	
2215	<p>Masimo Headband For use with TF-I Sensors 12 / box, non-sterile Not made with natural rubber latex</p>	Box	
21514	<p>LNC Replacement Cover For use with all LNC patient cables 25 covers and 50 pins</p>	Pack	
4507	<p>RD Reusable Sensor Replacement Latch For use with RD reusable RD sensor 25 latches / bag</p>	Bag	

SatShare® Interface Cables

Note: Do not interface the Radical-7 with any monitor or module not specifically listed in this chart.


PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1324	D001	<p>4 ft. SatShare Cable for use with: Datex-Ohmeda</p> <ul style="list-style-type: none"> AS/3 Patient Monitor (M-NSAT) System 5 Compact Monitor (M-NSAT) 	EA	
1326	D002	<p>4 ft. SatShare Cable for use with: Datex-Ohmeda</p> <ul style="list-style-type: none"> AS/3 Patient Monitor (M-NESTPR-00-00 / M-NESTR-00-00) CardioCap / 5 Monitor Light Monitor System 5 Compact Monitor (M-NESTPR-00-00 / M-NESTR-00-00) 	EA	
2223	D006	<p>4 ft. SatShare Cable for use with: Datex-Ohmeda</p> <ul style="list-style-type: none"> System 5 Compact Monitor (M-PRESTN..01) 	EA	
2969	D007	<p>4 ft. SatShare Cable for use with: Datex-Ohmeda</p> <ul style="list-style-type: none"> System 5 Anesthesia Monitor (E-PRESTN) 	EA	
1325	DS01	<p>4 ft. SatShare Cable for use with:</p> <p>Corometrics</p> <ul style="list-style-type: none"> Maternal / Fetal Monitor Model 118 Maternal / Fetal Monitor Model 120 <p>Critikon</p> <ul style="list-style-type: none"> Dinamap Monitor (117310) <p>Datascope</p> <ul style="list-style-type: none"> Expert Monitor with Nellcor SpO₂ (HS-500) Accutor Monitor with Nellcor SpO₂ Passport Monitor with Nellcor SpO₂ <p>North American Dräger</p> <ul style="list-style-type: none"> Narcomed 4 Anesthesia Machine (MP202) Vitalert 3200 Monitor (203) <p>Medical Data Electronics</p> <ul style="list-style-type: none"> Escort Monitor <p>Witt Biomedical</p> <ul style="list-style-type: none"> Frontend05 Patient Monitor 	EA	

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1324	D001	4 ft. SatShare Cable for use with: Datex-Ohmeda <ul style="list-style-type: none"> AS/3 Patient Monitor (M-NSAT) System 5 Compact Monitor (M-NSAT) 	EA	
1365	DS03	4 ft. SatShare Cable for use with: Dynascope <ul style="list-style-type: none"> Fukuda Denshi DS-5300W (HS-500) 	EA	
2979	DS04	4 ft. SatShare Cable for use with: Dynascope <ul style="list-style-type: none"> Fukuda Denshi DS-5400 (HS-500) 	EA	
1357	HP03	5 ft SatShare Cable for use with: Philips / Agilent / HP <ul style="list-style-type: none"> Omnicare / Viridia Model 24 / 26 Series Monitor System (M1204A) Series 50XM Fetal Monitor 	EA	
1321	MQ01	4 ft. SatShare Cable for use with: GE / Marquette <ul style="list-style-type: none"> Solar 7000 / 8000 Series Monitor (Tram 250SL, Tram 450SL, Tram 650SL, Tram 850A&SL, and SpO₂ Solar Parameter Module) Dash3000 Monitor 	EA	
1789	MS01	2.5 ft. SatShare Cable for use with: Datex-Ohmeda <ul style="list-style-type: none"> S5 (E-Masimo-00) Dräger <ul style="list-style-type: none"> Infinity Delta (SmartPOD with Masimo SET) GE / Marquette <ul style="list-style-type: none"> Solar 8000 (Tram 451M, Tram 851M, Masimo SET single parameter module) Philips <ul style="list-style-type: none"> IntelliVue MP40 Monitor (MMS Module [Masimo SET]) GE <ul style="list-style-type: none"> Carescape B850 (Patient Data Module [Masimo SET]) 	EA	

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1324	DO01	<p>4 ft. SatShare Cable for use with: Datex-Ohmeda</p> <ul style="list-style-type: none"> AS/3 Patient Monitor (M-NSAT) System 5 Compact Monitor (M-NSAT) 	EA	
1368	NK02	<p>3 ft. SatShare Cable for use with: Nihon Kohden</p> <ul style="list-style-type: none"> BSM-2301 (Type 9266) BSM-5105 (Type 9266) 	EA	
1840	PM01	<p>4 ft. SatShare Cable for use with: Philips Medical</p> <ul style="list-style-type: none"> M3/M4 (MMS Module) Intellivue Monitor (MMS Module) 	EA	
1323	SL01	<p>4 ft. SatShare Cable for use with: SpaceLabs Medical</p> <ul style="list-style-type: none"> PC Bedside Monitor (90303) [90465] PC2 Central Monitor (90315, 90316) [90465] <p>SpaceLabs Medical</p> <ul style="list-style-type: none"> Adult/Neonatal Monitor (90621) 	EA	
1362	SL02	<p>4 ft. SatShare Cable for use with: SpaceLabs Medical</p> <ul style="list-style-type: none"> Ultraview 1700 Bedside Monitor (90387) [90496] Ultraview UCW Monitor (90385) [90496] 	EA	
1322	SM01	<p>4 ft. SatShare Cable for use with: Siemens Medical Systems</p> <ul style="list-style-type: none"> Sirecust 1281 Monitor(99731472E2529) SC 6002XL, SC7000 	EA	

Note: Do not interface the Radical-7 with any monitor or module not specifically listed in this chart.

SatShare Extension Cable

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)*	UM	Host Connector
1605		SatShare Extension Cable, 4 ft.	EA	

Index

A

Accessories • 113

B

Bedside Device Accessories • 116

Blood Pressure Hoses and Adapters • 27

D

Demo Tools • 115

Direct Connect Sensors (for use with Masimo SET® Devices) • 110

Disposable Blood Pressure Cuffs • 22

Doctella™ • 39

E

EMMA™ • 14

G

GE® Healthcare Compatible Devices • 44

H

Handheld Device Accessories • 121

I

ISA™ CO2 • 13

iSpO2® Rx • 47

K

Kite® and UniView™ Supplemental Display • 40

L

LNC Patient Cables (for use with Masimo SET® Devices) • 98

LNCS Adapter Cables • 103

LNCS Red Patient Cables (for use with Masimo rainbow SET Devices) • 96

LNCS Reusable Sensors • 95

LNCS Single Patient Use Sensors with Replaceable Tapes • 89

LNCS Specialty Sensors • 91

LNCS Sterile Sensors • 93

LNCS® Adhesive Sensors • 86

M

Masimo Patient SafetyNet™ System • 32

Masimo rainbow SET™ Parameter Options • 51

Masimo rainbow® Pulse CO-Oximeter® Devices • 5

Masimo SET® Pulse Oximeter Devices • 29

MightySat™ Rx • 45

M-LNC Patient Cables • 104

M-LNCS Adapter Cables • 108

Monitoring Devices • 5

Mountable Patient Cable Extensions and Adapter • 112

N

NomoLine Sampling Lines • 15

NomoLine® Capnography • 13

Noninvasive Blood Pressure (NIBP) • 22

O

O3® Regional Oximeter Monitoring • 11

P

Parameter Options and Upgrade Tools • 51

Philips® Compatible Devices • 42

R

rainbow Acoustic Monitoring® Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRA®) • 71

rainbow® Adhesive Sensors (for rainbow SET Devices with SpO2, SpCO®, and SpMet) • 58

rainbow® Adhesive Sensors R1 series (for rainbow SET Devices with SpHb, SpO2, and SpMet) • 57

rainbow® Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet) • 62

rainbow® Patient Cables (for rainbow SET Devices) • 63

rainbow® Reusable Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet) • 61

RD rainbow Lite SET™ Adhesive Sensors (for rainbow SET Devices with SpO2, and ORi) • 56

RD rainbow SET™ Adhesive Sensors (for rainbow SET Devices with SpHb®, SpO2, and SpMet®) • 55

RD SET Adapter Cables • 85

RD SET Patient Cables • 81

RD SET Reusable Sensors • 79

RD SET Specialty Sensors • 78

RD SET™ Cabled Sensors • 75

RD SET™ Disposable Sensors • 74

Red Direct Connect Sensors (for Masimo rainbow SET Devices) • 110

Replica™ • 38

Reusable Blood Pressure Cuffs • 25

S

SatShare Extension Cable • 132

SatShare® Interface Cables • 129

SedLine® Brain Function Monitoring • 9

Sensor Accessories • 124

Sensors and Cables • 55

Software • 116

rainbow® Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO2, SpMet) • 59

T

Temperature • 28

Training Kits • 113



PLCO-003021/PLM-10062H 0719