

Liebel-Flarsheim Single Use Empty Syringes

Syringe Components Sheet

P/N: 600172
Size: 130 ml
Injectors: Liebel-Flarsheim Angiomat 6000
Medrad Mark II, III, IV
Medtron 82M/HP
Liebel-Flarsheim Optistat



Contents: Syringe with Handifill
Package: 50 units/box
Pressure Limit Syringe : 1200 PSI / 84 bar
Connectors: Linden Luer

Linden-Luer Features:
- Easy to use connection
- Reduced thread pitch
- Reduced area for contrast to collect

Materials: Polypropylene
Polycarbonate
Natural Rubber
Polyethylene
Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation

Use: Single Use – do not reuse

CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237

Country of Origin: USA



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany