

Liebel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

P/N: 601195
Lengths: 150 cm
Description: Spiral Line for CT



Connectors: LL f/m
Package: 100 units/box
Pressure Limit Syringe: 300 PSI / 21 bar

Materials: Polycarbonate
PVC
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer: Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237
Country of Origin: USA



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1-3
53761 Hennef
Germany