

Liebel-Flarsheim Single Use Empty Syringes

Syringe Components Sheet

P/N: 800099
Size: 200 ml
Injectors: Liebel-Flarsheim CT 9000
Liebel-Flarsheim CT 9000 ADV



Contents: Multipack with Syringe,
Spiral Tube and Handifill
Package: 50 units/box
Pressure Limit Syringe : 1200 PSI / 84 bar
Pressure Limit Tube : 300 PSI / 21 bar
Connectors: Linden Luer
Materials: Polypropylene
Polycarbonate
Synthetic Rubber
Polyvinyl Chloride
Latex free

Linden-Luer Features:
- Easy to use connection
- Reduced thread pitch
- Reduced area for
Contrast to collect

Front Loading Syringe
- for unsurpassed ease of
use in CT

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123
Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237
Country of Origin: USA



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

