

Liebel-Flarsheim Single Use Empty Syringes

Syringe Components Sheet

P/N: 801801
Size: 60 ml
Injectors: Liebel-Flarsheim Optistar LE



Contents: Syringe with Spike
Package: 50 units/box
Pressure Limit Syringe : 200 PSI / 14 bar

Materials: Polypropylene
Synthetic Rubber
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 2 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237
Country of Origin: USA

tyco
Healthcare

MALLINCKRODT

Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany