

Mallinckrodt Single Use Empty Syringes

Syringe Components Sheet

P/N: 844016
Size: 2 x 200 ml Syringes with Handifils and 150 cm coiled Y-Tube with Single Check Valve
Injectors: Mallinckrodt OptiVantage DH



Contents: 2 Syringes with Handifil and 150 cm coiled Y-tube with Single Check Valve

Package: 20 units/box

Pressure Limit: Syringes: 400 PSI / 28 bar
Tube: 400 PSI / 28 bar

Connectors: Linden Luer

Materials: Polypropylene
Polycarbonate
Synthetic Rubber, PVC
Silicone

Latex free

Sterilisation: Ethylene Oxide (EO)

Shelf Life: 3 years from Date of Sterilisation

Use: Single Use – do not reuse

CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237

Country of Origin: USA

Linden-Luer Features:

- Easy to use connection
- Reduced thread pitch
- Reduced area for Contrast to collect

Front Loading Syringe

- for unsurpassed ease of use in CT



Tyco Healthcare Deutschland GmbH
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

