

Mallinckrodt Single Use Empty Syringes

Syringe Components Sheet

P/N: 844020
Size: 200 ml Syringe with Handifil and 150 cm coiled Y-Tube
Injectors: Mallinckrodt OptiVantage DH




Contents: Syringe with Handifil and 150 cm coiled Y-tube
Package: 50 units/box
Pressure Limit: Syringe: 400 PSI / 28 bar
 Tube: 400 PSI / 28 bar

Connectors: Linden Luer
Materials: Polypropylene
 Polycarbonate
 Synthetic Rubber
 PVC

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123
Manufacturer : Liebel-Flarsheim Company
 2111 East Galbraith Road, Cincinnati, Ohio 45237
Country of Origin: USA

Linden-Luer Features:
 - Easy to use connection
 - Reduced thread pitch
 - Reduced area for Contrast to collect

Front Loading Syringe
 - for unsurpassed ease of use in CT



Tyco Healthcare Deutschland GmbH
 Josef-Dietzgen-Str. 1
 53761 Hennef
 Germany

