

Liebel-Flarsheim Single Use Empty Syringes

Syringe Components Sheet

P/N: 900102
Size: 150 ml
Injectors: Liebel-Flarsheim Angiomat Illumena



<u>Contents:</u>	Multipack with Syringe, Spiral Tube and Handfill	<u>Syringe Feature:</u> Ultra clear syringe barrel
<u>Package:</u>	50 units/box	
<u>Pressure Limit Syringe :</u>	1200 PSI / 84 bar	
<u>Pressure Limit Tube :</u>	300 PSI / 21 bar	
<u>Connectors:</u>	Quick Connection	
<u>Materials:</u>	Copolyester Polycarbonat Synthetic Rubber Latex free	
<u>Sterilisation:</u>	Ethylene Oxide (EO)	
<u>Shelf Life:</u>	3 years from Date of Sterilisation	
<u>Use:</u>	Single Use – do not reuse	
<u>CE:</u>	0123	
<u>Manufacturer :</u>	Liebel-Flarsheim Company 2111 East Galbraith Road, Cincinnati, Ohio 45237	
<u>Country of Origin:</u>	USA	



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

