

Liebel-Flarsheim Single Use Empty Syringes

Syringe Components Sheet

P/N: 900104
Size: 150 ml
Injectors: Liebel-Flarsheim Angiomat Illumena



Contents: Multipack with Syringe,
Spiral Tube and Handfill
Package: 50 units/box
Pressure Limit Syringe : 1200 PSI / 84 bar
Pressure Limit Tube : 300 PSI / 21 bar
Connectors: Linden Luer
Materials: Copolyester
Polycarbonat
Synthetic Rubber
Latex free

Linden-Luer Features:
- Easy to use connection
- Reduced thread pitch
- Reduced area for
contrast to collect

Syringe Feature:
Ultra clear syringe barrel

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237
Country of Origin: USA

tyco
Healthcare

MALLINCKRODT

Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

