

Liebel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

P/N: HPxxRA
Lengths: 50 cm / 75 cm / 120 cm
Description: Standard High Pressure Line

**with rotating adapter
for Angiography**



Connectors: Female Luer with Rotating Male Adapter
Package: 25 units/box
Pressure Limit Syringe : 1200 PSI / 84 bar

Materials: PVC
Polycarbonate
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 5 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer : Medtron AG,
Hauptstraße 255, 66128 Saarbrücken
Country of Origin: Germany

tyco
Healthcare

MALLINCKRODT

Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1-3
53761 Hennef
Germany

