

Finale[®]

RADIAL COMPRESSION DEVICE



*Compression at
the turn of a dial*

Finale[®]

Merit Medical - your partner for radial procedures



Clear dial provides easy monitoring of access site.



Ergonomic dial: one hand to increase compression and two hands to decrease compression, preventing un-intended pressure release.



Comfortable Velcro[®] band. Available in two lengths

Catalog Number	Band Length
FCD100	10" (25.4cm)
FCD200	12" (30.5cm)

Packaged 10 per box



PRELUDE[®] SHEATH INTRODUCERS

Excellent transitions facilitate smooth entry into patient and vessel. Resists kinking during insertion.



DIAGNOSTIC CARDIOLOGY CATHETERS

Performa[®] Ultimate diagnostic angiography catheters for procedures performed using the radial approach.



Velcro[®] is a registered trademark of Velcro Industries
US & Foreign Patents Pending.

Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24 • France 0800 91 60 30

Germany 0800 182 0871 • Ireland (republic) 1800 553 163 • Luxembourg 8002 25 22 • Netherlands 0800 022 81 84 • Sweden 020 792 445 • UK 0800 973 115