



**NAVIGATE.
CAPTURE.
SUCCEED.**



EN Snare®

Endovascular Snare System

ONE Snare®

Endovascular Snare System

EN Snare[®]

Endovascular Snare System

The only delivery catheter with **reinforced strain-resistant hub and peel-away insertion tool** on the market.

EN Snare[®] loop tips stay as much as **58% closer** to vessel walls than the competition.¹

Embedded Radiopaque Marker Band

identifies the snare delivery catheter's precise location under imaging

Interlaced Nitinol Loops

designed to resist prolapse and improve capture by allowing the snare to move as a single-unit through the vessel

Improved Vessel Navigation


facilitated by 15 degree curved tip on 6F and 7F snare delivery catheters

Enhanced Visualization

provided by platinum strands incorporated into the loops

High Performance Snare Delivery Catheter

designed for strength and to minimize the potential for kinking and buckling



REINFORCED CATHETER HUB

NEW flexible, reinforced, strain relief hub designed to reduce buckling.

ONE Snare[®]
2.4X more buckle force¹


EN Snare[®]
1.9X more buckle force²



ONE Snare®

Endovascular Snare System

ONE Snare® delivery catheter provides **2.4X more buckle force** than the competition.¹



Designed for Better Foreign Body Capture with 90 degree angle loop


Excellent Support from Nitinol core wire that provides shape-memory and superelasticity

Kink Resistant Catheter designed to improve strength and ensure procedural success

Nitinol & Gold-Plated Tungsten Loop retains snare's circular shape to support capture of foreign bodies, while allowing superior visibility.




High Visibility under fluoroscopy due to gold-plated tungsten loop



REDUCED PROCEDURAL TIMES

UNIQUE peel-away insertion tool for quick and efficient snare loading with EN Snare® and ONE Snare.®

ONE Snare
Peel-Away insertion tool



Competitor's non-peelable introducer

EN Snare

EN Snare®

Endovascular Snare System



EN SNARE ORDERING INFORMATION

Catalog Number	Loop Diameter	Snare Length	Catheter Size	Catheter Length
Standard Snare Kit				
EN2006010	6-10 mm	120 cm	6F	100 cm
EN2006015	9-15 mm	120 cm	6F	100 cm
EN2006020	12-20 mm	120 cm	6F	100 cm
EN2007030	18-30 mm	120 cm	7F	100 cm
EN2007045	27-45 mm	120 cm	7F	100 cm
Mini Snare Kit				
EN1003004	2-4 mm	175 cm	3.2F	150 cm
EN1003008	4-8 mm	175 cm	3.2F	150 cm

ONE Snare®

Endovascular Snare System



ONE SNARE ORDERING INFORMATION

Catalog Number	Loop Diameter	Snare Length	Catheter Size	Catheter Length
Standard Snare Kit				
ONE500	5 mm	120 cm	4F	100 cm
ONE1000	10 mm	120 cm	4F	100 cm
ONE1500	15 mm	120 cm	6F	100 cm
ONE2000	20 mm	120 cm	6F	100 cm
ONE2500	25 mm	120 cm	6F	100 cm
ONE3000	30 mm	120 cm	6F	100 cm
ONE3500	35 mm	120 cm	6F	100 cm
Petite Snare Kit				
ONE1001	10 mm	65 cm	4F	48 cm
ONE2501	25 mm	65 cm	6F	48 cm
Microsnare Kit				
ONE200	2 mm	175 cm	2.3-3F	150 cm
ONE201	2 mm	200 cm	2.3-3F	175 cm
ONE400	4 mm	175 cm	2.3-3F	150 cm
ONE401	4 mm	200 cm	2.3-3F	175 cm
ONE700	7 mm	175 cm	2.3-3F	150 cm
ONE701	7 mm	200 cm	2.3-3F	175 cm
Snare Replacement Catheters				
ONE4000	-	-	4F	100 cm
ONE6000	-	-	6F	100 cm

Footnotes:

¹ Data on file, Merit 2012.

² Data on file, Merit 2014.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115