



## Accessories and Consumables

CATALOGUE 2018.04

Mindray Building, Keji 12th Road South,  
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China  
Tel: +86 755 8188 8998 Fax: +86 755 26582680  
E-mail: PMLS-MarCom@mindray.com [www.mindray.com](http://www.mindray.com)

**mindray** is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.  
© 2017 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.  
P/N: ENG-Accessories Catalogue-210210x132P-20180417





## Welcome to the Mindray Accessories Catalogue

This catalogue will provide you with the parts and accessories that connect to your Mindray Patient Monitor, Electrocardiograph, Defibrillator. Each Mindray product is the product of a special brand of patient focused, clinician-friendly design. For this reason, you can expect the same service, focus and quality with our parts and accessories.

## Finding the Right Part

This catalog has been designed to make finding the right part easy. Chapters are organized by specific parameter categories. Simply locate the type of part you are looking for under the appropriate category.

### Note:

This catalog is not an Operating Instructions Manual. This catalog will assist you in identifying the correct parts and accessories to connect to your Mindray product, please refer to the Operating Instructions Manual. Warnings, Precautions and Notes can also be found in the Operating Instructions.

## CONTENTS

01

### ● Patient Monitor Accessories

|                                      |     |
|--------------------------------------|-----|
| ECG                                  | 2   |
| SpO <sub>2</sub>                     | 30  |
| NIBP                                 | 44  |
| TEMP                                 | 54  |
| CO <sub>2</sub> AG                   | 62  |
| IBP                                  | 77  |
| ICP                                  | 84  |
| RM                                   | 85  |
| C.O.                                 | 87  |
| ICG                                  | 90  |
| BIS                                  | 92  |
| PiCCO <sub>2</sub> ScvO <sub>2</sub> | 95  |
| EEG                                  | 100 |
| NMT                                  | 103 |
| rSO <sub>2</sub>                     | 105 |
| Mounting                             | 106 |
| Others                               | 113 |

02

### ● Electrocardiograph Accessories

116

03

### ● Defibrillator Accessories

120



ECG Accessories




Integrated ECG Cable Assembly

- Integrated design, convenient for use and maintenance
- Meeting the requirements of EC53
- Outstanding shielding property and anti-interference performance, protecting ECG signal from being interfered
- Excellent defibrillation-proof performance, well protecting the equipment
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

Integrated ECG Cable Assembly  
For BeneView, iPM, iMEC, uMEC series monitors, BeneHeart defibrillator

| Picture | Model   | Part No.      | Description                                                                                    | Purchasing Unit |
|---------|---------|---------------|------------------------------------------------------------------------------------------------|-----------------|
|         | EA6251B | 040-000961-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 5-Lead, Defib-Proof,<br>AHA, Snap, 3.6 m | Each            |
|         | EA6252B | 040-000963-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 5-Lead, Defib-Proof,<br>IEC, Snap, 3.6 m | Each            |
|         | EA6231B | 040-000965-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 3-Lead, Defib-Proof,<br>AHA, Snap, 3.6 m | Each            |
|         | EA6232B | 040-000967-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 3-Lead, Defib-Proof,<br>IEC, Snap, 3.6 m | Each            |

| Picture                                                                             | Model   | Part No.      | Description                                                                                    | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|------------------------------------------------------------------------------------------------|-----------------|
|    | EA6251A | 040-000960-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 5-Lead, Defib-Proof,<br>AHA, Clip, 3.6 m | Each            |
|    | EA6252A | 040-000962-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 5-Lead, Defib-Proof,<br>IEC, Clip, 3.6 m | Each            |
|   | EA6231A | 040-000964-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 3-Lead, Defib-Proof,<br>AHA, Clip, 3.6 m | Each            |
|  | EA6232A | 040-000966-00 | ECG cable and wires (integrative):<br>Adu/Ped, 12 Pin 3-Lead, Defib-Proof,<br>IEC, Clip, 3.6 m | Each            |

**For PM, MEC series monitors**


| Picture                                                                               | Model   | Part No.      | Description                                                                     | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|---------------------------------------------------------------------------------|-----------------|
|    | EA6151B | 0010-30-12244 | ECG cable and wires: Adu/Ped,<br>6 Pin 5-lead, Defib-Proof, AHA,<br>Snap, 3.6 m | Each            |
|    | EA6152B | 0010-30-12245 | ECG cable and wires: Adu/Ped,<br>6 Pin 5-lead, Defib-Proof, IEC,<br>Snap, 3.6 m | Each            |
|   | EA6131B | 0010-30-12246 | ECG cable and wires: Adu/Ped,<br>6 Pin 3-lead, Defib-Proof, AHA,<br>Snap, 3.6 m | Each            |
|  | EA6132B | 0010-30-12247 | ECG cable and wires: Adu/Ped,<br>6 Pin 3-lead, Defib-Proof, IEC,<br>Snap, 3.6 m | Each            |

Trunk Cables

- Easy to replace leadwires
- Meeting the requirements of EC53
- Outstanding shielding property and anti-interference performance, protecting ECG signal from being interfered
- Excellent defibrillation-proof performance, well protecting the equipment
- ESU-proof, ensuring ECG signals not interfered during operation
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, BeneHeart defibrillator

| Picture                                                                            | Model  | Part No.      | Description                                                           | Purchasing Unit |
|------------------------------------------------------------------------------------|--------|---------------|-----------------------------------------------------------------------|-----------------|
|   | EV6201 | 0010-30-42719 | ECG trunk cable: 3/5-lead, Adu/Ped, 12 Pin, Defib-Proof, AHA/IEC, 3 m | Each            |
|  | EV6211 | 0010-30-42723 | ECG trunk cable: 3/5-lead, Adu/Ped, 12 Pin, ESU-Proof, AHA/IEC, 3 m   | Each            |

| Picture                                                                               | Model  | Part No.      | Description                                                              | Purchasing Unit |
|---------------------------------------------------------------------------------------|--------|---------------|--------------------------------------------------------------------------|-----------------|
|    | EV6202 | 0010-30-42720 | ECG trunk cable: 3-lead, Ped/Neo, 12 Pin, Defib-Proof, AHA/IEC, 3 m      | Each            |
|    | EV6212 | 0010-30-42724 | ECG trunk cable: 3-lead, Ped/Neo, 12 Pin, ESU-Proof, AHA/IEC, 3 m        | Each            |
|   | EV6222 | 040-000754-00 | ECG Trunk Cable, 3-lead, Ped/Neo, 12 Pin, TPU, Defib-Proof, AHA/IEC, 3 m | Each            |
|  | EV6203 | 0010-30-42721 | ECG trunk cable: 12-lead, Adu/Ped, 12 Pin, Defib-Proof, AHA, 3 m         | Each            |



| Picture                                                                           | Model  | Part No.      | Description                                                      | Purchasing Unit |
|-----------------------------------------------------------------------------------|--------|---------------|------------------------------------------------------------------|-----------------|
|  | EV6204 | 0010-30-42722 | ECG trunk cable: 12-lead, Adu/Ped, 12 Pin, Defib-Proof, IEC, 3 m | Each            |

For PM, MEC series monitors

| Picture                                                                             | Model  | Part No.      | Description                                                               | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|---------------------------------------------------------------------------|-----------------|
|    | EV6101 | 0010-30-42782 | ECG trunk cable: 3/5-lead, Adu/Ped, 6 Pin, Defib-Proof, TPU, AHA/IEC, 3 m | Each            |
|   | EV6102 | 0010-30-12377 | ECG trunk cable: 3-lead, Neo, 6 Pin, Defib-Proof, TPU, AHA/IEC, 3 m       | Each            |
|  | EV6122 | 040-000753-00 | ECG Trunk Cable, 3-lead, Ped/Neo, 6 Pin, Defib-Proof, TPU, AHA/IEC, 3 m   | Each            |

ECG Leadwire – IEC

- Easy to replace trunk cables
- Meeting the requirements of EC53
- Outstanding shielding property and anti-interference performance, protecting ECG signal from being interfered
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

Match with 3/5-lead cables (0010-30-42719, 0010-30-42723, 0010-30-42782)

| Picture                                                                              | Model   | Part No.      | Description                                                | Purchasing Unit |
|--------------------------------------------------------------------------------------|---------|---------------|------------------------------------------------------------|-----------------|
|   | EL6502A | 0010-30-42728 | 5-Lead ECG wires, Clip, Adu, TPU, IEC, 0.6 m/1m            | Each            |
|  | EL6504A | 0010-30-42730 | 5-Lead ECG wires, Clip, Adu/Ped, TPU, IEC, long, 1 m/1.4 m | Each            |

| Picture                                                                             | Model   | Part No.      | Description                                     | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|-------------------------------------------------|-----------------|
|    | EL6502B | 0010-30-42736 | 5-Lead ECG wires, Snap, Adu, TPU, IEC, 1 m/1.4m | Each            |
|    | EL6302B | 0010-30-42733 | 3-Lead ECG wires, Snap, Adu/Ped, TPU, IEC, 1m   | Each            |
|   | EL6302A | 0010-30-42725 | 3-Lead ECG wires, Clip, Adu, TPU, IEC, 0.6 m    | Each            |
|  | EL6304A | 0010-30-42732 | 3-Lead ECG wires, Clip, Adu/Ped, TPU, IEC, 1 m  | Each            |

Match with 3-lead cables (0010-30-42720, 0010-30-42724, 0010-30-12377)

| Picture                                                                              | Model   | Part No.      | Description                                  | Purchasing Unit |
|--------------------------------------------------------------------------------------|---------|---------------|----------------------------------------------|-----------------|
|   | EL6306A | 0010-30-42897 | 3-Lead ECG wires, Clip, Neo, TPU, IEC, 1m    | Each            |
|   | EL6308B | 0010-30-42901 | 3-Lead ECG wires, Snap, Ped, TPU, IEC, 0.6 m | Each            |
|  | EL6308A | 0010-30-42899 | 3-Lead ECG wires, Clip, Ped, TPU, IEC, 0.6 m | Each            |

Match with 3-lead Neonatal cables (040-000754-00, 040-000753-00)

| Picture                                                                           | Model   | Part No.      | Description                                                        | Purchasing Unit |
|-----------------------------------------------------------------------------------|---------|---------------|--------------------------------------------------------------------|-----------------|
|  | EL6312B | 040-000147-00 | 3-lead ECG wires, Single-patient use, Ped/Neo, Snap, IEC, DIN, 1 m | Each            |
|  | EL6312A | 040-000149-00 | 3-lead ECG wires, Single-patient use, Ped/Neo, Clip, IEC, DIN, 1 m | Each            |


Match with 12-lead cables (0010-30-42721, 0010-30-42722)

| Picture                                                                               | Model   | Part No.      | Description                                                    | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|----------------------------------------------------------------|-----------------|
|    | EL6802B | 0010-30-42907 | 12-Lead ECG limb wires, Snap, Adu/Ped, Gray, TPU, IEC, 0.8 m   | Each            |
|    | EL6804B | 0010-30-42909 | 12-Lead ECG Chest Wires, Snap, Adu/Ped, White, TPU, IEC, 0.6 m | Each            |
|   | EL6802A | 0010-30-42903 | 12-Lead ECG limb wires, Clip, Adu/Ped, Gray, TPU, IEC, 0.8 m   | Each            |
|  | EL6804A | 0010-30-42905 | 12-Lead ECG Chest Wires, Clip, Adu/Ped, White, TPU, IEC, 0.6 m | Each            |





For Telemetry(TMS-6016)

| Picture                                                                             | Model   | Part No.      | Description                                            | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|--------------------------------------------------------|-----------------|
|    | EY6301B | 115-004868-00 | 3-Lead ECG wires, Telemetry, IEC, Snap, Adu/Ped, 0.8 m | Each            |
|    | EY6501B | 115-004870-00 | 5-Lead ECG wires, Telemetry, IEC, Snap, Adu/Ped, 0.8 m | Each            |
|   | EY6301A | 115-004872-00 | 3-Lead ECG wires, Telemetry, IEC, Clip, Adu/Ped, 0.8 m | Each            |
|  | EY6501A | 115-004874-00 | 5-Lead ECG wires, Telemetry, IEC, Clip, Adu/Ped, 0.8 m | Each            |

| Picture                                                                               | Model   | Part No.      | Description                                                         | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|---------------------------------------------------------------------|-----------------|
|    | EY6303B | 040-001138-00 | 3-Lead ECG wires, Telemetry, IEC, Snap, Adu/Ped, 0.8 m, Water-proof | Each            |
|    | EY6503B | 040-001140-00 | 5-Lead ECG wires, Telemetry, IEC, Snap, Adu/Ped, 0.8 m, Water-proof | Each            |
|   | EY6303A | 040-001142-00 | 3-Lead ECG wires, Telemetry, IEC, Clip, Adu/Ped, 0.8 m, Water-proof | Each            |
|  | EY6503A | 040-001144-00 | 5-Lead ECG wires, Telemetry, IEC, Clip, Adu/Ped, 0.8 m, Water-proof | Each            |

For Telemetry (TM-80)

| Picture                                                                             | Model | Part No.      | Description                                                | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------|-----------------|
|    |       | 009-004769-00 | 3-Lead, Adu/Ped, Snap, 36", reusable, TPU, IEC, BeneVision | Each            |
|                                                                                     |       | 009-004768-00 | 3-Lead, Adu/Ped, Snap, 24", reusable, TPU, IEC, BeneVision | Each            |
|    |       | 009-004775-00 | 3-Lead, Adu/Ped, Clip, 36"reusable, TPU, IEC, BeneVision   | Each            |
|                                                                                     |       | 009-004774-00 | 3-Lead, Adu/Ped, Clip, 24", reusable, TPU, IEC, BeneVision | Each            |
|   |       | 009-004785-00 | 5-Lead, Adu/Ped, Snap, 36", reusable, TPU, IEC, BeneVision | Each            |
|                                                                                     |       | 009-004784-00 | 5-Lead, Adu/Ped, Snap, 24", reusable, TPU, IEC, BeneVision | Each            |
|  |       | 009-004789-00 | 5-Lead, Adu/Ped, Clip, 36", reusable,TPU, IEC, BeneVision  | Each            |
|                                                                                     |       | 009-004788-00 | 5-Lead, Adu/Ped, Clip, 24", reusable, TPU, IEC, BeneVision | Each            |

| Picture                                                                               | Model   | Part No.      | Description                                                 | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|-------------------------------------------------------------|-----------------|
|    | EY6304B | 009-004797-00 | 6-Lead, Adu/Ped, Snap, 36", reusable, TPU, IEC, BeneVision  | Each            |
|                                                                                       |         | 009-004796-00 | 6-Lead, Adu/Ped, Snap, 24", reusable, TPU, IEC, BeneVision  | Each            |
|    | EY6504B | 009-004801-00 | 6-Lead, Adu/Ped, Clip,36", reusable, TPU, IEC, BeneVision   | Each            |
|                                                                                       |         | 009-004800-00 | 6-Lead, Adu/Ped, Clip, 24", reusable, TPU, IEC, BeneVision  | Each            |
|   | EY6304A | 009-004780-00 | 3-Lead, Adu/Ped, Snap, 24", single-patient, IEC, BeneVision | Each            |
|  | EY6504A | 009-004792-00 | 5-Lead, Adu/Ped, Snap, 24", single-patient, IEC, BeneVision | Each            |

ECG Leadwire – AHA

- Easy to replace trunk cables
- Meeting the requirements of EC53
- Outstanding shielding property and anti-interference performance, protecting ECG signal from being interfered
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

Match with 3/5-lead cables (0010-30-42719, 0010-30-42723, 0010-30-42782)

| Picture                                                                            | Model   | Part No.      | Description                                               | Purchasing Unit |
|------------------------------------------------------------------------------------|---------|---------------|-----------------------------------------------------------|-----------------|
|   | EL6501A | 0010-30-42727 | 5-Lead ECG wires, Clip, Adu, TPU, AHA, 0.6 m/1m           | Each            |
|  | EL6503A | 0010-30-42729 | 5-Lead ECG wires, Clip, Adu/Ped, TPU, AHA, long, 1m/1.4 m | Each            |

| Picture                                                                               | Model   | Part No.      | Description                                      | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|--------------------------------------------------|-----------------|
|    | EL6501B | 0010-30-42735 | 5-Lead ECG wires, Snap, Adu, TPU, AHA, 1m/1.4m   | Each            |
|    | EL6301B | 0010-30-42734 | 3-Lead ECG wires, Snap, Adu/Ped, TPU, AHA, 1m    | Each            |
|   | EL6301A | 0010-30-42726 | 3-Lead ECG wires, Clip, Adu, TPU, AHA, 0.6 m     | Each            |
|  | EL6303A | 0010-30-42731 | 3-Lead ECG wires, Clip, Adu/Ped, TPU, AHA, (1 m) | Each            |

Match with 3-lead cables (0010-30-42720, 0010-30-42724, 0010-30-12377)

| Picture                                                                            | Model   | Part No.      | Description                                  | Purchasing Unit |
|------------------------------------------------------------------------------------|---------|---------------|----------------------------------------------|-----------------|
|   | EL6305A | 0010-30-42896 | 3-Lead ECG wires, Clip, Neo, TPU, AHA, 1m    | Each            |
|   | EL6307B | 0010-30-42900 | 3-Lead ECG wires, Snap, Ped, TPU, AHA, 0.6 m | Each            |
|  | EL6307A | 0010-30-42898 | 3-Lead ECG wires, Clip, Ped, TPU, AHA, 0.6 m | Each            |

Match with 3-lead Neonatal cables (040-000754-00, 040-000753-00)


| Picture                                                                             | Model   | Part No.      | Description                                  | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|----------------------------------------------|-----------------|
|  | EL6311B | 040-000146-00 | 3-Lead ECG wires, Disposable, AHA, Snap, 1 m | Each            |
|  | EL6311A | 040-000148-00 | 3-Lead ECG wires, Disposable, AHA, Clip, 1 m | Each            |

Match with 12-lead cables (0010-30-42721, 0010-30-42722)

| Picture                                                                             | Model   | Part No.      | Description                                                    | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|----------------------------------------------------------------|-----------------|
|    | EL6801A | 0010-30-42902 | 12-Lead ECG limb wires, Clip, Adu/Ped, Gray, TPU, AHA, 0.8 m   | Each            |
|    | EL6803A | 0010-30-42904 | 12-Lead ECG Chest wires, Clip, Adu/Ped, Brown, TPU, AHA, 0.6 m | Each            |
|   | EL6801B | 0010-30-42906 | 12-Lead ECG limb wires, Snap, Adu/Ped, Gray, TPU, AHA, 0.8 m   | Each            |
|  | EL6803B | 0010-30-42908 | 12-Lead ECG Chest wires, Snap, Adu/Ped, Brown, TPU, AHA, 0.6 m | Each            |

For Telemetry(TMS-6016)

| Picture                                                                               | Model   | Part No.      | Description                                            | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|--------------------------------------------------------|-----------------|
|    | EY6302B | 115-004867-00 | 3-Lead ECG wires, Telemetry, AHA, Snap, Adu/Ped, 0.8 m | Each            |
|    | EY6502B | 115-004869-00 | 5-Lead ECG wires, Telemetry, AHA, Snap, Adu/Ped, 0.8 m | Each            |
|   | EY6302A | 115-004871-00 | 3-Lead ECG wires, Telemetry, AHA, Clip, Adu/Ped, 0.8 m | Each            |
|  | EY6502A | 115-004873-00 | 5-Lead ECG wires, Telemetry, AHA, Clip, Adu/Ped, 0.8 m | Each            |

| Picture                                                                             | Model   | Part No.      | Description                                                         | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|---------------------------------------------------------------------|-----------------|
|    | EY6304B | 040-001137-00 | 3-Lead ECG wires, Telemetry, AHA, Snap, Adu/Ped, 0.8 m, Water-proof | Each            |
|    | EY6504B | 040-001139-00 | 5-Lead ECG wires, Telemetry, AHA, Snap, Adu/Ped, 0.8 m, Water-proof | Each            |
|   | EY6304A | 040-001141-00 | 3-Lead ECG wires, Telemetry, AHA, Clip, Adu/Ped, 0.8 m, Water-proof | Each            |
|  | EY6504A | 040-001143-00 | 5-Lead ECG wires, Telemetry, AHA, Clip, Adu/Ped, 0.8 m, Water-proof | Each            |

For Telemetry (TM-80)

| Picture                                                                               | Model | Part No.      | Description                                                | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------|-----------------|
|    |       | 009-004766-00 | 3-Lead, Adu/Ped, Snap,36", reusable,TPU, AHA, BeneVision   | Each            |
|                                                                                       |       | 009-004765-00 | 3-Lead, Adu/Ped, Snap, 24", reusable, TPU, AHA, BeneVision | Each            |
|    |       | 009-004772-00 | 3-Lead, Adu/Ped, Clip,36", reusable, TPU, AHA, BeneVision  | Each            |
|                                                                                       |       | 009-004771-00 | 3-Lead, Adu/Ped, Clip, 24", reusable, TPU, AHA, BeneVision | Each            |
|   |       | 009-004783-00 | 5-Lead, Adu/Ped, Snap, 36", reusable, TPU, AHA, BeneVision | Each            |
|                                                                                       |       | 009-004782-00 | 5-Lead, Adu/Ped, Snap, 24", reusable, TPU, AHA,BeneVision  | Each            |
|  |       | 009-004787-00 | 5-Lead, Adu/Ped, Clip, 36", reusable, TPU, AHA, BeneVision | Each            |
|                                                                                       |       | 009-004786-00 | 5-Lead, Adu/Ped, Clip, 24", reusable, TPU, AHA, BeneVision | Each            |




| Picture                                                                             | Model | Part No.      | Description                                                 | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------------------|-----------------|
|    |       | 009-004795-00 | 6-Lead, Adu/Ped, Snap, 36", reusable, TPU, AHA, BeneVision  | Each            |
|                                                                                     |       | 009-004794-00 | 6-Lead, Adu/Ped, Snap, 24", reusable, TPU, AHA, BeneVision  | Each            |
|                                                                                     |       | 009-004799-00 | 6-Lead, Adu/Ped, Clip, 36", reusable, TPU, AHA, BeneVision  | Each            |
|                                                                                     |       | 009-004798-00 | 6-Lead, Adu/Ped, Clip, 24", reusable, TPU, AHA, BeneVision  | Each            |
|   |       | 009-004777-00 | 3-Lead, Adu/Ped, Snap, 24", single-patient, AHA, BeneVision | Each            |
|  |       | 009-004790-00 | 5-Lead, Adu/Ped, Snap, 24", single-patient, AHA, BeneVision | Each            |

Electrode

- Latex free
- DEHP free
- Good biocompatibility, avoiding allergic reactions to patient

| Picture                                                                             | Model    | Part No.      | Description                                   | Purchasing Unit |
|-------------------------------------------------------------------------------------|----------|---------------|-----------------------------------------------|-----------------|
|  | 31499224 | 0010-10-12304 | Adult ECG Electrode (Kendall, Medi Trace 210) | 10 pcs/pouch    |
|                                                                                     | H124SG   | 900E-10-04880 | Neonatal ECG Electrode (Kendall, H124SG)      | 50pcs/pouch     |
|  | 2249     | 0509-10-00094 | Adult ECG electrode (3M, 2249)                | 25 pcs/pouch    |
|                                                                                     | 2245     | 9000-10-07469 | Pediatric ECG electrode (3M, 2245)            | 50 pcs/pouch    |
|                                                                                     |          | 040-002711-00 | Adult ECG electrode (INTCO)                   | 5 pcs/pouch     |
|                                                                                     |          | 040-002833-00 | Pediatric/Neonatal ECG electrode (INTCO)      | 30 pcs/pouch    |

Match with 3-lead Neonatal cables (040-000754-00, 040-000753-00)

| Picture                                                                             | Model                | Part No.        | Description                                                                | Purchasing Unit              |
|-------------------------------------------------------------------------------------|----------------------|-----------------|----------------------------------------------------------------------------|------------------------------|
|    | 1050NPSM<br>KittyCat | 0681-00-0098-01 | Neonatal 3-Lead Pre-wired electrodes, Radio Opaque, AAMI, 18" (45 cm)      | 100pouch/box<br>(3pcs/pouch) |
|    | 1051NPSM<br>KittyCat | 0681-00-0098-02 | Neonatal 3-Lead Pre-wired electrodes, Radio Translucent, AAMI, 18" (45 cm) | 100pouch/box<br>(3pcs/pouch) |
|   | 0406060              | 040-003251-00   | Disposable neonatal 3-lead pre-wired electrode, radio opaque, AHA, 60cm    | 50pouch/box<br>(3pcs/pouch)  |
|  | 0406059              | 040-003252-00   | Disposable neonatal 3-lead pre-wired electrode, radio opaque, IEC, 60cm    | 50pouch/box<br>(3pcs/pouch)  |

| Picture                                                                             | Model   | Part No.      | Description                                                                  | Purchasing Unit             |
|-------------------------------------------------------------------------------------|---------|---------------|------------------------------------------------------------------------------|-----------------------------|
|  | 0406062 | 040-003254-00 | Disposable neonatal 3-lead pre-wired electrode, radio translucent, AHA, 60cm | 50pouch/box<br>(3pcs/pouch) |
|  | 0406061 | 040-003255-00 | Disposable neonatal 3-lead pre-wired electrode, radio translucent, IEC, 60cm | 50pouch/box<br>(3pcs/pouch) |

Veterinary Accessory (Reusable)

| Picture                                                                              | Model | Part No.      | Description    | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------|---------------|----------------|-----------------|
|  |       | 9101-20-58104 | Crocodile clip | 5pcs/box        |

SpO<sub>2</sub> Accessories

Mindray SpO<sub>2</sub> Accessories

Mindray SpO<sub>2</sub> Cable

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Flexible and durable cables
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Easy to change sensor, meeting clinical requirements for patient use
- Latex free

For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, VS-600/900, BeneHeart defibrillators

| Picture                                                                            | Model | Part No.      | Description                                                     | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|-----------------------------------------------------------------|-----------------|
|  | 562A  | 009-004600-00 | Mindray SpO <sub>2</sub> extension cable, 7 Pin, 2.5m(high-end) | Each            |
|                                                                                    | 562B  | 040-001443-00 | Mindray SpO <sub>2</sub> extension cable, 7 Pin, 1.2 m          | Each            |

For PM, MEC series monitors, VS-800 vital signs monitor

| Picture                                                                             | Model | Part No.      | Description                                            | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------|-----------------|
|  | 561A  | 0010-20-42594 | Mindray SpO <sub>2</sub> extension cable, 6 Pin, 2.5 m | Each            |
|                                                                                     |       | 040-001274-00 | Mindray SpO <sub>2</sub> extension cable, 6 Pin, 1.2 m | Each            |

For TMS-6016 telemetry monitoring system

| Picture                                                                             | Model | Part No.      | Description                                                                  | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------------|-----------------|
|  |       | 0152-30-39689 | Mindray SpO <sub>2</sub> telemetry module with extension cable, 4 Pin, 0.5 m | Each            |
|                                                                                     |       | 0152-21-39686 | Telemetry adapter cable (metal connector to plastic connector)               | Each            |

For TM80 telemetry monitoring system

| Picture                                                                               | Model | Part No.      | Description                                                      | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------|-----------------|
|  |       | 115-029488-00 | Mindray SpO <sub>2</sub> module for BeneVision TM80, 6 Pin, 0.5m | Each            |

Mindray SpO<sub>2</sub> Sensor

Finger-Clip Sensor (Reusable)

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- High quality photoelectric element, ensuring precise measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Perfect performance against light interference, can be used in environment of strong light
- ESU-proof, ensuring SpO<sub>2</sub> signals not interfered during operation
- Strict electric safety specification, guaranteeing safety for use
- Few pit structure, not easily staining, convenient for cleaning
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient

For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                            | Model | Part No.      | Description                                              | Purchasing Unit |
|--------------------------------------------------------------------|-------|---------------|----------------------------------------------------------|-----------------|
| A white, reusable finger-clip sensor with a cable and a connector. | 512F  | 512F-30-28263 | Reusable sensor, adult, finger-clip, 1.1 m, >30 kg       | Each            |
| A white, reusable finger-clip sensor with a cable and a connector. | 512H  | 512H-30-79061 | Reusable sensor, pediatric, finger-clip, 1.1 m, 10-30 kg | Each            |

Finger-Tip Sensor (Reusable)

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- High quality photoelectric element, ensuring precise measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Perfect performance against light interference, can be used in environment of strong light
- ESU-proof, ensuring SpO<sub>2</sub> signals not interfered during operation
- Strict electric safety specification, guaranteeing safety for use
- Silicone rubber sheath, not likely to break in case of drop, hardly sensor off
- Few pit structure, not likely staining, convenient for cleaning
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient

For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                           | Model | Part No.      | Description                                             | Purchasing Unit |
|-------------------------------------------------------------------|-------|---------------|---------------------------------------------------------|-----------------|
| A white, reusable finger-tip sensor with a cable and a connector. | 512E  | 512E-30-90390 | Reusable sensor, adult, finger-tip, 1.1 m, >30 kg       | Each            |
| A white, reusable finger-tip sensor with a cable and a connector. | 512G  | 512G-30-90607 | Reusable sensor, pediatric, finger-tip, 1.1 m, 10-30 kg | Each            |

Bandage Sensor (Reusable)

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- High quality photoelectric element, ensuring precise measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- ESU-proof, ensuring SpO<sub>2</sub> signals not interfered during operation
- Strict electric safety specification, guaranteeing safety for use
- Compact and lightweight, single buckle design, easy to handle, not likely to loose
- Silicone rubber bandage (518B), no sharp edges, no hazard to hurt the neonates, no breakage from drop
- Foam bandage (518C), comfortable, permeable to air
- Few pit structure, not likely staining, convenient for cleaning
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient.

For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                                             | Model | Part No.      | Description                                                          | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|----------------------------------------------------------------------|-----------------|
|   | 518B  | 518B-30-72107 | Reusable sensor, Neo, Foot (adult/pediatric, finger), 1.2 m, <5 kg   | Each            |
|  | 518C  | 040-000407-00 | Reusable sensor, with foam wrap replacement, Neo, Foot, 1.2 m, <5 kg | Each            |

| Picture                                                                             | Model | Part No.      | Description                            | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|----------------------------------------|-----------------|
|  |       | 115-004895-00 | Disposable wrap for 518C sensor, <5 kg | 12 pcs/pouch    |

Ear Sensor (Reusable)

- Suitable for SpO<sub>2</sub> monitoring for the case sensor cannot be applied to fingers or toes
- Ergonomic design, comfortable to use
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient.

For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                                              | Model | Part No.      | Description                                                       | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------------------------|-----------------|
|  | 513A  | 115-034527-00 | Reusable SpO <sub>2</sub> sensor, adu/ped, ear-clip, 1.6m, >30 kg | Each            |

Veterinary Sensor (Reusable)

For all SpO<sub>2</sub> Cables of Veterinaty Monitor

| Picture                                                                           | Model | Part No.      | Description                                                | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------|-----------------|
|  | 551B  | 115-036222-00 | Veterinary SpO <sub>2</sub> Sensor, 2-clip, reusable, 1.2m | Each            |

For iMEC, uMEC series patient monitor

| Picture                                                                            | Model | Part No.      | Description                                                                  | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------------|-----------------|
|   | 512F  | 115-012807-00 | 512F integrative reusable SpO <sub>2</sub> sensor, Adult, Finger, >30 kg, 3m | Each            |
|  | 518B  | 115-020887-00 | Integrative reusable SpO <sub>2</sub> s ensor, (Adu/Neo , bandage, 7 pin, 3m | Each            |

Non-Adhesive Sensor (Disposable)

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- High quality photoelectric element, ensuring precise measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- ESU-proof, ensuring SpO<sub>2</sub> signals not interfered during operation
- Strict electric safety specification, guaranteeing safety for use
- Velcro® bandage, convenient for changing application site
- Latex free
- Good biocompatibility, avoiding allergic reactions tc patient.

For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                                               | Model | Part No.      | Description                                                                  | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------------|-----------------|
|    | 520A  | 009-005087-00 | 520A disposable SpO <sub>2</sub> sensor, non-adhesive, Adu, >30kg, 0.5m      | 20 pcs/box      |
|   | 520P  | 009-005088-00 | 520P disposable SpO <sub>2</sub> sensor, non-adhesive, Ped, 10 to 50kg, 0.5m | 20 pcs/box      |
|  | 520I  | 009-005089-00 | 520I disposable SpO <sub>2</sub> sensor, non-adhesive, Inf, 3 to 20kg, 1m    | 20 pcs/box      |



| Picture                                                                           | Model | Part No.      | Description                                                                          | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------------------------|-----------------|
|  | 520N  | 009-005090-00 | 520N disposable SpO <sub>2</sub> sensor, non-adhesive, Adult (>30kg ) /Neo(<3kg), 1m | 20 pcs/box      |

### Adhesive Sensor (Disposable)

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- High quality photoelectric element, ensuring precise measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- ESU-proof, ensuring SpO<sub>2</sub> signals not interfered during operation
- Strict electric safety specification, guaranteeing safety for use
- Medical adhesive tape , great performance for fixation
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient.

#### For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                                             | Model | Part No.      | Description                                                              | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------------|-----------------|
|   | 521A  | 009-005091-00 | 521A disposable SpO <sub>2</sub> sensor, adhesive, Adu, >30kg, 0.5m      | 20pcs/box       |
|  | 521P  | 009-005092-00 | 521P disposable SpO <sub>2</sub> sensor, adhesive, Ped, 10 to 50kg, 0.5m | 20 pcs/box      |

#### For all Mindray SpO<sub>2</sub> Cables and PM-50/60 pulse oximeter

| Picture                                                                             | Model | Part No.      | Description                                                                      | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|----------------------------------------------------------------------------------|-----------------|
|  | 521I  | 009-005093-00 | 521I disposable SpO <sub>2</sub> sensor, adhesive, Inf, 3 to 20kg, 1m            | 20 pcs/box      |
|  | 521N  | 009-005094-00 | 521N disposable SpO <sub>2</sub> sensor, adhesive, Adult (>30kg ) /Neo(<3kg), 1m | 20 pcs/box      |

Nellcor SpO<sub>2</sub> Accessories

Nellcor SpO<sub>2</sub> Cables

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Flexible and durable cables
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Easy to change sensor, meeting clinical requirements for patient use
- Latex free

For BeneVision, BeneView, iPM, iMEC series monitors, VS-600/900, BeneHeart defibrillators

| Picture                                                                           | Model | Part No.      | Description                                            | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------|-----------------|
|  | 572A  | 0010-20-42712 | Nellcor SpO <sub>2</sub> extension cable, 8 Pin, 2.5 m | Each            |

For VS-800, PM series monitors

| Picture                                                                            | Model | Part No.      | Description                                            | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------|-----------------|
|  | 571B  | 0010-20-42595 | Nellcor SpO <sub>2</sub> extension cable, 6 Pin, 2.5 m | Each            |

Nellcor SpO<sub>2</sub> Sensors

Reusable Nellcor SpO<sub>2</sub> Sensors

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Strict electric safety specification, guaranteeing safety for use
- Outstanding cable jacket, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient

For all the Nellcor Cables (0010-20-42712, 0010-20-42595)




| Picture                                                                               | Model   | Part No.      | Description                                                | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|------------------------------------------------------------|-----------------|
|    | DS-100A | 9000-10-05161 | Nellcor reusable sensor, adult, 1 m, >40 kg                | Each            |
|   | OXI-P/I | 9000-10-07308 | Nellcor reusable sensor, pediatric/Inf, 1 m, 3-40 kg       | Each            |
|  | OXI-A/N | 9000-10-07336 | Nellcor reusable sensor, adult/Neo, 1 m, 1-<3 kg or >40 kg | Each            |

| Picture                                                                           | Model | Part No.      | Description                                                                                    | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------------------------------|-----------------|
|  | D-YS  | 0010-10-12476 | Nellcor multisite sensor, adult/pediatric/Inf/Neo, 1 m, 1-<3 kg or >40 kg (A/N), 3-40 kg (P/I) | Each            |
|                                                                                   | V-SAT | 9101-10-58134 | Nellcor veterinary SpO2 sensor                                                                 | Each            |

Nellcor SpO<sub>2</sub> Sensor (Disposable)

- Ergonomic design, precise engineering and clinical testing guaranteeing reliable measurement
- Well anti-electromagnetic interference, suitable for complex electrical environment
- Strict electric safety specification, guaranteeing safety for use
- Various sensors meeting requirements of different patients
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient

For all the Nellcor Cables (0010-20-42712, 0010-20-42595)

| Picture                                                                               | Model | Part No.      | Description                                                  | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------|-----------------|
|    | MAX-A | 0010-10-12202 | Nellcor disposable sensor, adult, >30 kg, 0.5 m              | 24 pcs/box      |
|   | MAX-P | 0010-10-12203 | Nellcor disposable sensor, pediatric, 10 to 50 kg, 0.5 m     | 24 pcs/box      |
|  | MAX-I | 0010-10-12204 | Nellcor disposable sensor, Inf, 3 to 20 kg, 0.5 m            | 24 pcs/box      |
|  | MAX-N | 0010-10-12205 | Nellcor disposable sensor, adult/Neo, >40 kg or <3 kg, 0.5 m | 24 pcs/box      |

Non-invasive Blood Pressure Accessories



Non-invasive Blood Pressure Tubing

- Not deform even if for long-term use
- Reasonable hardness avoids severe bending, ensures the safety of deflation
- Reasonable hardness ensures well transmission of blood pressure signal
- TPU ensures good air tightness and long life
- Latex free, PVC free
- Good biocompatibility, free from biological hazard to skin

For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, VS-600/900,BeneHeart defibrillator

| Picture                                                                              | Part No.      | Description                                                                            | Purchasing Unit |
|--------------------------------------------------------------------------------------|---------------|----------------------------------------------------------------------------------------|-----------------|
|   | 6200-30-09688 | NIBP Tubing, Adu/Ped/Inf, with connectors, 3 m                                         | Each            |
|   | 115-012522-00 | NIBP Tubing, Adu/Ped/Inf, with connectors, 3 m (For iMEC, uMEC series patient monitor) | Each            |
|  | 6200-30-11560 | NIBP Tubing, Neo, with connectors, 3 m                                                 | Each            |

For VS-800, MEC, PM series monitors

| Picture                                                                           | Part No.      | Description                   | Purchasing Unit |
|-----------------------------------------------------------------------------------|---------------|-------------------------------|-----------------|
|  | 509B-30-06259 | NIBP tubing, Adu/Ped/Inf, 3 m | Each            |
|  | 509B-30-06260 | NIBP tubing, Neo, 3 m         | Each            |

Adapted with the tubing (6200-30-09688 and 509B-30-06259) and the neonatal cuffs

| Picture                                                                            | Part No.      | Description                                             | Purchasing Unit |
|------------------------------------------------------------------------------------|---------------|---------------------------------------------------------|-----------------|
|  | 040-000688-00 | NIBP Cuff Tubing Adapter (Adult tubing to Neonate cuff) | Each            |

Non-invasive Blood Pressure Cuff (Reusable)

CM1200 Series

- Soft and comfortable. Low hazard to skin even if a long-term use
- Easy to clean. The cuff wrap can not be damped or stained by liquid if duly cleaned
- Pilling-proof. Not deform even if for long-term use
- TPU bladder ensures good air tightness and long life
- Latex free, PVC free
- Good biocompatibility, free from biological hazard to skin



Connected with the tubing 6200-30-09688, 509B-30-06259 and 115-012522-00

| Model  | Part No.      | Description                                                                                    | Purchasing Unit |
|--------|---------------|------------------------------------------------------------------------------------------------|-----------------|
| CM1200 | 115-002480-00 | Reusable cuff, Small Inf, 7-13 cm<br>Connected with the tubing 509B-30-06260 and 6200-30-11560 | Each            |
| CM1201 | 0010-30-12157 | Reusable cuff, Inf, 10-19 cm, with connector                                                   | Each            |
| CM1202 | 0010-30-12158 | Reusable cuff, Child, 18-26 cm, with connector                                                 | Each            |
| CM1203 | 0010-30-12159 | Reusable cuff, Adu, 25-35 cm, with connector                                                   | Each            |
| CM1204 | 0010-30-12160 | Reusable cuff, Large Adu, 33-47 cm, with connector                                             | Each            |
| CM1205 | 0010-30-12161 | Reusable cuff, Thigh, 46-66 cm, with connector                                                 | Each            |

CM1300 Series

- Easy to clean. No bladder, can be cleaned and disinfected directly
- Color-coded by size for different patient
- Avoid inaccurate measurement. No-bladder design avoids inaccurate measurements due to bladder twisting
- Pilling-proof. Not deform even if for long-term use
- Round angle design improves safety
- Easy to use. The index line for proper size and placement
- Latex free, PVC free
- Good biocompatibility



Connected with the tubing 6200-30-09688, 509B-30-06259 and 115-012522-00

| Model  | Part No.      | Description                                                                                                | Purchasing Unit |
|--------|---------------|------------------------------------------------------------------------------------------------------------|-----------------|
| CM1300 | 040-000968-00 | Reusable bladderless cuff, Small Inf, 7-13 cm<br>Connected with the tubing 509B-30-06260 and 6200-30-11560 | Each            |
| CM1301 | 040-000973-00 | Reusable bladderless cuff, Inf, 10-19 cm, with connector                                                   | Each            |
| CM1302 | 040-000978-00 | Reusable bladderless cuff, Child, 18-26 cm, with connector                                                 | Each            |
| CM1303 | 040-000983-00 | Reusable bladderless cuff, Adu, 25-35 cm, with connector                                                   | Each            |
| CM1304 | 040-000988-00 | Reusable bladderless cuff, Large Adu, 33-47 cm, with connector                                             | Each            |
| CM1305 | 040-000993-00 | Reusable bladderless cuff, Thigh, 46-66 cm, with connector                                                 | Each            |
| CM1306 | 115-015930-00 | Reusable bladderless cuff, Adu, long wrap, 25-35 cm, with connector                                        | Each            |
| CM1307 | 115-015931-00 | Reusable bladderless cuff, Large Adu, long wrap, 33-47 cm, with connector                                  | Each            |

Non-invasive Blood Pressure Cuff (Single Use)

CM1500 Series (Adult and Child)

- Soft and comfortable. Low hazard to skin even if a long-term use
- Avoid inaccurate measurement. No bladder design avoids inaccurate measurements due to bladder twisting
- For single patient only, avoiding cross infection
- Round angle design improves safety
- Easy to use. The index line for proper size and placement
- Latex free, PVC free
- Good biocompatibility, free from biological hazard to skin



Connected with the tubing 6200-30-09688, 509B-30-06259 and 115-012522-00

| Model  | Part No.      | Description                                             | Purchasing Unit |
|--------|---------------|---------------------------------------------------------|-----------------|
| CM1501 | 001B-30-70697 | Single-patient use cuff, Inf, 10-19 cm                  | 10 pcs/box      |
| CM1502 | 001B-30-70698 | Single-patient use cuff, Child, 18-26 cm                | 10 pcs/box      |
| CM1503 | 001B-30-70699 | Single-patient use cuff, Adu, 25-35 cm                  | 10 pcs/box      |
| CM1504 | 001B-30-70700 | Single-patient use cuff, Large Adu, 33-47 cm            | 10 pcs/box      |
| CM1505 | 001B-30-70701 | Single-patient use cuff, Thigh, 46-66 cm                | 10 pcs/box      |
| CM1506 | 115-016969-00 | Single-patient use cuff, Adu, long wrap, 25-35 cm       | 10 pcs/box      |
| CM1507 | 115-016970-00 | Single-patient use cuff, Large Adu, long wrap, 33-47 cm | 10 pcs/box      |

CM1500 Series (Neonatal)

- Soft and comfortable. Low hazard to skin even if a long-term use
- Easy to clean. The stain can be cleaned by wet cloth
- Avoid inaccurate measurement
  - No-bladder design avoids inaccurate measurements due to bladder twisting
- For single patient only, avoiding cross infection
- Transparent design for observing the patient skin condition conveniently
- Easy to use. The index line for proper size and placement
- Latex free, PVC free
- Good biocompatibility, free from biological hazard to skin



Connected with the tubing 6200-30-11560 and 509B-30-06260

| Model   | Part No.      | Description                               | Purchasing Unit |
|---------|---------------|-------------------------------------------|-----------------|
| CM1500A | 001B-30-70692 | Single-patient use cuff, Neo, 3.1-5.7 cm  | 20 pcs/box      |
| CM1500B | 001B-30-70693 | Single-patient use cuff, Neo, 4.3-8.0 cm  | 20 pcs/box      |
| CM1500C | 001B-30-70694 | Single-patient use cuff, Neo, 5.8-10.9 cm | 20 pcs/box      |
| CM1500D | 001B-30-70695 | Single-patient use cuff, Neo, 7.1-13.1 cm | 20 pcs/box      |
| CM1500E | 001B-30-70696 | Single-patient use cuff, Neo, 8-15 cm     | 20 pcs/box      |

For Telemetry BP10

| Picture                                                                             | Part No.      | Description                                                       | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------------|-------------------------------------------------------------------|-----------------|
|  | 040-002379-00 | NIBP tubing                                                       | Each            |
|  | 115-020374-00 | CM1102 reusable cuff, Child,18-26cm,with connector                | Each            |
|                                                                                     | 115-020375-00 | CM1103 reusable cuff, Adu,25-35cm,with connector                  | Each            |
|                                                                                     | 115-020376-00 | CM1104 reusable cuff, Large Adu,33-47cm,with connector            | Each            |
|  | 115-020377-00 | CM1602 Single-patient cuff liner, Ped, 18-26cm (for CM1102)       | Each            |
|                                                                                     | 115-020378-00 | CM1603 Single-patient cuff liner, Adu, 25-35cm (for CM1103)       | Each            |
|                                                                                     | 115-020379-00 | CM1604 Single-patient cuff liner, Large Adu, 33-47cm (for CM1104) | Each            |



NIBP Accessory Kit

For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, VS-600/900, BeneHeart defibrillator

| Part No.      | Description                                                                                                                                                                                               | Purchasing Unit |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 0651-30-77024 | Tubing + Adu cuff (25-35 cm)<br>Including:<br>6200-30-09688 NIBP Tubing, Adu/Ped/Inf, with connectors, 3 m 1 pcs<br>CM1203 Reusable cuff, Adu, 25-35 cm, with connector 1 pcs                             | Set             |
| 0651-30-77025 | Tubing + Ped cuff (18-26 cm)<br>Including:<br>6200-30-09688 NIBP Tubing, Adu/Ped/Inf, with connectors, 3 m 1 pcs<br>CM1202 Reusable cuff, Ped, 18-26 cm, with connector 1 pcs                             | Set             |
| 0651-30-77026 | Tubing + Neo disposable cuff (5.8-10.9 cm)<br>Including:<br>6200-30-11560 NIBP Tubing, Neo, with connectors, 3 m 1 pcs<br>CM1500C (001B-30-70679) Disposable cuff, Neo, 5.8-10.9 cm, with connector 1 pcs | Set             |

For VS-800, MEC, PM series monitors

| Part No.      | Description                                                                                                                                                  | Purchasing Unit |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 115-003557-00 | Tubing + Adu cuff (25-35 cm)<br>Including:<br>509B-30-06259 NIBP tubing, Adu/Ped/Inf, 3 m 1 pcs<br>CM1203 Reusable cuff, Adu, 25-35 cm, with connector 1 pcs | Set             |
| 115-003558-00 | Tubing + Ped cuff (18-26 cm)<br>Including:<br>509B-30-06259 NIBP tubing, Adu/Ped/Inf, 3 m 1 pcs<br>CM1202 Reusable cuff, Ped, 18-26 cm, with connector 1 pcs | Set             |

| Part No.      | Description                                                                                                                                                                              | Purchasing Unit |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 115-003559-00 | Tubing + Neo disposable cuff (5.8-10.9 cm)<br>Including:<br>509B-30-06260 NIBP tubing, Neo, 3 m 1 pcs<br>CM1500C (001B-30-70679) Disposable cuff, Neo, 5.8-10.9 cm, with connector 1 pcs | Set             |

Temperature Accessories



Reusable Temperature Probes

- Available in Rectal/Esophageal and Skin Surface Styles
- Accuracy ± 0.1°C
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient

For BeneVision, BeneView, iPM, iMEC, uMECseries monitors, BeneHeart defibrillator

| Picture                                                                               | Model  | Part No.      | Description                                                 | Purchasing Unit |
|---------------------------------------------------------------------------------------|--------|---------------|-------------------------------------------------------------|-----------------|
|    | MR401B | 0011-30-37392 | Reusable Temp Probe, Adu, Esophageal/Rectal, 2 Pin, 3 m     | Each            |
|    | MR402B | 0011-30-37394 | Reusable Temp Probe, Ped/Neo, Esophageal/Rectal, 2 Pin, 3 m | Each            |
|   | MR403B | 0011-30-37393 | Reusable Temp Probe, Adu, Skin, 2 Pin, 3.6 m                | Each            |
|  | MR404B | 0011-30-37395 | Reusable Temp Probe, Ped/Neo, Skin, 2 Pin, 3.6 m            | Each            |

For MEC, PM series monitors

| Picture                                                                             | Model | Part No.      | Description                                                 | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------------------|-----------------|
|    | MR401 | 0011-30-90440 | Reusable Temp Probe, Adu, Esophageal/Rectal, Audio, 3 m     | Each            |
|    | MR402 | 0011-30-90441 | Reusable Temp Probe, Ped/Neo, Esophageal/Rectal, Audio, 3 m | Each            |
|   | MR403 | 0011-30-90442 | Reusable Temp Probe, Adu, Skin, Audio, 3.6 m                | Each            |
|  | MR404 | 0011-30-90443 | Reusable Temp Probe, Ped/Neo, Skin, Audio, 3.6 m            | Each            |

Temperature cables

- Extension cables are connected with the disposable temperature probes
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

Extension Cables: For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, BeneHeart defibrillator

| Picture                                                                             | Model  | Part No.      | Description                                                                    | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|--------------------------------------------------------------------------------|-----------------|
|  | MR420B | 0011-30-37391 | Extension Cable For Disposable Probes, 2 Pin, 3 m, round connector             | Each            |
|  |        | 040-001235-00 | Extension Cable For Disposable Probes, 2 Pin, 3 m, Rectangle Connector (Molex) | Each            |

Extension Cables: For MEC, PM series monitors

| Picture                                                                               | Model | Part No.      | Description                                                        | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------|-----------------|
|  | MR420 | 0011-30-90444 | Extension cable for disposable probes, Audio, 3 m, Round Connector | Each            |

| Picture                                                                           | Model | Part No.      | Description                                                                    | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------------------|-----------------|
|  |       | 040-001264-00 | Extension cable for disposable probes, Audio, 3 m, Rectangle Connector (Molex) | Each            |

**Adapter Cable: For BeneVision, BeneView, iPM, iMEC, uMEC series monitors, BeneHeart defibrillator**

| Picture                                                                           | Model | Part No.      | Description                                                                                                                                                     | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|  |       | 0010-30-43056 | TEMP Adapter Cable, 2 pin To Audio, 0,3 m<br>Note: With this adapter, Audio Temperature probe can be used for BeneView series monitors, BeneHeart defibrillator | Each            |

**Disposable Temperature Probes**

- Available in Rectal/Esophageal and Skin Surface Styles
- Accuracy ± 0.1°C
- Flexible and durable cables
- Latex free
- Good biocompatibility, avoiding allergic reactions to patient

**Match with the extension cables 0011-30-37391 and 0011-30-90444**

| Picture                                                                             | Model | Part No.      | Description                                                                         | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------------------------------------------|-----------------|
|  | MR411 | 0011-30-90446 | MR411 disposable temp probe, Adu/Ped/Neo, Esophageal/Rectal, round connector, 0,6 m | Each            |
|  | MR412 | 0011-30-90447 | MR411 disposable temp probe, Adu/Ped/Neo, Skin, round connector, 0,6 m              | Each            |

**Match with the extension cables 040-001235-00 and 040-001264-00**

| Picture                                                                               | Model | Part No.      | Description                                                               | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|---------------------------------------------------------------------------|-----------------|
|   | MR411 | 040-003294-00 | Disposable general purpose temperature probe, 9FR, Molex connector, 0,76m | 20 pcs/box      |
|  | MR412 | 040-003295-00 | Disposable temperature probe, skin surface, Molex connector, 0,76m        | 20 pcs/box      |

SmarTemp™ Probes

- Available in oral/axillary and rectal styles
- The probe cover is disposable, avoiding cross-infection
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free
- Good biocompatibility, avoiding allergicreactions to patient

For VS600/800/900

| Picture                                                                             | Part No.      | Description                          | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------------|--------------------------------------|-----------------|
|    | 6006-30-39598 | Oral/Axillary probe, 0.4~4 m         | Each            |
|   | 6006-30-39599 | Rectal probe, 0.4~4 m                | Each            |
|  | M09A-20-62062 | Blue sheath for Oral/Auxiliary probe | Each            |

| Picture                                                                             | Part No.         | Description                 | Purchasing Unit              |
|-------------------------------------------------------------------------------------|------------------|-----------------------------|------------------------------|
|  | M09A-20-62062-51 | Red sheath for Rectal probe | Each                         |
|  | M09A-20-62124    | Probe covers                | 20 pcs/box                   |
|                                                                                     | M09A-30-62128    | Probe covers                | 2000 pcs/case<br>(100 boxes) |

Tympanic Thermometer

For VS900

| Picture                                                                               | Part No.      | Description                   | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------------|-------------------------------|-----------------|
|   | 115-038397-00 | Genius 2 tympanic thermometer | Each            |
|  | 100-000200-00 | Probe covers for Genius 2     | 96pcs/box       |

CO<sub>2</sub> Accessories

Side-stream CO<sub>2</sub> Accessories (Mindray)

- Sampling line for intubated patients, Nasal sample cannula for non-intubated patients
- Sampling lines are disposable
- Water trap can be reused (Should be replaced every month)
- Latex free
- Good biocompatibility

For BeneView, iPM, iMEC, PM, MEC-2000 series monitors, BeneHeart defibrillator

| Picture                                                                             | Model       | Part No.                         | Description                                         | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------------|----------------------------------|-----------------------------------------------------|-----------------|
|   | 60-13100-00 | 115-043022-00<br>(9200-10-10530) | Dryline water trap, Adu/Ped<br>for dual-slot module | 10 pcs/box      |
|  | 60-13200-00 | 115-043023-00<br>(9200-10-10574) | Dryline water trap, Neo<br>for dual-slot module     | 10 pcs/box      |

For BeneVision, BeneView series monitors

| Picture                                                                             | Model | Part No.                         | Description                                                  | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|----------------------------------|--------------------------------------------------------------|-----------------|
|  |       | 115-043024-00<br>(100-000080-00) | M02C DRYLINE II water trap Adu/Ped<br>for single-slot module | 10 pcs/box      |
|  |       | 115-043025-00<br>(100-000081-00) | M02C DRYLINE II water trap Neo<br>for single-slot module     | 10 pcs/box      |

For BeneVision, BeneView, iPM, iMEC, PM, MEC-2000 series monitors, BeneHeart defibrillator

| Picture                                                                               | Model       | Part No.                         | Description                   | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------------|----------------------------------|-------------------------------|-----------------|
|   | 60-15200-00 | 115-043017-00<br>(9200-10-10533) | Sampling line, Adu/Ped, 2.5 m | 25 pcs/box      |
|  | 60-15300-00 | 115-043018-00<br>(9200-10-10555) | Sampling line, Neo, 2.5 m     | 25 pcs/box      |

| Picture                                                                             | Model       | Part No.                         | Description                                      | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------------|----------------------------------|--------------------------------------------------|-----------------|
|    | 4000        | 115-043001-00<br>(M02A-10-25937) | CO <sub>2</sub> Nasal sample cannula, Adu, 2.1 m | 25 pcs/box      |
|                                                                                     | 4100        | 115-043002-00<br>(M02A-10-25938) | CO <sub>2</sub> Nasal sample cannula, Pad, 2.1 m | 25 pcs/box      |
|                                                                                     | 4200        | 115-043003-00<br>(M02B-10-64509) | CO <sub>2</sub> Nasal sample cannula, Inf, 2.1 m | 25 pcs/box      |
|    | 60-14100-00 | 115-043020-00<br>(9000-10-07486) | Dryline airway adapter, straight                 | 10 pcs/box      |
|   | 60-14200-00 | 115-043021-00<br>(9000-10-07487) | Dryline airway adapter, elbow                    | 10 pcs/box      |
|  |             | 115-043019-00<br>(040-001187-00) | Dryline airway adapter, straight, neonate        | 10 pcs/box      |

Side-stream CO<sub>2</sub> Accessory Kit  
For BeneView, iPM, iMEC, uMEC, PM, MEC-2000 series monitors, BeneHeart defibrillator

| Part No.      | Description                                                                                                                                                                                                                                                                                                                | Purchasing Unit |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 6800-30-50618 | Side-stream CO <sub>2</sub> kit, Adu/Ped<br>Including:<br>60-13100-00 Dryline water trap, Adu, 2 pcs<br>4000 CO <sub>2</sub> Nasal sample cannula, Adu, 2 pcs<br>4100 CO <sub>2</sub> Nasal sample cannula, Ped, 2 pcs<br>60-15200-00 Sampling line, Adu, 2 pcs<br>60-14100-00 Dryline air way adapter, straight, 2 pcs    | Set             |
| 6800-30-50467 | Sides-tream CO <sub>2</sub> kit, Neo<br>Including:<br>60-13200-00 Dryline water trap, Neo, 2 pcs<br>4200 CO <sub>2</sub> Nasal sample cannula, Inf, 2 pcs<br>60-15300-00 Sampling line, Neo, 2 pcs<br>60-14100-00 Dryline airway adapter, straight, 2 pcs                                                                  | Set             |
| 115-021054-00 | M02C Sidestream CO <sub>2</sub> kit, Adu/Ped<br>Including:<br>M02C DRYLINE II water trap Adu/Ped, 2 pcs<br>4000 CO <sub>2</sub> Nasal sample cannula, Adu, 2 pcs<br>4100 CO <sub>2</sub> Nasal sample cannula, Ped, 2 pcs<br>60-15200-00 Sampling line, Adu, 2 pcs<br>60-14100-00 Dryline air way adapter, straight, 2 pcs | Set             |
| 115-021055-00 | M02C Sidestream CO <sub>2</sub> kit, Neo<br>Including:<br>M02C DRYLINE II water trap Neo, 2 pcs<br>4200 CO <sub>2</sub> Nasal sample cannula, Inf, 2 pcs<br>60-15300-00 Sampling line, Neo, 2 pcs<br>040-001187-00 Dryline air way adapter, straight, 2 pcs                                                                | Set             |





Micro-stream CO<sub>2</sub> Accessories (Oridion)

- Sampling line (FilterLine) for intubated patients, Nasal sample cannula and oral nasal cannula (Smart CapnoLine, CapnoLine and NIV-Line) for non-intubated patients
- Sampling lines are disposable
- Easy to use, without watertrap
- Latex free
- Good biocompatibility

For BeneVision, BeneView, iPM, iMEC, PM series monitors, BeneHeart defibrillator

| Picture                                                                             | Model   | Part No.      | Description                                | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|--------------------------------------------|-----------------|
|    | XS04620 | 0010-10-42560 | FilterLine Set Adult/Pediatric, 2 m        | Each            |
|                                                                                     | 007768  | 0010-10-42563 | FilterLine Set Adult/Pediatric Long, 4 m   | Each            |
|   | 008174  | 0010-10-42577 | NIV-Line Adult, 2 m                        | Each            |
|                                                                                     | 008175  | 0010-10-42578 | NIV-Line Pediatric, 2 m                    | Each            |
|  | XS04624 | 0010-10-42561 | FilterLine H Set Adult/Pediatric, 2 m      | Each            |
|                                                                                     | 007737  | 0010-10-42564 | FilterLine H Set Adult/Pediatric Long, 4 m | Each            |

| Picture                                                                               | Model  | Part No.      | Description                                    | Purchasing Unit |
|---------------------------------------------------------------------------------------|--------|---------------|------------------------------------------------|-----------------|
|    | 006324 | 0010-10-42562 | FilterLine H Set Infant/Neonatal, 2 m          | Each            |
|                                                                                       | 007738 | 0010-10-42565 | FilterLine H Set Infant/Neonatal Long, 4 m     | Each            |
|    | 009818 | 0010-10-42566 | Smart CapnoLine Plus Adult, 2 m                | Each            |
|   | 007266 | 0010-10-42567 | Smart CapnoLine Pediatric, 2 m                 | Each            |
|  | 009822 | 0010-10-42568 | Smart CapnoLine Plus O <sub>2</sub> Adult, 2 m | Each            |

| Picture                                                                             | Model  | Part No.      | Description                                          | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|------------------------------------------------------|-----------------|
|    | 009826 | 0010-10-42570 | Smart CapnoLine Plus O <sub>2</sub> Adult, Long, 4 m | Each            |
|                                                                                     |        |               |                                                      |                 |
|    | 007269 | 0010-10-42569 | Smart CapnoLine O <sub>2</sub> Pediatric, 2 m        | Each            |
|                                                                                     | 007743 | 0010-10-42571 | Smart CapnoLine O <sub>2</sub> Pediatric Long, 4 m   | Each            |
|   | 008177 | 0010-10-42572 | CapnoLine H Adult, 2 m                               | Each            |
|                                                                                     | 008179 | 0010-10-42574 | CapnoLine H Infant/Neonatal, 2 m                     | Each            |
|  | 008180 | 0010-10-42575 | CapnoLine H O <sub>2</sub> Adult, 2 m                | Each            |
|                                                                                     | 008181 | 0010-10-42576 | CapnoLine H O <sub>2</sub> Pediatric, 2 m            | Each            |

For BeneVision N1, VS-900, BeneHeart D3

| Picture                                                                             | Model | Part No.      | Description                        | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|------------------------------------|-----------------|
|  |       | 045-003134-00 | CO <sub>2</sub> adapter (reusable) | Each            |
|                                                                                     |       |               |                                    |                 |


Micro-stream CO<sub>2</sub> Accessory Kit

| Part No.                                                                     | Description                                                                                                                                     | Purchasing Unit |
|------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| For BeneVision, BeneView, iPM, iMEC series monitors, BeneHeart defibrillator |                                                                                                                                                 |                 |
| 040-000027-00                                                                | Microstream CO <sub>2</sub> kit, Ped<br>Including:<br>XS04620 FilterLine Set Adult/Pediatric 1 pcs<br>008175 NIV-Line Pediatric 1 pcs           | Set             |
| 040-000026-00                                                                | Microstream CO <sub>2</sub> kit, Neo<br>Including:<br>006324 FilterLine H Set Infant/Neonatal 1 pcs<br>008179 CapnoLine H Infant/Neonatal 1 pcs | Set             |
| 6800-30-50619                                                                | Microstream CO <sub>2</sub> kit, Adu<br>Including:<br>XS04620 FilterLine Set Adult/Pediatric, 1 pcs<br>008174 NIV-Line Adult, 1 pcs             | Set             |

Main-stream CO<sub>2</sub> Accessories (Capnostat)

- Designed for intubated patients. Non-intubated patient need to wear mask to use mainstream CO2 module
- Short response time
- Adapter and mask are disposable
- Latex free
- Good biocompatibility

For BeneVision, BeneView, iPM, iMEC series monitors

| Picture                                                                             | Model   | Part No.      | Description                                        | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------|---------------|----------------------------------------------------|-----------------|
|    | /       | 6800-30-50760 | Capnostat CO <sub>2</sub> sensor<br>(Cable: 2.4 m) | Each            |
|   | 9960LGE | 0010-10-42669 | Capnostat CO <sub>2</sub> mask, large Adu          | Each            |
|                                                                                     | 9960STD | 0010-10-42670 | Capnostat CO <sub>2</sub> mask, Adu                | Each            |
|                                                                                     | 9960PED | 0010-10-42671 | Capnostat CO <sub>2</sub> mask, Ped                | Each            |
|  | 6063    | 0010-10-42662 | Single-Patient Use Air way Adapter,<br>Adu/Ped     | Each            |

| Picture                                                                               | Model | Part No.      | Description                                                     | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|-----------------------------------------------------------------|-----------------|
|    | 6312  | 0010-10-42664 | Single-Patient Use Air way Adapter,<br>Neo                      | Each            |
|    | 6412  | 0010-10-42663 | Single-Patient Use Air way Adapter,<br>Adu/Ped, with mouthpiece | Each            |
|   | 8751  | 0010-10-42668 | Capnostat CO <sub>2</sub> sensor holding clips                  | 5 pcs/pack      |
|  | 6934  | 0010-10-42667 | Cable management straps                                         | 5 pcs/pack      |

Main-stream CO<sub>2</sub> Accessory Kit  
For BeneVision, BeneView, iPM, iMEC series monitors

| Part No.      | Description                                                                                                                                                                                                                                                  | Purchasing Unit |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 6800-30-50613 | Mainstream CO <sub>2</sub> kit, Adu/Ped/Neo<br>Including:<br>Capnostat CO <sub>2</sub> sensor. 1 pcs<br>6063 Single-Patient Use Airway Adapter, Adu/Ped, 1 pcs<br>6312 Single-Patient Use Airway Adapter, Neo, 1 pcs<br>8751 Cable management straps, 1 pack | Set             |

Anesthesia Gas Accessories

- Only for intubated patients
- Sampling lines are disposable
- Water trap can be reused (Should be replaced every month)
- Latex free
- Good biocompatibility

For BeneVision, BeneView series monitors

| Picture                                                                               | Model | Part No.                         | Description                                               | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|----------------------------------|-----------------------------------------------------------|-----------------|
|   |       | 115-043024-00<br>(100-000080-00) | M02C DRYLINE II water trap Adu/Ped for single-slot module | 10 pcs/box      |
|  |       | 115-043025-00<br>(100-000081-00) | M02C DRYLINE II water trap Neo for single-slot module     | 10 pcs/box      |

For BeneView, iPM, PM series, iMEC 15 monitors

| Picture                                                                           | Model       | Part No.                         | Description                                         | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------------|----------------------------------|-----------------------------------------------------|-----------------|
|  | 60-13100-00 | 115-043022-00<br>(9200-10-10530) | Dryline water trap, Adu/Ped<br>for dual-slot module | 10 pcs/box      |
|  | 60-13200-00 | 115-043023-00<br>(9200-10-10574) | Dryline water trap, Neo<br>for dual-slot module     | 10 pcs/box      |

For BeneVision, BeneView, iPM, PM series, iMEC 15 monitors

| Picture                                                                             | Model       | Part No.                         | Description                   | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------------|----------------------------------|-------------------------------|-----------------|
|   | 60-15300-00 | 115-043018-00<br>(9200-10-10555) | Sampling line, Neo, 2.5 m     | 25 pcs/box      |
|  | 60-15200-00 | 115-043017-00<br>(9200-10-10533) | Sampling line, Adu/Ped, 2.5 m | 25 pcs/box      |

| Picture                                                                              | Model       | Part No.                         | Description                                  | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------------|----------------------------------|----------------------------------------------|-----------------|
|   | 60-14100-00 | 115-043020-00<br>(9000-10-07486) | Dryline airway adapter, straight             | 10 pcs/box      |
|   | 60-14200-00 | 115-043021-00<br>(9000-10-07487) | Dryline airway adapter, elbow                | 10 pcs/box      |
|  |             | 115-043019-00<br>(040-001187-00) | Dryline airway adapter, straight,<br>neonate | 10 pcs/box      |

Anesthesia Gas Accessory Kit

| Part No.                                                 | Description                                                                                                                                                                                                                                                                                                                            | Purchasing Unit |
|----------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| <b>For BeneView series, iPM series, iMEC 15 monitors</b> |                                                                                                                                                                                                                                                                                                                                        |                 |
| 6800-30-50621                                            | Multi-Gas accessory kit<br>Including:<br>60-13100-00 Dryline water trap, Adu/Ped, 5 pcs<br>60-13200-00 Dryline water trap, Neo, 5 pcs<br>60-14100-00 Dryline airway adapter, straight, 2 pcs<br>60-14200-00 Dryline airway adapter, elbow, 2 pcs<br>60-15200-00 Sampling line, Adu/Ped, 4 pcs<br>60-15300-00 Sampling line, Neo, 4 pcs | Set             |
| <b>For PM series monitors</b>                            |                                                                                                                                                                                                                                                                                                                                        |                 |
| 9200-30-10575                                            | Multi-Gas accessory kit<br>Including:<br>60-13100-00 Dryline water trap, Adu/Ped, 5 pcs<br>60-13200-00 Dryline water trap, Neo, 5 pcs<br>60-14100-00 Dryline airway adapter, straight, 2 pcs<br>60-14200-00 Dryline airway adapter, elbow, 2 pcs<br>60-15200-00 Sampling line, Adu/Ped, 4 pcs<br>60-15300-00 Sampling line, Neo, 4 pcs | Set             |
| <b>For BeneVision, BeneView series monitors</b>          |                                                                                                                                                                                                                                                                                                                                        |                 |
| 115-034234-00                                            | Dual-slot Multi-Gat accessory kit<br>Including:<br>M02C Dryline II watertrap, Adu/Ped, 5pcs<br>60-14100-00 Dryline airway adapter, straight, 2pcs<br>60-14200-00 Dryline airway adapter, elbow, 2pcs<br>60-15200-00 Sampling line, Adu/Ped, 2pcs                                                                                       | Set             |

Invasive Blood Pressure (IBP) Accessories

Invasive Blood Pressure Cables

- Compatible solution with major monitor IBP module interface and disposable pressure transducer brands in the market
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

For BeneVision, BeneView, iPM, iMEC series monitors, BeneHeart defibrillator

| Picture                                                                              | Model  | Part No.      | Description                            | Purchasing Unit |
|--------------------------------------------------------------------------------------|--------|---------------|----------------------------------------|-----------------|
|  | IM2201 | 001C-30-70759 | 12 Pin IBP Cable (for ICU Medical), 4m | Each            |

| Picture                                                                             | Model  | Part No.      | Description                                                | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|------------------------------------------------------------|-----------------|
|    | IM2202 | 001C-30-70757 | 12 Pin IBP Cable (for BD), 4m                              | Each            |
|    | IM2207 | 0010-21-43082 | 12 Pin IBP Cable (for Memscap, SP844 82031 transducer), 4m | Each            |
|   | IM2211 | 0010-21-12179 | 12 Pin IBP Cable (for Edwards), 4m                         | Each            |
|  | IM2206 | 115-017849-00 | 12 Pin IBP cable (for Utah), 4m                            | Each            |

For MEC-2000, PM series monitors

| Picture                                                                               | Model  | Part No.      | Description                                               | Purchasing Unit |
|---------------------------------------------------------------------------------------|--------|---------------|-----------------------------------------------------------|-----------------|
|    | IM2101 | 001C-30-70760 | 6 Pin IBP Cable (for ICU Medical), 4m                     | Each            |
|    | IM2102 | 001C-30-70758 | 6 Pin IBP Cable (for BD), 4m                              | Each            |
|   | IM2103 | 0010-21-43094 | 6 Pin IBP Cable (for Edwards), 4m                         | Each            |
|  | IM2107 | 0010-21-43081 | 6 Pin IBP Cable (for Memscap, SP844 82031 transducer), 4m | Each            |



| Picture                                                                           | Model  | Part No.      | Description                    | Purchasing Unit |
|-----------------------------------------------------------------------------------|--------|---------------|--------------------------------|-----------------|
|  | IM2106 | 115-017848-00 | 6 Pin IBP cable (for Utah), 4m | Each            |

Adapter Cable: For BeneVision, BeneView, iPM, iMEC series monitors, BeneHeart defibrillator

| Picture                                                                           | Model | Part No.      | Description                         | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|-------------------------------------|-----------------|
|  |       | 0010-20-42795 | IBP cable adapter (12 pin to 6 pin) | Each            |

Y-type IBP cable: For BeneView, iPM series patient monitor

| Picture                                                                             | Model  | Part No.      | Description                                               | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|-----------------------------------------------------------|-----------------|
|  | IM2204 | 040-001029-00 | Y-type IBP cable (switch one connector to two connectors) | Each            |

Bracket and Mounting Clamp

| Picture                                                                              | Model  | Part No.      | Description                                               | Purchasing Unit |
|--------------------------------------------------------------------------------------|--------|---------------|-----------------------------------------------------------|-----------------|
|   | 42602  | M90-000133--- | Steady rest for IBP transducer (ICU Medical)              | Each            |
|   | 42394  | M90-000134--- | Mounting clamp for IBP transducer (ICU Medical)           | Each            |
|  | 682275 | 0010-10-12156 | Mounting background and accessory for IBP transducer (BD) | Each            |

Connected with Abbott (ICU Medical) Cables (001C-30-70759, 001C-30-70760, 115-020496-00)

| Picture                                                                           | Model | Part No.      | Description                         | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|-------------------------------------|-----------------|
|  | 42584 | 0010-10-42638 | Disposable transducer (ICU Medical) | 1 pc/box        |

Connected with BD cables (001C-30-70757, 001C-30-70758, 115-020492-00)

| Picture                                                                            | Model   | Part No.      | Description                | Purchasing Unit |
|------------------------------------------------------------------------------------|---------|---------------|----------------------------|-----------------|
|  | DT-4812 | 6000-10-02107 | Disposable transducer (BD) | 5 pcs/box       |

Invasive Blood Pressure Accessory Kit

For BeneVision, BeneView, iPM, iMEC series monitors, BeneHeart defibrillator


| Part No.      | Description                                                                                                                                                                                                                                                                                            | Purchasing Unit |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 6800-30-50876 | IBP accessory kit (ICU Medical, 12 Pin)<br>Including:<br>IM2201 12 Pin IBP Cable (for ICU Medical) 1 pcs<br>42584 Disposable transducer (ICU Medical) 5 pcs<br>M90-000133--- Steady rest for IBP transducer (ICU Medical) 1 pcs<br>M90-000134--- Mounting clamp for IBP transducer (ICU Medical) 1 pcs | Set             |
| 6800-30-50877 | IBP accessory kit (BD, 12 Pin)<br>Including:<br>IM2202 12 Pin IBP Cable (for BD) 1 pcs<br>DT-4812 Disposable transducer (BD) 5 pcs<br>0010-10-12156 Mounting background and accessory for IBP transducer (BD) 1 pcs                                                                                    | Set             |

For MEC-2000, PM series monitors

| Part No.      | Description                                                                                                                                                                                                                                                                                          | Purchasing Unit |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 6000-30-02164 | IBP accessory kit (ICU Medical, 6 Pin)<br>Including:<br>IM2101 6 Pin IBP Cable (for ICU Medical) 1 pcs<br>42584 Disposable transducer (ICU Medical) 5 pcs<br>M90-000133--- Steady rest for IBP transducer (ICU Medical) 1 pcs<br>M90-000134--- Mounting clamp for IBP transducer (ICU Medical) 1 pcs | Set             |
| 0010-30-42925 | IBP accessory kit (BD, 6 Pin)<br>Including:<br>IM2102 6 Pin IBP Cable (for BD) 1 pcs<br>DT-4812 Disposable transducer (BD) 5 pcs<br>0010-10-12156 Mounting background and accessory for IBP transducer (BD) 1 pcs                                                                                    | Set             |

Invasive Blood Pressure (ICP) Accessories

For BeneVision, BeneView, iPM, iMEC series monitors, BeneHeart defibrillator

| Picture                                                                            | Model | Part No.      | Description      | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|------------------|-----------------|
|   |       | 009-005460-00 | 12 Pin ICP cable | Each            |
|  |       | 040-002336-00 | ICP sensor kit   | Each            |

Respiratory Mechanics (RM) Accessories

- Precise volume, flow, and pressure data are obtained by RM module
- Respiratory mechanics monitoring is applicable to neonate

For BeneView series monitors

| Picture                                                                               | Model    | Part No.      | Description                                   | Purchasing Unit |
|---------------------------------------------------------------------------------------|----------|---------------|-----------------------------------------------|-----------------|
|    | PN155362 | 0010-30-42678 | Reusable flow sensor, adult/pediatric, 1.8 m  | Each            |
|   | PN155500 | 0010-30-42680 | Single patient use flow sensor, infant, 1.8 m | Each            |
|  | CP12901  | 6800-20-50328 | RM connecting plug                            | Each            |

For BeneVision, BeneView series monitors

| Picture                                                                           | Model  | Part No.      | Description                                                   | Purchasing Unit |
|-----------------------------------------------------------------------------------|--------|---------------|---------------------------------------------------------------|-----------------|
|  | MR4414 | 040-001948-00 | SPIRIT Single patient use flow sensor, neonate, 1.8m          | Each            |
|                                                                                   | MR4415 | 040-001950-00 | SPIRIT Single patient use flow sensor, neonate, 3.3m          | Each            |
|  | MR4412 | 040-001947-00 | SPIRIT Single patient use flow sensor, adult/pediatric, 1.8 m | Each            |
|                                                                                   | MR4413 | 040-001949-00 | SPIRIT Single patient use flow sensor, adult/pediatric, 3.3 m | Each            |

Respiratory Mechanics (RM) Accessory Kit

For BeneView series monitors

| Part No.      | Description                                                                                                                           | Purchasing Unit |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 6800-30-50612 | RM accessory kit, Adu/Ped<br>Including:<br>PN155362 Reusable flow sensor, adult/pediatric, 1 pcs<br>CP12901 RM connecting plug, 1 pcs | Set             |

Cardiac Output Accessories

Cardiac Output Cables

- Secure connection to ensure clear transmission signal
- Outstanding cable material, enduring repeated cleaning and disinfection
- Convenient for cleaning
- Latex-free

For BeneVision, BeneView, iPM series monitors

| Picture                                                                              | Model  | Part No.      | Description       | Purchasing Unit |
|--------------------------------------------------------------------------------------|--------|---------------|-------------------|-----------------|
|  | CO7702 | 0010-30-42743 | 12 Pin C.O. Cable | Each            |

For MEC-2000, PM series monitors

| Picture                                                                               | Model  | Part No.      | Description      | Purchasing Unit |
|---------------------------------------------------------------------------------------|--------|---------------|------------------|-----------------|
|  | CO7703 | 900E-30-04952 | 6 Pin C.O. Cable | Each            |

Cardiac Output Sensors

- Secure connection to ensure clear transmission signal
- Latex-free
- Good biocompatibility, avoiding allergic reactions to patient

Connected with all the CO cables (CO7702 and CO7703)

| Picture                                                                             | Model  | Part No.      | Description                                      | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|--------------------------------------------------|-----------------|
|    | SP4042 | 6000-10-02079 | Injection liquid temperature sensor (BD)         | Each            |
|    | SP5045 | 6000-10-02080 | Injection liquid temperature sensor housing (BD) | Each            |
|  | MX387  | 6000-10-02081 | 12CC control syringe w/rotator (Medex)           | Each            |

Cardiac Output Accessory Kit

For BeneVision, BeneView, iPM series monitors

| Picture                                                                             | Part No.      | Description                                                                                                                                                                                                                 | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|  | 6800-30-50617 | C.O. accessory kit<br>Including:<br>CO7702 12 Pin C.O. Cable,<br>SP4042 Injection liquid temperature sensor (BD)<br>SP5045 Injection liquid temperature sensor housing (BD)<br>MX387 12CC control syringe w/rotator (Medex) | Set             |

For MEC-2000, PM series monitors

| Picture                                                                             | Part No.      | Description                                                                                                                                                                                                               | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|  | 6000-30-02165 | C.O. accessory kit<br>Including:<br>CO7703 6 Pin C.O. Cable<br>SP4042 Injection liquid temperature sensor (BD)<br>SP5045 Injection liquid temperature sensor housing (BD)<br>MX387 12CC control syringe w/rotator (Medex) | Set             |

Impedance Cardiography (ICG) Accessories

For BeneVision, BeneView series monitors

| Picture                                                                            | Model   | Part No.      | Description              | Purchasing Unit |
|------------------------------------------------------------------------------------|---------|---------------|--------------------------|-----------------|
|   | N1301-3 | 100-000149-00 | ICG patient cable-normal | Each            |
|  | N1201   | 100-000148-00 | ICG sensor               | 4pcs/pouch      |

Impedance Cardiography (ICG) Accessory Kit

For BeneVision, BeneView series monitors

| Part No.      | Description                                                                                                 | Purchasing Unit |
|---------------|-------------------------------------------------------------------------------------------------------------|-----------------|
| 115-034123-00 | ICG accessory kit (normal)<br>Including:<br>N1201 ICG sensor 4 pcs<br>N1301-3 ICG patient cable-normal 1pcs | Set             |

Bispectral Index (BIS) Accessories

Bispectral Index (BIS) Cables

- Consisting of BIS module and patient cable, connected to monitors
- Outstanding cable material, enduring repeated cleaning and disinfection
- Convenient for cleaning
- Latex-free

For BeneVision, BeneView series monitors

| Picture                                                                             | Model       | Part No.      | Description             | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------------|---------------|-------------------------|-----------------|
|   | 186-0195-MR | 6800-30-50761 | BIS module and cables   | Each            |
|  | 186-0224-MR | 115-005707-00 | BISx4 module and cables | Each            |

Bispectral Index (BIS) Sensors

- Noninvasive, suitable for ICU, OR and clinical study
- Flexible design adjusts to different head sizes
- Connector provides secure click-in connection with push button release
- Single patient use, avoiding cross-infection
- Latex-free
- Good biocompatibility, avoiding allergic reactions to patient

For BeneVision, BeneView series monitors

| Picture                                                                               | Model    | Part No.      | Description                            | Purchasing Unit |
|---------------------------------------------------------------------------------------|----------|---------------|----------------------------------------|-----------------|
|    | 186-0106 | 0010-10-42672 | BIS sensor, Adu (QUATRO, 186-0106)     | Each            |
|   | 186-0200 | 0010-10-42673 | BIS sensor, Ped (QUATRO, 186-0200)     | Each            |
|  | 186-0212 | 040-000392-00 | BISx4 bilateral sensor, Adu (186-0212) | Each            |

Bispectral Index (BIS) Accessory Kits

For BeneVision, BeneView series monitors

| Part No.      | Description                                                                                                                       | Purchasing Unit |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 6800-30-50878 | BIS accessory kit, Adu<br>Including:<br>186-0195-MR BIS module and cables<br>186-0150-MR BIS sensor kit (5 pcs/pkg)               | Set             |
| 6800-30-50144 | BIS accessory kit, Ped<br>Including:<br>186-0195-MR BIS module and cables<br>186-0206-MR BIS sensor kit (5 pcs/pkg)               | Set             |
| 115-005614-00 | BISx4 accessory kit, Adu<br>Including:<br>186-0224-MR BISx4 module and cables<br>186-0212 BISx4 bilateral sensor, Adu (5 pcs/pkg) | Set             |

PiCCO<sub>2</sub> Accessories

PiCCO<sub>2</sub> Cables

- Secure connection to ensure clear transmission signal
- Outstanding cable material, enduring repeated cleaning and disinfection
- Convenient for cleaning
- Latex-free

For BeneVision, BeneView series monitors

| Picture                                                                               | Model  | Part No.      | Description                                                                           | Purchasing Unit |
|---------------------------------------------------------------------------------------|--------|---------------|---------------------------------------------------------------------------------------|-----------------|
|   | CO7701 | 040-000816-00 | PiCCO <sub>2</sub> 12 Pin Cable (for injecta temperature and blood temperature)       | Each            |
|  | IM2203 | 040-000815-00 | PiCCO <sub>2</sub> Y-Type IBP Cable (for artery pressure and central venuos pressure) | Each            |





For BeneVision, BeneView series monitors

| Part No.      | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Purchasing Unit |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| 115-041293-00 | PiCCO <sub>2</sub> Accessory kit, Pediatric, 2pcs IBP Cable (Utah+BD/Argon)<br>Including:<br>12 Pin IBP Y Cable (040-002827-00)<br>12 Pin Utah IBP Cable (115-017849-00)<br>12 Pin BD IBP Cable (001C-30-70757)<br>CO7701 12 Pin PiCCO <sub>2</sub> Cable<br>PC80105 2 Pin Injecta Temperature Sensor Cable<br>PV2013L07 PiCCO <sub>2</sub> catheter Pediatric<br>PV8215 PiCCO <sub>2</sub> Monitoring kit<br>PiCCO <sub>2</sub> Mounting clamp for PV8215 (040-002903-00)             | Set             |
| 115-041296-00 | PiCCO <sub>2</sub> Accessory kit, Adult, 2pcs IBP Cable (Utah+ICU Medical)<br>Including:<br>12 Pin IBP Y Cable (040-002827-00)<br>12 Pin Utah IBP Cable (115-017849-00)<br>12 Pin ICU Medical IBP Cable (001C-30-70759)<br>CO7701 12 Pin PiCCO <sub>2</sub> Cable<br>PC80105 2 Pin Injecta Temperature Sensor Cable<br>2015L20N PiCCO <sub>2</sub> catheter Adult<br>PV8215 PiCCO <sub>2</sub> Monitoring kit<br>PiCCO <sub>2</sub> Mounting clamp for PV8215 (040-002903-00)          | Set             |
| 115-041294-00 | PiCCO <sub>2</sub> Accessory kit, Pediatric, 2pcs IBP Cable (Utah+ICU Medical)<br>Including:<br>12 Pin IBP Y Cable (040-002827-00)<br>12 Pin Utah IBP Cable (115-017849-00)<br>12 Pin ICU Medical IBP Cable (001C-30-70759)<br>CO7701 12 Pin PiCCO <sub>2</sub> Cable<br>PC80105 2 Pin Injecta Temperature Sensor Cable<br>PV2013L07 PiCCO <sub>2</sub> catheter Pediatric<br>PV8215 PiCCO <sub>2</sub> Monitoring kit<br>PiCCO <sub>2</sub> Mounting clamp for PV8215 (040-002903-00) | Set             |

ScvO<sub>2</sub> Accessories

For BeneVision, BeneView series monitors



| Picture                                                                             | Model     | Part No.      | Description                     | Purchasing Unit |
|-------------------------------------------------------------------------------------|-----------|---------------|---------------------------------|-----------------|
|  | PV2022-35 | 040-000920-00 | CeVOX fiber optic probe (35 cm) | Each            |
|                                                                                     | PV2013L07 | 040-000919-00 | CeVOX fiber optic probe (37 cm) | Each            |


ScvO<sub>2</sub> Accessory Kit

| Part No.      | Description                                        | Purchasing Unit |
|---------------|----------------------------------------------------|-----------------|
| 115-008293-00 | CeVOX Optical Module + Cable + 35 cm Optical Probe | Set             |
| 115-008294-00 | CeVOX Optical Module + Cable + 37 cm Optical Probe | Set             |

EEG Accessories


For BeneVision, BeneView series monitors

| Picture                                                                             | Model | Part No.      | Description                                                            | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------|-----------------|
|    |       | 040-001594-00 | EEG cable, 14 PIN                                                      | Each            |
|   |       | 040-001596-00 | EEG subdermal needle electrodes (disposable, L=10 mm, adu/ped, 10 pcs) | Each            |
|  |       | 040-001598-00 | EEG cup electrodes (reusable, Ag/AgCl, D=10 mm, adu/ped, 10 pcs)       | Each            |

| Picture                                                                              | Model | Part No.      | Description                                                      | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------|-----------------|
|   |       | 040-001602-00 | EEG cup electrodes (reusable, Ag/AgCl, D=10 mm, ped/neo, 10 pcs) | Each            |
|   |       | 040-001595-00 | EEG scrub cream                                                  | Each            |
|  |       | 040-001597-00 | EEG conductive cream                                             | Each            |


EEG Accessory Kit

For BeneVision, BeneView series monitors

| Picture                                                                            | Part No.      | Description                                                                                                                                                                                                                                  | Purchasing Unit |
|------------------------------------------------------------------------------------|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|   | 115-018153-00 | EEG accessory kit (needle electrodes, adu/ped)<br>Including:<br>040-001594-00 EEG cable<br>040-001596-00 EEG subdermal needle electrodes (disposable, L=10 mm, adu/ped, 10 pcs)                                                              | set             |
|   | 115-018154-00 | EEG accessory kit (cup eletrodes, adu/ped)<br>Including:<br>040-001594-00 EEG cable<br>040-001598-00 EEG cup electrodes (reusable, Ag/AgCl, D=10 mm, adu/ped, 10 pcs)<br>040-001595-00 EEG scrub cream<br>040-001597-00 EEG conductive cream | set             |
|  | 115-018155-00 | EEG accessory kit (cup eletrodes, ped/neo)<br>Including:<br>040-001594-00 EEG cable<br>040-001602-00 EEG cup electrodes (reusable, Ag/AgCl, D=10 mm, ped/neo, 10 pcs)<br>040-001595-00 EEG scrub cream<br>040-001597-00 EEG conductive cream | set             |

NMT Accessories

For BeneVision, BeneView series monitors

| Picture                                                                               | Model   | Part No.      | Description           | Purchasing Unit |
|---------------------------------------------------------------------------------------|---------|---------------|-----------------------|-----------------|
|    | NM13101 | 040-001462-00 | NMT main cable        | Each            |
|   | NM13401 | 040-001463-00 | NMT transducer cable  | Each            |
|  | NM13701 | 040-001464-00 | NMT stimulation cable | Each            |

For BeneVision, BeneView series monitors

| Picture                                                                           | Part No.      | Description                              | Purchasing Unit |
|-----------------------------------------------------------------------------------|---------------|------------------------------------------|-----------------|
|  | 9000-10-07469 | Pediatric ECG electrode (3M, 2245)       | 50 pcs/pouch    |
|                                                                                   | 040-002833-00 | Pediatric/neonatat ECG electrode (INTCO) | 30 pcs/pouch    |
|  | 040-002258-00 | Bandage for NMT transducer, disposable   | 20 pcs/pouch    |

NMT Accessory Kit

For BeneVision, BeneView series monitors

| Picture                                                                             | Part No.      | Description                                                                                                                                                                                                                                                                   | Purchasing Unit |
|-------------------------------------------------------------------------------------|---------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|  | 115-018586-00 | NMT accessory kit<br>Including:<br>040-001462-00 NMT main cable<br>040-001463-00 NMT transducer cable<br>040-001464-00 NMT stimulation cable<br>040-002833-00 Pediatric/neonatat ECG electrode (INTCO), 30 pcs<br>040-002258-00 Bandage for NMT transducer, disposable, 20pcs | set             |

INVOS rSO<sub>2</sub> Accessories

For BeneVision series monitors

| Picture                                                                               | Model | Part No.      | Description                                                                                                                                                                                                                                                                           | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|    |       | 115-033947-00 | INVOS rSO <sub>2</sub> accessory kit, Adu<br>Including:<br>INVOS 5100C cable, channel 1&2, with Pre-ampliter A 1pcs<br>INVOS 5100C reusable sensor, channel 1 (blue) 1pcs<br>INVOS 5100C reusable sensor, channel 2(brown) 1pcs<br>INVOS SomaSensor disposable sensor, Adu>40kg, 2pcs | Each            |
|   |       | 115-033948-00 | INVOS rSO <sub>2</sub> accessory kit, Ped<br>Including:<br>INVOS 5100C cable, channel 1&2, with Pre-ampliter A 1pcs<br>INVOS 5100C reusable sensor, channel 1 (blue) 1pcs<br>INVOS 5100C reusable sensor, channel 2(brown) 1pcs<br>INVOS SomaSensor disposable sensor, Ped<40kg, 2pcs | Each            |
|  |       | 115-033949-00 | INVOS rSO <sub>2</sub> accessory kit, Neo<br>Including:<br>INVOS 5100C cable, channel 1&2, with Pre-ampliter A 1pcs<br>INVOS cable + disposable sensor, Neo<5kg, brain/body, 2pcs                                                                                                     | Each            |

Wall Mount Bracket

| Picture                                                                             | Model | Part No.      | Description                                            | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------|-----------------|
|    |       | 045-000926-00 | T5 wall mount bracket                                  | Each            |
|    |       | 045-000929-00 | T8 wall mount bracket                                  | Each            |
|   |       | 045-000936-00 | T8 display wall mount bracket                          | Each            |
|  |       | 045-000931-00 | iPM/iMEC/uMEC wall mount bracket<br>(iMEC 15 excluded) | Each            |
|                                                                                     |       | 045-000934-00 | T8 keyboard wall mount bracket                         | Each            |

| Picture                                                                               | Model | Part No.      | Description                     | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|---------------------------------|-----------------|
|    |       | 045-000942-00 | PM-7000/8000 wall mount bracket | Each            |
|    |       | 045-000947-00 | PM-9000 wall mount bracket      | Each            |
|   |       | 045-000949-00 | VS-800 wall mount bracket       | Each            |
|  |       | 045-001059-00 | VS-600/900 wall mount bracket   | Each            |

| Picture                                                                            | Model | Part No.      | Description                                                                              | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------------------------|-----------------|
|   |       | 045-000933-00 | SMR wall mount bracket                                                                   | Each            |
|   |       | 045-003054-00 | VS-600/900 basic wall mount bracket                                                      | Each            |
|  |       | 045-003055-00 | iPM/iMEC/uMEC basic wall mount bracket (iMEC 15 excluded)                                | Each            |
|                                                                                    |       | 045-003312-00 | iMEC 15 wall mount bracket, GCX wall bracket + adapter for iMEC15 (Connected plate 6308) | Each            |

Trolley

| Picture                                                                              | Model | Part No.      | Description                                                                                                                          | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|   |       | 045-000915-00 | N17/N15/T5 rolling stand                                                                                                             | Each            |
|   |       | 045-000921-00 | T8 rolling stand                                                                                                                     | Each            |
|  |       | 045-000924-00 | N12/iPM/iMEC/uMEC rolling stand (iMEC 15/uMEC 15 excluded)                                                                           | Each            |
|                                                                                      |       | 045-003255-00 | N12/iPM/iMEC/uMEC rolling stand, ≤6kg, fixed-angle, quick-lock, with two brackets and barrel fix mounting (iMEC 15/uMEC 15 excluded) | Each            |

| Picture                                                                             | Model | Part No.      | Description                                                                                                | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------------------------------------------------|-----------------|
|    |       | 045-003052-00 | VS-600/900 basic rolling stand                                                                             | Each            |
|    |       | 045-003053-00 | iPM/iMEC/uMEC basic rolling stand (iMEC 15 excluded)                                                       | Each            |
|                                                                                     |       | 045-003311-00 | iMEC 15 rolling stand, ≤8kg, angle adjustable, with two baskets+ adapter for iMEC15( Connected plate 6308) | Each            |
|  |       | 045-000905-00 | PM-7000/8000 rolling stand                                                                                 | Each            |

| Picture                                                                               | Model | Part No.      | Description                         | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|-------------------------------------|-----------------|
|    |       | 045-000908-00 | PM-9000 rolling stand               | Each            |
|    |       | 045-000911-00 | VS-800 rolling stand                | Each            |
|   |       | 045-001057-00 | VS-600/900 rolling stand            | Each            |
|  |       | 045-002960-00 | VS-600/900 rolling stand (enhanced) | Each            |





| Picture                                                                           | Model | Part No.      | Description                                         | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|-----------------------------------------------------|-----------------|
|  |       | 8000-30-90169 | Bedrail hook                                        | Each            |
|                                                                                   |       | 8000-30-90170 | Bedrail clamp                                       | Each            |
|                                                                                   |       | 115-044101-00 | Transition plate for bedrail mounting(for uMEC6/10) | Each            |
|                                                                                   |       | 115-044102-00 | Transition plate for bedrail mounting(for uMEC7/12) | Each            |

Others

| Picture                                                                              | Model | Part No.      | Description                             | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------|---------------|-----------------------------------------|-----------------|
|   |       | 8000-21-10361 | Nurse call cable for T6/T8/T9           | Each            |
|                                                                                      |       | 009-003116-00 | Nurse call cable for iPM / iMEC series  | Each            |
|                                                                                      |       | 009-003118-00 | Defib. sync cable for iPM / iMEC series | Each            |
|  |       | 009-000027-00 | Defib. sync cable for T6/T8/T9          | Each            |
|                                                                                      |       | 009-003249-00 | Defib. sync cable for T1                | Each            |

| Picture                                                                           | Model | Part No.      | Description                                  | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|----------------------------------------------|-----------------|
|  |       | 009-003117-00 | Analog output cable<br>for iPM / iMEC series | Each            |
|  |       | 6800-20-50146 | Analog output cable<br>for T6/T8/T9          | Each            |
|                                                                                   |       | 009-002944-00 | Analog output cable<br>for T1                | Each            |
|                                                                                   |       | A30-000001--- | Thermal paper(50mmx20m)                      | Each            |

| Picture                                                                             | Model | Part No.      | Description                        | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|------------------------------------|-----------------|
|  |       | 009-003650-00 | cable protecting tube, 20cm        | 5pcs/pkg        |
|                                                                                     |       | 009-003649-00 | cable protecting tube, 40cm        | 5pcs/pkg        |
|  |       | 009-003903-00 | Accessories management tape , 62cm | 2pcs/pkg        |




Electrocardiograph accessories

- Apply for 12-lead ECG acquisition on electrocardiograph
- Both ECG cable with banana plugs, Limb Clamps, Chest Bulbs and ECG cable with Electrode clips are available
- Both IEC and AHA accessories are available

For BeneHeart R3/R12

| Picture                                                                            | Model | Part No.      | Description                             | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|-----------------------------------------|-----------------|
|   |       | 040-001579-00 | 12-lead ECG cable, IEC, Φ4 banana plugs | Each            |
|  |       | 040-001582-00 | 12-lead ECG cable, AHA, Φ4 banana plugs | Each            |

| Picture                                                                               | Model | Part No.      | Description                                                   | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|---------------------------------------------------------------|-----------------|
|    |       | 040-001587-00 | Limb clamps, 4/set, IEC                                       | 4 pcs/set       |
|    |       | 040-001586-00 | Limb clamps, 4/set, AHA                                       | 4 pcs/set       |
|    |       | 040-001585-00 | Chest bulbs, 6/set                                            | 6 pcs/set       |
|   |       | 040-001645-00 | 12-lead ECG cable, IEC, Clip                                  | Each            |
|                                                                                       |       | 040-001643-00 | 12-lead ECG cable, AHA, Clip                                  | Each            |
|  |       | 115-024339-00 | 12-lead ECG cable + tab adapters +tab electrode, Adu/ped, IEC | Each            |
|                                                                                       |       | 115-024341-00 | 12-lead ECG cable + tab adapters +tab electrode, Adu/ped, AHA | Each            |
|                                                                                       |       | 040-001646-00 | Tab adapters (use between the cable and tab electrodes)       | 10 pcs/set      |

| Picture                                                                            | Model | Part No.      | Description                              | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|------------------------------------------|-----------------|
|   |       | 040-001908-00 | Tab electrode                            | 100 pcs/pkg     |
|   |       | 115-024338-00 | Carrying bag for R3                      | Each            |
|   |       | 045-001385-00 | Trolley for R3, with case and cable arm  | Each            |
|  |       | 045-001384-00 | Trolley for R12, with case and cable arm | Each            |

| Picture | Model | Part No.      | Description                                  | Purchasing Unit |
|---------|-------|---------------|----------------------------------------------|-----------------|
|         |       | M002-10-69954 | Rolling thermal paper (80mmx20m)             | Each            |
|         |       | 095-002708-00 | Folding thermal paper (80mmx14m)             | Each            |
|         |       | 095-002775-00 | Thermal paper, A4 size, 210mmx295mm<br>100p  | Each            |
|         |       | 095-002773-00 | Thermal paper, A4 size, 210mmx295mm<br>150p  |                 |
|         |       | 095-002776-00 | Thermal paper, US letter 215mmx280mm<br>100p | Each            |
|         |       | 095-002774-00 | Thermal paper, US letter 215mmx280mm<br>150p |                 |



Defibrillator Accessories



Reusable Internal paddles

- Apply for exposed heart defibrillation in OR during the heart surgery
- Integrate both cable and paddles in one cable
- Offer three paddle sizes to be chosen based on different patient types (1", 2" and 3")
- Auto-identify the paddle connection and reduce the energy level under 50J (default setting is 10J) for internal defibrillation use by BeneHeart D6 and D3
- Autoclavable

For BeneHeart defibrillator D6/D3

| Picture                                                                               | Model | Part No.      | Description              | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|--------------------------|-----------------|
|    |       | 115-018366-00 | Internal paddles, 1 inch | Each            |
|   |       | 115-018367-00 | Internal paddles, 2 inch | Each            |
|  |       | 115-018368-00 | Internal paddles, 3 inch | Each            |

External Paddles and Cables

- Applicable for both adults and pediatric patients, and easy to switch
- Safe for defibrillation energy delivery
- Space-saving spiral cable
- Flexible and durable cables
- Outstanding cable material, enduring repeated cleaning and disinfection
- Latex free

For BeneHeart defibrillator D6/D3

| Picture                                                                           | Model | Part No.      | Description                                                         | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|---------------------------------------------------------------------|-----------------|
|  |       | 0651-30-76994 | External paddles kit, Adu/Ped, English, with conductive gel (250 g) | Each            |

Pads Cable

- High voltage resistance and high safety
- Flexible and durable cables
- Outstanding cable material enduring repeated cleaning and disinfection

For BeneHeart defibrillator D6/D3

| Picture                                                                             | Model  | Part No.      | Description                                     | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|-------------------------------------------------|-----------------|
|  | MR6702 | 040-000545-00 | Cable of electrode pads with test load (50 ohm) | Each            |

Multifunctional Pads

- Two models for either adult or pediatric
- Multifunctional for defibrillation, pacing, and ECG monitoring
- Wide applicable temperature scope
- Highly adhesive
- Outstanding cable material enduring repeated cleaning and disinfection
- Disposable after use, and no pollution to the environment

For BeneHeart defibrillator D6/D3

| Picture                                                                              | Model | Part No.      | Description                                                                                      | Purchasing Unit |
|--------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------------------------------------|-----------------|
|   | MR60  | 0651-30-77007 | Adult multifunction electrode pads, cable length 1.2 m preconnectable (0.45m cable out-of-pouch) | 5pcs/box        |
|  | MR61  | 0651-30-77008 | Pediatric multifunction electrode pads, cable length 1.2 m (no cable out-of-pouch)               | 5pcs/box        |

For BeneHeart defibrillator D1

| Picture                                                                           | Model | Part No.      | Description                                                      | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------------|-----------------|
|  | MR62  | 115-035426-00 | Adult multifunction electrode pads, with auto-identification     | 5pcs/box        |
|  | MR63  | 115-035427-00 | Pediatric multifunction electrode pads, with auto-identification | 5pcs/box        |

Test Load

- Used to test the performance of the main unit and multifunctional cable
- High reliability and safety

For BeneHeart defibrillator D6/D3

| Picture                                                                             | Model  | Part No.      | Description                               | Purchasing Unit |
|-------------------------------------------------------------------------------------|--------|---------------|-------------------------------------------|-----------------|
|  | MR6905 | 040-000413-00 | Test load<br>(for use with 0651-30-77010) | Each            |

CPR Sensor

For BeneHeart defibrillator D6/D3

| Picture                                                                             | Model | Part No.      | Description                                                                                                                                      | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
|  |       | 115-044836-00 | CPR sensor kit (with battery)<br>Including:<br>CPR sensor with battery<br>CPR sensor cable<br>Disposable CPR sensor adhesive tape,<br>3pcs       | Each            |
|  |       | 115-044871-00 | CPR sensor kit (without battery)<br>Including:<br>CPR sensor without battery<br>CPR sensor cable<br>Disposable CPR sensor adhesive tape,<br>3pcs | Each            |



Carrying Case

- Carrying a defibrillator for pre-hospital or field emergency treatment
- High safety coefficient in load bearing
- High applicability in extreme environments
- Outstanding material enduring repeated cleaning and disinfection

| Picture                                                                            | Model | Part No.      | Description                                            | Purchasing Unit |
|------------------------------------------------------------------------------------|-------|---------------|--------------------------------------------------------|-----------------|
|   |       | 115-018610-00 | BeneHeart D3 package & protection cover                | Each            |
|   |       | 115-008361-00 | BeneHeart D6 luxury carrying case and shield cover kit | Each            |
|  |       | 115-023421-00 | D1 carrying case with rescue kit                       | Each            |

Conductive Gel

| Picture                                                                             | Model | Part No.      | Description           | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|-----------------------|-----------------|
|  | 15-25 | 0000-10-10775 | Conductive gel, 250 g | Each            |

Table Mounting

| Picture                                                                               | Model | Part No.      | Description                                                                         | Purchasing Unit |
|---------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------------------------------------------|-----------------|
|    |       | 115-013412-00 | Vehicle/ table mounting for BeneHeart D3                                            | Each            |
|    |       | 115-005061-00 | Vehicle/ table mounting for BeneHeart D6                                            | Each            |
|    |       | 115-030528-00 | D6 vehicle dock station CE, DC power input connector, with cigarette DC power cable | Each            |
|  |       | 115-030529-00 | D6 vehicle dock station CE, AC power input connector, without AC power cable        | Each            |



Battery Charger

| Picture                                                                           | Model | Part No.      | Description                                                | Purchasing Unit |
|-----------------------------------------------------------------------------------|-------|---------------|------------------------------------------------------------|-----------------|
|  |       | 115-009187-00 | External Li-ion battery charger + 1 power cord (GB)        | Each            |
|                                                                                   |       | 115-009188-00 | External Li-ion battery charger + 1 power cord (US)        | Each            |
|                                                                                   |       | 115-009189-00 | External Li-ion battery charger + 1 power cord (India)     | Each            |
|                                                                                   |       | 115-009190-00 | External Li-ion battery charger + 1 power cord (EU)        | Each            |
|                                                                                   |       | 115-009191-00 | External Li-ion battery charger + 1 power cord (Brazil)    | Each            |
|                                                                                   |       | 115-009192-00 | External Li-ion battery charger + 1 power cord (UK)        | Each            |
|                                                                                   |       | 115-025630-00 | External Li-ion battery charger + 1 power cord (Australia) | Each            |
|                                                                                   |       | 115-033660-00 | External Li-ion battery charger + 1 power cord (Swiss)     | Each            |

Thermal Paper

| Picture | Model | Part No.      | Description             | Purchasing Unit |
|---------|-------|---------------|-------------------------|-----------------|
|         |       | A30-000001--- | Thermal paper(50mmx20m) | Each            |
|         |       | M002-10-69954 | Thermal paper(80mmx20m) | Each            |

Mobile Trolley

| Picture                                                                             | Model | Part No.      | Description                                                                                                 | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------------------------------------------------------------------|-----------------|
|  |       | 115-015823-00 | Mobile trolley kit for D6<br>Including:<br>115-005061-00 mounting for D6<br>115-013418-00 trolley for D3/D6 | Each            |
|                                                                                     |       | 115-015825-00 | Mobile trolley kit for D3<br>Including:<br>115-013412-00 mounting for D3<br>115-013418-00 trolley for D3/D6 | Each            |

Cabinet

| Picture                                                                             | Model | Part No.      | Description                                     | Purchasing Unit |
|-------------------------------------------------------------------------------------|-------|---------------|-------------------------------------------------|-----------------|
|  |       | 045-001140-00 | AED cabinet (with mounting kit, lock and alarm) | Each            |