

Liebel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

P/N: 801107
Lengths: 230 cm
Injectors: Liebel-Flarsheim Optistar MR



Contents: Y-connector with
230 cm patient line and single check valve
Package: 100 units/box
Pressure Limit Syringe: 200 PSI / 14 bar

Materials: Polycarbonate
Polyvinyl Chloride
Silicone
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer: Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237

Country of Origin: USA



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

