

Liabel-Flarsheim

Single Use Empty Syringes and Connecting Lines

Syringe and Tube Components Sheet

P/N: 801800
Size: 2 x 60 ml Syringe and
230 cm Connecting Line
Injectors: Liabel-Flarsheim Optistar LE



Contents: 2 x 60 ml Syringe with spikes and Y-connector with 230 cm patient line

Package: 50 units/box

Pressure Limit Syringe : 200 PSI / 14 bar

Materials: Polypropylene
Synthetic Rubber
Polycarbonate
Polyvinyl Chloride
Silicone
Latex Free

Sterilisation: Ethylene Oxide (EO)

Shelf Life: 3 years from Date of Sterilisation

Use: Single Use – do not reuse

CE: 0123

Manufacturer : Liabel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237

Country of Origin: USA



Tyco Healthcare Europe
Liabel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany

