

Mallinckrodt Single Use Connection Lines

Tubes Components Sheet

P/N: 844010
Lengths: 150 cm
Description: Coiled Y-Tube with 13 cm Legs




Connectors: LL f/m
Package: 50 units/box
Pressure Limit: 400 PSI / 28 bar

Materials: Polycarbonate
PVC
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 3 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237
Country of Origin: USA



Tyco Healthcare Deutschland GmbH
Josef-Dietzgen-Str. 1
53761 Hennef
Germany