

Liebel-Flarsheim Single Use Empty Syringes

Syringe Components Sheet

P/N: 900105
Size: 200 ml
Injectors: Liebel-Flarsheim Angiomat Illumena



Contents: Syringe with Handfill

Package: 50 units/box

Pressure Limit Syringe : 900 PSI / 62 bar

Pressure Limit Tube : 200 PSI / 21 bar

Connectors: Linden Luer

Materials:
Polypropylene
Polycarbonate
Synthetic Rubber
Polyvinyl Chloride
Polyethylen

Sterilisation: Ethylene Oxide (EO)

Shelf Life: 3 years from Date of Sterilisation

Use: Single Use – do not reuse

CE: 0123

Manufacturer : Liebel-Flarsheim Company
2111 East Galbraith Road, Cincinnati, Ohio 45237

Country of Origin: USA

Linden-Luer Features:

- Easy to use connection
- Reduced thread pitch
- Reduced area for contrast to collect

tyco
Healthcare

MALLINCKRODT

Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1
53761 Hennef
Germany