

Liabel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

P/N: CTLxxx
Lengths: 75 cm / 120 cm / 150 cm
Description: Pressure Line for CT



Connectors: LL f/m
Package: 50 units/box
Pressure Limit Syringe : 300 PSI / 21 bar

Materials: PVC
ABS
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 5 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123

Manufacturer : Medtron AG,
Hauptstraße 255, 66128 Saarbrücken
Country of Origin: Germany



Tyco Healthcare Europe
Liabel-Flarsheim Products
Josef-Dietzgen-Str. 1-3
53761 Hennef
Germany