

Liebel-Flarsheim Single Use Connection Lines

Tubes Components Sheet

P/N: MRTSW
Lengths: 85 cm Filling Tube,
5 cm Connecting Tube
Description: Automatic filling tube
with drip chamber for MR



Connectors: Male Luer Lock
Package: 50 units/box
Pressure Limit Syringe : 300 PSI / 21 bar
Materials: PVC, ABS, Polycarbonat,
Silicon, PE,
Versapor
Latex Free

Sterilisation: Ethylene Oxide (EO)
Shelf Life: 5 years from Date of Sterilisation
Use: Single Use – do not reuse
CE: 0123
Manufacturer : Medtron AG,
Hauptstraße 255, 66128 Saarbrücken
Country of Origin: Germany



Tyco Healthcare Europe
Liebel-Flarsheim Products
Josef-Dietzgen-Str. 1-3
53761 Hennef
Germany

