


DRAINAGE SOLUTIONS

Everything you
Need for Drainage
Procedures



MERIT MAK-NV™ MINI ACCESS KIT NON-VASCULAR

Excellent transitions for a smooth introduction.

- Radiopaque marker tip clearly shows tip of introducer.
- Superior tip design decreases buckling and introducer peel back during insertion.
- Echo-enhanced needle increases visualization for ultrasound placement.
- 4F introducer ID allows 0.018" guide wire to be left in place for safety during larger wire placement.

ORDERING INFORMATION

Catalog Number	Introducer	Dilator	Guide Wire 0.018" (0.46 mm) 60 cm	Guide Wire Heavy Duty/Double Ended 0.038" (0.97 mm) 3mm J	Needle
MAK-NV-001	6F x 20 cm	4F	SS/PT	150 cm	21G x 15 cm Trocar
MAK-NV-002	6F x 20 cm	4F	NT/PT	150 cm	21G x 15 cm Trocar
MAK-NV-003	6F x 20 cm	4F	SS/PT	-	21G x 15 cm Trocar
MAK-NV-004	6F x 20 cm	4F	NT/PT	-	21G x 15 cm Trocar
MAK-NV-005*	-	-	-	-	21G x 15 cm Trocar
MAK-NV-006*	6F x 20 cm	4F	-	-	21G x 15 cm Trocar
MAK-NV-007*	6F x 20 cm	4F	-	-	-
MAK-NV-008	6F x 20 cm	4F	SS/PT	80 cm	21G x 15 cm Trocar
MAK-NV-009	6F x 20 cm	4F	NT/PT	80 cm	21G x 15 cm Trocar
MAK-NV-010	6F x 20 cm	4F	SS/PT	150 cm	21G x 15 cm Chiba
MAK-NV-011	6F x 20 cm	4F	NT/PT	150 cm	21G x 15 cm Chiba
MAK-NV-012	6F x 20 cm	4F	SS/PT	-	21G x 15 cm Chiba
MAK-NV-013	6F x 20 cm	4F	NT/PT	-	21G x 15 cm Chiba
MAK-NV-014*	-	-	-	-	21G x 15 cm Chiba
MAK-NV-015*	6F x 20 cm	4F	-	-	21G x 15 cm Chiba
MAK-NV-016	6F x 20 cm	4F	SS/PT	80 cm	21G x 15 cm Chiba
MAK-NV-017	6F x 20 cm	4F	NT/PT	80 cm	21G x 15 cm Chiba
MAK-NV-018	6F x 20 cm	4F	SS/PT	150 cm	21G x 20 cm Trocar
MAK-NV-019	6F x 20 cm	4F	NT/PT	150 cm	21G x 20 cm Trocar
MAK-NV-020*	-	-	-	-	21G x 20 cm Trocar
MAK-NV-021*	-	-	-	-	21G x 20 cm Chiba
MAK005-60	-	-	SS/PT	-	-
MAK001N-60	-	-	NT/PT	-	-

SS = Stainless Steel PT = Platinum Tip NT = Nitinol

Packaged 10 per box.

*Packaged 5 per box.



RESOLVE® BILIARY LOCKING DRAINAGE CATHETER

Our ReSolve® Biliary Locking Drainage Catheter is designed for fast, unimpeded drainage and offers up to 64% more drainage area* than other drainage catheters, with multiple large holes and different size configurations.

- Atraumatic internal pigtail reduces tissue irritation
- Depth markers and radiopaque marker band for fast and accurate positioning

ORDERING INFORMATION

Catalog Number	French Size	Working Length	Drainage Location	Pigtail to Distal Shaft Drainage Hole	Drainage Holes	Securement Device
RBC-8-038	8.5	40 cm	14 cm	5 cm	17	–
RBC-10-038	10	40 cm	14 cm	5 cm	17	–
RBC-12-038	12	40 cm	14 cm	5 cm	18	–
RBC-14-038	14	40 cm	14 cm	5 cm	18	–
RBDC-8-038	8.5	40 cm	11 cm	2 cm	17	–
RBDC-10-038	10	40 cm	11 cm	2 cm	17	–
RBDC-12-038	12	40 cm	11 cm	2 cm	18	–
RBDC-14-038	14	40 cm	11 cm	2 cm	18	–
RBC-8-REV	8.5	40 cm	14 cm	5 cm	17	Revolution™
RBC-10-REV	10	40 cm	14 cm	5 cm	17	Revolution™
RBC-12-REV	12	40 cm	14 cm	5 cm	18	Revolution™
RBC-14-REV	14	40 cm	14 cm	5 cm	18	Revolution™
RBDC-8-REV	8.5	40 cm	11 cm	2 cm	17	Revolution™
RBDC-10-REV	10	40 cm	11 cm	2 cm	17	Revolution™
RBDC-12-REV	12	40 cm	11 cm	2 cm	18	Revolution™
RBDC-14-REV	14	40 cm	11 cm	2 cm	18	Revolution™
RBC-8-SFX	8.5	40 cm	14 cm	5 cm	17	StayFix®
RBC-10-SFX	10	40 cm	14 cm	5 cm	17	StayFix®
RBC-12-SFX	12	40 cm	14 cm	5 cm	18	StayFix®
RBC-14-SFX	14	40 cm	14 cm	5 cm	18	StayFix®
RBDC-8-SFX	8.5	40 cm	11 cm	2 cm	17	StayFix®
RBDC-10-SFX	10	40 cm	11 cm	2 cm	17	StayFix®
RBDC-12-SFX	12	40 cm	11 cm	2 cm	18	StayFix®
RBDC-14-SFX	14	40 cm	11 cm	2 cm	18	StayFix®

*Data on File.



RESOLVE® LOCKING DRAINAGE CATHETER

Designed for abscess, nephrostomy and miscellaneous drainage procedures, the ReSolve® Locking Drainage Catheter comes packaged with a metal stiffening cannula, a flexible stiffening cannula*, and a trocar styllet.

- Catheter tip positioned within pigtail curve to reduce tissue irritation.
- Depth markers located on the catheter shaft aid accurate placement and act as a reference for correct positioning.
- Large drainage holes to maximize drainage.
- Tapered tip and hydrophilic coating for a smooth introduction to reduce patient discomfort.
- Flexible polyurethane material to prevent kinking and ensure excellent catheter rebound.
- Repositioning Tool (included) for unlocking and repositioning the catheter at any time.

ORDERING INFORMATION

Catalog Number	French Size	Working Length	Drainage Holes	Marker Band
RLC-6-038	6.5	25 cm	8	-
RLC-7-038	7.5	25 cm	8	-
RLC-8-038	8.5	25 cm	8	-
RLC-10-038	10	25 cm	8	-
RLC-12-038	12	25 cm	9	-
RLC-14-038	14	25 cm	9	-
RLC15-6-038MB	6.5	15 cm	8	Yes
RLC15-7-038MB	7.5	15 cm	8	Yes
RLC15-8-038MB	8.5	15 cm	8	Yes
RLC15-10-038MB	10	15 cm	8	Yes
RLC15-12-038MB	12	15 cm	9	Yes
RLC15-14-038MB	14	15 cm	9	Yes
RLC-6-038MB	6.5	25 cm	8	Yes
RLC-7-038MB	7.5	25 cm	8	Yes
RLC-8-038MB	8.5	25 cm	8	Yes
RLC-10-038MB	10	25 cm	8	Yes
RLC-12-038MB	12	25 cm	9	Yes
RLC-14-038MB	14	25 cm	9	Yes
RLC40-6-038MB	6.5	40 cm	8	Yes
RLC40-7-038MB	7.5	40 cm	8	Yes
RLC40-8-038MB	8.5	40 cm	8	Yes
RLC40-10-038MB	10	40 cm	8	Yes
RLC40-12-038MB	12	40 cm	9	Yes
RLC40-14-038MB	14	40 cm	9	Yes

Packaged 5 per box

Max guide wire: 0.038" (0.97mm)

RC-RTOOL Repositioning Tool Packaged 15 per box

* Flexible stiffening cannula not included with 6.5F and 7.5F catheters.



DRAINAGE CATHETER TRAYS

Conveniently containing the ReSolve® Locking Drainage Catheter, MAK-NV™ Introducer System, Merit Drainage Depot™ Drainage Bag, and dilators, these trays offer all you need from stick to fixation.



ORDERING INFORMATION


Catalog Number	French Size	Working Length	20 cm Dilator French Size	Marker Band	Securement Device
RLC-NV-6	6.5	25 cm	-	No	-
RLC-NV-8	8.5	25 cm	8	No	-
RLC-NV-10	10	25 cm	8, 10	No	-
RLC-NV-12	12	25 cm	8, 10, 12	No	-
RLC-NV-14	14	25 cm	8, 10, 12, 14	No	-
RLC-NV-6-SFX	6.5	25 cm	-	No	StayFIX (small)
RLC-NV-8-SFX	8.5	25 cm	8	No	StayFIX (small)
RLC-NV-10-SFX	10	25 cm	8, 10	No	StayFIX (small)
RLC-NV-12-SFX	12	25 cm	8, 10, 12	No	StayFIX (small)
RLC-NV-14-SFX	14	25 cm	8, 10, 12, 14	No	StayFIX (large)
RLC-NV-6-REV	6.5	25 cm	-	No	Revolution
RLC-NV-8-REV	8.5	25 cm	8	No	Revolution
RLC-NV-10-REV	10	25 cm	8, 10	No	Revolution
RLC-NV-12-REV	12	25 cm	8, 10, 12	No	Revolution
RLC-NV-14-REV	14	25 cm	8, 10, 12, 14	No	Revolution
RLCMB-NV-6	6.5	25 cm	-	Yes	-
RLCMB-NV-8	8.5	25 cm	8	Yes	-
RLCMB-NV-10	10	25 cm	8, 10	Yes	-
RLCMB-NV-12	12	25 cm	8, 10, 12	Yes	-
RLCMB-NV-14	14	25 cm	8, 10, 12, 14	Yes	-
RLCMB-NV-6-SFX	6.5	25 cm	-	Yes	StayFIX (small)
RLCMB-NV-8-SFX	8.5	25 cm	8	Yes	StayFIX (small)
RLCMB-NV-10-SFX	10	25 cm	8, 10	Yes	StayFIX (small)
RLCMB-NV-12-SFX	12	25 cm	8, 10, 12	Yes	StayFIX (small)
RLCMB-NV-14-SFX	14	25 cm	8, 10, 12, 14	Yes	StayFIX (large)
RLCMB-NV-6-REV	6.5	25 cm	-	Yes	Revolution
RLCMB-NV-8-REV	8.5	25 cm	8	Yes	Revolution
RLCMB-NV-10-REV	10	25 cm	8, 10	Yes	Revolution
RLCMB-NV-12-REV	12	25 cm	8, 10, 12	Yes	Revolution
RLCMB-NV-14-REV	14	25 cm	8, 10, 12, 14	Yes	Revolution

Packaged as individual units

MAK-NV includes: 21G x 15 cm Trocar needle, 6F x 20 cm introducer, 0.018" x 60 cm SS/PT guide wire and 0.038" x 150 cm Guide Wire.

ReSolve+™ LOCKING DRAINAGE CATHETER WITH SECUREMENT DEVICE

Having your drainage catheter and securement device conveniently packaged together saves both time and inventory space.



ORDERING INFORMATION

Catalog Number	French Size	Working Length	Marker Band	Securement Device
RLC-6-SFX	6.5	25 cm	-	StayFIX® (small)
RLC-8-SFX	8.5	25 cm	-	StayFIX® (small)
RLC-10-SFX	10	25 cm	-	StayFIX® (small)
RLC-12-SFX	12	25 cm	-	StayFIX® (small)
RLC-14-SFX	14	25 cm	-	StayFIX® (large)
RLCMB-6-SFX	6.5	25 cm	Yes	StayFIX® (small)
RLCMB-8-SFX	8.5	25 cm	Yes	StayFIX® (small)
RLCMB-10-SFX	10	25 cm	Yes	StayFIX® (small)
RLCMB-12-SFX	12	25 cm	Yes	StayFIX® (small)
RLCMB-14-SFX	14	25 cm	Yes	StayFIX® (large)
RLC-6-REV	6.5	25 cm	-	Revolution™
RLC-8-REV	8.5	25 cm	-	Revolution™
RLC-10-REV	10	25 cm	-	Revolution™
RLC-12-REV	12	25 cm	-	Revolution™
RLC-14-REV	14	25 cm	-	Revolution™
RLCMB-6-REV	6.5	25 cm	Yes	Revolution™
RLCMB-8-REV	8.5	25 cm	Yes	Revolution™
RLCMB-10-REV	10	25 cm	Yes	Revolution™
RLCMB-12-REV	12	25 cm	Yes	Revolution™
RLCMB-14-REV	14	25 cm	Yes	Revolution™

Packaged 5 per box

RESOLVE® NON-LOCKING DRAINAGE CATHETER

Designed for drainage of abscess, nephrostomy, and miscellaneous drainage procedures, the ReSolve® comes in a J or straight tip and is packaged with a metal stiffening cannula and a trocar stylet.

- Large drainage holes to maximize drainage.
- Tapered tip and hydrophilic coating for a smooth introduction to reduce patient discomfort.
- Color-coded hubs that match French sizing charts for easy identification.
- Flexible polyurethane material to prevent kinking and to ensure excellent catheter rebound.

ORDERING INFORMATION

Catalog Number	French Size	Working Length	Tip Configuration	Drainage Holes
RNL-6-038S	6.5	20 cm	Straight	4
RNL-7-038S	7.5	25 cm	Straight	4
RNL-8-038S	8.5	25 cm	Straight	4
RNL-10-038S	10	25 cm	Straight	4
RNL-12-038S	12	25 cm	Straight	5
RNL-6-038J	6.5	20 cm	J Curve	4
RNL-7-038J	7.5	25 cm	J Curve	4
RNL-8-038J	8.5	25 cm	J Curve	4
RNL-10-038J	10	25 cm	J Curve	4
RNL-12-038J	12	25 cm	J Curve	5

Packaged 5 per box
Max guide wire: 0.038" (0.97 mm)



VALVED ONE-STEP™ CENTESIS CATHETER

The self-closing valve prevents air infiltration and fluid leakage. When activated with a luer, the valve fully opens to allow maximum flow.

- Self-closing valve prevents air infiltration and fluid leakage.
- Four large holes to maximize drainage.
- Tapered tip for a smooth introduction to reduce patient discomfort.
- Echo-enhanced needle for increased needle visibility.
- Bevel indicator gives a tactile and visual reference for bevel orientation.
- Clear hub increases flashback visibility with ergonomic design for easy gripping and handling.



ORDERING INFORMATION

Catalog Number	French Size	Working Length	Drainage Holes	Needle Luer
PIGTAIL				
OSC-4F-12PLT	4	12	4	Fixed Luer
OSC-4F-17PLT	4	17	4	Fixed Luer
OSC-5F-12PLT	5	12	4	Fixed Luer
OSC-5F-17PLT	5	17	4	Fixed Luer
OSC-6F-12PLT	6	12	4	Fixed Luer
OSC-6F-17PLT	6	17	4	Fixed Luer
OSC-8F-12PLT	8	12	4	Fixed Luer
OSC-8F-17PLT	8	17	4	Fixed Luer
OSC-4F-12PT	4	12	4	Slip Fit
OSC-4F-17PT	4	17	4	Slip Fit
OSC-5F-12PT	5	12	4	Slip Fit
OSC-5F-17PT	5	17	4	Slip Fit
OSC-6F-12PT	6	12	4	Slip Fit
OSC-6F-17PT	6	17	4	Slip Fit
OSC-8F-12PT	8	12	4	Slip Fit
OSC-8F-17PT	8	17	4	Slip Fit
STRAIGHT				
OSC-4F-7LT	4	7	4	Fixed Luer
OSC-4F-10LT	4	10	4	Fixed Luer
OSC-4F-15LT	4	15	4	Fixed Luer
OSC-5F-7LT	5	7	4	Fixed Luer
OSC-5F-10LT	5	10	4	Fixed Luer
OSC-5F-15LT	5	15	4	Fixed Luer
OSC-6F-7LT	6	7	4	Fixed Luer
OSC-6F-10LT	6	10	4	Fixed Luer
OSC-6F-15LT	6	15	4	Fixed Luer
OSC-8F-7LT	8	7	4	Fixed Luer
OSC-8F-10LT	8	10	4	Fixed Luer
OSC-8F-15LT	8	15	4	Fixed Luer
OSC-4F-7T	4	7	4	Slip Fit
OSC-4F-10T	4	10	4	Slip Fit
OSC-4F-15T	4	15	4	Slip Fit
OSC-5F-7T	5	7	4	Slip Fit
OSC-5F-10T	5	10	4	Slip Fit
OSC-5F-15T	5	15	4	Slip Fit
OSC-6F-7T	6	7	4	Slip Fit
OSC-6F-10T	6	10	4	Slip Fit
OSC-6F-15T	6	15	4	Slip Fit
OSC-8F-7T	8	7	4	Slip Fit
OSC-8F-10T	8	10	4	Slip Fit
OSC-8F-15T	8	15	4	Slip Fit

Packaged 5 per box

ONE-STEP™ CENTESIS CATHETER

For ease of insertion, the One-Step™ Centesis Catheter has a smooth transition for decreased resistance and an echo-enhanced needle for more accurate placement.

- Available with luer locking or slip fit needle.
- Four large holes to maximize drainage.
- Tapered tip for a smooth introduction to reduce patient discomfort.
- Unique hub design for easy gripping and handling during procedure.
- Echo-enhanced needle for increased needle visibility.



ORDERING INFORMATION

Catalog Number	French Size	Working Length	Drainage Holes	Needle Luer
PIGTAIL				
OSC-4F-12PL	4	12	4	Fixed Luer
OSC-4F-17PL	4	17	4	Fixed Luer
OSC-5F-12PL	5	12	4	Fixed Luer
OSC-5F-17PL	5	17	4	Fixed Luer
OSC-6F-12PL	6	12	4	Fixed Luer
OSC-6F-17PL	6	17	4	Fixed Luer
OSC-8F-12PL	8	12	4	Fixed Luer
OSC-8F-17PL	8	17	4	Fixed Luer
OSC-4F-12PS	4	12	4	Slip Fit
OSC-4F-17PS	4	17	4	Slip Fit
OSC-5F-12PS	5	12	4	Slip Fit
OSC-5F-17PS	5	17	4	Slip Fit
OSC-6F-12PS	6	12	4	Slip Fit
OSC-6F-17PS	6	17	4	Slip Fit
OSC-8F-12PS	8	12	4	Slip Fit
OSC-8F-17PS	8	17	4	Slip Fit
STRAIGHT				
OSC-4F-7L	4	7	4	Fixed Luer
OSC-4F-10L	4	10	4	Fixed Luer
OSC-4F-15L	4	15	4	Fixed Luer
OSC-5F-7L	5	7	4	Fixed Luer
OSC-5F-10L	5	10	4	Fixed Luer
OSC-5F-15L	5	15	4	Fixed Luer
OSC-6F-7L	6	7	4	Fixed Luer
OSC-6F-10L	6	10	4	Fixed Luer
OSC-6F-15L	6	15	4	Fixed Luer
OSC-8F-7L	8	7	4	Fixed Luer
OSC-8F-10L	8	10	4	Fixed Luer
OSC-8F-15L	8	15	4	Fixed Luer
OSC-4F-7S	4	7	4	Slip Fit
OSC-4F-10S	4	10	4	Slip Fit
OSC-4F-15S	4	15	4	Slip Fit
OSC-5F-7S	5	7	4	Slip Fit
OSC-5F-10S	5	10	4	Slip Fit
OSC-5F-15S	5	15	4	Slip Fit
OSC-6F-7S	6	7	4	Slip Fit
OSC-6F-10S	6	10	4	Slip Fit
OSC-6F-15S	6	15	4	Slip Fit
OSC-8F-7S	8	7	4	Slip Fit
OSC-8F-10S	8	10	4	Slip Fit
OSC-8F-15S	8	15	4	Slip Fit

Packaged 5 per box

THORACENTESIS AND PARACENTESIS SET*

Includes necessary products to complete a thoracentesis and paracentesis procedure, including three safety needles to ensure a safe working environment.

* Not available in countries where CE mark is required.



TRAY MAY INCLUDE

- 1 Each Valved One-Step Centesis Catheter
- 1 Each ChloraPrep® 3 mL tinted*
- 1 Each Drape tri-layer
- 1 Each Fenestrated drape
- 1 Each Lidocaine 20 mL
- 1 Each Lidocaine label
- 1 Each 18G x 1" shielded needle
- 1 Each 22G x 1.5" shielded needle
- 1 Each 25G x 1" shielded needle
- 1 Each Medallion® Syringe 10 mL
- 1 Each 10 mL Luer lock syringe
- 1 Each 60 mL Luer lock syringe
- 1 Each #11 Futura® Safety Scalpel
- 1 Each ShortStop® Temporary Sharps Holder
- 10 Each Gauze 4" x 4" (10 cm x 10 cm)
- 1 Each Non-vented female cap
- 3 Each Specimen vials (10 mL) with caps
- 1 Each 3-Way large bore stopcock
- 1 Each 48" (122 cm) tubing with spike for evacuated container
- 1 Each 96" (244 cm) wall suction tubing
- 1 Each Stepped adapter
- 1 Each Bandage dressing
- 1 Each CSR wrap

ORDERING INFORMATION

Catalog Number	Catheter French Size	Catheter Length	Tip	Luer
TAPS-100	-	-	-	-
TAPS-5F-12PLT	5	12	PIGTAIL	FIXED
TAPS-5F-17PLT	5	17	PIGTAIL	FIXED
TAPS-6F-12PLT	6	12	PIGTAIL	FIXED
TAPS-6F-17PLT	6	17	PIGTAIL	FIXED
TAPS-8F-12PLT	8	12	PIGTAIL	FIXED
TAPS-8F-17PLT	8	17	PIGTAIL	FIXED
TAPS-5F-12PT	5	12	PIGTAIL	SLIP
TAPS-5F-17PT	5	17	PIGTAIL	SLIP
TAPS-6F-12PT	6	12	PIGTAIL	SLIP
TAPS-6F-17PT	6	17	PIGTAIL	SLIP
TAPS-8F-12PT	8	12	PIGTAIL	SLIP
TAPS-8F-17PT	8	17	PIGTAIL	SLIP
TAPS-5F-7LT	5	7	STRAIGHT	FIXED
TAPS-5F-10LT	5	10	STRAIGHT	FIXED
TAPS-6F-7LT	6	7	STRAIGHT	FIXED
TAPS-6F-10LT	6	10	STRAIGHT	FIXED
TAPS-8F-7LT	8	7	STRAIGHT	FIXED
TAPS-8F-10LT	8	10	STRAIGHT	FIXED
TAPS-5F-7T	5	7	STRAIGHT	SLIP
TAPS-5F-10T	5	10	STRAIGHT	SLIP
TAPS-6F-7T	6	7	STRAIGHT	SLIP
TAPS-6F-10T	6	10	STRAIGHT	SLIP
TAPS-8F-7T	8	7	STRAIGHT	SLIP
TAPS-8F-10T	8	10	STRAIGHT	SLIP

Packaged 5 per box

ChloroPrep® is a registered trademark of CareFusion 2200, Inc.

SAFETY PARACENTESIS PROCEDURE TRAY™ *

Includes necessary products to complete a paracentesis procedure, including three safety needles to ensure a safe working environment.

* Not available in countries where CE mark is required.

ORDERING INFORMATION

Catalog Number	Catheter French Size	Catheter Length	Tip	Luer
SPPT-100	-	-	-	-
SPPT-5F-12PL	5	12	PIGTAIL	FIXED
SPPT-5F-17PL	5	17	PIGTAIL	FIXED
SPPT-6F-12PL	6	12	PIGTAIL	FIXED
SPPT-6F-17PL	6	17	PIGTAIL	FIXED
SPPT-8F-12PL	8	12	PIGTAIL	FIXED
SPPT-8F-17PL	8	17	PIGTAIL	FIXED
SPPT-5F-12PS	5	12	PIGTAIL	SLIP
SPPT-5F-17PS	5	17	PIGTAIL	SLIP
SPPT-6F-12PS	6	12	PIGTAIL	SLIP
SPPT-6F-17PS	6	17	PIGTAIL	SLIP
SPPT-8F-12PS	8	12	PIGTAIL	SLIP
SPPT-8F-17PS	8	17	PIGTAIL	SLIP
SPPT-5F-7L	5	7	STRAIGHT	FIXED
SPPT-5F-10L	5	10	STRAIGHT	FIXED
SPPT-6F-7L	6	7	STRAIGHT	FIXED
SPPT-6F-10L	6	10	STRAIGHT	FIXED
SPPT-8F-7L	8	7	STRAIGHT	FIXED
SPPT-8F-10L	8	10	STRAIGHT	FIXED
SPPT-5F-7S	5	7	STRAIGHT	SLIP
SPPT-5F-10S	5	10	STRAIGHT	SLIP
SPPT-6F-7S	6	7	STRAIGHT	SLIP
SPPT-6F-10S	6	10	STRAIGHT	SLIP
SPPT-8F-7S	8	7	STRAIGHT	SLIP
SPPT-8F-10S	8	10	STRAIGHT	SLIP

Packaged 5 per box

Chloraprep® is a registered trademark of CareFusion 2200, Inc.



TRAY MAY INCLUDE

- 1 Each One-Step™ centesis catheter
- 1 Each Chloraprep® 3 mL tinted*
- 1 Each Drape tri-layer
- 1 Each Fenestrated drape
- 1 Each Lidocaine 20 mL
- 1 Each Lidocaine label
- 1 Each 18G x 1" shielded needle
- 1 Each 22G x 1.5" shielded needle
- 1 Each 25G x 1" shielded needle
- 1 Each Medallion® Syringe 10 mL
- 1 Each 10 mL Luer lock syringe
- 1 Each 60 mL Luer lock syringe
- 1 Each #11 Futura® Safety Scalpel
- 1 Each ShortStop® Temporary Sharps Holder
- 10 Each Gauze 4" x 4" (10 cm x 10 cm)
- 1 Each Non-vented female cap
- 3 Each Specimen vials (10 mL) with caps
- 1 Each 3-Way large bore stopcock
- 1 Each 48" (122 cm) tubing with spike for evacuated container
- 1 Each 96" (244 cm) wall suction tubing
- 1 Each Stepped adapter
- 1 Each Bandage dressing
- 1 Each CSR wrap

CATHETER SECUREMENT

The Revolution™ Catheter Securement Device saves time, provides convenience and enhances compliance with safety regulations for securing catheters with a needleless procedure. The StayFIX® Fixation Device design secures catheters in place without the use of sutures, reducing the exposure to sharps and sharps-related accidents.



Revolution™

Catheter Securement Device

- Secures non vascular catheters from 6F-24F.
- Hydrocolloid is absorbent around puncture site, and conforms well to body anatomy.
- Provides safe and secure fixation for up to 7 days.
- Included shower cap covers insertion site for added protection while showering.

ORDERING INFORMATION

CFG624 - Packaged 25 per box

StayFIX®

Fixation Device

- Secures non-vascular catheters from 5F - 22F.
- Provides safe and secure fixation for up to 7 days.

ORDERING INFORMATION

680ME 5F-14F - Packaged 25 per box

685ME 12F-22F - Packaged 20 per box

DRAINAGE BAGS AND TUBING

Merit offers a selection of drainage bags and tubing to meet your specific drainage needs.

- Adjustable length tubing to fit a variety of patients
- Twist valve for easy disposal of fluids
- Internal anti-reflux valve prevents backflow
- Cloth backing for patient comfort
- Easy-to-read markings to identify fluid volume

ORDERING INFORMATION

MDD600	600 mL bag
MDD600PML	600 mL bag with polypropylene Luer
FZ624	600 mL cloth-backed bag

Packaged 20 per box

K08-01907	Jackson-Pratt 100 mL bulb
K10-03450	MDD600 drainage bag with Jackson-Pratt® 100 mL bulb designed for vacuum drainage

Packaged 10 per box

ABSCESS DRAINAGE KIT

K10-00230B	OSHA-compliant closed waste system with 1400 mL MDD bag, tubing and CRV Valve Designed for quick and safe drainage elimination.
-------------------	--

Packaged 15 per box

DRAINAGE CONNECTOR TUBING

K10-04129	30 cm connecting tubing
K10-04130	42 cm connecting tubing
K10-04131	30 cm connecting tubing with one-way stopcock
K10-04132	42 cm connecting tubing with one-way stopcock

Packaged 25 per box

DRAINAGE TUBING

FAS5248	48" large bore tubing with pinch clamp and non-vented spike. Packaged 25 units per box.
K08-02027	Tubing kit with barbed male luer. Packaged 25 units per box.
K10-04750	48" Large bore tubing with roller clamp and 16G x 1.5" (38.1 mm) needle. Packaged 15 per box



Jackson-Pratt® is a registered trademark of Cardinal Health.

PI

DRAINAGE ACCESSORIES

With a breadth of products to complete all of your drainage procedure needs, Merit is your One-Stop-Shop for drainage procedure products.

VALVE ADAPTER

ORDERING INFORMATION

VA-40 Designed to be back loaded over a guide wire and attached to the MAK-NV to facilitate aspiration or infusion. Accepts up to a 4F device

Packaged 10 per box

DILATORS

ORDERING INFORMATION

K12-08322 8F X 20 cm

K12-08323 10F X 20 cm

K12-08324 12F X 20 cm

K12-08325 14F X 20 cm

Packaged 10 per box

LARGE BORE STOPCOCKS

Clear body with large lumen designed to increase flow rate.

ORDERING INFORMATION

L3SNC 3-way large bore stopcock

L4SNC 4-way large bore stopcock

Packaged 50 per box

VacLok®

Vacuum Pressure Syringe

- Clear body design with multiple locking positions.
- Maintains negative pressure while reducing physician strain.

ORDERING INFORMATION

VAC120 20 mL VacLok Syringe

VAC130 30 mL VacLok Syringe

Packaged 20 per box

VAC110 10 mL VacLok Syringe

VAC160 60 mL VacLok Syringe

Packaged 25 per box

ACCESS NEEDLE

ORDERING INFORMATION

K12-08326 18G X 15 cm Chiba

Packaged 10 per box

InQwire®

Guide Wires

InQwire® guide wires are precoated, creating a smooth, consistent surface for drainage catheters to track over. 80 cm length reduces excess wire and risk for contamination.

ORDERING INFORMATION

IQ35F80B Bentson 0.035"(0.89 mm)

6680-61 Super Stiff 0.035"(0.89 mm)

IQ35F80J1O5RS Rosen (1.5 mm J) 0.035"(0.89 mm)

IQ35F80J3 Rosen (3 mm J-tip) 0.035"(0.89 mm)



VacLok®
Vacuum
Pressure Syringe

InQwire®
Diagnostic Guide Wire

WASTE BASINS

Merit's innovative waste basins provide quick and safe fluid waste elimination by providing a seal-tight lid to prevent fluid leakage after a procedure. Our waste basins can be used in all surgical environments as they are designed with particulate-free materials.

BackStop®

Disposal System

Designed for quick and safe drainage disposal.

ORDERING INFORMATION

MBS100 500 mL waste basin

MBSP100 500 mL waste basin with Temporary Sharps Holder

Packaged 10 per box

MiniStop®

Disposal System

Designed for smaller drainage procedures.

ORDERING INFORMATION

MMS100 100 mL waste basin

MMSP100 100 mL waste basin with Temporary Sharps Holder

Packaged 20 per box

DugOut®

Closed Waste Basin

Designed for quick and safe disposal of the largest volume drainage procedures.

ORDERING INFORMATION

DUG100 1,000 mL waste basin

DUG1SS 1,000 mL waste basin with ShortStop
Temporary Sharps Holder

Packaged 6 per box



BackStop+
Disposal System

RLS™

Resolve Label Set

Created to help clinicians to standardize and improve patient care.

ORDERING INFORMATION

RLS100 Orange, Includes Pen, 5 labels per sheet

Packaged 25 per box



RLS™
Resolve Label Set



Corporate Headquarters, South Jordan, Utah

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
+1 801-208-4300
+1 800 35 MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30
(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115