



BLUperc® PDT kit with BLUselect® tube



BLUgriggs® PDT kits with BLUselect® tube

BLUperc® and BLUgriggs® Percutaneous Dilation Tracheostomy Kits

Components

BLUperc® and BLUgriggs® Percutaneous Dilation Tracheostomy Kits Components

| BLUperc® PDT kits with BLUselect® include | | | |
|--|---|--|--|
| All Trays and Kits Include (101/562/---) | Additional components in 101/561/--- includes a cuffed tube 101/563/--- includes a Suctionaid® tube | Additional components in 101/595/--- includes a cuffed tube 101/563/--- includes a Suctionaid® tube | Additional components in 101/573/--- |
| <ul style="list-style-type: none"> Lubricating jelly (2) Gauze sponges (10) Split tracheostomy dressing Safety scalpel 5 mL syringe(s) 10 mL syringe Introducer cannula and needle assembly Introducer needle(s) J-tip guidewire pre-loaded in single handed feeder 14 FR short dilator Single stage dilator pre-loaded with guiding catheter Curved hemostat forceps Tracheostomy tube introducer(s) Needle safety device Tracheostomy tube holder with brush Disconnection wedge | <ul style="list-style-type: none"> BLUselect® tracheostomy tube Inner cannulas (2) If Suctionaid tube, vacuum control valve Obturator Patient labels Dedicated introducer | <ul style="list-style-type: none"> BLUselect® tracheostomy tube Inner cannulas (2) Needle driver Polypropylene sutures PEEP keep swivel adaptor If Suctionaid® tube, vacuum control valve Obturator Patient labels ChlorPrep applicators (2) 5 mL Lidocaine and Epinephrine 1.5% (2) Filter straw Fenestrated drape Hypodermic needles Dedicated introducer Polypropylene sutures | <ul style="list-style-type: none"> Needle driver Polypropylene sutures PEEP peep swivel adaptor ChlorPrep applicators (2) 5 mL Lidocaine and Epinephrine 1.5% (2) Filter straw Fenestrated drape Hypodermic needles Polypropylene sutures |

| BLUgriggs® PDT kits include | | | |
|--|--|--|---|
| Trays and Kits include 101/543/--- includes a cuffed tube 101/893/--- includes a Suctionaid® tube | | Additional components in 101/541/--- includes a cuffed tube 101/891/--- includes a Suctionaid® tube | Additional components in 101/540/--- includes a cuffed tube 101/892/--- includes a Suctionaid® tube |
| <ul style="list-style-type: none"> Lubricating jelly (2) Gauze sponges (10) Split tracheostomy dressing Safety scalpel 5 mL syringe(s) 10 mL syringe Introducer cannula and needle assembly Introducer needle(s) J-tip guidewire pre-loaded in single handed feeder | <ul style="list-style-type: none"> 14 FR short dilator BLUselect® tracheostomy tube Inner cannulas (2) Needle safety device Tracheostomy tube holder with brush Disconnection wedge If Suctionaid tube, vacuum control valve Obturator Patient labels | <ul style="list-style-type: none"> Guidewire dilating forceps | <ul style="list-style-type: none"> Guidewire dilating forceps Needle driver Polypropylene sutures PEEP keep swivel adaptor ChlorPrep applicators 5 mL Lidocaine and Epinephrine 1.5% Filter straw Fenestrated drape Hypodermic needles |

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN ALL GEOGRAPHIES

MMSPCA-0570

Smiths Medical ASD, Inc.
6000 Nathan Lane North
Minneapolis, MN 55442, USA
Tel: 1-614-210-7300
Toll-Free USA: 1-800-258-5361
www.smiths-medical.com

Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. BLUgriggs, BLUperc, Portex and the Smiths Medical design marks are trademarks of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries ©2019 Smiths Medical. All rights reserved. TR194567EN-102019

**Rx
ONLY**

smiths medical