

# Clinical Accessories

Trusted. Reliable. Simple.



# Contents

<b>PATIENT MONITORING</b>	<b>4</b>	<b>Flow Sensors – Respiratory</b>	<b>49</b>	<b>Lullaby™ Warmer</b>	<b>76</b>
<b>ECG</b>	<b>6</b>	Respiratory Flow Sensors and		Temperature Probes	76
Monitoring Recording Paper	6	Expiratory Valve Assemblies	49	Temperature Probe Patches	76
Trunk Cables	7	<b>Flow Sensors – Anesthesia</b>	<b>50</b>	<b>Shuttle™</b>	<b>76</b>
Leadwires	8-10	<b>Nebulizer Accessories</b>	<b>51</b>	Patient Restraint	76
Trunk Cable Accessories	11	<b>Anesthesia O2 Cells</b>	<b>51</b>		
<b>Pulse Oximetry</b>	<b>12</b>	<b>Respiratory Safety Guard</b>	<b>52</b>		
Sensors & Cables	12-18	<b>Pressure Catheters</b>	<b>52</b>		
<b>Temperature</b>	<b>19</b>	<b>DIAGNOSTIC CARDIOLOGY</b>	<b>53</b>		
Cables & Probes	19-20	<b>MAC™, CASE™ and CardioSoft™</b>	<b>55</b>		
<b>Invasive Blood Pressure</b>	<b>21</b>	<b>ECG</b>	<b>55</b>		
Cables	21-22	Cables	55		
<b>Cuffs</b>	<b>23</b>	Leadwires	55-56		
Hoses	23	Accessories	57		
SOFT-CUF	24	Recording Paper	58		
DURA-CUF	25	Electrodes	59		
CLASSIC-CUF	26	<b>SEER® Holter Monitors</b>	<b>60</b>		
RADIAL-CUF	27	<b>ECG – CM 4000</b>	<b>60</b>		
Neonatal	28	Leadwire Sets	60		
<b>Cardiac Output</b>	<b>29</b>	Accessories	61		
Cables	29-30	<b>ECG – SEER® 12</b>	<b>61</b>		
Other Accessories	30	Leadwire Sets	61		
<b>Adequacy of Anesthesia</b>	<b>31</b>	Accessories	61		
Cable	31-32	<b>ECG – SEER® 1000</b>	<b>62</b>		
Sensors	31-32	Leadwire Sets	62		
<b>Neurology</b>	<b>33</b>	Accessories	63		
EEG Leadsets	33	<b>ECG – SEER® Light</b>	<b>64</b>		
EEG Electrodes	34	Leadwire and Cable Kits	64		
AEP Accessories	34	Accessories	65		
<b>Patient Spirometry</b>	<b>35</b>	<b>MATERNAL INFANT CARE</b>	<b>66</b>		
Sensors	35	<b>Maternal Fetal Care</b>	<b>68</b>		
Tubes	36	<b>Mini Telemetry</b>	<b>68</b>		
Kits	36	Recording Thermal Paper	68		
<b>Gas Monitoring</b>	<b>37</b>	Remote Event Marker	68		
Airway Adapter Sets	37	<b>Novii™</b>	<b>69</b>		
Gas Sampling Lines	38	Novii™ Patch	69		
Airway Adapters	38	Cables	69		
Water Traps	38	<b>Corometrics™</b>	<b>70</b>		
<b>ANESTHESIA &amp; RESPIRATORY CARE</b>	<b>39</b>	Cables	70		
<b>Absorbents</b>	<b>41</b>	Transducer Belts	70		
Absorbers	41	Ultrasound Transducers	70		
Canisters	41	<b>Internal Monitoring</b>	<b>71</b>		
Dust Filter	41	Cables	71		
<b>Vaporizer Bottle Adapters</b>	<b>41</b>	Electrode	71		
<b>Filters &amp; Breathing Masks</b>	<b>42</b>	Pad	71		
<b>Reservoir Bags &amp; Gas Scavenging</b>	<b>43</b>	<b>Infant Care</b>	<b>72</b>		
Anesthesia Breathing Bags	43	<b>Giraffe™/Panda™</b>	<b>72-73</b>		
Scavenging Bags, Tubes, & Connectors	43	Air filters	72		
<b>Patient Breathing Circuits</b>	<b>44</b>	Clip	72		
Circuit Startup Kits	44	Cover	73		
Circuit Kits	45	<b>Phototherapy Accessories</b>	<b>74</b>		
Circuits	46-47	Pad Covers	74		
Circuit Components	48	<b>Resuscitation</b>	<b>75</b>		
		T-Piece	75		

## Welcome

GE Healthcare offers you a simple, reliable one-stop solution for high-quality GE Healthcare-compatible clinical accessories. We provide an expanded inventory of products that have been verified by our engineers. The compatible devices listed may not be available for sale and/or cleared/approved in all countries. Check with your representative for product availability.

# PATIENT MONITORING



## **PATIENT MONITORING**

### **ECG**

- Leadwires
- Recording Paper
- Trunk Cable
- Trunk Cables Accessories

### **Pulse Oximetry**

- Cables
- Clip
- Other Accessories
- Probes
- Sensors
- Tape & Wrap

### **Temperature**

- Cables
- Probes

### **Invasive Blood Pressure**

- Cables

### **Cuffs**

- Hoses
- SOFT-CUFF
- DURA-CUF
- CLASSIC-CUF
- RADIAL-CUF
- Neonatal

### **Cardiac Output**

- Cables
- Catheter Kits
- Other Accessories

### **Adequacy of Anesthesia**

- Cable
- Electrodes
- Sensors

### **Neurology**

- AEP Accessories
- EEG Electrodes
- EEG Leadsets

### **Patient Spirometry**

- Kits
- Sensors
- Tubes

### **Gas Monitoring**

- Airway Adapter Sets
- Airway Adapters
- Gas Sampling Lines
- Water Traps

# ECG



## MONITORING RECORDING PAPER



### MS Recording Thermal Paper

Length	CARESCAPE™ B450/B650 Monitors, CARESCAPE™ B850 Monitor with PRN 50-M+ printer	B20(i)/B40(i) Monitors	B1x5 Monitors	CARESCAPE™ V100 and VC150 Monitors	Part Number	Quantity
57.2 mm x 24.4 m				X	2106825-001	10
50 mm x 25.6 m	X	X	X		2106824-001	20
50 mm x 27.4 m	X	X	X		2106823-001	48

For ECG paper for Diagnostic cardiology and Maternal infant care equipment, please go to pages 58 and 68.

*Not available in all countries*





2106305-003



2106310-002

ECG Trunk Cables – IEC			
Description	Length	Part Number	Quantity
Neonatal DIN 3-lead	3.6 m	2106306-003	1
	1.2 m	2106306-004	1
3/5-lead	3.6 m	2106305-003	1
	1.2 m	2106305-004	1
	6 m	2106311-002	1
3/5-lead w/ESU filter	3.6 m	2106310-002	1
6-lead	3.6 m	2106307-003	1
	1.2 m	2106307-004	1
12-lead	3.6 m	2106308-003	1
	1.2 m	2106308-004	1



2106307-003



2106308-004

## Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module E-PSM
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ Monitor with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry

Specifications available on request



ECG Trunk Cables – IEC			
Description	Length	Part Number	Quantity
Neonatal trunk cable	3.6 m	2080557-003	1
	1.2 m	2080557-004	1

Not available in all countries



## ECG Leadwires Grabber – IEC Compliant

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm*	2106390-003	1
		130 cm	2106390-004	
5-lead	Molded	74 cm*	2106389-003	1
		130 cm	2106389-004	
		mix 74 cm, 130 cm	2106389-006	
	Non-molded	74 cm *	2106391-003	
		130 cm	2106391-004	
		mix 74 cm, 130 cm	2106391-006	
5-lead C2-6	Non-molded	130 cm	2106393-004	1
		74 cm	2106393-002	
6-lead	Non-molded	74 cm*	2106398-001	1
		130 cm	2106398-002	
		mix 74 cm, 130 cm	2106398-003	

### Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

\*74cm 3-, 5-, and 6-leadwire sets are compatible with ApexPro™

*Not available in all countries*







2106385-004



2106383-006



2106381-004



2106394-002



2106398-001

## ECG Leadwires Snap – IEC Compliant

Description	Connection	Length	Part Number	Quantity
3-lead	Molded	74 cm*	2106385-003	1
		130 cm	2106385-004	
5-lead	Molded	74 cm*	2106383-003	1
		130 cm	2106383-004	
		mix 74 c m, 130 cm	2106383-006	
	Non-molded	74 cm*	2106381-003	
130 cm		2106381-004		
5-lead C2-6	Non-molded	130 cm	2106394-004	1
		74 cm	2106394-002	
6-lead	Non-molded	74 cm*	2106396-001	1
		130 cm	2106396-002	

## Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor

Specifications available on request

\*74cm 3-, 5-, and 6-leadwire sets are compatible with ApexPro™

*Not available in all countries*



ECG Leadwires – IEC Compliant				
Description	Connection	Length	Part Number	Quantity
4-lead	Snap	74 cm/29 in	2066468-002	1
		130 cm/51 in	2066468-006	
	Grabber	74 cm/29 in	2066468-004	1
		130 cm/51 in	2066468-008	

### Product compatibility

- CARESCAPE™ ONE

Specifications available on request

*Not available in all countries*



2106403-003



2106403-001



2106399-001



2106399-002

## ECG Trunk Cable Accessories

Description		Part Number	Quantity
DIN Converter	3-Lead	2106403-003	1
	5-Lead	2106403-001	1
Leadwire separator		2106399-001	3
Bedsheet clip		2106399-002	1

## Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry

Specifications available on request

*Not available in all countries*



# Pulse Oximetry



## TRUSIGNAL™

GE Healthcare provides compatibility with three major brands of Pulse Oximetry (SpO2) measurement to meet all of your needs:

- GE Healthcare TruSignal™
- MASIMO
- Nellcor



TS-F2-GE



TS-SA4-GE



TS-E2-GE

TruSignal™ Integrated Reusable Sensors						
Description	Connector	Patient	Patient Weight	Length	Part Number	Quantity
Finger	GE	Adult, Pediatric	≥ 20 kg	2 m	TS-F2-GE	1
				4 m	TS-F4-GE	1
	Ohmeda			1 m	TS-F1-H	1
				Datex	4 m	TS-F4-N
FingerTip	GE	Adult	> 30 kg		4 m	TS-SA4-GE
PediTip		Pediatric	15-30 kg	3 m	TS-SP3-GE	1
Ear	GE	Adult, Pediatric	≥ 10 kg	2 m	TS-E2-GE	1
				Ohmeda	4 m	TS-E4-GE
	TS-E4-H				1	
	Datex				TS-E4-N	1

*Not available in all countries*





TS-F-D



TS-SA-D/TS-SP-D




TS-E-D



TS-W-D



TS-SE-3



TruSignal™ Reusable Sensors

Description	Patient	Patient Weight	Length	Part Number	Quantity
Finger	Adult, Pediatric	≥ 20 kg	1 m	TS-F-D	1
PediTip	Pediatric	15-30 kg		TS-SP-D	1
Fingertip	Adult	> 30 kg		TS-SA-D	
Ear	Adult, Pediatric	≥ 10 kg		TS-E-D	1
Wrap	Adult, Pediatric, Infant	≥ 3 kg		TS-W-D	1
Sensitive Skin	Adult, Pediatric, Infant, Neonate	No limit		TS-SE-3	3

Requires Compatible Interconnect Cable



TS-G3

TruSignal™ Interconnect Cables

Connector	Length	Part Number	Quantity
GE	3 m	TS-G3	1
Ohmeda		TS-H3	
Datex		TS-N3	



TruSignal™ Disposable Sensors					
Description	Patient	Patient Weight	Length	Part Number	Quantity
Adult/Pediatric	Pediatric, infant	≥ 20 kg	0.3 m	TS-AP-10	10
				TS-AP-25	25
AllFit	Adult, pediatric, infant, neonate	No limit	0.5 m	TS-AF-10	10
				TS-AF-25	25
Adhesive Wrap	Adult	> 20 kg	0.3 m	TS-AAW-10	10
				TS-AAW-25	25
	Pediatric, infant	3-20 kg		TS-PAW-10	10
				TS-PAW-25	25

Requires Compatible Interconnect Cable

TruSignal™ Sensor Accessories				
Description	Use With	Size	Part Number	Quantity
Replacement wide adhesive tape	TruSignal™ Sensitive Skin Sensor, TruSignal™ Wrap Sensor	-	OXY-RTW	100
Replacement adhesive tape (blue)	TruSignal™ AllFit Sensor	-	OXY-RT	100
Replacement adhesive tape (teddy bears)			OXY-RTB	
Replacement foam wrap	TruSignal™ Sensitive Skin Sensor	Small	OXY-RWS	24
		Medium	OXY-RWM	
		Large	OXY-RWL	
Infant foam sandal	TruSignal™ Sensitive Skin Sensor	-	OXY-SND	3
Replacement Headband	TruSignal™ Ear Sensor	-	OXY-HB	3



### Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM) or Patient Side Module (E-PSM)
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 2000/2500/3000/4000/5000 patient monitors
- Solar™ monitors with CARESCAPE™ Patient Data Module (PDM)
- ECG trunk cables are also compatible with GE Healthcare ECG leadwire sets
- ApexPro™ CH telemetry
- CARESCAPE™ ONE Monitor
- Giraffe Warmer
- Panda Warmer

Specifications available on request

*Not available in all countries*



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Masimo™ SpO2 technology.



LNCS Digit Sensor



DCI-P Digit sensor



LNCS Tip-Clip Ear Sensor

**Masimo™ Reusable LNCS and RD SET Sensors**

Sensors	Masimo™ Reference	Type	Patient	Patient Weight	Part Number	Quantity
LNCS Digit Sensor	1863	DC-I	Adult	> 30 kg	2027258-001	1
	1864	DCI-P	Pediatric	10-50 kg	2027259-001	
LNCS Tip-Clip Ear Sensor	1895	TC-I	Adult	> 30 kg	2027261-001	1
RD SET SpO2 Sensor	4050	DCI	Adult	> 30 kg	MASIMO-RD-4050	1
	4051	DCI-P	Pediatric	10-50 kg	MASIMO-RD-4051	

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.



2027263-002

**Masimo™ LNCS & RD SET Interconnect Cables**

Description	Size	Masimo™ Reference	Part Number	Quantity
LNC-10 Interconnect Cable	For GE DINAMAP	1814	2027263-001	1
	For GE	2016	2027263-002	1
RED Interconnect Cable	LNC-04	2055	2059005-001	1
	LNC-10	2056	2059006-001	1
	LNC-14	2057	2059007-001	1
RD SET Interconnect Cable	GE-05 (For GE)	4084	MASIMO-RD-4084	1
	GE-12 (For GE)	4085	MASIMO-RD-4085	1
RD SET Interconnect Cable	MD14-05, 14-pin	4080	MASIMO-RD-4080	1
	MD14-12, 14-pin	4081	MASIMO-RD-4081	1
	MD20-05, 20-pin	4103	MASIMO-RD-4103	1
	MD20-12, 20-pin	4104	MASIMO-RD-4104	1

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

*Not available in all countries*



2051070-001

**Masimo™ LNCS & RD SET Disposable Sensors**

Description	Patient Weight	Masimo™ Reference	Part Number	Quantity	
LNCS SpO2	Adult	>30 kg	1859	2027253-001	20
	Pediatric	10-50 kg	1860	2027254-001	20
	Infant	3-20 kg	2328	2051068-001	20
	Neonatal PT	< 1 kg	2330	2051070-001	20
	Adult / Neonatal	> 40 kg / < 3 kg	2329	2051069-001	20
RD SET® SpO2	Adult	>30 kg	4000	MASIMO-RD-4000	20
	Pediatric	10-50 kg	4001	MASIMO-RD-4001	20
	Infant	3-20 kg	4002	MASIMO-RD-4002	20
	Adult / Neonatal	> 40 kg / < 3 kg	4003	MASIMO-RD-4003	20
	Premature Baby	< 1 kg	4004	MASIMO-RD-4004	20

Requires Compatible Interconnect Cable

These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

**Masimo™ LNCS & RD SET Sensor Accessories**

Description	Masimo™ Reference	Part Number	Quantity	
LNCS Replacement Tape	Infant	2307	2051071-001	100
	Neonatal	2308	2051073-001	100
	Neonatal PT	2309	2051074-001	100
RD Set Replacement Tape	Infant	4046	MASIMO-RD-4046	102
	Neonatal	4047	MASIMO-RD-4047	102
RD Set Replacement Wrap	Neonatal PT	4048	MASIMO-RD-4048	10


These accessories may only be available from GE Healthcare when purchased with compatible GE Healthcare medical devices. Please consult with your local GE Healthcare representative if purchasing replacement accessories.

*Not available in all countries*





For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



**Nellcor™ Interconnect Cables**

Description	Connector	Length	Part Number	Quantity
OxiMax™ SpO2 Interconnect Cable	GE	1.2 m	2021406-002	1
		3 m	2021406-001	1
	GE for Corometrics™ only	2.9 m	2025350-001	1
OxiSmart SpO2 Interconnect Cable	GE	2.9 m	2006644-001	1
Interconnect Cables*	DOC-4	1.2 m	2059000-001	1
	DOC-8	2.4 m	2059001-001	1
	DOC-10	3 m	2008773-001	1

\* These accessories are sold ONLY with the equipment

414248-001



70124033



70124021



70124035

**Nellcor™ Reusable Sensors**

Sensors	Patient	Patient Weight	Nellcor™ Reference	Part Number	Quantity
DuraSensor Finger Probe	Adult	< 30 kg	DS100A-1	70124021	1
Multisite SpO2 Sensor	Adult	> 1 kg	D-YS	70124033	1
Two-piece SpO2 sensor	Neonatal / Adult	< 3 kg / > 40 kg	OXI-A/N	70124035	1 (delivered with 24 adhesive strips)
	Pediatric, infant	3-40 kg	OXI-P/I	414248-001	1 (delivered with 24 adhesive strips)

Requires Compatible Interconnect Cable

These accessories are sold ONLY with the equipment

*Not available in all countries*



For CARESCAPE™ B450 / B650 / B850 monitors with a CARESCAPE™ Patient Data Module (PDM) with Nellcor™ SpO2 technology.



70124027



70124032

### Nellcor™ Disposable Sensors

Sensors	Nellcor™ Reference	Patient	Patient Weight	Part Number	Quantity
Adhesive SpO2 sensor	MAX-A	Adult	> 30 kg	70124027	24
Long Finger Adhesive sensor	MAXAL			2028117-001	24
Adhesive SpO2 sensor	MAXP	Pediatric	10-50 kg	70124022	24
	MAXI	Infant	3-20 kg	70124026	24
	MAX-N	Adult	> 40 kg	70124032	24
	Neonate	< 3 kg			

Requires Compatible Interconnect Cable  
These accessories are sold ONLY with the equipment



### Nellcor™ Fixation Clip

Description	Nellcor™ Reference	For Use With	Part Number	Quantity
SpO2 Ear-Clip	D-YSE	D-YS	70124034	1

These accessories are sold ONLY with the equipment




### Nellcor™ Sensor Adhesive Tape

Product	Patient	Nellcor™ Reference	Use With	Part Number	Quantity
Sensor Wrap Replacement Tape	Adult/Neonatal	ADH-A/N	OXI-A/N, D-YS	2016130-001	100
	Pediatric/Infant	ADH-P/I	OXI-P/I, D-YS	2016131-001	100

These accessories are sold ONLY with the equipment

*Not available in all countries*

# Temperature



2104178-001



2104178-002



2104179-001

## Temperature Cables

Description	Connection	Length	Part Number	Quantity
400 Series for Reusable Probes	Single	0.5 m	2104178-001	1
	Dual	0.5 m	2104178-002	1
400 Series for disposable Probes	Single	3.6 m	2104179-001	1

## Product compatibility

- CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors

*Not available in all countries*





2107176-032



2107176-033

Temperature Probes Reusable					
Description	Connection	Matching Code on the Label	Length	Part Number	Quantity
Esophageal/rectal	Adult	1010132011GE	280 cm	2107176-031	1
	Pediatric	1010132010GE	280 cm	2107176-032	1
Skin	Adult	1010132023GE	280 cm	2107176-033	1

### Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP and E-PT modules
- B105/125 patient monitors
- B20/B40 and B20i/B40i patient monitors
- Dash™ 3000/4000/5000 patient monitors
- Solar™ 8000M with TRAM module

Specifications available on request

*Not available in all countries*



# Invasive Blood Pressure



Argon™ - Single



Edwards™ LS TruWave™ - Single



Argon™ - Dual



Edwards™ LS TruWave™ - Dual

## Invasive Blood Pressure Cables

Description	Connection	Length	Part Number	Quantity
Argon™	Single	3.6 m	2104166-001	1
		1.2 m	2104166-002	
	Dual	3.6 m	2104166-003	1
		1.2 m	2104166-004	
Edwards™ LS TruWave™	Single	3.6 m	2104162-001	1
		1.2 m	2104162-002	
	Dual	3.6 m	2104162-003	1
		1.2 m	2104162-004	


## Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-SPMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

*Not available in all countries*

# INVASIVE BLOOD PRESSURE CABLES



ICU Medical Transpac-IV® - Dual



Utah/Biotrans™ - Single



Utah/Biotrans™ - Dual



Dual Adapter Cable

Invasive Blood Pressure Cables				
Description	Connection	Length	Part Number	Quantity
ICU Medical Transpac-IV®	Single	3.6 m	2104158-001	1
		1.2 m	2104158-002	
	Dual	3.6 m	2104158-003	1
		1.2 m	2104158-004	
Utah/Biotrans™	Single	3.6 m	2104107-001	1
		1.2 m	2104107-002	
	Dual	3.6 m	2104107-003	1
		1.2 m	2104107-004	
Dual Adapter Cable	-	30 cm	2103966-001	1

## Product compatibility

- CARESCAPE™ Monitors B450/B650/B850 with CARESCAPE™ Patient Data Module (PDM), Patient Side Module E-PSMP, E-PP, E-PT, E-COP(SV) and E-PiCCO™ modules
- B105/B125 patient monitors
- B20/B40 and B20i/B40i patient monitors

Specifications available on request

*Not available in all countries*

# Cuffs



Hoses					
Description	Patient	Compatibility	Length	Part Number	Quantity
 <p>2058203-001 2058203-002 2058203-003</p>	Adult/ Pediatric	V100/VC150, DINAMAP Pro Series, DINAMAP Procure	3.6 m	2058205-001	1
			7.2 m	2058205-002	
		Block style monitor connector for Coro, PDM and some older legacy monitors	1.2 m	2058203-001	
			3.6 m	2058203-002	
			7.2 m	2058203-003	
		CARESCAPE ONE	1.2 m	2087389-001	
			3.6 m	2087389-002	
			7.2 m	2087389-003	
		 <p>2058204-002 2058204-003</p>	Neonate	CARESCAPE ONE	
3.6 m	2089791-003				
Block style monitor connector for Coro, PDM and some older legacy monitors	2.4 m			2058204-002	
	3.6 m			2058204-003	
V100/VC150, DINAMAP Pro Series, DINAMAP Procure	3.6 m			2058206-001	

Not available in all countries

# SOFT-CUF



Our most popular cuff, the CRITIKON™ SOFT-CUF™ is designed with soft, absorbent material with rounded corners for maximizing patient comfort during extended stays at the hospital. These cuffs are ideal for patients with longer term stays or for protocols where the cuff travels with the patient through their total length of time in the hospital. It is economically priced to help support your infection control efforts.

20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Infant	Orange	8-13 cm	SFT-P1-2A
Child*	Green	12-19 cm	SFT-P2-2A-INT
Child Long	Green	12-19 cm	SFT-P2-2A-L
Small adult*	Royal Blue	17-25 cm	SFT-A1-2A-INT
Small adult Long	Royal Blue	17-25 cm	SFT-A1-2A-L
Adult*	Navy	23-33 cm	SFT-A2-2A-INT
Adult Long*	Navy	23-33 cm	SFT-A2-2A-L-INT
Large Adult*	Pink	31-40 cm	SFT-A3-2A-INT
Large Adult Long	Pink	31-40 cm	SFT-A3-2A-L
Thigh	Brown	38-50 cm	SFT-T1-2A

Easy-to-use, fast **DINACLICK connector** enables nurses to link the cuff to the NIBP hose quicker than a pair of traditional tube connectors. The audible click lets users know there is a proper connection.



\*Has a respective part number with an "-INT" suffix, denoting that the Instructions For Use (IFU) is available in multiple languages

## ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult	2059303-001
20 Cuffs : 6 Infant, 10 Child, 4 Small Adult	2059303-002
20 Cuffs : 3 Small Adult, 7 Adult, 3 Adult Long, 3 Large Adult, 4 RADIAL-CUF (forearm)	2059303-003
20 Cuffs : 2 Infant, 2 Child, 2 Small Adult, 6 Adult, 2 Adult Long, 2 Large Adult, 2 RADIAL-CUF (forearm), 2 Thigh	2059303-004

Not available in all countries







CRITIKON™ DURA-CUF™ is gentle and designed for durability. Constructed with materials not made with PVC or rubber latex.

5 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Infant	Orange	8-13 cm	DUR-P1-2A
Child	Green	12-19 cm	DUR-P2-2A
Child long	Green	12-19 cm	DUR-P2-2A-L
Small Adult	Royal Blue	17-25 cm	DUR-A1-2A
Small Adult Long	Royal Blue	17-25 cm	DUR-A1-2A-L
Adult w/ Hanger	Navy	23-33 cm	DUR-A2-2A
Adult long	Navy	23-33 cm	DUR-A2-2A-L
Large Adult	Wine	31-40 cm	DUR-A3-2A
Large Adult Long	Wine	31-40 cm	DUR-A3-2A-L
Thigh	Brown	38-50 cm	DUR-T1-2A



**ASSORTMENTS DINACLICK CONNECTOR**

PRODUCT	2-TUBE DINACLICK
3 cuffs: 1 Small Adult, 1 Adult, 1 Large Adult	2059301-001
6 cuffs : 2 Infant, 3 Child, 1 Small Adult	2059301-002
6 cuffs : 1 Small Adult, 3 Adult, 1 Adult Long, 1 Large Adult	2059301-003
6 cuffs : 1 Infant, 1 Child, 1 Small Adult, 1 Adult, 1 Large Adult, 1 Thigh	2059301-004

*Not available in all countries*

# CLASSIC-CUF

CRITIKON™ CLASSIC-CUF™ BP cuffs are designed and clinically validated for limited reuse on multiple patients. Made with durable vinyl construction, this time-tested CRITIKON blood pressure cuff is both easy-to-clean and strong enough to withstand repeated inflations – yet economical enough for single patient- use.



5 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
Infant	Orange	8-13 cm	CLA-P1-2A
Child	Green	12-19 cm	CLA-P2-2A
Small Adult	Light Blue	17-25 cm	CLA-A1-2A
Adult	Navy	23-33 cm	CLA-A2-2A
Adult long	Navy	23-33 cm	CLA-A2-2A-L
Large Adult	Pink	31-40 cm	CLA-A3-2A
Large Adult Long	Pink	31-40 cm	CLA-A3-2A-L
Thigh	Brown	38-50 cm	CLA-T1-2A



## ASSORTMENTS DINACLICK CONNECTOR

PRODUCT	2-TUBE DINACLICK
3 Cuffs : 1 Small Adult, 1 Adult, 1 Large Adult	2059303-001
20 Cuffs : 5 Small Adult, 10 Adult, 5 Large Adult	2059303-002

Not available in all countries





GE Healthcare offers the first BP forearm cuff clinically validated for use on obese or overly muscular adult patients. It provides proper fit and comfort to patients, while giving reliable blood pressure readings that support quality clinical care.

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE DINACLICK
OBESE ADULT	PURPLE	26-36 CM (FOREARM)	SFT-F1-2A



*Not available in all countries*

# NEONATAL



Providing big benefits for your tiniest patients, GE Healthcare's neonatal cuffs are designed with the accuracy you trust and the gentle comfort you expect. We offer a full range of neonatal sizes with the choice of extra soft, absorbent material or easy-to-clean vinyl.

The playful, color-coded animal designs help you find the right size at a glance and they are economically priced for single-patient-use.

20 cuffs/box

SIZE	COLOR	LIMB CIRCUMFERENCE	2-TUBE NEO-SNAP	
			CLASSIC-CUF	SOFT-CUF
1	Orange	3-6 cm	CLA-N1-2B	SFT-N1-2B
2	LT. Blue	4-8 cm	CLA-N2-2B	SFT-N2-2B
3	Green	6-11 cm	CLA-N3-2B	SFT-N3-2B
4	Navy	7-13 cm	CLA-N4-2B	SFT-N4-2B
5	Pink	8-15 cm	CLA-N5-2B	SFT-N5-2B


## ASSORTMENTS NEO-SNAP CONNECTOR

PRODUCT	2-TUBE NEO-SNAP
20 cuffs SOFT-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5	2059305-002
20 cuffs CLASSIC-CUF : 2 Neo #1, 3 Neo #2, 5 Neo #3, 5 Neo #4, 5 Neo #5	2059306-002

*Not available in all countries*



# Cardiac Output



ICU-Bed 1  
Medical

10:23

HR /min 50, J150

**60**

ECG PVC 0 /min

ST mm

V1	0.4	II	0.4	V5	0.4
V2	0.4	III	0.4	V6	0.4
V3	0.4	aVF	0.4	I	0.4
V4	0.4	aVL	0.0		

SpO2 % 90, J OFF

**97**

P1 Art 1 mmHg Syst 80, J160

**120/82** (97)

P2 CVP mmHg Contrl Mean 0, J15

**(13)**

P4 PA mmHg Contrl Syst 10, J40


**25/10** (16)

C.O. & C.I.

C.O. 4.5 C.I. 2.2

l/min l/min/m2

## CARDIAC OUTPUT



### Cardiac Output Cables

Description	Length	For use with	Part Number	Quantity
Cardiac Output Cable	3.6 m	CARESCAPE™ Monitors B450, B650 and B850 with the CARESCAPE™ Patient Data Module (PDM)	2104180-001	1
	2.4 m		2104180-002	1

Specifications available on request

Not available in all countries

## OTHER ACCESSORIES

Other Accessories				
Description	For use with	Manufacturer Reference	Part Number	Quantity
CO Set for Room Temperature Injectate, disposable	E-COP, E-COPsv	93610	16574VS	1
Catheter Connecting Cable	E-COP, E-COPsv	COM2-CC	16590-HEL	1
Thermodilution Catheter, disposable	E-COP, E-COPsv	131HF7	M1238258	5

These accessories are sold ONLY with the equipment

## PiCCO™ CABLES

PiCCO™ Cardiac Output Cables				
Description	Length	For use with	Part Number	Quantity
Continuous Cardiac Output Cable	4.1m	E-PiCCO module	M1154585	1
Injectate Sensor Cable	0.9m		M1182448	1

These accessories are sold ONLY with the equipment



M1182448

GE Healthcare compatible Cardiac Output probes are available from multiple manufacturers. You can purchase the following product portfolios from our partners directly:

- Vyair Medical
- Bath Probes
- Argon Medical
- Spectramed Temperature Probes
- Edwards Life Sciences
- Injectate Bath Probes
- Flow-through Injectate Temperature Probes
- Getinge
- PiCCO Accessories

*Not available in all countries*



# Adequacy of Anesthesia



## ENTROPY

### Sensors

Product	Use	Part Number	Quantity
Entropy™ EasyFit Sensor	Disposable	M1174413	25
Entropy™ Sensor		M1038681	

### Cable

Product	Length	Part Number	Quantity
Entropy™ Cable	3.5 m	M1050784	1

### Product compatibility

E-ENTROPY module


The GE Healthcare Entropy™ module, E-ENTROPY, and accessories are indicated for adult and pediatric patients older than 2 years within a hospital for monitoring the state of the brain by data acquisition of electroencephalograph (EEG) and frontal electromyograph (FEMG) signals during general anesthesia.

*Not available in all countries*



M1038681

M1174413





Electromyography (EMG)  
Requires 5 electrodes



Kinemyography (KMG)  
Requires 2 electrodes



888418/897439



888416

Sensors				
Product	Use	Patient	Part Number	Quantity
NMT MechanoSensor	Reusable	Adult	888418	1
		Pediatric	897439	
NMT ElectroSensor		Adult, Pediatric	888416	1
NMT Regional Block Adapter		Adult, Pediatric	888417	1

Accessories				
Product	Use	Patient	Part Number	Quantity
NMT Electrodes	Disposable	Adult, Pediatric	57268-HEL	30

Cable			
Product	Length	Part Number	Quantity
NMT Sensor Cable	3 m	888414	1


### Product compatibility

E-NMT module

*Not available in all countries*



# Neurology



898050

Reusable Preconfigured EEG Leadsets			
Product	Length	Part Number	Quantity
Basic montage	0.6 m	898050	1
AEP montage		898052	1

*Not available in all countries*

# EEG ELECTRODES



2083319-001



2083320-001

Disposable EEG Electrodes		
Product	Part Number	Quantity
EEG Clip electrodes Ag/AgCl	M1240794	40 pouches of 25 (1,000 total)
	572686	500
EEG Prewired cup electrodes	2083319-001	46 pouches of 10 (460 total)
EEG Prewired wet gel electrodes	2083320-001	42 pouches of 10 (420 total)

All electrodes within a set-up should be similar

## Product compatibility

E-EEG(X) module and N-EEG(X) headbox

AEP Accessories		
Product	Part Number	Quantity
Disposable ear tips for use with earphone set	8001233	20

To use with 898052

## Product compatibility

E-EEG(X) module and N-EEG(X) headbox

*Not available in all countries*



# Patient Spirometry



896952



8001948



733910-HEL



73393

Sensors				
Description	Type	Approved with Modules	Part Number	Quantity
D-lite™ +	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE <sup>1</sup> , E-sCAiOVX	896952	50
Pedi-lite+			8001948	
D-lite™	Re-usable		733910-HEL	1
Pedi-lite			73393	

Specifications available on request

(1) Not available in the US

Not available in all countries

# SPIROMETRY



890031

## Spirometry Tubes

Description	Length	Type	Approved with Modules	Part Number	Quantity
Spirometry Tubes	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE <sup>1</sup> , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	890031	5
	3m			884101	
	6m			891191-HEL	



D-lite™ ++ kit

## Kits

(include gas sample line, sensor and spirometry tube)

Description	Length	Type	Approved with Modules	Part Number	Quantity
D-lite™ ++ kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE <sup>1</sup> , E-sCAiOVX	2104297-001	20
	3m			2104297-002	
D-lite™ + kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE <sup>1</sup> , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8004381	20
	3m			8004382	
Pedi-lite+ kit	2m	Single Patient Disposable	E-sCOV, E-sCOVX, E-sCAiOV, E-sCAiOVE <sup>1</sup> , E-sCAiOVX, E-COV, E-COVX, E-CAiOV, E-CAiOVX	8002718	20
	3m			M1032634	

Specifications available on request



Pedi-lite+ kit

(1) Not available in the US

Not available in all countries



# Gas Monitoring



## ICU AIRWAY ADAPTER SETS



### Airway Gas Sampling Sets

Description	Type	Approved with Modules	Length	Part Number	Quantity
Adult/Pediatric	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX	3 m	8090121531	20 pcs/box
			2 m	8090121521	
Infant		E-sCO, E-sCOV, E-sCOVX	2 m	8090121522	20 pcs/box

*Not available in all countries*

# GAS SAMPLING LINES AND ADAPTERS



Gas Sampling Lines					
Description	Type	Approved with Modules	Length	Part Number	Quantity
PVC/PE gas sampling line	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE <sup>1</sup> , E-sCAiOV, E-sCAiOVX, E-sCAiOVE <sup>1</sup> , N-CAiO	2 m	2097307-001	10
			3 m	2097307-002	
			6 m	2097307-003	5



3090121002

Gas Monitoring Airway Adapter				
Description	Type	Approved with Modules	Part Number	Quantity
Airway Adapter with Luer connector	Single Patient Disposable	E-miniC, E-sCO, E-sCOV, E-sCOVX, E-sCAiO, E-sCAiOE <sup>1</sup> , E-sCAiOV, E-sCAiOVX, E-sCAiOVE <sup>1</sup> , N-CAiO	3090121002	25

# WATER TRAPS



881319-HEL



M1200227



876446-HEL



M1182629



8002174

Water Traps – ICU				
Description	Type	Approved with Modules	Part Number	Quantity
D-fend™+	Single Patient Disposable	E-CO, E-COV, E-COVX	881319-HEL	10
D-fend™ Pro+		E-sCO, E-sCOV, E-sCOVX	M1200227	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

Water Traps – Anesthesia				
Description	Type	Approved with Modules	Part Number	Quantity
D-fend™	Multi-Patient	E-CAiO, E-CAiOV, E-CAiOVX	876446-HEL	10
D-fend™ Pro		E-sCAiO, E-sCAiOE <sup>1</sup> , E-sCAiOV, E-sCAiOVX, E-sCAiOVE <sup>1</sup>	M1182629	
Mini D-fend™		E-miniC	8002174	

Specifications available on request

*Not available in all countries*



ANESTHESIA & RESPIRATORY CARE



## **ANESTHESIA & RESPIRATORY CARE**

### **Absorbents**

- Absorbers
- Canisters
- Dust Filter

### **Vaporizer Bottle Adapters**

### **Filters & Breathing Masks**

### **Reservoir Bags & Gas Scavenging**

- Anesthesia Breathing Bags
- Scavenging Bags, Tubes, & Connectors

### **Patient Breathing Circuits**

- Circuits
- Circuit Components
- Circuit Kits
- Circuit Startup Kits

### **Flow Sensors - Respiratory**

- Cables
- Heater
- Respiratory Flow Sensors and Expiratory Valve Assemblies

### **Flow Sensors - Anesthesia**

### **Nebulizer Accessories**

### **Anesthesia O2 Cells**

### **Respiratory Safety Guard**

### **Pressure Catheters**





# Anesthesia & Respiratory care



Disposable CO2 Absorbers			
Description	Compatibility	Part Number	Quantity
AMSORB® Plus G-Can Absorber, 1.0L	All Aisys™, Avance™, Aespire™, 9100c and ADU variants	2105489-003	8
AMSORB® Plus Prefilled Cartridge, 1.2L	All Aestiva™ variants	2105489-004	12
AMSORB® Plus CareCan Absorber, 1.4L	Carestation™ 600 series & 750 series	2105489-006	8

Reusable CO2 Absorbent Components			
Description	Compatibility	Part Number	Quantity
AMSORB® Plus Jerican Loosefill Absorbent, 5.0L	GE Healthcare Reusable Canisters	2105489-007	2
ABS Reusable CO <sub>2</sub> Canister	All Aisys™, Avance™, Aespire™ & 9100c variants	1407-7004-000	1
CBS Reusable CO <sub>2</sub> Canister	Carestation™ 600 series & 750 series	2071165-001-S	1
Foam Dust Filter	ABS Reusable Canister	2105489-001	40

Easy-Fil Vaporizer Bottle Adapters			
Vaporizer Compatibility	Color Indicator	Part Number	Quantity
Isoflurane Easy-Fil	Purple	2106572-001	1
Halothane Easy-Fil	Red	2106572-002	
Sevoflurane Easy-Fil	Yellow	2106572-003	

# FILTERS AND BREATHING MASKS



HMEFs and Filters – Disposable				
Description	Size	Sampling Port	Part Number	Quantity
HMEF	Neonatal	Included	2106570-006	50
	Pediatric		2106570-009	
	Adult		2106570-010	
Bacterial and Viral Filter	Adult	Included	2106570-007	50
HEPA Filter	Adult		2106570-008	



Anesthesia Breathing Masks – Disposable			
Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-012	50
	1 – Small Child	2106570-013	
22mm	2 – Child	2106570-014	50
	3 – Large Child	2106570-015	
	4 – Adult	2106570-016	50
	5 – Large Adult	2106570-017	
	6 – Extra Large Adult	2106570-018	

Anesthesia Breathing Masks – Reusable			
Mask Connection	Mask Size	Part Number	Quantity
15mm	0 – Neonatal	2106570-011	1
		2106570-029	10
	1 – Small Child	2106570-019	1
		2106570-030	10
22mm	2 – Child	2106570-020	1
		2106570-031	10
	3 – Large Child	2106570-021	1
		2106570-032	10
	4 – Adult	2106570-022	1
		2106570-033	10
	5 – Large Adult	2106570-023	1
		2106570-034	10

*Not available in all countries*



# RESERVOIR BAGS AND GAS SCAVENGING

## Anesthesia Manual Breathing Bags – Reusable



Bag Connection Size	Material	Bag Size	Part Number	Quantity
22mm	Clear Silicone	0.5L	2106570-001	1
			2106570-025	5
		1.0 L	2106570-002	1
			2106570-026	5
		2.0 L	2106570-003	1
			2106570-027	5
		3.0 L	2106570-004	1
			2106570-028	5

## Passive Scavenging Tubes – Disposable



Connection	Length	Type	Part Number	Quantity
30M-30M	5 m	Tube	2105488-006	1
			2105488-022	10

## Gas Scavenging Bag – Reusable



Size	Compatibility	Connection	Part Number	Quantity
3.0L	Aisys, Avance, Aespire™ (ABS Breathing Systems)	30M	2106570-024	1

## Straight Connectors – Scavenging

Connection	Reprocessable	Part Number	Quantity
30M-19M	Disposable	M1003947	1
		1500-3376-000	5
30M-309M		M1003134	1

*Not available in all countries*



## PATIENT BREATHING CIRCUIT KITS – START UP

All startup kits include patient breathing circuit, reservoir bag tube, reservoir bag, gas sampling lines and gas sampling exhaust lines. These are primarily intended for initial anesthesia machine setup.



### Patient Breathing Circuit Kits – Start Up – Adult

Type	Gas Sampling Line	Gas Sampling Exhaust Line	Limb Circuit	Breathing Circuit	Breathing Bag		Part Number	Quantity
				Length	Manual Tube 1.5 M	Size		
Disposable	3m PVC/PE	0.18m & 1m	Dual	Expandable up to 3 m	Straight	3L	2096533-001	1
						2L	2096533-003	1
			Single Split Lumen	1.8 m	Angled*	2L	2096536-002	2
					Straight	2L	2096536-003	2
					Angled*	2L	2105488-019	1
Reusable			Dual	1.5 m	Straight	2L	2096534-009	1
					Angled*	2L	2096536-001	1

\*For use with bag support arm

Not available in all countries



# PATIENT BREATHING CIRCUITS

Our wide range of patient breathing circuits and components offer robust and high quality options for both disposable and reusable applications.

Patient Breathing Circuit Kits									
Type	Size	3m PVC/PE Gas Sampling Line	Circuit Type	Breathing Circuit	Breathing Bag		Includes	Part Number	Qty
				Length	Manual Tube	Size			
Disposable	Adult	Yes	Dual Limb	Expandable up to 3 m	Not included	3L	-	2096535-006	20
				1.8 m	Angled*, 1 M		-	2105488-013	20
					Straight, 1.5 M		-	2105488-014	20
		No	Single Split Lumen	1.8 m	Straight	3L	-	2105488-010	20
				Dual Limb	Expandable up to 1.8 m	Not included	-	2105488-011	20
			1.5 m		Not included	-	2x Water Traps	2096535-001	20
	Straight, 0.8 M			-	1x Water Trap	2105488-017	20		
	Not Included				Angled*, 2 M	2L		2105488-012	20
	Neonatal	Yes	Dual Limb	1.5 m	Not included	1L	HME & Filter	2105488-001	20
		No						2105488-021	1
Reusable	Adult	No	Dual Limb	1.5 m	Not included	2x Water Traps	2096534-005	1	
	Pediatric						2096534-006	1	

\*For use with bag support arm



Not available in all countries

# PATIENT BREATHING CIRCUITS



2096534-007  
Reusable, Pediatric



2105488-020  
Disposable, Pediatric



2105488-015  
Disposable, Adult



2105488-009  
Disposable, Adult, Inteliflo® Single Split Lumen

## Patient Breathing Circuits

Type	Size	Circuit Type	Breathing Circuit Length	Patient Elbow	Gas Sampling Port Location	Part Number	Quantity
Reusable	Adult	Dual Limb	1.5 m	Double Swivel, Removable	Elbow	2096534-007	1
	Pediatric					2096534-008	1
Disposable	Adult	Single Split Lumen	1.8 m	Fixed, Removable	Elbow	2105488-009	20
		Dual Limb	1.5 m	Double Swivel, Removable	Patient Wye	2096535-002	20
	Fixed, Removable			Elbow	2105488-015	20	
	Pediatric	Dual Limb	1.5 m	Double Swivel, Removable	Patient Wye	2105488-020	20



## Non-Recirculating Breathing Circuits – Jackson Rees

Type	Breathing Circuit Length	Bag Size	Gas Sampling Port Location	Part Number	Quantity
Disposable	0.8 m	1 L	Elbow	2105488-008	20

Not available in all countries



Reusable Patient Breathing Circuit Tubes					
Size	Material	Length	Connections	Part Number	Quantity
Adult	Hytrel®	0.35 m	22F-22F (90 Deg)	2096534-001	1
		0.75 m		2096534-002	1
		0.9 m		2096534-003	1
		1.5 m		2096534-004	1
	Silicone	1 m	22F-22F	2096534-101	1
		1.5 m		2096534-102	1
Pediatric	Silicone	1.5 m	22F-10F	2096534-103	1
	Hytrel®	1.1 m	22F-15F	2105488-007	1
				2105488-023	10



Not available in all countries

# PATIENT BREATHING CIRCUITS



Patient Breathing Circuit Components				
Description	Size	Reprocessable	Part Number	Quantity
Y-Piece Connector Kit, 15M-10F	Neonatal	Reusable	2096534-104	10
Patient Angled Y-Piece			2105488-002	10
Patient Straight Y-Piece - Swivel			2105488-004	10
Straight Connector - CBS Bag Arm, 22M-22M	Adult	Reusable	2096533-004	1
ABS Bag Arm Connector, 22M-22M	All		2105488-005	10
Patient Elbow - Swivel	Adult	Reusable	2105488-003	10
Patient Elbow		Disposable	2105488-018	50
Gas Return Adapter, 22M-22F	Adult	Disposable	2105488-016	10

*Not available in all countries*







Respiratory Flow Sensors and Expiratory Valve Assemblies (EVA)				
Description	Reprocessable	Patient Type	Part Number	Quantity
Flow Sensor Only	Reusable	All	1505-3231-000	1
Whole Assembly (housing and flow sensor)			1505-3848-000	1
Housing Only (includes water trap)			1505-8568-000	1
Flow Sensor Only (Neonatal)		Neonatal	1505-3272-000	1
Flow Sensor Only	Single Patient Use	All	1505-3231-001	1
Whole Assembly (housing and flow sensor)			1505-3848-001	1
Housing Only (includes water trap)			1505-8568-001	1
Flow Sensor Only			1505-3231-010	10
Whole Assembly (housing and flow sensor)			1505-3848-010	10
Housing Only (includes water trap)			1505-8568-010	10

## Product compatibility

CARESCAPE™ R860 and Engström Ventilators

*Not available in all countries*

## FLOW SENSORS – RESPIRATORY



M1200693



M1188723

### Respiratory Flow Sensor and Expiratory Valve Components

Description	Patient Type	Part Number	Quantity
Flow Sensor Cable	Neonatal	1505-5604-000	1
Exhalation Valve Heater	All	M1200693	1
Exhalation Valve Heater Cable		M1188723	1

### Product compatibility

CARESCAPE™ R860 and Engström Ventilators

## FLOW SENSORS – ANESTHESIA



2087640-001-S



2089610-001-S



2069358-001-S

### Anesthesia Flow Sensors

Machine Compatibility	Autoclavable	Part Number	Quantity
Aisys™, Avance™, Aespire™, Aestiva™	Yes	2087640-001-S	1
	No	2089610-001-S	
Carestation™ 600 series & 750 series	Yes	2069358-001-S	1

*Not available in all countries*



## NEBULIZER ACCESSORIES



Aerogen® Nebulizer Accessories			
Description	Patient Type	Part Number	Quantity
Aerogen® Pro Nebulizer kit with adult and ped tee and cable	Adult and Pediatric	1505-3846-000	1
Nebulizer head with filler cap	All	AG-AP1000	
Nebulizer Cable		1505-5602-000	
Nebulizer tee with plug	Adult	AG-AP1010	5
	Pediatric	AG-AP1020	
	Neonatal	AG-AP1035	
Nebulizer replacement kit		AG-AP1100	1
Aerogen® Solo Nebulizer	Adult	AG-AS3100-S	5
		AG-AS3200-S	10
Nebulizer tee	Adult	AG-AS3010	10
	Pediatric	AG-AS3020	
	Neonatal	AG-AS3035	

### Product compatibility

CARESCAPE™ R860 and Engström Ventilators

These accessories are sold ONLY with the equipment.

## ANESTHESIA O2 CELLS

Anesthesia O2 Cells			
Description	Machine Compatibility	Part Number	Quantity
O2 Sensor ONLY	Aisys, Avance, Aespire™	M-10	1
O2 Sensor and cable		2082006-001	
O2 Sensor with Adapter and O-Ring	Carestation™ 600 series & 750 serie	2071018-001-S	



2071018-001-S

*Not available in all countries*



## RESPIRATORY INSPIRATORY SAFETY GUARD



Respiratory Inspiratory Safety Guard			
Machine Compatability	Patient Type	Part Number	Quantity
CARESCAPE™ R860	All	5830446	1
		5514015	10

## PRESSURE CATHETERS




Respiratory Pressure Catheters			
Description	Machine Compatability	Part Number	Quantity
Intratracheal Pressure Catheter	CARESCAPE™ R860 and Engström Ventilators	M1045564	10

*Not available in all countries*



# DIAGNOSTIC CARDIOLOGY



## DIAGNOSTIC CARDIOLOGY

### MAC™, CASE™ and CardioSoft™

#### ECG

- Accessories
- Cables
- Electrodes
- Leadwires
- Recording Paper

### SEER® Holter Monitors

#### ECG – CM 4000

- Leadwire Sets
- Accessories

#### ECG – SEER® 12

- Leadwire Sets
- Accessories

#### ECG – SEER® 1000

- Leadwire Sets
- Accessories

#### ECG – SEER® Light

- Leadwire and Cable Kits
- Accessories

# MAC™, CASE™ and CardioSoft™



## CABLES AND LEADWIRES – IEC

### ECG Trunk Cables – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
Trunk Cable 10-lead IEC	220 cm	2104722-001			X	X	X	X		



### Integrated Cable and Leadwire Set – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
ECG Value Cable 10-lead Cable / Leadwire IEC	Chest leads: 80 cm Limb leads: 125 cm Trunk cable: 220cm	2104726-001			X	X	X	X		



*Not available in all countries*



# LEADWIRES – IEC

## ECG Leadwire Set – Banana – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
Leadwire Set, Base 10, banana, IEC	Limb leads: 125 cm	2104724-002	X	X	X	X	X	X	X	X
Leadwire Set, 4 Add-on, banana, IEC-HEIM	Chest leads: 80 cm	2104750-001	X	X					X	X



2104724-002



2104750-001

## ECG Leadwire Set – Grabber – IEC

Description	Length	Part Number	MAC™						CASE™	CardioSoft™
			VU360	5500	3500	2000	800	600		
Base 10, IEC	Limb lead: 125 cm	2104748-001	X	X	X	X	X	X	X	X
4 Add-on, IEC-HEIM	Chest lead: 80 cm	2104725-002	X	X					X	X



2104748-001



2104725-002

Not available in all countries





# CABLES AND LEADWIRES – IEC



2104729-001



2104767-001

## Accessories

Description	Part Number	MAC™						CASE™	CardioSoft™	Quantity
		VU360	5500	3500	2000	800	600			
Adapter Clip, Wide Mouth Clip	2104729-001	X	X	X	X	X	X	X	X	14
ECG Leadwire Port Plug	2104767-001	X	X					X	X	4

# CARDIOSOFT™



## CardioSoft™ Accessories

Description	Part Number	Quantity
Supply BP Bag Belt TONOPORT	2104824-001	1

*Not available in all countries*

# ECG RECORDING PAPER



## DCAR ECG Recording Thermal Paper

Folding	White patient data area	Paper size	Red grid wide	Queue	MAC™						CASE™	Part Number	Quantity
					VU360	5500	3500	2000	800	600			
z-fold	No	A4 *	195 mm	Hole		X	X				X	2104769-001	300 sheets/ pack, 8 packs
			205 mm	Block	X			X				2104772-001	150 sheets/ pack, 10 packs
		21.5 cm x 28 cm	200 mm	Hole		X	X				X	2104768-001	300 sheets/ pack, 8 packs
		80 mm x 90 mm	75 mm	Block						X		2104773-001	280 sheets/ pack, 10 packs
		110 mm x 140 mm	100 mm	Block					X			2104775-001	200 sheets/ pack, 10 packs
		Yes	21.5 cm x 28 cm	155 mm	Hole	X	X	X					2104770-001
	Hole				X		X	X				2104771-001	150 sheets/ pack, 16 packs

For ECG paper for Patient monitoring, please go to page 6.



Not available in all countries

## REUSABLE ECG ELECTRODES

GE Healthcare reusable bulb electrodes with silver-silver chloride (Ag/AgCl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.



Reusable Bulb Electrode		
Description	Part Number	Case Quantity
Reusable Bulb Electrode	2104783-001	6



Bulb Electrode Accessories			
Description	Compatibility	Part Number	Case Quantity
Clamp Limb AHA Value	MAC™ 600, MAC™ 800, MAC™ 2000 ECG System	2104784-001	4/set
Clamp Limb IEC Value		2104785-001	

Specifications available on request

*Not available in all countries*



# SEER<sup>®</sup> Holter Monitors



## CM 4000 LEADWIRE SETS



2105803-005

### CM 4000 Leadwire Set - IEC

Description	Length	Part Number	Quantity
3-leads	80 cm	2105803-001	1
5-leads		2105803-003	1
7-leads		2105803-005	1

Specifications available on request

*Not available in all countries*



CM 4000 Accessories		
Description	Part Number	Case Quantity
Carry Pouch Reusable	2054933-014	1
Carry Bag Disposable	2054933-015	1
USB Cable	2054933-013	1

Specifications available on request



2104808-001

SEER® 12 – IEC				
Description	Length	Type	Part Number	Quantity
12-Channel Holter ECG Cable	130 cm	Round	2104808-001	1
	105 cm	Flat	2104809-001	

SEER® 12 ECG cables are compatible with SEER® 12 monitors and CardioMem™ CM 3000-12 BT monitors



2104803-001

SEER® 12 Holter Monitor Accessories		
Description	Part Number	Quantity
SEER® 12 lead portable pouch, waist belt, shoulder strap, 1 battery compartment	2104803-001	1
SEER® 12 Carrying pouch, 2 battery compartment	2014606-297	1

2104803-001 is compatible with SEER® 12, CM 3000, CM 3000-12  
 2014606-297 must be purchased with service  
 Specification available on request

*Not available in all countries*



2104811-001

SEER® 1000 Leadwire Set – IEC

Description	Length	Part Number	Quantity
5 leadwire, 2 CH	75 cm	2104810-001	1
7 leadwire, 3 CH		2104811-001	
5 leadwire, 2 CH	105 cm	2104812-001	
7 leadwire, 3 CH		2104813-001	
3 leadwire, 3 CH	75 cm	2104818-001	
	105 cm	2104819-001	

**Product compatibility**

SEER® 1000 ECG cables are compatible with SEER® 1000 24-hour Holter monitor, SEER® 1000 48-hour Holter monitor and SEER® 1000 7-day Holter monitors

Specifications available on request

Not available in all countries



2067634-073

## SEER® 1000 Holter Monitor Accessories

Description	Part Number	Quantity
Pouch with waist belt and shoulder strap	2067634-073	1
Pouch Belt Clip	2067634-166	1

Specification available on request



2067634-166

*Not available in all countries*

SEER® Light Leadwire Set Kits – IEC			
Description	Length	Part Number	Quantity
Carry Pouch and 2 Channel Patient cable	140 cm	2042686-003	1
Carry Pouch and 3 Channel Patient cable		2042686-004	
Carry Pouch and 2 Channel Patient cable	60 cm	2042686-007	1
Carry Pouch and 3 Channel Patient cable		2042686-008	

## Product compatibility

SEER® Light and SEER® Light Extend  
 Specifications available on request

*Not available in all countries*



SEER® Light Accessories		
Description	Part Number	Quantity
Compact Flash	2008750-003	1
Connect Kit	2020044-292	1
Controller Kit	2042685-007	1

### Product compatibility

SEER® Light and SEER® Light Extend  
Specifications available on request



SEER® Light Holter Monitor Accessories		
Description	Part Number	Quantity
SEER® Light, shoulder straps, waist belt, portable pouch	2104827-001*	1

### Product compatibility

2104827-001 is compatible with SEER® Light and SEER® Light Extend  
Specification available on request

*Not available in all countries*



MATERNAL INFANT CARE



## MATERNAL INFANT CARE

### Maternal Fetal Care

#### Mini Telemetry

- Recording Thermal Paper
- Remote Event Marker

#### Novii™

- Cables
- Novii™ Patch
- Tape

#### Corometrics™ and Mini Telemetry

- Cables
- Transducer Belts
- Transducer Button and Loop

#### Internal Monitoring

- Cables
- Electrode
- Pad

### Infant Care

#### Giraffe™/Panda™

- Air filters
- Clip
- Cover

#### Phototherapy accessories

- Pad Covers

#### Resuscitation

- T Piece

#### Lullaby™ Warmer


- Temperature Probes
- Temperature Probe Patches

#### Shuttle™

- Patient Restraint

# Maternal Fetal Care

- Mini Telemetry
- Novii™
- Corometrics™



## MIC Fetal Monitoring Recording Thermal Paper

Scaling	Corometrics™ 170 series	Corometrics™ 250cx series	Part Number	Quantity
30-240 bpm scaling	X	X	2104907-001	40 packs
50-210 bpm scaling	X	X	2104908-001	

For ECG paper for Patient monitoring and Diagnostic cardiology, go to pages 6 and 58.

Specifications available on request

## Accessories

Description	Part Number	Quantity
Remote Event Marker	2107349-001	1



2107349-001

*Not available in all countries*





Novii™ Accessories			
Description		Part Number	Quantity
Novii™ Patch		107-PT-004-10	10
		107-PT-004-50	50
Novii™ Interface Cable	GE FECG	105-PT-102	1
	GE MECG	105-PT-104	
	GE UA	105-PT-106	
3M™ Red Dot™ 2236 Skin Prep Tape		100-PT-007	1
FECG/MECG Y Adapter Cable, 120 Series		1442AAO	1


*Not available in all countries*





MECG Accessories - Vyaire			
Description		Part Number	Quantity
Gray MECG Cable	3 Lead Set, Model 118	1554AAO	1
	3 Lead Set IEC, Model 118	1554BAO	1
120 Green MECG Cable	3 Lead Set	1553AAO	1
	3 Lead Set IEC	1553BAO	1

Sold with equipment only.



2015827-001

Transducer Belts			
Description		Part Number	Quantity
Single patient use	Button-style, blue	2015827-001	10



2108347-001

Ultrasound Transducers		
Description	Part Number	Quantity
Ultrasound Transducer Button and Loop	2108346-001	1
TOCO Transducer Button and Loop	2108347-001	1

*Not available in all countries*





2109002-001



2109006-001

## Fetal Spiral Electrode Accessories

Description	Type	Part Number	Quantity
Cable	Rectangular Connector	2109005-001	1
	Round Connector	2109004-001	1
ATTACHMENT PAD	Large, Single Patient Use	2109006-001	50
SPIRAL ELECTRODE	Single Helix, Single Patient Use	2109002-001	50

For ECG Electrodes for Patient Monitoring, go to pages 14 and 15.

*Not available in all countries*



# Infant Care

- Giraffe™/Panda™
- Phototherapy accessories
- Resuscitation



## GIRAFFE™/PANDA™



2107356-001



2107357-001

Accessories		
Description	Part Number	Quantity
Air Filters	2107356-001	10
Bag and mask holder clip	2107357-001	1

*Not available in all countries*





#### Giraffe™ Hood Covers

Description	Part Number	Quantity
OmniBed™ Cover	6600-0825-800	1
Incubator Cover	6600-0846-800	1

Sold with equipment only

*Not available in all countries*





Flat cushioned cover with straps  
*The flat, cushioned BiliSoft™ Cover lets you swaddle baby and phototherapy pad together.*



Cushioned cover with developmental positioning foot roll and straps  
*The BiliSoft™ Nest offers sick babies the boundaries and support they need.*

BiliSoft™ Disposable Pad Covers					
Description	BiliSoft™ 1.0	BiliSoft™ 2.0	Size	Part Number	Quantity
Disposable Pad Covers, Flat		X	Small	2105422-001	20
				2105418-001	50
	X		Small	2107201-001	20
				2107200-001	50
Disposable Nest Pad Covers		X	Small	2105420-001	15
	X			2107202-001	
	X	X	Large	2105421-001	

**Product compatibility**

The BiliSoft™ Pad Covers and Nest are designed for use with the Bilisoft™ 1.0 and 2.0 Phototherapy System, which is used to treat both premature and full-term infants for hyperbilirubinemia.

Available in small and large sizes to match the size of the Light Pad\*

\* Note: the BiliSoft™ 2.0 Disposable Pad Cover (Flat and Nest) in Size Large also fits the BiliSoft™ Phototherapy System (Legacy) large pad.

Specifications available on request



BiliSoft™ 1.0



BiliSoft™ 2.0

Not available in all countries

**T-Piece Circuit:** Disposable resuscitator designed for neonates, provides accurate reliable ventilation and pulmonary resuscitation.



2109061-001  
Disposable T Piece  
with Mask

**Infant Resuscitator T-Piece Circuit Kit**

Description	Mask Size	Part Number	Quantity
Disposable T-Piece	-	2109059-001	20
Disposable T-Piece with Mask	0	2109060-001	10
	1	2109061-001	10

**Product compatibility**

With a universal connector, it fits with all commonly used manual resuscitator units and does not require an adapter. Specifications available on request



Panda™ Warmer



StAr

Not available in all countries

## LULLABY™ WARMER

Temperature Probes - Vyaire		
Description	Part Number	Quantity
Disposable Temperature Probe	6600-0196-700	50
	6600-0208-700	10
Reusable Temperature probe – Lullaby™ Warmer	2044384-001	1

Sold with equipment only in most countries.

Temperature Probe Patches - Vyaire		
Description	Part Number	Quantity
Patch Heat Reflecting ROHS	0203-1980-300	50
Heat Reflecting Patches	2044385-001	50

Sold with equipment only in most countries.

## SHUTTLE™

Transport		
Description	Part Number	Quantity
Patient restraint	M1169540	5

*Not available in all countries*







Product may not be available in all countries and regions.  
Full product technical specifications are available upon request.  
Contact a GE Healthcare representative for more information.

**Please visit: [www.gehealthcare.com/promotional-locations](http://www.gehealthcare.com/promotional-locations)**

Data subject to change.

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges.

From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

© 2020 General Electric Company – All rights reserved.

GE, the GE monogram, MAC<sup>™</sup>, CASE<sup>™</sup>, CardioSoft<sup>™</sup>, CardioMem<sup>™</sup>, SEER<sup>®</sup>, Aespire<sup>™</sup>, Aestiva<sup>™</sup>, Aisys<sup>™</sup>, Carestation<sup>™</sup>, CARESCAPE<sup>™</sup>, Solar<sup>™</sup>, Dash<sup>™</sup>, Entropy<sup>™</sup>, D-lite<sup>™</sup>, D-fend<sup>™</sup>, Corometrics<sup>™</sup>, Shuttle<sup>™</sup>, Panda<sup>™</sup>, Lullaby<sup>™</sup> and TruSignal<sup>™</sup> are trademarks of General Electric Company. ICU Medical Transpac is a registered trademark of ICU Medical, Inc.; Utah/Biotrans is a registered trademark of Utah Medical Products, Inc. Edwards Lifesciences TruWave is a trademark of Edwards Lifesciences Corporation. Argon is a trademark of Argon Medical Devices. AMSORB and G-Can are trademarks of Armstrong Medical. Hytrel is a registered trademark of E. I. du Pont de Nemours and Company or its affiliates. Inteliflo is a registered trademark of Flexicare Medical Ltd. Nellcor and Oximax are trademarks of a Medtronic Company. Masimo and SET are trademarks of Masimo Corporation. All other trademarks are the property of their respective owners.

Reproduction in any form is forbidden without prior written permission from GE. Nothing in this material should be used to diagnose or treat any disease or condition. Readers must consult a healthcare professional.

JB77537XE(1)