

CardioSoft Accessories

Part Number	Description	Case Quantity
2104824-001	Supply BP Bag Belt TONOPORT V	1

ECG Cables and Leadwires - MAC, CASE and Cardiosoft

ECG cables and leadwires for the MAC[™], CASE[™], Cardiosoft[™] systems help you reliably make the correct connection to GE Healthcare devices. These high-quality reusable cables complement your equipment, enabling accurate measurements. GE Healthcare cables and leadwires undergo strict testing to help ensure quality, reliability and safety.



ECG Cables and Leadwires - AHA Compliant

Part Number	Description	MAC VU360	MAC 5500	MAC 3500	MAC 2000	MAC 800	MAC 600	CASE	Cardio- Soft
ECG Trunk Ca	ables – AHA								
2104723-001	Trunk Cable 10-lead AHA			×	×	×	X		
Integrated C	able and Leadwire Set - AHA								
2104727-001	ECG Value Cable 10-lead Cable / Leadwire AHA				×	×	×		
ECG Leadwires - Banana - AHA									
2104724-001	Leadwire Set, Base 10, banana, AHA	×	х	×	×	×	×	×	×
2104751-001	Leadwire Set, 4 Add-on, banana, AHA-HEIM	x	х					×	x
2104752-001	Leadwire Set, 3 Add-on, banana, AHA-PED (V3R, V4R, V7)	x	х					×	
ECG Leadwir	es – Grabber – AHA								
2104749-001	Leadwire Set, Base 10, grabber, AHA	×	х	х	×			х	x
2104725-001	Leadwire Set, 4 Add-on, grabber, AHA-HEIM	x	×					×	x

Specification available on the website



ECG Cables and Leadwires - IEC Compliant

Part Number	Description	MAC VU360	MAC 5500	MAC 3500	MAC 2000	MAC 800	MAC 600	CASE	Cardio- Soft
ECG Trunk C	ables – IEC								
2104722-001	Trunk Cable 10-lead IEC			×	×	х	х		
Integrated C	able and Leadwire Set – IEC								
2104726-001	ECG Value Cable 10-lead Cable / Leadwire IEC				×	×	x		
ECG Leadwin	res - Banana - IEC								
2104724-002	Leadwire Set, Base 10, banana, IEC	×	×	×	×	×	×	×	х
2104750-001	Leadwire Set, 4 Add-on, banana, IEC-HEIM	x	x					×	х
ECG Leadwin	res - Grabber - IEC								
2104748-001	Leadwire Set, Base 10, grabber, IEC	×	×	×	×			×	х
2104725-002	Leadwire Set, 4 Add-on, grabber, IEC-HEIM	х	x					×	Х
Accessories									
2104729-001	Adapter Clip, Wide Mouth Clip, 14/pack	×	×	×	×	×	×	×	х
2104767-001	ECG Leadwire port plug 4/pack	×	×					X	X

ECG Recording Paper

Part Number	Description	Case Quantity
Size A4 Paper ¹		
2104769-001	Thermal paper A4, red grid 195mm wide, z-fold, hole queue	300 sheets, 8 packs
2104772-001	Thermal paper A4, red grid 205mm wide, z-fold, block queue	150 sheets, 10 packs
Size 8.5 x 11 in	ch Red Grid Paper	
2104768-001	Thermal paper 8.5in x 11in, red grid 200mm wide, z-fold, hole queue	300 sheets, 8 packs
2104773-001	Thermal paper 80mm wide, red grid 75mm wide, z-fold, block only, patient monitoring paper queue	280 sheets, 10 packs
2104775-001	Thermal paper 110mm wide, red grid 100mm wide, z-fold, block queue	200 sheets, 10 packs, 40 packs /case
Size 8.5 x 11 inc	ch Red Grid/White Data Area	
2104770-001	Thermal paper 8.5in x 11in, white patient data area, red grid 155mm wide, z-fold, hole queue	300 sheets, 8 packs
2104771-001	Thermal paper 8.5in x 11in, white patient data area, red grid 155mm wide, z-fold, hole queue	150 sheets, 16 packs



ECG Recording Paper (continued) Product compatibility

Part Number	Equipment						
Part number	MAC VU360	MAC 5500	MAC 3500	MAC 2000	MAC 800	MAC 600	CASE
2104769-001		×	×				Х
2104772-001	×			×			
2104768-001		×	×				×
2104773-001						×	
2104775-001					×		
2104770-001	×	×	×				
2104771-001	×		×	×			

Reusable Bulb Electrode

GE Healthcare reusable bulb electrodes with silver-silver chloride (Ag/Agcl) plating are designed to ensure positive contact with the skin and provide maximum signal quality.



Part Number	Description	Case Quantity
2104783-001	Reusable Bulb Electrode	6

Bulb Electrode Accessories

Part Number	Description	Case Quantity
2104784-001	Clamp Limb AHA Value 4/ set	4
2104785-001	Clamp Limb IEC Value 4/ set	4

Product compatibility

MAC 600, MAC 800, MAC 2000 ECG Systems

Specification available on the website

ECG Cables and Leadwires – SEER Holter Monitor

Quality ECG cables and leadwires for your Holter devices play an important role in ensuring that they perform at the high level you have come to expect from GE Healthcare equipment.



ECG Cables and Leadwires - AHA Compliant

Part Number	Description	Case Quantity
SEER 12		
2107228-001	12-Channel Holter ECG Cable, Round, AHA, 130 cm (51.2 in)	1
2107229-001	12-Channel Holter ECG Cable, Flat AHA, 105 cm (41.3 in)	1



DIAGNOSTIC CARDIOLOGY | ECG Cables and Leadwires – SEER Holter Monitor

ECG Cables and Leadwires – AHA Compliant (continued)

Part Number	Description	Case Quantity
SEER 1000		
2104821-001	SEER [™] 1000 Leadwire set, 3 leadwire, 3 CH, 75 cm (29.5 in), AHA	1
2104822-001	SEER 1000 Leadwire set, 3 leadwire, 3 CH, 105 cm (41 in), AHA	1
2104814-001	SEER 1000 Leadwire set, 5 leadwire, 2 CH, 75 cm (29.5 in), AHA	1
2104816-001	SEER 1000 Leadwire set, 5 leadwire, 2 CH, 105 cm (41 in), AHA	1
2104815-001	SEER 1000 Leadwire set, 7 leadwire, 3 CH, 75 cm (29.5 in), AHA	1
2104817-001	SEER 1000 Leadwire set, 7 leadwire, 3 CH, 105 cm (41 in), AHA	1



ECG Cables and Leadwires - IEC Compliant

Part Number	Description	Case Quantity
SEER 12		
2104808-001	12-Channel Holter ECG Cable, Round, IEC, 130 cm	1
2104809-001	12-Channel Holter ECG Cable, Flat, IEC, 105 cm	1
SEER 1000		
2104810-001	SEER 1000 Leadwire set, 5 leadwire, 2 CH, 75 cm (29.5 in), IEC	1
2104811-001	SEER 1000 Leadwire set, 7 leadwire, 3 CH, 75 cm (29.5 in), IEC	1
2104812-001	SEER 1000 Leadwire set, 5 leadwire, 2 CH, 105 cm (41 in), IEC	1
2104813-001	SEER 1000 Leadwire set, 7 leadwire, 3 CH, 105 cm (41 in), IEC	1
2104818-001	SEER 1000 Leadwire set, 3 leadwire, 3 CH, 75 cm (29.5 in), IEC	1
2104819-001	SEER 1000 Leadwire set, 3 leadwire, 3 CH, 105 cm (41 in), IEC	1

Product compatibility

SEER 12 ECG cables are compatible with SEER 12 monitors and CardioMem[™] CM 3000-12 BT monitors

SEER 1000 ECG cables are compatible with SEER 1000 24-hour Holter monitor, SEER 1000 48-hour Holter monitor and SEER 1000 7-day Holter monitors

Specification available on the website

SEER Holter Monitor Accessories

Part Number	Description	Case Quantity
2104827-001	SEER Light Portable Pouch, shoulder straps, waist belt	1
2104803-001	12 Lead Portable Pouch, shoulder straps, waist belt	1

Product compatibility

- 2104827-001 is compatible with SEER Light and SEER Light Extend
- 2104803-001 is compatible with SEER 12, cm 3000, cm 3000-12

Specification available on website

