



Patient monitoring and telemetry supplies

PHILIPS

Disposable ECG Electrodes

Disposable Adult ECG Monitoring Electrodes

Product #	Description	Packaging	Use With:
REVISED 989803148801 CE	Small Adult, foam, solid gel snap ECG electrode, shelf life: 6 months Size: 1.6" x 1.26"	30/pack 1500/case	Adult ECG Monitoring
REVISED 989803148821 CE	Small Adult, radiolucent foam, solid gel ECG electrode, shelf life: 6 months Size: 1.4" x 1.6"	30/pack 600/case	Adult ECG Monitoring
40489E CE	Adult paper tape ECG electrode Recommended for diaphoretic patients Dimensions: 2-3/8" (59mm) diameter Shelf life: 6 months	30/pack 300/case	Adult ECG Monitoring
40493D CE	Adult foam, wet gel ECG electrode Dimensions: 2-1/8" (54mm) diameter Shelf life: 6 months	5/pack 300/case	Adult ECG Monitoring, M2735A, M2738A, M2720A
40493E CE	Adult foam, wet gel ECG electrode Dimensions: 2-1/8" (54mm) diameter Shelf life: 6 months	30/pack 300/case	Adult ECG Monitoring, M2735A, M2738A, M2720A



989803148801



989803148821



40489E



40493D / 40493E

Disposable Adult ECG Monitoring Electrodes (continued)


Product #	Description	Packaging	Use With:
REVISÉD M4612A CE	Adult foam, solid gel ECG Electrode Teardrop shape, 2" (50mm) diameter Shelf life: 6 months	5 electrodes per pouch/60 pouches per case	Adult ECG Monitoring
REVISÉD M4613A CE	Adult foam, solid gel ECG Electrode Teardrop shape, 2" (50mm) diameter Shelf life: 6 months	30 electrodes per pouch/10 pouches per case	Adult ECG Monitoring
13941E CE	Adult cloth, wet gel ECG electrode. Contains non-allergenic skin balance gel. Dimensions: 2-1/8" (54mm) diameter Shelf life: 6 months	30/pack 300/case	Adult ECG Monitoring
M2202A CE	Adult Radiotranslucent Foam Electrode Shelf life: 6 months	30/pack 300/case	Adult ECG Monitoring
13942E CE	Adult Plastic Tape ECG Electrode, Disp. Shelf life: 6 months	30/pack 300/case	Adult ECG Monitoring
M4606A CE	Skin prep paper	100 per bag	use to prep skin for electrodes



M4612A / M4613A



13941E / 13942E




M4606A

Disposable Adult Metallic ECG Lead Wire Electrodes - AAMI Color Coding

Product #	Description	Packaging	Use With:
989803156201 CE	Preattached Metallic 3-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2 inches x 1.8 inches Shelf life: 24 months unopened	1 set per pack, 100 sets/box	Adult ECG Monitoring
989803156241 CE	Preattached Metallic 5-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2 inches x 1.8 inches Shelf life: 24 months unopened	1 set per pack, 60 sets/box	Adult ECG Monitoring



989803156201



989803156241

Disposable Adult Metallic ECG Lead Wire Electrodes - IEC Color Coding

Product #	Description	Packaging	Use With:
989803156211 CE	Preattached Metallic 3-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm Shelf life: 24 months unopened	1 set per pack, 100 sets/box	Adult ECG Monitoring
989803156251 CE	Preattached Metallic 5-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm Shelf life: 24 months unopened	1 set per pack, 60 sets/box	Adult ECG Monitoring

Disposable Adult Radiolucent ECG Lead Wire Electrode - AAMI Color Coding

Product #	Description	Packaging	Use With:
989803156221 CE	Preattached radiolucent 3-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2" x 1.8" Shelf life: 24 months unopened	1 set per pack, 100 sets/box	Adult ECG Monitoring
989803156261 CE	Preattached radiolucent 5-leadwire (36") electrode set, non-shielded, solid gel. Dimensions: 1.2" x 1.8" Shelf life: 24 months unopened	1 set per pack, 60 sets/box	Adult ECG Monitoring

Disposable Adult Radiolucent ECG Lead Wire Electrode - IEC Color Coding

Product #	Description	Packaging	Use With:
989803156231 CE	Preattached radiolucent 3-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm Shelf life: 24 months unopened	1 set per pack, 100 sets/box	Adult ECG Monitoring
989803156271 CE	Preattached radiolucent 5-leadwire (1m) electrode set, non-shielded, solid gel. Dimensions: 30mm x 45mm Shelf life: 24 months unopened	1 set per pack, 60 sets/box	Adult ECG Monitoring

Disposable Pediatric ECG Monitoring Electrodes

Product #	Description	Packaging	Use With:
13950B CE	Disposable cloth electrode with KCl foam gel retaining ring, snap. (Individually wrapped) Diameter: 1.25" (32mm) 6 month shelf life	30 per box 10 boxes/case	Pediatric/Neonatal Monitors
13951C CE	Disposable, repositionable solid gel electrode, snap. (3 electrodes mounted on plastic card) Dimensions: 7/8" x 1-1/4" (22mm x 33mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors
13955C CE	Disposable, repositionable solid gel electrode, snap. (3 electrodes mounted on plastic card) Dimensions: 7/8" x 7/8" (22mm x 22mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors



13950B



13951C



13955C

Disposable Infant/Pediatric Metallic Lead Wire Electrodes: AAMI Color Coding

Product #	Description	Packaging	Use With:
13952A	Preattached metallic (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 1-1/4" (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors
13953A CE	Preattached metallic (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 7/8" (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors



13952A

Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: AAMI Color Coding

Product #	Description	Packaging	Use With:
13952B CE	Preattached radiolucent (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 1-1/4" (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors
13953B CE	Preattached radiolucent (18") 3-leadwire electrode, non-shielded, solid gel. Dimension: 7/8" x 7/8" (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors



13952B

Disposable Infant/Pediatric Metallic Lead Wire Electrodes: IEC Color Coding

Product #	Description	Packaging	Use With:
13952D CE	Preattached metallic (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 33mm, (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors
13953D CE	Preattached metallic (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 22mm, (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors



13953D

Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: IEC Color Coding

Product #	Description	Packaging	Use With:
13952E CE	Preattached radiolucent (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 33mm, (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors
13953E CE	Preattached radiolucent (46cm) 3-leadwire electrode, non-shielded, solid gel. Dimension: 22mm x 22mm, (socket size: 1.5mm) 6 month shelf life	3 per pack, 300/case	Pediatric/Neonatal Monitors

EEG Electrode (electroencephalogram)

Product #	Description	Packaging	Use With:
M1935A CE	Disposable EEG snap electrode	25/pack 100/case	Adult EEG Monitoring



M1935A