

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



1 A



1 B



2



3

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Lead Sets AAMI Color Coding (continued)

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*		
(Also for use with Efficia CM monitors)									
4	5 limb	M1968A	ICU, grabber	Reusable	Limb 1m (3.3ft) chest 1,6m (5.3 ft)	1 each	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)		
		M1644A	ICU, snap						
		989803192181	ICU, DIN to Tab						
5		M1973A	OR, grabber		Limb 70 cm (28") chest 1,3m (2.3ft)		20 each	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)	
6	5 limb	M1647A	Miniclip grabber, unshielded						
7		989803173131	ICU, grabber	Single patient use					
	5 chest	M1976A	ICU, grabber	Reusable	1m (3.3ft)	1 each			Efficia CM12 (863303), CM120 (863302), CM150 (863304)
		M1602A	ICU, snap Reusable Length: 1m (3.3ft)						
		989803192201	ICU, DIN to Tab				5+5 Trunk cable M1949A, M1663A		
8		M1979A	OR, grabber				Efficia CM12 (863303), CM120 (863302), CM150 (863304)		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Trunk cables

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:* (Also for use with Efficia CM monitors)
1	A	M1669A CE	ECG Trunk Cable, 12 pin connector, AAMI and IEC labels included	Reusable	2,7 m (9.9ft)	1 cable	Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 863088, 863317; Efficia Monitors 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322; Defibrillators M3534A, M3535A, M3536A (Heart Start MRx), M4735A (HeartStart XL), 861290 HeartStartXL+, DFM 100, Efficia CM12 (863303), CM120 (863302), CM150 (863304)
		989803170171 (Orange) CE	OR ECG Trunk Cable, 12-pin connector. AAMI and IEC labels included. For use in operating room			1 cable	Multi-Measurement Modules: M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863064, 863065, 863066, 863068, 863317 Efficia Monitors 863300, 863301, 863302, 863303, 863304 Defibrillator: DFM 100, Efficia CM12 (863303), CM120 (863302), CM150 (863304)
2		M1668A CE	ECG Trunk Cable, 12 pin connector, AAMI and IEC labels included			1 cable	Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863064, 863065, 863066, 863068, 863086, 863087, 863088, 863317; Efficia Monitors 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322; Defibrillators: M3534A, M3535A, M3536A (Heart Start MRx), M4735A (HeartStart XL), 861290 HeartStartXL+, DFM 100, Xper Flex Cardio: 860338, 850335, Efficia CM12 (863303), CM120 (863302), CM150 (863304)
		989803170181 CE	OR ECG Trunk Cable, 12-pin connector. AAMI and IEC labels included. For use in operating room			1 cable	Multi-Measurement Modules: M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863065, 863066, 863068, 863317 Efficia Monitors 863300, 863301, 863302, 863303, 863304, 863320, 863321, 863322 Defibrillator: DFM 100, Efficia CM12 (863303), CM120 (863302), CM150 (863304)
3	6	M1667A CE	ECG Trunk Cable, 12 pin connector, AAMI and IEC labels included			1 cable	Multi Measurement Server: M3000A, M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18, MP2: M8102A Opt C12, MP5: M8105A Opt C12, Defibrillator HeartStartXL+ 861290

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



1



2



3



4

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Trunk cables (continued)

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:* (Also for use with Efficia CM monitors)
4		M1663A CE	ECG Trunk Cable, 12 pin connector, AAMI and IEC labels included	Reusable	2m (6.5ft)	1 cable	Multi-Measurement Modules: M3000A, M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18; Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG Modules: M1001B, M1002B, Efficia Monitors: 863302, 862202, 863304, 863321, 863322 Defibrillators: M3535A (HeartStart MRx), M4735 (HeartStart XL), Defibrillator Heart Start XL+ 861290, Efficia CM12 (863303), CM120 (863302), CM150 (863304)
	10 (5+5)	M1949A CE	ECG Trunk Cable, 12 pin connector, AAMI and IEC labels included		2,7 m (9.9ft)	1 cable	Multi-Measurement Modules: M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18; Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG, Efficia Monitors: 863302, 863303, 863304, 863321, 863322 Defibrillators: M3535A, M3536A (HeartStart MRx), M4735A (HeartStart XL), Defibrillator Heart Start XL+ 861290 XPer Flex Cardio 860335, 860338, Efficia CM12 (863303), CM120 (863302), CM150 (863304)
5	10 (6+4)	M1665A CE	ECG Trunk Cable, 12 pin connector, AAMI and IEC labels included			1 cable	Multi-Measurement Modules: M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18; Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family




ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Trunk cable adapter for non Philips monitors (connects non Philips trunk cables to IntelliVue Lead sets)


	# of Leads	Product #	Description	Reusable/ Single patient use	Packaging	Use with:*
1	5 (can also be used for 3 lead ECG)	989803174501 CE	ECG Trunk Cable Adapter- GE, AAMI/IEC	Reusable	1 adapter	Connects to GE 5-lead monitoring trunk cables with MultiLink-style connectors (such as product numbers 41-2931-XXX, 411910-0011)
2		989803174511 CE	ECG Trunk Cable Adapter - Spacelabs, AAMI/IEC			Connects to Spacelabs 5-lead monitoring trunk cables with TruLink-style connectors (such as Spacelabs 700-0008-06 shielded trunk)
3		989803174521 CE	ECG Trunk Cable Adapter - Datascope AAMI			Connects to Datascope 5-lead monitoring AAMI trunk cables with DIN-style connector (0012-00-0620-02)
		989803174531 CE	ECG Trunk Cable Adapter - Mindray, AAMI/IEC			Connects to Mindray 5-lead monitoring trunk cables (such as 0012-00-1255-02)
		989803177881 CE	ECG Trunk Cable Adapter - Datascope IEC			Connects to Datascope 5-lead monitoring IEC trunk cables with DIN-style connector
		989803182171 CE	ECG Trunk Cable Adapter - Draeger/Fukuda Denshi, AAMI/IEC			Connect Philips 3- and 5-lead leadsets to Draeger monitoring trunk cables (such as 3368391, 5950196, 5191221)
4		989803182181 CE	ECG Trunk Cable Adapter - Nihon Kohden, AAMI /IEC			Connect any Philips 3- and 5-lead leadsets to Nihon Kohden monitoring trunk cables (such as JC-906P)

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



4



5



6



7



8

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Lead Sets AAMI Color Coding


	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
							(Also for use with Efficia CM monitors)
1	3	M1671A	ICU, grabber	Reusable	1m (3.3ft)	1 each	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
		M1673A	ICU, snap				
2		M1675A	OR, grabber				
		M1622A	Miniclip, grabber, unshielded		45cm (18")		
3		M1624A	Miniclip, grabber, unshielded		70cm (28")		
		989803173121	ICU grabber	Single Patient Use		20 each	
	4 chest	M1532A	ICU, grabber	Reusable	1m (3.3ft)	1 each	6+4 trunk cable M1665A
		M1537A	ICU, snap				
		989803192241	ICU, DIN to Tab				
		M1557A	OR, grabber				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



9



10


ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Lead Sets AAMI Color Coding (continued)

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
							(Also for use with Efficia CM monitors)
	5 chest	989803192141	ICU, grabber	Single patient use	1m (3.3ft)	20 each	5+5 Trunk cable M1949A,M1663A
9	6	M1680A	ICU, grabber	Reusable	Limb 1m (3.3ft) chest 1,6m (5.3 ft)	1 each	6 lead Trunk cable M1667A or 6+4 Trunk cable M1665A
		M1682A	ICU, snap				
		989803192221	ICU, DIN to Tab				
10		M1684A	OR, grabber				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family





ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Lead Sets IEC Color Coding

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:* (Also for use with Efficia CM monitors)
1	3	M1672A CE	ICU, grabber	Reusable	1m (3.3ft)	1 each	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
		M1674A CE	ICU, snap				
2		M1678A CE	OR, grabber		70cm (28")		
3		M1626A CE	Miniclip, grabber, unshielded				
		989803174201 CE	ICU, grabber	Single patient use	1m (3.3ft)	20 each	6+4 trunk cable M1665A
4	4 chest	M1533A CE	ICU, grabber	Reusable			
		M1538A CE	ICU, snap				
		989803192251 CE	ICU, DIN to Tab				
5		M1558A CE	OR, grabber				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Lead Sets IEC Color Coding (continued)


	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:* (Also for use with Efficia CM monitors)
4	5 limb	M1971A CE	ICU, grabber	Reusable	Limb 1m (3.3ft) chest 1,6m (5.3 ft)	1 each	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
		M1645A CE	ICU, snap				
		989803192191 CE	ICU, DIN to Tab				5 lead Trunk cable M1668A, OR trunk cable 989803170181 or 5+5 trunk cable M1949A, M1663A
		M1974A CE	OR, grabber				
5		M1648A CE	Miniclip grabber, unshielded		Limb 70 cm (28in) chest 1,3m (2.3ft)		Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
		989803174211 CE	ICU, grabber	Single patient use		20 each	5 lead Trunk cable M1668A, OR trunk cable 989803170181 or 5+5 trunk cable M1949A, M1663A
6	5 chest	M1978A CE	ICU, grabber	Reusable	1 m (3.3ft)	1 each	5+5 Trunk cable M1949A, M1663A
		M1604A CE	ICU, snap				
		989803192211 CE	ICU, DIN to Tab				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



7



8



9

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Lead Sets IEC Color Coding (continued)

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
							(Also for use with Efficia CM monitors)
7	5 chest	M1984A CE	OR, grabber	Reusable	1m (3.3ft)	1 each	Efficia CM12 (863303), CM120 (863302), CM150 (863304)
		989803192151 CE	ICU, grabber	Single patient use		20 each	5+5 Trunk cable M1949A, M1663A
8	6	M1681A CE	ICU, grabber	Reusable	Limb 1m (3.3ft) chest 1,6m (5.3 ft)	1 each	6 lead Trunk cable M1667A or 6+4 Trunk cable M1665A
		M1683A CE	ICU, snap				
		989803192231 CE	ICU, DIN to Tab				
9		M1685A CE	OR, grabber				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

Radiolucent ECG Lead Sets – AAMI/IEC Color Coding

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	5	M1649A CE	5 unshielded radiolucent grabber leads with Safety DIN connectors, AAMI/IEC color markers included	Reusable	91cm (36")	5 individual leads	5 lead trunk cable M1668A, OR trunk cable 989803170181 or 5+5 trunk cable M1949A, M1663A

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.


ECG Cables and Leads – Bedside ECG IntelliVue Cable Family



1



2



3



4




5



6



7



8

ECG Cables and Leads – Bedside ECG IntelliVue Cable Family

IntelliVue Cable family - Leads Accessories

Product #	Description	Packaging	Use with:*
989803193821	DIN to Tab Adapters	10 adapters	All DIN to Tab lead sets
1 989803148841 C€	IntelliVue Cable Management Kit	5 organizers	All cables
2 989803148861 C€	IntelliVue Trunk Cable Covers (5 leads)	10 covers	M1663A or M1949 5+5 trunk cables
3 M1636A	Organizer for unshielded 3-lead sets	10 organizers	All unshielded 3 lead set
4 M1638A	Organizer for unshielded 5-lead sets		All unshielded 5 lead sets
M1503A	Organizer for shielded 3-lead sets		All shielded 3-lead sets
5 M1664A	Organizer for shielded 4-lead sets		All shielded 4-lead sets
M1504A	Organizer for shielded 5-lead sets		All shielded 5-lead sets
6 M1679A	Organizer for shielded 6-lead sets		All shielded 6-lead sets
7 M1509A	Bedsheet clip for trunk cables	3 clips	All trunk cables
8 M2281A	Patient cable organizer (hook)	8 hooks	All cables

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Sure Signs one piece trunk/leadset cable



Sure Signs one piece trunk/leadset cable

One Piece Trunk/Leadset cables AAMI and IEC Color Coding

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	3	989803143171	Grabber IEC C€	Reusable	2,5m (8.2ft)	1 cable	SureSigns monitors 863063 86064,86065, 863066, 86068, 863077, 863085, 863086, 863087, 863088, 863300, 863301, 863302, 863303, 863304, 863317, 863320, 863321, 863322, Measurement Server M3001A M3001AL, MP2 M8102A, MP5 M8105A, X2 M3002A, Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
2		989803143181	Grabber AAMI				
3	5	989803143191	Grabber, IEC C€				SureSigns monitors 86064,86065, 863066, 86068, 863086, 863087, 863088, 863300, 863301, 863302, 860303, 863304, 863317, 863320, 863321, 863322, Measurement Server M3001A M3001AL, MP2 M8102A, MP5 M8105A, X2 M3002A, Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
	6	989803143201	Grabber, AAMI				

ECG Cables and Leads- Hemodynamic Monitoring

Leads for Hemodynamic Monitoring

	# of Leads	Product #	Description	Reusable/ Single patient use	Length		Use with:*
4	4	989803174231	Grabber AAMI for 16 Lead Hemo C€	Reusable	1,6m (3.3ft)	1 each	Xper Flex Cardio physiomonitring system FC2010 and FC2020. Use with M1949A 5+5 ECG trunk cable
5		989803174241	Grabber IEC for 16 Lead Hemo C€				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Efficia ECG Cables



1



2



3



4



5



6

Efficia ECG Cables

One Piece Sets

	Product #	Description	IEC Part number	Length	Use with:*
1	989803160731	3-Lead grabber	989803160741	2.5 m	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
2	989803160751	3-Lead snap	989803160761		
3	989803160771	5-Lead grabber	989803160781		
4	989803160791	5-Lead snap	989803160801		

3-Lead Sets

	Product #	Description	IEC Part number	Use with:*
5	989803160651	Grabber	989803160661	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
6	989803160671	Snap	989803160681	

5-Lead Sets

	Product #	Description	IEC Part number	Use with:*
	989803160691	Grabber	989803160701	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
	989803160711	Snap	989803160721	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue MX 40 Wearable Patient Monitoring Supplies – Patient Cables




IntelliVue MX 40 Wearable Patient Monitoring Supplies – Patient Cables

Patient Cables AAMI color coding

	# of Leads	Product #	Parameter	Lead Termination	Reusable/ Single patient use	Length	Packaging	Use with:*
1	3	989803171801	ECG	Grabbers	Reusable	85cm (34")	1 each	MX40 1.4 GHz - SH Patient Worn Monitor 865350 (US only), MX40 2.4 GHz - SH Patient Worn Monitor 865351
2		989803171811	ECG+ SpO ₂					
3	5	989803171821	ECG	Snaps				
4		989803171831	ECG	Grabbers				
5		989803171841	ECG+ SpO ₂	Snaps				
6		989803171851	ECG+ SpO ₂	Grabbers				
7		989803172031	ECG					
8		989803172051	ECG+ SpO ₂					
9	6	989803171861	ECG		Reusable			
10		989803171871	ECG+ SpO ₂					

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue MX 40 Wearable Patient Monitoring Supplies – Patient Cables



IntelliVue MX 40 Wearable Patient Monitoring Supplies – Patient Cables

Patient Cables IEC Color coded

# of Leads	Product #	Parameter	Lead Termination	Reusable/ Single patient use	Length	Packaging	Use with:*
3	989803171901	ECG CE	Grabbers	Reusable	85cm (34")	1 each	MX40 1.4 GHz - SH Patient Worn Monitor 865350 (US only), MX40 2.4 GHz - SH Patient Worn Monitor 865351
	989803171911	ECG+ SpO ₂ CE					
5	989803171931	ECG CE		Single patient use			
	989803171951	ECG+ SpO ₂ CE					
	989803172131	ECG CE					
	989803172151	ECG+ SpO ₂ CE					
6	989803171961	ECG CE		Reusable		1 each	
	989803171971	ECG+ SpO ₂ CE					

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue MX 40 Wearable Patient Monitoring Supplies – Patient Cables



IntelliVue MX 40 Wearable Patient Monitoring Supplies – Patient Cables

Miscellaneous MX40 Cable/Supplies

Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
989803172221 CE	Patient Monitor ECG Trunk Cable AAMI/IEC	Reusable	2,7m (9 ft)	1 each	MX40 1.4 GHz – SH Patient Worn Monitor 865350 (US only), MX40 2.4 GHz – SH Patient Worn Monitor 865351
989803172241 CE	Extender cable		1m (3.2 ft)		
989803172211 CE	IntelliVue Adapter Cable Bedside/Tele ECG+SpO ₂		30cm (1 ft)		
989803176491	Protective Cap for reusable Leads			10 caps	All MX40 reusable leads
989803176501	Protective Cap for adapter cable 989803172211				SpO ₂ socket of adapter cable 989803172211

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue Telemetry ECG Cables and Leads – Patient Cables





IntelliVue Telemetry ECG Cables and Leads – Patient Cables

IntelliVue Telemetry Leadsets AAMI

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	3	989803151971	Grabber and shower shield	Reusable	85cm (34")	1 each	Telemetry 862231 None U.S. TRx4851A 2.4 Ghz IntelliVue Tele TRX, 862439: U.S. only TRx4841A 1.4. Ghz IntelliVue Tele TRX
2		989803151991	Snap and shower shield				
3		989803173141	Grabber and shower shield	Single patient use		20 each	
4	5	989803152051	Grabber and shower shield	Reusable		1 each	
5		989803152071	Snap and shower shield				
6		989803173151	Grabber and shower shield	Single patient use		20 each	
7	6	989803152131	Grabber and shower shield	Reusable		1 each	
8		989803152151	Snap and shower shield				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue Telemetry ECG Cables and Leads - Patient Cables



IntelliVue Telemetry ECG Cables and Leads – Patient Cables

IntelliVue Telemetry Leadsets IEC

	# of Leads	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	3	989803151981	Grabber and shower shield CE	Reusable	85cm (34")	1 each	Telemetry 862231 None U.S. TRx4851A 2.4 Ghz IntelliVue Tele TRX, 862439: U.S. only TRx4841A 1.4. Ghz IntelliVue Tele TRX
2		989803152001	Snap and shower shield CE				
3	5	989803152061	Grabber and shower shield CE				
4		989803152081	Snap and shower shield CE				
5	6	989803152141	Grabber and shower shield CE				
6		989803152161	Snap and shower shield CE				

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue Telemetry ECG Cables and Leads - Patient Cables



IntelliVue Telemetry ECG Cables and Leads – Patient Cables

IntelliVue Telemetry Shields and Alignment Guides

	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	989803143481 C€	Tether cable connects PWD to TeleMon tether cable or MP5 tether cable	Reusable	51 cm (20")	1 cable	862231 TRX4851A 2.4 GHz IntelliVue Tele TRX, 862439 1.4 GHZ IntelliVue Tele TRX, M8105A IntelliVue MP5
2	989803143491 C€	TeleMon tether cable connects the M2636C to the PWD tether cable		2 m (79")		862231 TRX4851A 2.4 GHz IntelliVue Tele TRX, 862439 1.4 GHZ IntelliVue Tele TRX
3	989803146911 C€	MP5 tether cable connects the MP5 to the PWD tether cable		862231 TRX4851A 2.4 GHz IntelliVue Tele TRX, 862439 1.4 GHZ IntelliVue Tele TRX, M8105A IntelliVue MP5		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue Telemetry ECG Cables and Leads - Patient Cables



IntelliVue Telemetry ECG Cables and Leads – Patient Cables

IntelliVue Telemetry Shields and Alignment Guides

	Product #	Description	Packaging	Use with:*
1	989803153031	3 lead detachable shield replacements	10 shields	All IntelliVue Telemetry lead sets
2	989803153041	5 lead detachable shield replacements		
3	989803153051	6 lead detachable shield replacements		
	989803140401	Single ECG alignment guide protective insert, untethered		IntelliVue Telemetry: 862231:TRx4851A Transceiver; 862439: TRx4841A Transceiver; 862108: M2601B Transmitter
	989803140411	Single ECG alignment guide protective insert, tethered		
	989803140421	Double ECG alignment guide protective insert, untethered		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.


EEG Supplies (Electroencephalogram)




1



2



3




4



5



6



7



8

EEG Supplies (Electroencephalogram)

EEG Supplies (Electroencephalogram)

	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	M2268A C€	5 lead trunk cable	Reusable	2,7m (8.8ft)	1 cable	M1027A EEG Module
2	989803180541 C€	4 channel trunk cable				
3	M1931A C€	5 lead EEG adult cup electrode set		80cm (31.5")	5 leads	
4	M1932A C€	5 lead EEG pediatric cup electrode set				
5	M1934A C€	5 lead EEG miniclip lead set				
6	989803180561 C€	9 lead EEG adult cup electrode set			9 leads	Cup electrodes
	989803180531 C€	9 lead EEG miniclip lead set				
7	989803180521 C€	9 lead EEG infant and pediatric cup electrode set				
	M1937A C€	EEG conductive paste			1 tube	
8	M1935A C€	EEG Snap electrodes			100 each	Miniclip EEG lead sets

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

IntelliVue NMT (Neuro Muscular Transmission) Patient Cable



1

IntelliVue NMT (Neuro Muscular Transmission) Patient Cable

IntelliVue NMT (Neuro Muscular Transmission) Patient Cable

	Product #	Description	Reusable/ Single patient use	Length	Packaging	Use with:*
1	989803174581 C€	Patient cable with 12 pin instrument connector	Reusable	3m (9.8ft)	1 each	IntelliVue Neuro muscular Transmission Module 865383

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes



1



2



3



4




5



6



7



8

Disposable ECG Electrodes

Disposable Adult ECG Monitoring Electrodes


Product #	Description	Packaging	Use with:*
1 989803148801	Small Adult, foam, solid gel snap ECG electrode, shelf life: 6 months Size: 1.6" x 1.26"	30/pack 1500/case	Adult AEG monitoring Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304)
2 989803148821	Small Adult, radiolucent foam, solid gel ECG electrode, shelf life: 6 months Size: 1.4" x 1.6"	30/pack 600/case	
3 40489E	Adult paper tape ECG electrode, shelf life: 6 months Recommended for diaphoretic patients Dimensions: 2 3/8" (59mm) diameter	30/pack 300/case	Adult ECG Monitoring
40493D	Adult foam, wet gel ECG electrode Dimensions: 2 1/8" (54mm) diameter Shelf life: 6 months	30/pack 60 packs/case	Efficia CM10 (863301), CM12 (863303), CM100 (863300), CM120 (863302), CM150 (863304) Adult ECG Monitoring: M2735A, M2738A, M2720A; Cableless Toco+ MP Transducer: 866075; Cableless ECG/IUP Transducer 866077 to Use with column for 40493D and 40493E
4 40493E	Adult foam, wet gel ECG electrode Dimensions: 2 1/8" (54mm) diameter Shelf life: 6 months		
M4612A	Adult foam, solid gel ECG electrode Teardrop shape, 2" (50mm) diameter Shelf life: 6 months		
5 M4613A	Adult foam, solid gel ECG electrode Teardrop shape, 2" (50mm) diameter Shelf life: 6 months		
6 13941E	Adult cloth, wet gel ECG electrode Contains non-allergenic skin balance gel Dimensions: 54mm (2 1/8") diameter Shelf life: 6 months		
M2202A	Adult radiotranslucent foam electrode Shelf life: 6 months		
7 13942E	Adult plastic tape ECG electrode, Disposable Shelf life: 6 months		
8 M4606A	Skin prep paper	100 per bag	Use to prep skin for electrodes

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes



1



2

Disposable ECG Electrodes

Disposable Adult Metallic ECG Lead Wire Electrodes - AAMI Color Coding

Product #	Description	Packaging	Use with:*
1 989803156201	Preattached metallic 3-leadwire (36") electrode set, non-shielded, solid gel Dimensions: 1.2" x 1.8" Shelf life: 24 months, unopened	1 set/pack 100 sets/box	Adult ECG monitoring
2 989803156241	Preattached metallic 5-leadwire (36") electrode set, non-shielded, solid gel Dimensions: 1.2" x 1.8" Shelf life: 24 months, unopened	1 set/pack 60 sets/box	

Disposable Adult Metallic ECG Lead Wire Electrodes - IEC Color Coding

Product #	Description	Packaging	Use with:*
989803156211	Preattached metallic 3-leadwire (1m) electrode set, non-shielded, solid gel Dimensions: 30mm x 45mm, Shelf life: 24 months, unopened	1 set/pack 100 sets/box	Adult ECG monitoring
989803156251	Preattached metallic 5-leadwire (1m) electrode set, non-shielded, solid gel Dimensions: 30mm x 45mm Shelf life: 24 months, unopened	1 set/pack 60 sets/box	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.



Disposable ECG Electrodes

Disposable Adult Radiolucent ECG Lead Wire Electrode - IEC Color Coding

Product #	Description	Packaging	Use with:*
989803156231	Preattached radiolucent 3-leadwire (1m) electrode set, non-shielded, solid gel Dimensions: 30mm x 45mm Shelf life: 24 months, unopened	1 set/pack 100 sets/box	Adult ECG monitoring
989803156271	Preattached radiolucent 5-leadwire (1m) electrode set, non-shielded, solid gel Dimensions: 30mm x 45mm Shelf life: 24 months, unopened	1 set/pack 60 sets/box	

Disposable Adult Radiolucent ECG Lead Wire Electrode - AAMI Color Coding


Product #	Description	Packaging	Use with:*
989803156221	Preattached radiolucent 3-leadwire (36in) electrode set, non-shielded, solid gel Dimensions: 1.2" x 1.8" Shelf life: 24 months, unopened	1 set/pack 100 sets/box	Adult ECG monitoring
989803156261	Preattached radiolucent 5-leadwire (36in) electrode set, non-shielded, solid gel Dimensions: 1.2" x 1.8" Shelf life: 24 months, unopened	1 set/pack 60 sets/box	

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes



1



2



3


Disposable ECG Electrodes

Disposable Pediatric ECG Monitoring Electrodes

Product #	Description	Packaging	Use with:*
1 13950B	Disposable cloth electrode with KCl foam gel retaining ring, snap (Individually wrapped) Diameter: 1.25" (32mm) Shelf life: 6 months	30/box 10 boxes/case	Pediatric/neonatal monitors
2 13951C	Disposable, repositionable solid gel electrode, snap Dimensions: 7/8" x 1 1/4" (22mm x 33mm) 3 electrodes mounted on plastic card Shelf life: 6 months	3/pack 100 packs/case	
3 13955C	Disposable, repositionable solid gel electrode, snap Dimensions: 7/8" x 7/8" (22mm x 22mm) 3 electrodes mounted on plastic card Shelf life: 6 months		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes





Disposable ECG Electrodes

Disposable Infant/Pediatric Metallic Lead Wire Electrodes: AAMI Color Coding

Product #	Description	Packaging	Use with:*
1 13952A	Pre-attached metallic 46 cm (18") leadwire electrode, nonshielded, AAMI solid gel. Dimension: 22mm x 33mm (7/8" x 1 1/4"); socket size 1.5mm. Shelf life: 6 months.	3 electrodes/pack; 300/case	IntelliVue MP2 - M8102A, IntelliVue MP5 -M8105A, MMS - M3000A, M3001A, M3002A, M1001A, M1001B, M1002A, M1002B, 78260A, 78261A, 78801A, 78801B, 78831A, 78832A, 78833A, 78833B, 78834A, 78834C, Monitors (Neonatal/Pediatric), 863063, 863064, 863065, 863066, 863068, 863077
13953A	Pre-attached metallic 45cm (18") leadwire electrode, nonshielded, AAMI solid gel. Dimension: 22mm x 22mm (7/8" x 7/8"); socket size 1.5mm; Shelf life: 6 months	3 electrodes/pack; 300/case	
2 989803183101	NeoLead, metallic, pre-attached leadwire (24") pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. AAMI, Standard DIN connector. 1" round.	3 electrodes/pack 50 packs/box	ECG modules: M1001B, M1002B; Multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors: M8105A, M8102A, 863064, 863065, 863066, 863068, 860338, 860335; Requires ECG trunk cable: M1668A or 989803170181
3 989803183081	Mini NeoLead, metallic pre-attached leadwire (24") pediatric/neonatal ECG electrodes. Hydrogel/Hydrocolloid adhesive. AAMI, Standard DIN connector, 3/4" round		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes





Disposable ECG Electrodes

Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: AAMI Color Coding

Product #	Description	Packaging	Use with:*
1 13952B	Preattached radiolucent (18") 3-leadwire electrode, non-shielded, solid gel Dimension: 7/8" x 1-1/4" (socket size: 1.5mm) Shelf life: 6 months	3 electrodes/pack 300/case	MP2 - M8102A, MP5 -M8105A, MMS - M3000A, M3001A, M3002A, ECG modules -M1001A, M1001B, M1002A, M1002B, 78260A, 78261A, 78801A, 78801B, 78831A, 78832A, 78833A, 78833B, 78834A, 78834C, Monitors (Neonatal/Pediatric) 863063, 863064, 863065, 863066, 863068, 863077; Requires ECG trunk cable: M1668A or 989803170181
13953B	Preattached radiolucent (18") 3-leadwire electrode, non-shielded, solid gel Dimension: 7/8" x 7/8" (socket size: 1.5mm) Shelf life: 6 months		
2 989803183121	NeoLead, radiolucent pre-attached leadwire (24") pediatric/neonatal ECG electrodes. AAMI; Standard DIN connector 1" round.	3 electrodes per pouch 50 per box	ECG modules: M1001B, M1002B; multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors MP5 -M8105A, MP2-M8102A, 863064, 863065, 833066, 863068, 860338, 860335; Requires ECG trunk cable: M1668A or 989803170181
3 989803183141	Micro NeoLead, radiolucent pre-attached leadwire (24") pediatric/neonatal ECG electrodes. AAMI, standard DIN connector 1" x 3/8".		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes





Disposable ECG Electrodes

Disposable Infant/Pediatric Metallic Lead Wire Electrodes: IEC Color Coding

Product #	Description	Packaging	Use with:*
1 13952D	Preattached metallic (46cm) 3-leadwire electrode, non-shielded, solid gel Dimension: 22mm x 33mm, (socket size: 1.5mm) Shelf life: 6 months	3 electrodes /pack; 300/case	MP2 - M8102A, MP5 -M8105A, MMS - M3000A, M3001A, M3002A, ECG modules -M1001A, M1001B, M1002A, M1002B, 78260A, 78261A, 78801A, 78801B, 78831A, 78832A, 78833A, 78833B, 78834A, 78834C, Monitors (Neonatal/Pediatric) 863063, 863064, 863065, 863066, 863068, 863077; Requires ECG trunk cable: M1668A or 989803170181
13953D	Preattached metallic (46cm) 3-leadwire electrode, non-shielded, solid gel Dimension: 22mm x 22mm, (socket size: 1.5mm) Shelf life: 6 months		
2 989803183111	NeoLead, metallic pre-attached leadwire (60 cm) pediatric/neonatal ECG electrodes; IEC; standard DIN connector, 25 mm round	3 electrodes /pack; 50 per box	ECG modules: M1001B, M1002B; multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors MP5 -M8105A, MP2-M8102A, 863064, 863065, 833066, 863068, 860338, 860335; Requires ECG trunk cable: M1668A or 989803170181
3 989803183091	Mini NeoLead, metallic pre-attached leadwire (60 cm) pediatric/neonatal ECG electrodes; AAMI; standard DIN connector. 19 mm round		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.

Disposable ECG Electrodes



Disposable ECG Electrodes

Disposable Infant/Pediatric Radiolucent Lead Wire Electrodes: IEC Color Coding

Product #	Description	Packaging	Use with:*
1 13952E	Preattached radiolucent (46cm) 3-leadwire electrode, non-shielded, solid gel Dimension: 22mm x 33mm, (socket size: 1.5mm), Shelf life: 6 months	3 electrodes/pack; 300/case	MP2 - M8102A, MP5 -M8105A, MMS - M3000A, M3001A, M3002A, ECG modules -M1001A, M1001B, M1002A, M1002B, 78260A, 78261A, 78801A, 78801B, 78831A, 78832A, 78833A, 78833B, 78834A, 78834C, Monitors (Neonatal/Pediatric)- 863063, 863064, 863065, 863066, 863068, 863077; Requires ECG trunk cable: M1668A or 989803170181
13953E	Preattached radiolucent (46cm) 3-leadwire electrode, non-shielded, solid gel Dimension: 22mm x 22mm, (socket size: 1.5mm) Shelf life: 6 months		
2 989803183131	NNeoLead, radiolucent pre-attached leadwire (60 cm) pediatric/neonatal ECG electrodes; IEC, standard DIN controller, 25 mm round	3/pack 50 packs/box	ECG modules: M1001B, M1002B; multi-measurement servers: M3000A, M3001A, M3002A; Patient monitors MP5 -M8105A, MP2-M8102A, 863064, 863065, 833066, 863068, 860338, 860335; Requires ECG trunk cable: M1668A or 989803170181
3 989803183151	Micro NeoLead, radiolucent pre-attached leadwire (60 cm) pediatric/neonatal ECG electrodes; IEC, standard DIN connector, 25mm x 9.5 mm		

* Please reference the [MySupplies](#) or the [e-store](#) website for more details on hardware compatibility.