



Fluid Management

CONTROL | CONVENIENCE | CONFIDENCE



Fluid Management

Merit Medical offers an integrated suite of fluid management products designed to effectively manage patient fluid levels, contrast media, and waste in diagnostic and interventional procedures. From tubing to transducers, fluid spikes to drainage bags, Merit Medical provides high quality angiographic accessories you can rely on to deliver safe and efficient fluid management.

MANIFOLDS & STOPCOCKS

Merit® Manifolds

- Up to 5 ports
- Right & left-handed orientation
- 200 & 500 PSI (14 & 34 Bar) pressure ratings
- 6" extension tubing & rotating Luer
- Block body & half body styles
- TRAM® and TRAM-P™ with integrated transducers
- Devos™ manifold with integrated check relief valve

Marquis® Series Stopcock

- Large, easy-grip handle for effortless grip and turn
- Pressure ratings: 200, 500, 1050 & 1200 PSI (14, 34, 72 & 83 Bar)
- Luer options include rotator, swivel nut, fixed and slip
- 1-way, 3-way & 4-way styles

CONTRAST MANAGEMENT

Miser® Management System

- Ball and burette design to prevent air from entering fluid lines
- Dual one-way check valves' microbial barriers block cross contamination
- Automatic refill during each contrast withdrawal for quick and easy preparation

CT Contrast Transfer Sets


- Large and small-bore configurations
- Multiple spike options
- Halkey-Roberts® swabbable valve

WASTE MANAGEMENT

Merit Disposal Depot™ Waste Bag

This closed system fully contains fluid waste for easy disposal. The bags come in many configurations for connecting directly to a manifold or for convenient placement on the back table.

- Large and small-bore tubing available
- Fully customizable
- With and without reflux valves



CONTROL

with a broad selection of fluid management products to match your preferences.

CONFIDENCE

with fluid and contrast management solutions that may help limit exposure for both patients and health care workers.

CONVENIENCE

having all the items you need delivered together, streamlining clinical set up time.



Merit®
Pressure Monitoring
Tubing

Meritrans®
Pressure Transducer

Merit Disposal Depot™ Waste Bag
Closed Waste System

TUBING

Merit® Tubing

- High pressure contrast injection tubing rated at 500, 900 & 1200 PSI (34, 62 & 83 Bar)
- Pressure monitoring tubing
- Fluid administration sets

TRANSDUCERS & ACCESSORIES

Meritrans® Disposable Pressure Transducer

With reliable precision, customized configurations and a clear fluid path, the Meritrans can be tailored to your clinical needs.

- Hand-held & pole-mounted set-ups
- Cable lengths: 12" & 48" (30.5 & 122 cm)
- Assembled with or without bonded stopcocks
- Mounting plate and pole clamp accessories



SYRINGES

CCS™, Inject8™ & Inject10ⁿ™

Coronary Control Syringes

Merit Syringes feature a clear polycarbonate barrel, solid plunger body and are available in a variety of sizes and configurations.

- Rotating and fixed Luers
- Available locking ring
- Optional smart tip
- Inject8 and Inject10ⁿ generate more pressure with less force



Medallion® Syringes

The Medallion Syringe features a clear polycarbonate barrel, colored plunger and fixed male Luer.

- 7 color options
- Optional sword handle
- Custom printing or frosted field available



Fluid Management Kits, Packs & Trays

OPTIMIZE EFFICIENCY, REDUCE COST

Merit Medical offers pre-assembled fluid management kits and complete procedure packs to improve efficiency and deliver cost-savings. Our skilled representatives leverage years of experience and Merit's broad portfolio of Cath Lab products to design kits that meet your requirements. While each design is unique, all kits and packs feature quality products, manufacturer-direct service, and our value-driven expertise in designing cost-effective Cath Lab solutions.

Choose from our comprehensive portfolio of Cath Lab products:

- Pen & Label Sets
- InQwire® Guidewires
- Angiography Needles
- BackStop® Disposal Systems
- ShortStop® Temporary Sharps Holders
- Prelude® Introducer Sheaths
- Futura® Safety Scalpels
- Pressure Infusor Bags
- Inflation Devices
- Torque Devices
- Hemostasis Valves
- Equipment & Table Covers
- Patient Drapes
- OR & Hand Towels
- Gowns & Gloves
- Bowls, Basins & Cups
- Dressings, Gauze & Sutures
- Surgical Instruments



Contact your Merit Representative for a complimentary Cath Lab utilization assessment.

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

Australia
+61 (0) 3 8373 4943

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Brazil
+55 11 4561 - 0062

Canada
1.800.364.4370

China
+86-10 8561 0788

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Hong Kong
+852-2207 6888

India
+91-80-41223376

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Singapore
+65 67536311

South Korea
+8210 2205 0707

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 22 518 02 30

(Français)
+41 22 518 02 52

(Italiano)
+41 22 518 00 35

UK
0800 973 115